

Gibco选择性抗生素 让选择变得不再困难



Gibco™多种高质量的选择试剂，为您提供满足您研究需要的独特解决方案，您可以轻松选择适合您实验需求的选择性抗生素。

- 快速筛选稳定的真核细胞系
- 双重筛选
- 经济有效的选择
- 筛选适用于多种生物体

选择性抗生素	最常见的选择用途	规格	货号
杀稻瘟菌素 (Blasticidin)	真核细胞和细菌	20 mL	A1113902
		10 x 1 mL	A1113903
		50 mg	R21001
Geneticin™ (G-418)	真核细胞	20 mL	10131035
		100 mL	10131027
		1 g	11811023
		5 g	11811031
		25 g	11811098
潮霉素B (Hygromycin B)	双重筛选和真核生物	20 mL	10687010
霉酚酸 (Mycophenolic Acid)	哺乳动物和细菌	500 mg	11814019
嘌呤霉素 (Puromycin)	真核细胞和细菌	20 mL	A1113802
		10 x 1 mL	A1113803
博莱霉素 (Zeocin)	哺乳动物、昆虫、酵母、细菌和植物	50 mL	R25005
		8 x 1.25 mL	R25001
羧苄青霉素二钠盐 (Carbenicillin, Disodium Salt)	土壤杆菌和大肠杆菌	5 g	10177012

选择指南: Gibco选择性抗生素及其对应的携带抗性标记的载体

选择性抗生素	细胞类型	选择性抗生素的浓度(g/ml)*	携带抗性标记的载体
Blasticidin S	HeLa	1-3	<ul style="list-style-type: none"> • pcDNA6 vectors (Constitutive mammalian expression, 推荐货号V22020) • BsdCassette™ vectors (Constructing customized Blasticidinresistant vectors, 推荐货号V51020) • pIB/V5-His-TOPO™ vectors (Stable insect expression, 推荐货号K89020) • pMIB/V5-His vectors (Secreted insect expression, 推荐货号V803001) • pCoBlast (Selection vector for DEST™, 推荐货号K512001, K513001) • pPIC6 vector (Constitutive expression in Pichia pastoris, 推荐货号V21020) • pPIC6α vector (Secreted expression in Pichia pastoris, 推荐货号V21520)
	NIH3T3	5-10	
	CHO	5-10	
	COS-1	3-10	
	293 HEK	5-10	
	S2 Drosophila	~5	
	S. cerevisiae	~25	
Zeocin	HeLa	~150	<ul style="list-style-type: none"> • pcDNA4 vectors (Constitutive mammalian expression, 推荐货号V86420) • T-REx™ System, pcDNA4/TO vectors (Inducible mammalian expression, 推荐货号V102020) • pSecTag2 vector (Secreted mammalian expression, 推荐货号V90020) • ZeoCassette™ vectors (Constructing Zeocin™-resistant vectors, 推荐货号V50120) • InsectSelect™ System, pIZ/V5-His vectors (Stable expression in insect cells, 推荐货号V800001) • pPICZ and pICZ vectors (Inducible expression in Pichia pastoris, 推荐货号V19520, V19020) • pGAPZ and pGAPZ vectors (Constitutive expression in Pichia pastoris, 推荐货号V20020, V20520)
	NIH3T3	~400	
	CHO	~250	
	COS-1	~400	
	293 HEK	200-400	
	Jurkat T cell	~200	
	Sf9 insect	~250	
	S2 Drosophila	~75	
	Pichia pastoris	~100	
S. cerevisiae	200-300		
Geneticin	HeLa	200-400	<ul style="list-style-type: none"> • pcDNA3.4™ vectors (Constitutive mammalian expression, 推荐货号A14697) • pShooter™ vectors (Intracellular protein targeting, 推荐货号V82320, V89120, V82220, V82120, V82020) • pDisplay™ vectors (Protein display, 推荐货号V66020) • pBlue-TOPO™, pGlow-TOPO™ vectors (Assessing promoter activity, 推荐货号K483101, K483001)
	NIH3T3	600-1,000	
	CHO	~400	
	293 HEK	600-800	
	Jurkat T cell	600-700	
Hygromycin B	HeLa	~550	<ul style="list-style-type: none"> • pcDNA5 vectors (Constitutive mammalian expression, 推荐货号V601020) • pSecTag2/Hygro vector (Secreted mammalian expression, 推荐货号V91020) • pREP4/pCEP4 vectors (Episomal mammalian expression, 推荐货号V00450, V04450) • pCoHygro (Selection vector for DEST™, 推荐货号K412001, K413001)
	CHO	~250	
	Jurkat T cell	~1,000	
	S2 Drosophila	200-300	

* 应使用杀伤曲线来确定您的细胞系所需的选择性抗生素的最佳浓度



赛默飞
官方微信



赛默飞
生命科学官方微信

免费服务电话: 800 820 8982/400 820 8982
 销售服务信箱: sales.china@thermofisher.com
 技术咨询信箱: LifeScience-CNTS@thermofisher.com

ThermoFisher
SCIENTIFIC