

life lab

PRODUCTS, INFORMATION, AND SCIENTAINMENT
ISSUE 18 | 2018

Inspiring discovery

Gene expression
solutions that give you
extra confidence

page 6

Scientists performing
inspiring science

page 14

Imagine the beauty
in simplified cell analysis

page 18

thermo
scientific

applied
biosystems

invitrogen

gibco

iontorrent

ThermoFisher
SCIENTIFIC

invitrogen

Flow cytometry

Let's get to the science



Stimulation
antibody



Conjugated
antibodies



Cell health



Flow RNA



Flow cytometry
instrument

Analyze single cells like never before

Within a single cell are insights waiting to be translated, answers waiting to be found, and discoveries waiting to be made. Get to them in new ways, in this new era of flow cytometry.

Harness the power of flow at
[thermofisher.com/flowcytometry](https://www.thermofisher.com/flowcytometry)

ThermoFisher
SCIENTIFIC

TABLE OF CONTENTS

- 4 Party inspiration
- 5 5 tips for successful CRISPR editing
- 6 Gene expression solutions that give you extra confidence
- 8 With high-quality resources, your research can change the world




- 17 Finding inspiration
- 18 Imagine the beauty in simplified cell analysis



- 20 The ingenuity behind impactful assays
- 23 Designed with the customer in mind

- 24 The fast and the inspired
- 26 Promising solutions for precision oncology research



PARTY INSPIRATION

MAKE YOUR HOLIDAY, NEW YEAR'S, OR BIRTHDAY PARTY MORE FUN WITH SCIENCE

Mixology

Gel spheres—Have you ever seen cocktails topped with floating gel spheres? These spheres are created by mixing a liquid with sodium alginate in a syringe. The liquid needs to have a pH range of 4–10 for the liquid to gel. Fruit juices are a popular choice and also add a pop of flavor to champagne or sparkling water.

Cocktail popsicles—For hot days, there is no better way to cool off than with a popsicle. Make it an adult-friendly treat by combining your favorite cocktail with a tasty fruit puree and freezing. The mixers from the cocktail and the fruit puree will blend and help freeze even higher-proofed alcohols.

Nitrogen-chilled cocktails—These are not for the faint of heart. Instantly chill cocktails using liquid nitrogen instead of ice for a great conversation starter at parties. Remember, liquid nitrogen should only be handled when wearing safety glasses, protective clothing, and gloves.



Appetizers

Trail mix test tubes—Fill test tubes (unused, of course) with trail mix or mixed nuts to create a fun snack for your guests.

Cheese molecules—Let partygoers make their own molecules from a plate of cheese cubes, olives, fruit, and veggies.



5 TIPS

Innovative CRISPR technology offers numerous possibilities. Our CRISPR trainer and technical support scientist, Yashashree Joshi, PhD, would like to share the tips she's picked up along the way.

FOR SUCCESSFUL CRISPR EDITING

1. Aim for 30–70% confluency.

Optimize cell confluency—seeding density, passage number, and the type of media can impact editing efficiency. Ideal starting confluency when using a Cas9 protein for lipid transfection is 30–70% and 70–90% for electroporation.

2. Test 2–3 guide RNAs.

Target the early exon within the gene to disrupt the reading frame. This will identify the gRNA with the highest editing efficiency.

[thermofisher.com/trueguide](https://www.thermofisher.com/trueguide)

3. Start with lipid-mediated delivery.

Know what works for your cells, and progress to other methods like electroporation as needed.

[thermofisher.com/crisprtransfection](https://www.thermofisher.com/crisprtransfection)

4. Optimize use of controls.

Positive and negative controls help determine gRNA amount and transfection conditions that give optimal efficiency with the highest cell viability.

[thermofisher.com/crisprcontrols](https://www.thermofisher.com/crisprcontrols)

5. Confirm gene-editing efficiency.

For robust future screening experiments, verify the control target and select the conditions for optimal editing efficiency.

[thermofisher.com/geneeditdetect](https://www.thermofisher.com/geneeditdetect)

[thermofisher.com/ngs](https://www.thermofisher.com/ngs)

[thermofisher.com/sangersequencing](https://www.thermofisher.com/sangersequencing)

Yashashree Joshi, PhD

Technical Applications Scientist
CRISPR Workshop Trainer
Thermo Fisher Scientific



What is your favorite thing about leading this course?

Helping people! Scientists come into this course with their own projects and I get to help them understand the entire workflow. We held our first course in Winter 2016 and I am still in touch with many of the attendees—they tell me about how well their experiments are working!


What is one thing you want everyone to know?

The coolest thing is we have products for the whole genome editing workflow—something none of the other suppliers can say.

Why are you passionate about science?

It is the basis of all advancement and technology that everyone is benefiting from. Basic science research helps unravel the fundamental processes, which can have huge applications. For instance, CRISPR was discovered as the bacterial immune system, but who knew it would have such vast implications in genome editing. It now holds the promise of correcting disease-causing mutations and more.

[thermofisher.com/crisprworkshop](https://www.thermofisher.com/crisprworkshop)



GENE EXPRESSION SOLUTIONS THAT GIVE YOU **EXTRA** CONFIDENCE

Tools that give you peace of mind from the start

TAQMAN ADVANCED miRNA ASSAYS—GET GREATER FLEXIBILITY WITH FORMATS INCLUDING TAQMAN PLATES, ARRAY CARDS, AND OPENARRAY PANELS

Applied Biosystems™ TaqMan® Advanced miRNA assays achieve high-sensitivity and high-specificity microRNA quantification using qPCR. This solution provides a productive workflow (using a universal RT step) that is optimized for the analysis of multiple miRNA targets from challenging low-level RNA research samples in serum; plasma; formalin-fixed, paraffin-embedded (FFPE) tissue; and blood.

thermofisher.com/advancedmiRNA



APPLIED BIOSYSTEMS™ CLARIOM™ ASSAYS FOR FAST TRANSCRIPTOME PROFILING—DON'T LET DATA ANALYSIS BE YOUR BOTTLENECK

Take control of your RNA biomarker discovery projects. Free, fast, intuitive, and powerful analysis software designed for biologists—allowing you to gain insight into the transcriptome in just three days.

thermofisher.com/fasttranscriptome



SEED YOUR INSPIRATION WITH GENETIC ANALYSIS WEBINARS

Learn how genetic analysis techniques, including sequencing, NGS, microarrays, and quantitative and digital PCR, are being used by researchers to advance knowledge and innovations in challenging disease and research areas.

thermofisher.com/gawebinars

GENERATE SUPERIOR SANGER SEQUENCING DATA

Choosing the right workflow from the start is the key to success

Proven through decades of results, Sanger sequencing is the gold-standard, simple technology that allows you to:


- Study diseases with clearly defined phenotypes
- Sequence 1–2 genes or up to 96 targets
- Sequence 1–96 samples at a time without barcoding
- Confirm next-generation sequencing (NGS) variants with up to 99.99% accuracy
- Get longer read lengths (up to 1,000 bp)

For many applications, Sanger sequencing remains the technology of choice due to the simple workflow, quick turnaround time, and cost-effectiveness when running a limited number of samples. It is also an important confirmatory method for NGS results. However, poor quality data can result from not using the right products or the right protocol. Poor-quality data increases the rerun rate, which adds to the overall cost and turnaround time. We have developed a proven protocol that provides a simple and fast Sanger sequencing workflow that can be completed in less than one workday,

from sample to answer. The protocol enables high-quality sequencing data using the Applied Biosystems™ BigDye™ Terminator v3.1 or BigDye™ Direct Cycle Sequencing Kits. Visit our page to download our latest protocols for the Applied Biosystems™ products that support every step of our recommended workflow.

Ready to move to better results?

thermofisher.com/sangerworkflow



Capillary electrophoresis workflow



**WITH HIGH-QUALITY RESOURCES,
YOUR RESEARCH CAN
CHANGE
THE WORLD**

Change begins in the lab,
and accuracy is essential to
yielding meaningful results

Equip your lab with:

- Access to innovative products developed for accurate performance
- Preprogrammed methods tailored to fit the nuances of everyday applications
- Equipment you can depend on for consistent quality of outcomes

MAKE THE MOST OUT OF EVERY DROP

In a world facing water scarcity, forward-thinking scientists are utilizing water quality technology to pioneer new applications.

Your lab can work toward improving lives when one of the most advanced water purification and analysis tools drives your research.

Make an impact with your work, whether you're:

- Developing solutions for clean drinking water
- Bolstering food production locally or globally
- Providing fundamental resources for accelerating biopharmaceutical advancements

Thermo Scientific™ high-performance measurement technologies provide the most accurate sample readings possible, whether you're checking pH, ion concentration (ISE), conductivity, dissolved oxygen, temperature, or parameters unique to your application.

Put the most versatile water quality products to work for you and your evolving needs.

thermofisher.com/waterquality



ESSENTIAL ACCURACY, MEANINGFUL RESULTS

Meet a comprehensive water quality portfolio

Precise beginnings lead to profound advancements. Our comprehensive water quality portfolio was developed by scientists, for scientists, with accuracy, flexibility and longevity in mind, so you can spend more time working and less time worrying.

We understand the importance of consistent outcomes. They are the difference between breakthroughs and redos. You need to be able to focus on moving your research forward without worrying about inconsistencies. With a legacy of proven expertise and a quality portfolio of pH, electrochemistry, and water purification products, we bring cohesive workflows and greater precision to your research so you have more bandwidth for what's important.

thermofisher.com/waterquality

CREATED SO THAT YOU CAN EXPAND YOUR POSSIBILITIES



SIMPLIFY DNA ELECTROPHORESIS WITH THE ONLY INTEGRATED GEL RUNNING AND IMAGING SYSTEM

The Invitrogen™ E-Gel™ Power Snap Electrophoresis System combines rapid, real-time nucleic acid analysis with high-resolution image capturing for superior convenience. The integrated design helps reduce workflow time and accelerate discovery.

- **Fast**—go from sample loading to image capture in as little as 15 minutes
- **Simple**—intuitive user interface with large touch screen and integrated operating system
- **Safe**—minimize handling of hazardous chemicals when used with Invitrogen™ E-Gel™ precast gel cassettes

Watch the video and request a demo at thermofisher.com/powersnap or thermofisher.com/egel





BE UNSTOPPABLE WITH FASTER, MORE ACCURATE ONE-STEP RT-PCR

The Invitrogen™ SuperScript™ IV One-Step RT-PCR System combines high-processivity Invitrogen™ SuperScript™ IV Reverse Transcriptase (RT) and high-fidelity Invitrogen™ Platinum™ SuperFi™ DNA Polymerase to provide unmatched one-step RT-PCR performance, even with challenging RNA samples.

- Two-phase hot-start activation mechanism for high specificity, increased yield, and room-temperature setup
- Superior performance with sensitivity down to 0.01 pg of RNA, target length up to 13.8 kb, and the fastest one-step RT-PCR protocol available
- Reliable target detection, even in RNA samples with suboptimal purity
- Fast and easy gDNA removal for superior accuracy and confidence in your results

Request a sample at
[thermofisher.com/ssiv-onestep](https://www.thermofisher.com/ssiv-onestep)



Your cloud-enabled Applied Biosystems™ ProFlex™ and SimpliAmp™ Thermal Cyclers can be securely accessed anytime, anywhere, with any mobile device or desktop computer, allowing you to:

- Design and share protocols
- Schedule time to use an instrument
- Start or stop a run
- Check run status

[thermofisher.com/proflex](https://www.thermofisher.com/proflex)
[thermofisher.com/simpliamp](https://www.thermofisher.com/simpliamp)

Powered by **Thermo Fisher Cloud** 

WHY ARE YOU PASSIONATE ABOUT SCIENCE?

I have been passionate about science ever since I was young. Exploring, learning, and understanding the world around me sparked my aspiration for studying science. I was especially fascinated by organism diversity and how their diversity is controlled through the regulation of genetic information during development. I am continuing my passion for science and experience in development and molecular biology to Thermo Fisher Scientific, where I am researching and developing molecular tools and solutions for customers.



Vicki Hurless, PhD, R&D Scientist, Molecular Biology
 Thermo Fisher Scientific

CRISPR SOLUTIONS DESIGNED TO ELEVATE YOUR RESEARCH

TRULY INSPIRATIONAL

We are passionate about helping you make your next scientific breakthrough, and our R&D team is constantly looking for ways to develop high-performing, reliable, and innovative products that will help save you precious time so you can accelerate your rate of discovery. See what researchers are saying about these innovations.

New, innovative solutions designed to take your research to the next level



Invitrogen™ LentiArray™
CRISPR Libraries



Invitrogen™ TrueCut™
Cas9 Protein v2

Invitrogen™ TrueGuide™
Synthetic gRNA

“If I hadn’t done this functional study, where, knocking out the *PLK4* gene, resulted in the inhibition of cell proliferation I wouldn’t have discovered that *PLK4* is essential for *AT/RT* growth. *PLK4* is only slightly elevated in these tumors. However, this gene is tightly regulated and slight increases in its expression result in an aggressive tumor phenotype. This overexpression can be targeted by inhibitors opening a new therapeutic prospective for children with *AT/RT*. Significantly, we also found it elevated in other embryonal tumors of the brain, what may have a larger impact in patient care ... the fact that we used the Invitrogen™ LentiArray™ Human Kinase CRISPR Library for the functional assay was key.”

– Simone T. Sredni, MD, PhD, Associate Professor of Pediatric Neurosurgery at Ann and Robert H. Lurie Children’s Hospital of Chicago / Northwestern University Feinberg School of Medicine
thermofisher.com/crisprlibraries



CURIOSITY WITH CRISPR

Olivier Humbert, PhD,
Staff Scientist in
Hans-Peter Kiem’s
laboratory, Fred
Hutchinson Cancer
Research Center



What inspired you to get into science and research?

What inspires you to get up every day?

Curiosity about how my surroundings and living things work; I love the excitement of carrying out experiments, the anticipation for results, and the fulfillment of knowing that I am contributing to the development of new medical treatments.

What drives your passion for science?

Knowing that I am working with cutting-edge technologies to develop therapies that will become available to patients in the next 5–10 years.

Can you provide an overview of your research?

Our goal is to treat hereditary diseases by correcting the underlying mutation in bone marrow stem cells so that these cells will be curative after they are transplanted back to the affected patient. The correction of mutations can be done by means of non-infectious viral vectors or more recently with nucleases that target

specific sites in DNA, also known as molecular scissors (such as CRISPR-Cas9). Our laboratory covers a wide range of research, spanning basic science all the way to clinical studies; examples of diseases we work on are blood disorders (sickle cell anemia), severe combined immunodeficiency, and Fanconi anemia (defects in DNA repair).

What are the current challenges you face?

One challenge is to correct mutations and engineer DNA in bone marrow stem cells with high efficiency and without affecting the identity of these cells so that they will differentiate into all blood cell types and reconstitute the hematopoietic system of a patient after transplantation. Another big challenge is in the scale-up, where we are looking to engineer the DNA of several hundred millions of cells obtained from patients.

What are the potential solutions?

One solution is to find a nuclease platform (such as CRISPR-Cas9) and nuclease delivery method (such as electroporation) that can treat large numbers of cells with minimum toxicity and high efficiency so that they will be therapeutic after they repopulate the patient. Another solution is to refine our definition of true, long-term hematopoietic stem cells, which will decrease the number of these cells that need to be treated.

Describe your results, if you have any.

We have so far done two transplant experiments in a preclinical animal model aimed at treating hemoglobinopathies (blood disorders). The results are very encouraging, but we want to follow up on our treatment for over a year to really understand how effective and safe our therapy is.

How are the advances in genome editing enabling you to achieve your research goals?

Current gene editing technologies allow us to target virtually any site in the genome in pretty much any cell type. In addition, the improved efficacy of TrueCut Cas9 Protein v2 allows us to reduce the amount used and minimize toxicity.

How are the tools and technologies today allowing you to lead to discovery?

Improved efficacy and safety of technology equals more potential to bring an effective therapy to the patient. These reagents are also available in large-scale [quantities] and in cGMP grade for easy translation to the clinic.

What thought do you want to leave us with?

I am lucky to be working in such an exciting field with lots of potential to cure a variety of diseases. Nevertheless, researchers need to be vigilant to take all necessary steps to not move this approach too quickly to the patient to avoid any major setbacks.

SCIENTISTS PERFORMING INSPIRING SCIENCE

Image courtesy of
Ekaterina Turlova,
University of Toronto

Whether using cells as a model of disease, leveraging cells to make protein, or using cells as therapy, cell biologists inspire us. Scientists like Samantha Yammine and Caitlin Vander Weele are going beyond the bench for public outreach and using science as art.

Samantha Yammine

**PhD candidate,
Neuroscience and
Stem Cell Biology,
University of Toronto**



What research do you do and do you have any challenges?

I am a PhD candidate in Dr. Derek van der Kooy's neurobiology research lab at the University of Toronto researching stem cell hierarchies in the developing and adult mammalian brain. Our lab studies a variety of different stem and precursor cell populations, including those of

the retina, corneal limbus, pancreas, and neural crest, but my thesis has revolved around a very rare population of neural stem cells. And by rare, I mean really rare—after I microdissect the thin stem cell niche surrounding the lateral ventricles in the mouse brain, my stem cells of interest comprise only 0.1% of the total cells dissected at most, and at some ages, they are only 0.01% of the population. But they are phenotypically really interesting cells—we recently published in *Stem Cells* that they act as a reserve pool for the more prominent neural stem cell population, so it's been worth troubleshooting new, sensitive techniques to try to learn more about them.

Have you found a solution? Since our cells are so rare in bulk samples, we knew we'd have to continue improving our purification methods and experiment with new assays that are sensitive enough to detect signals from single cells. Fortunately, single-cell analyses have become more and more popular and I've found a lot of fantastic collaborators in Toronto to help me implement these techniques. After many workshops, literature searches, and planning meetings, we are currently having success studying embryonic precursors with several single-cell transcriptomic platforms thanks to help from technical experts from several other local labs.

What are your next steps? I am excited to combine these new transcriptomics data with previous functional data I've collected on these rare neural stem cells and put forth some new additions to current neural stem and progenitor cell hierarchies. By better understanding the lineage relationships between the earliest cells of the mammalian brain, we will better appreciate how the diversity of the cells of the brain are created during development and homeostasis. Given that the brain has over 170 billion cells and that all of these cells come together to give us the ability to think and do, I find its creation to be one of the most fascinating biological concepts.

What do you want people to know about you? What motivates you to share your story? I created the Instagram account @science.sam so that I could bring more transparency to the process of science research. My goal was to show that there are so many fascinating aspects of science, including the basic biology behind experiments we do in the lab everyday. I share daily updates of research life and new science news through a personal lens to challenge my audience to change their perceptions of scientists and science for the better. I strive for outreach that is interesting and inclusive to everyone, so I use enough analogies and metaphors for the casual science enthusiast, but just enough details to keep fellow scientists interested, too.

A passionate advocate for equity, diversity, and inclusivity, Sam also posts inspiring interviews with people from underrepresented groups in the STEM fields to invite more people into the exciting world of science research. Thank you for inspiring us, Sam!

Instagram: science.sam
Twitter: SamanthaZY
www.samanthayamine.com

Caitlin Vander Weele

**PhD candidate,
Neuroscience,
Tye Lab,
Massachusetts
Institute of Technology**



What research do you do and do you have any challenges?

What happens to all the beautiful and interesting images generated from failed science experiments? Typically nothing. Scientists use them to help generate better tools and hypotheses; however, because they will never be published, they get stored away on hard drives for no one to see. I started posting some of the pretty images I gathered from my failed experiments (I'm a PhD student studying the neural circuitry underlying motivated behaviors) on Twitter and other social media outlets. Microscopy images were some of my most viewed and shared media and it struck me how effective they were at communicating science content. I wanted to create a home for these images to be seen and appreciated. On a whim of procrastination, Interstellate was born in May 2016.

Have you found a solution? To start, I made a Twitter account for the project (@interstellate_) and sent over 100 cold emails to neuroscience friends, colleagues, and principal investigators asking for image donations for the project. To be quite honest, I was surprised how enthusiastic and supportive the community was and within six months, I collected over 100 images from ~80 scientists in nine different countries! I tweeted images I collected and started assembling them into an 86-page full gloss magazine. Each page of the magazine features a stunning image generated by scientific research and briefly explains a


neuroscience-related concept. For example, the first couple of pages explain the different types of brain cells, how neurons communicate with one another, and how we study them in the laboratory. In October 2016, the digital copy of Interstellate Volume 1 went live (http://pub.lucidpress.com/Interstellate_Volume1/) and with help from our generous corporate sponsors, over 1,000 copies have been printed and distributed for free! I think Interstellate provides a platform to celebrate important but often overlooked steps of research. Interstellate's goals are science celebration and neuroscience awareness through art.

What are your next steps? Interstellate is still a really new initiative, so it will be exciting to see what it evolves into! Right now, I'm working on a second volume (and trying to finish my PhD!), which will debut in November 2017, and expanding the social media and outreach presence (you can follow us on Tumblr and Instagram @interstellate_). The whole project is one massive collaborative effort so the future of Interstellate really depends on the network of people who support it. I would love to see Interstellate used as an outreach tool for the public and to recruit the next generation of brain explorers.

A brain explorer that collects images of pretty neurons along the way, Caitlin is showing us the beauty in science. Thank you, Caitlin!

Instagram: interstellate_
Twitter: interstellate_
www.caitlinvanderweele.com
www.interstellate.com

REVOLUTIONIZING CANCER THERAPY



A HUGE STEP FORWARD FOR IMMUNOTHERAPY: FIRST EVER FDA-APPROVED CAR-T CELL THERAPY

August 30, 2017 marked the start of an exciting new era in personalized medicine. On that day, the US Food and Drug Administration (FDA) approved a therapy that genetically alters a patient's own cells to fight leukemia. Novartis™ KYMRIAH™ therapy, also known as tisagenlecleucel, is a treatment for pediatric acute lymphoblastic leukemia (ALL), and is now the first FDA-approved chimeric antigen receptor T (CAR-T) cell therapy to be commercially available. This historic approval brings awareness to emerging immunotherapies

and their potential to transform cancer treatment, and even cure specific cancers.

Emily Whitehead, 12, was the first child ever to receive this type of "living drug." She was one of 63 patients whose cases were used as evidence showing that the treatment had an 83% remission rate at the 90-day mark.

Read the full blog post at
thermofisher.com/ctsblog

Read the press release and find out how
Gibco™ CTS™ products played a role at
thermofisher.com/ctspr

ENABLING CAR-T RESEARCH FROM THE BENCH TO THE CLINIC

Gibco™ Cell Therapy Systems (CTS™) media and reagents are designed to help you translate your cell therapy to clinical applications. They are manufactured in accordance with cGMP for medical devices (21 CFR Part 820), and backed by extensive safety testing and traceability documentation to facilitate regulatory approval, so you can transition your cell therapy to the clinic with confidence.

thermofisher.com/ctsimmunotherapy



For several decades, scientists have been studying ways to use the immune system to help fight cancer. Exciting new discoveries have shown that immuno-oncology (I-O) research can provide potential anticancer therapies to patients who previously had very few treatment options available to them. I-O therapies represent a breakthrough in cancer treatment and have the potential to revolutionize the way we treat many forms of cancer.



Harness the power of lentiviral vectors for CAR-T applications—they have a safer integration site profile than the commonly used gamma retroviral vectors to deliver CARs into T cells. Discover the first optimized, high-titer lentiviral production system for suspension culture at any scale.

thermofisher.com/lentiviral

FINDING INSPIRATION

Words to spark the imagination

What's your favorite inspiring science quote?

Learn the art of science and the science of art and you'll find everything is connected.
– Leonardo da Vinci



Science and everyday life cannot and should not be separated. – Rosalind Franklin



What's the best lab or science career advice you've ever received?

Career advice: Never stop asking questions.
Lab advice: Never forget the controls.



The only dumb question is the one you didn't ask.



To achieve greatness in the lab you must have an unlimited imagination to question everything.



WHAT BOOKS DO YOU READ TO INSPIRE YOU?

THINKING, FAST AND SLOW

by Daniel Kahneman

WHEN BREATH BECOMES AIR

by Paul Kalanithi

BIG MAGIC

by Elizabeth Gilbert

WHAT IF?: SERIOUS SCIENTIFIC ANSWERS TO ABSURD HYPOTHETICAL QUESTIONS

by Randall Munroe

LAB GIRL

by Hope Jahren

IMAGINE THE BEAUTY IN SIMPLIFIED CELL ANALYSIS

SEE IT TO BELIEVE IT

Utilize high-content analysis (HCA) to address the diverse set of experiments and cell types that are emerging in cell-based assays

We have developed the Thermo Scientific™ CellInsight™ CX7 LZR HCA Platform that uses laser technology and its On-Stage Incubator to give you both the sensitivity and speed you need for imaging your live cells, along with the quantitative tools you need to measure the phenotypic changes in your sample in real time. New live-cell imaging capabilities are now available across all of our HCA instruments to enable you to better understand the dynamics of your cells in context.

thermofisher.com/hca



COUNT ON EASE


What if cell counting could be made more enjoyable and better than manual cell counting? We've added new features to the Invitrogen™ Countess™ II and Countess™ II FL Automated Cell Counters to do just that. These cell counters now feature:

- **Dilution calculator**—displays the volumes of cell sample and buffer needed
- **Report generator**—generates a report of results, graph(s), and images in a PDF file
- **Simplified workflow**—save protocols easily via the adjustments screen
- **Modified profile indicator**—shows that the gating settings have been changed

Learn more about these features or download the software to upgrade your existing Countess II or Countess II FL Automated Cell Counter at

thermofisher.com/countessupdate

thermofisher.com/countessii



5 STEPS TO PUBLICATION-QUALITY IMAGES

Want to get the most out of your fixed-cell imaging experiments? From preparing your cells for labeling to imaging, this brochure will help you obtain publication-quality images in just 5 steps. Download the brochure at thermofisher.com/5steps

Check out our 5-steps promotion on the back cover!

COUNT ON EFFICIENCY

Learn how you can use the Countess II FL Automated Cell Counter with your CRISPR experiments to optimize measurements of transduction efficiency. Download the application note at thermofisher.com/crisprcountess



THE INGENUITY BEHIND IMPACTFUL ASSAYS



Gerald BW Wertheim, MD, PhD
Assistant Professor
Children's Hospital of Philadelphia

DEVELOPING AN ASSAY THAT MEASURES DNA METHYLATION, AN INDEPENDENT PREDICTOR OF LEUKEMIA

Background: Acute myelogenous leukemia (AML) is generally a disease of older people, although it is also the second most common blood cancer in children. In AML and other leukemias, alterations in DNA methylation are well recognized, and the aggressiveness of AML tumors can be determined by studying methylation. Mutations in the genes that regulate DNA methylation have been shown to cause different forms of leukemia. Studies have also demonstrated that the prognosis of patients with AML can be predicted by this altered DNA methylation pattern with

as few as 17 loci. However, tests that directly measure multiple-locus DNA methylation are typically expensive and technically challenging, making them difficult to perform in a routine setting. Gerald Wertheim and his colleagues have developed a novel approach to simultaneously analyze DNA methylation patterns at the 17 target loci by integrating multiplexed branched DNA technology from Thermo Fisher Scientific and fluorescent microsphere technology from Luminex. The method uses techniques that are inexpensive and can be easily performed in a routine

setting. The technique is called expedited microsphere HpaII tiny fragment enrichment by ligation-mediated PCR, or xMELP. It accurately reflects the methylation levels at each analyzed locus and enables segregation of individuals with acute myeloid leukemia into prognostic subgroups. Dr. Wertheim and his colleagues have published papers in *The Journal of Molecular Diagnostics* and the journal *Clinical Chemistry* [1,2].

Tell us about your studies of DNA methylation.

My primary interest has been on developing methods to better diagnose leukemia and lymphoma and to assess their prognosis in patients. Ultimately, we hope that clinicians can use our findings to help guide their therapy choices. We have known for some time that mutations are present in leukemia tumor cells, but we have not had the information to risk-stratify patients optimally. So I have been focusing on DNA methylation, which is thought to control transcription in certain numbers of genes. There are a number of pieces of evidence demonstrating that differential DNA methylation is important for AML prognosis. Many researchers have shown that genes that regulate DNA

methylation are mutated in different types of leukemias, and some have directly shown that looking at DNA methylation alone can predict the outcome in patients with AML. Our goals are to find a set of genetic loci that vary in their DNA methylation patterns among patients, which can predict who will do well and who will do poorly with therapy, and to then develop a clinically relevant assay to detect levels of methylation at these loci.

How did you develop your analysis method?

With an ultimate goal of taking our assay to the clinic, we wanted to find a technique that was relatively inexpensive and could be multiplexed. We wanted to be able to look at multiple methylated regions to improve robustness, and we needed an assay that was highly reproducible. Because we use microsphere technology in the clinical lab, I initially tried a microsphere assay without branched DNA, and it did not have the appropriate analytical sensitivity. One of my collaborators was aware of the Invitrogen™ QuantiGene™ Plex Assay from Thermo Fisher Scientific and suggested that we look at it for our work. The QuantiGene Plex Assay uses branched DNA (bDNA)

to amplify the signal rather than amplifying the target, and it sufficiently improved our sensitivity. Having the branched DNA amplification signal has been critical. We found that this QuantiGene Plex Assay was remarkably quantitative, with the results obtained by the QuantiGene Plex Assay and qPCR being virtually identical. Thus, this looked like the robust, multiplex assay program that could ultimately be used in high-throughput clinical lab settings. Evaluation of methylation is done through the measurement of microspheres and does not require custom-made, solid-phase oligonucleotide microarrays or high-throughput sequencing. The examination of methylation levels is performed by flow cytometric analysis of fluorescent microspheres, alleviating the need for high-throughput sequencing. This assay format can be easily expanded for the evaluation of 80 loci without an increase in reactions; the simultaneous assessment offers advantages, including a streamlined workflow and relatively fast turnaround time. The technical staff from Thermo Fisher taught us how to set up the instruments and run the assay. They were readily available by phone to help when I had questions. We were also able to have Thermo Fisher run some of our samples so we could compare runs and rapidly troubleshoot problems. I think that going forward, the diagnostic field is going to have to look at some of these other aspects, such as RNA expression analysis or epigenetic events. Hopefully, our assay will give clinicians a more accurate indication of prognosis and will help guide their decisions.

HIGH-SENSITIVITY IMMUNOASSAY INNOVATION

Many tools have been developed to advance our understanding of cellular mechanisms. However, a simple, scalable, and affordable method for low-level quantitation in limited-serum samples has seemed out of reach. Our latest innovation—Invitrogen™ ProQuantum™ immunoassays—is a set of ready-to-use kits for quantifying low-abundance proteins in small sample volumes. The assays utilize matched antibody pairs for specificity and Applied Biosystems™ TaqMan® qPCR technology for amplification and detection. Cloud-enabled ProQuantum™ software is available for the quantitative analysis of protein concentrations interpolated from a standard curve.

thermofisher.com/proquantum

Powered by

Thermo Fisher Cloud



References

1. Wertheim GB et al. (2015) Validation of DNA methylation to predict outcome in acute myeloid leukemia by use of xMELP. *Clinical Chemistry* 61(1):249–258.
2. Wertheim GB et al. (2014) Microsphere-based multiplex analysis of DNA methylation in acute myeloid leukemia. *Journal of Molecular Diagnostics* 16(2):207–215.



The SeqStudio Genetic Analyzer is a compact, white and blue laboratory instrument. It features a large, color touchscreen display on the front panel, which shows a user interface with several green circular icons and a digital readout. The machine is positioned on a reflective surface, and a large green circular graphic highlights the device.

-  All-in-one cartridge with four-month stability
-  Sequencing and fragment analysis on the same plate
-  Fast: runs in as little as 30 minutes
-  Secondary analysis software with system purchase
-  Access your data anywhere with Thermo Fisher Cloud*
-  SmartStart™ orientation to get you up and running fast
-  One-year instrument warranty included

* Internet access required.

SeqStudio Genetic Analyzer— say hello to your new best lab partner

Same workflow, same trusted technology—now with an innovative all-in-one cartridge that reduces setup time from hours to minutes. The new Applied Biosystems™ SeqStudio™ Genetic Analyzer delivers gold-standard Sanger sequencing and fragment analysis with just a simple click. Get the same data quality, service, and support you've come to expect from Applied Biosystems™ technology, with a modernized experience at an affordable price.

System comparison



**SeqStudio
Genetic Analyzer**



**310
Genetic Analyzer**



**3130
Genetic Analyzer**

Ease of use	Integrated all-in-one cartridge installs in minutes with just one click.	Clean block and syringe and install on instrument. Connect capillary to gel block. Tape in place and fill with polymer. Dilute buffers and place on instrument.	Unscrew and install polymer bottle. Dilute running buffer and fill reservoirs. Screw array in detection window and snap clips and array in place. Run array install wizard.
Calibration	Auto calibration	Autosampler must be calibrated when capillary is replaced. Syringe calibration required when replaced.	Spatial calibration required when you install or replace array. Spectral calibration required when changing array length or polymer type.
Number of capillaries	4	1	4
Interactive touch screen with easy-to-use interface	Yes	No	No
Consumable tracking (RFID)	Yes	No	No
Capillary/array length	28 cm (in click-in cartridge)	47 cm, 61 cm	22 cm, 36 cm, 50 cm, 80 cm
Polymer type	POP-1 universal polymer	POP CAP, POP-4, POP-6	POP CAP, POP-4, POP-6, POP-7
Instrument size	49.5 x 64.8 x 44.2 cm (W x D x H)	61.0 x 55.9 x 86.4 cm (W x D x H)	149 x 55 x 81 cm (W x D x H)
Minimum sequencing run time	30 min	38 min	35 min
Minimum analysis run time	25 min	30 min	20 min
Sample capacity	96-well plate and 8-tube strips	1–96 sample tubes	2 sample plates (96- or 384-well)
Number of dyes	6	5	5
Applications	Sequencing and fragment analysis on same run	Sequencing, fragment analysis	Sequencing, fragment analysis
Sequencing read length	Up to 800 bp	Up to 600 bp	Up to 950 bp
Maximum sequencing throughput (base pair reads/day)	67K	15K	78K
Maximum fragment throughput (samples/day)	230	48	288
Connectivity	Wi-Fi and RJ-45 Ethernet ports	No	No
Remote monitoring and data storage	Yes	No	No
Configuration	Stand-alone, optional desktop or laptop computer	Computer required	Computer required
SmartStart orientation	Yes	No	No

Trade in your existing Applied Biosystems™ genetic analyzer and simplify your workflow.

For full trade-in details, contact your sales representative or go to thermofisher.com/seqstudiotradein



THE

FAST

AND THE INSPIRED

HOW KIERAN MULRONEY SPEEDS UP BACTERIAL PROFILE TESTING

Drug-resistant bacteria are spreading globally and previously treatable infections are inflicting death. With the need to stem inappropriate antibiotic use, Kieran Mulroney, PhD candidate at the Harry Perkins Institute of Medical Research, is developing a faster and more reliable diagnostic test to profile drug susceptibility in bacteria. Learn how the Invitrogen™ Attune™ NxT Flow Cytometer enables him to rapidly and accurately profile the effect of drugs on individual bacterial cells directly from clinical samples.

2018 ATTUNE NXT FLOW CYTOMETER: WHAT YOU NEED TO KNOW

The Attune NxT cytometer is designed for faster experimental run times with acoustic focusing technology. Get accurate analysis with minimal clogging and sample preparation for a broad range of sample types.

- **Efficient**—10x faster than traditional flow cytometric analyzers
- **Flexible**—up to 14 colors; field upgradeable
- **Evolutionary**—run difficult samples with a system that is less sensitive to clogging

Find out more at thermofisher.com/attune



Kieran Mulrone, BSc



Biomedical Science Honors, PhD candidate at the Harry Perkins Institute of Medical Research, Faculty of Health and Medical Sciences, The University of Western Australia

Describe the importance of your work.

My research focuses on rapid diagnostics in severe bacterial infections, specifically looking at the antibiotic resistance profile. Using the correct antibiotic is essential for a patient to survive.

Tells us how the Attune NxT cytometer helps your research.


We are looking for the effects of the antibiotic drug on individual cells in a very rapid manner. Traditional antibacterial sensitivity testing can take 1–3 days

depending on how difficult the bacteria are to grow, isolate, and then test for sensitivity. The Attune NxT cytometer allows us to tell whether or not the drug is effective within 1 hour, and gives inhibitory concentrations within 3 hours.

What do you hope to achieve with the Attune NxT cytometer?

We aim to characterize and test a full suite of drugs and bacterial combinations in order to move into the clinical space. We want to create a portable diagnostic and bring it to the areas of the world that need it the most.

Learn more by watching Kieran's interview at <http://bit.ly/2vEKVtn>



PROMISING SOLUTIONS

FOR PRECISION ONCOLOGY RESEARCH

Moving toward the development of a universal companion diagnostic in non-small cell lung cancer (NSCLC)

Accelerate results, from sample to report, and guide oncology treatment with the first *in vitro* diagnostic (IVD) next-generation sequencing (NGS)-based test. The Oncomine™ Dx Target Test* is the first test of its kind for detection of genetic mutations in NSCLC.

* For *In Vitro* Diagnostic Use.

CAN A SINGLE TEST HELP GUIDE TREATMENT FOR NSCLC?

A single CDx report in 4 days

The Oncomine Dx Target Test* is an IVD NGS-based test covering 23 genes and 3 CDx markers to help guide NSCLC treatment.

thermofisher.com/oncomine-dxtarget



ENABLE MORE

The expanded portfolio of Ion Torrent™ Cell-Free Research Assays now covers cfDNA and cfRNA, SNVs, indels, CNVs, fusions and/or expanded gene coverage, enabling multibiomarker analysis for more actionable insights.

thermofisher.com/cfdna-assays

HARNESS THE POWER OF THE IMMUNE REPERTOIRE

With our new Ion AmpliSeq™ Immune Repertoire Assay Plus, TCR β , analyze the immune repertoire to capture the diversity of T-cell receptor rearrangements with the first long read assay on the Ion Torrent™ NGS platform:

- Discover immune-response biomarkers
- Investigate markers for immune-mediated adverse events
- Optimize the function and manufacture of potentially therapeutic T cells

thermofisher.com/immune-repertoire

LOW FREQUENCY IS NO LONGER A PROBLEM

The Applied Biosystems™ TaqMan® Digital PCR Liquid Biopsy Assays, combined with OncoPrint™ cfDNA solutions offer a comprehensive set of tools for liquid biopsy applications.

thermofisher.com/dpcr-raremutation

NGS CHANGES COURSE FROM TRADITIONAL GENOMICS TO WATER TREATMENT

While NGS is commonly used in hospitals for cancer and genetics research, and in labs focused on forensic analysis, a group of researchers in Australia is using a simple water sample to determine which organisms, including vertebrates, native fish, and bacteria, have been in contact with that water sample. The end result is an improved process and improved outcomes for customers in the form of cleaner water.

The innovators at the Australian Water Quality Centre (AWQC) are the first in the Australian water industry to apply NGS for more efficient and targeted water treatment.

“It’s very useful to know the diversity and abundance of aquatic life in water sources, as they can be indicators of a healthy water system and the Ion Chef™ and Ion S5™ technology are helping us do just this,” said Karen Simpson, senior manager of laboratory services at the AWQC.

Conventional methods for detecting fecal contamination require days for culturing, and are limited because these methods target specific bacteria; parasites can be missed even after hours under microscopy. Use of NGS technology cuts down on the time to result for these scientists, enabling them to implement a rapid, targeted response to contaminants, and notify the public when applicable.

“Knowing more about what’s in water sources enables better informed decisions on how to treat the water before it’s supplied to customers,” said Ian Hunter, minister for water and the River Murray.

The AWQC has been able to verify the data it is generating against curated databases and is developing its own unique curated database to match new samples against, to ensure accuracy and consistency of testing. By enhancing traditional testing with NGS DNA technology, one intelligent sample will provide more information, enabling better water quality-management directions and risk assessments.

INTRODUCING THE ION GENESTUDIO S5 SERIES



A LINE OF HIGHLY VERSATILE NEXT GENERATION SEQUENCERS

The Ion GeneStudio™ S5 Series, a new line of benchtop next-generation sequencing (NGS) instruments, provides unmatched flexibility and scalability enabled by five Ion S5 chips, including the new Ion 550 chip, to facilitate wide-ranging experiments on a single platform. The new series also provides low cost per sample for small and large projects across multiple research applications. The Ion GeneStudio S5 Series delivers the highest level of dynamic range and read lengths of up to 600 base pairs to drive more powerful and cost-effective experiments in multiple areas, including cancer research, inherited disease, and microbial and infectious disease applications. With the flexible chip format, the Ion GeneStudio S5 platform can accommodate two sequencing runs per day with data analysis, delivering between 2-260 million reads and a total output between 0.5 GB to 50 GB.

Flexibility on the Ion GeneStudio S5 Series is achieved simply by switching between Ion S5 chips, not another sequencer, to match the appropriate throughput requirement

for each project. The interchangeable chip format enables users to cost-effectively process samples in increments, eliminating the need to delay experiments by waiting to batch samples. The Ion GeneStudio S5 Series includes the Ion GeneStudio S5, Ion GeneStudio S5 Plus and the Ion GeneStudio S5 Prime systems. Formatted for ease of upgradability, each delivers increasing levels of speed, throughput and scalability and is designed to easily integrate with Ion AmpliSeq technology and the Ion Chef System for automated library preparation and amplification, as well as Ion Torrent downstream bioinformatics and reporting tools for a seamless sequencing workflow.

When combined with the expanded Oncomine portfolio of assays for liquid biopsy and immuno-oncology announced today, the new series of instruments offers the industry's most comprehensive, end-to-end oncology solution based on a single platform for applications ranging from research to preclinical.

Learn more at
thermofisher.com/genestudio

“In my institute, many customers outsource for exome and whole genome studies. They come to me for speed and flexibility for different gene panels tailored to disease areas of interest,” said Morten Dunoe, laboratory director, Molecular Genetic Laboratory at the Juliane Marie Centre, Rigshospitalet, Copenhagen University Hospital. “This is the most elegant workflow on the market today. I can now offer a range of applications and deliver the data in the fastest turnaround time at costs that are competitive.”

Ion GeneStudio S5 Series systems and Ion 550 Chip are for research use only. Not intended for diagnostic procedures.


SCIENCE ANYWHERE

Science is advancing at a rapid pace. The way you work—how you collect, manage, analyze, share, and protect your data—has to keep up. Thermo Fisher Cloud helps take your lab to the next level with a wide variety of secure cloud-based tools.

Did you know? Every new Thermo Fisher Cloud member gets 10 GB of free storage. Based on the instruments you use, your free storage amount can be adjusted up to 1 TB.

Sign up today at thermofisher.com/cloud

Powered by **Thermo Fisher Cloud** 



Mobile lab apps put digital science tools at your fingertips wherever you are. Explore our lab management apps at thermofisher.com/apps



EQUIPPED AND READY

Having a Supply Center on-site is like having our entire warehouse in your lab. Customize your on-site product mix from our portfolio of trusted brands to best suit your research needs, and we provide the storage solutions. Our specialists work with you to review your current product usage and research to propose a product-stocking plan to meet those needs, helping improve your efficiency today and in the future.

Read what some of our customers have to say:

“The Supply Center lets us order the products and reagents we need in a quick and efficient manner. In stock items can get here within minutes following our order and anything not stocked in our facilities is usually processed and shipped out overnight/fairly quickly. It helps us save time and keeps the research happening.”

– Holly Tran, Brown University

“Throughout our busy days, we often do not have the time to wait for or to order new supplies. Having the Supply Center stocked downstairs saves me with easy access to critical supplies on hectic days.”

– Scott Mansfield, Regeneron Pharmaceuticals

“I like placing orders via the on-site Supply Center in our Nelson Building because I do not need to wait for products. I just grab them down the hall with no charge for shipping and handling. It saved us time and money, so that we can push our projects quickly!”

– Weihuan Cao, Rutgers University

“I can find the quality products I need right here (in-house). Lots of them stocked so I don’t have to wait to get the supplies I need for my experiments. Products are easy to find, in-house store well organized and run by a super guy ... he should wear a cape!”

– Kay Nicholson, Medical College of Wisconsin

Learn more about on-site stocking solutions for your lab at thermofisher.com/supplycenter

DESIGNED WITH THE CUSTOMER IN MIND

Modern western blot imaging is making a real difference with our customers



Rebecca Sinnott-Devaux,
Postdoctoral Associate,
University at Albany, SUNY, New York

What are you researching and how has the iBright system helped your lab?

In Dr. Herschkowitz's lab at the Cancer Research Center at the University at Albany, I am focused on the study of long noncoding RNAs and how they function in breast cancer progression. Over the past 30 years, advances in imaging technology and increased emphasis on early detection have led to a dramatic increase in the diagnosis of ductal carcinoma *in situ* (DCIS), a precancerous growth in the breast. Although only 40% of diagnosed DCIS cases are predicted to progress into breast cancer, we currently cannot distinguish between those DCIS cases that would remain harmless, and those that require immediate attention. Therefore, a DCIS diagnosis results in a cancer diagnosis, surgery, potentially radiation therapy, and potentially hormone therapy as well. This has created a current state of overdiagnosis and overtreatment, impacting tens of thousands of women every year. The only way to truly

combat overtreatment is to gain a better understanding of the biology underlying breast cancer progression to better inform patient treatment, or nontreatment, options.

Long noncoding RNAs (lncRNAs) are relatively new players in both normal biology and in impacting disease. Our studies focus on the expression of lncRNAs that are misregulated in breast cancer progression. Recently, we purchased the Invitrogen™ iBright™ CL1000 Imaging System from Thermo Fisher Scientific. This instrument was critical to our lab environment as it not only replaced older, less efficient technology, but also combined white light, UV, and chemiluminescence imaging in one unit. For my studies, this instrument has been wonderful for developing southern and western blots to look at RNA:protein interactions, but in reality the entire department uses this instrument, taking advantage of its versatility. It has been a great addition to our platform.

AN ELEGANT, USER-CENTERED SOLUTION

The iBright CL1000 and FL1000 Imaging Systems are high-performance instruments that enhance the western blot and gel imaging experience through advanced automated features and a streamlined, easy-to-use interface. Both models feature chemiluminescent blot, nucleic acid gel, and protein gel imaging modes. The FL1000 model also offers fluorescent blot imaging capability. Complete with a powerful 9.1 megapixel camera, proprietary Invitrogen™ Smart Exposure™ acquisition technology, and Thermo Fisher Cloud connectivity, the iBright systems enable customers to capture and analyze western blots faster and easier than before.

“We continue on the path of developing innovative products that streamline and modernize western blotting,” said Paul Haney, senior product manager at Thermo Fisher Scientific. “With the iBright Imaging Systems, we aim to improve the convenience, speed, and quality of western blot experimentation and circumvent the limitations associated with the capture of western blot images using x-ray film. Our iBright FL1000 model is designed to work seamlessly with our Thermo Scientific™ SuperSignal™ chemiluminescent substrates or our Invitrogen™ Alexa Fluor™ Plus reagents for maximum multiplexing abilities. To put it simply, the iBright FL1000 model is the most advanced western blot imaging system available in both design and capability.”

To learn more, or to request a live demo to see for yourself how the iBright Imaging Systems can streamline your research, go to thermofisher.com/ibright





#IAmAScientist

If you're reading this, a scientist somewhere at some moment of time has impacted your life. Scientists working hard toward advancing discovery have a shared passion for the work that they do, but who exactly are these scientists? We've created a space for scientists, like you, to share their stories and explore personal journeys of other scientists around the world to learn more about the people behind the science.

SHARING A MISSION

"My customers are in the San Francisco Bay Area and the Pacific Northwest. They constantly amaze and inspire me with their dedication to their research. We are connected by our common love for science and the possibilities that we can achieve because of it. Through science, we have the technology and the power to really make a difference. I care so much about what happens with my customers' experiments because we are on the journey together.



However, I know that it is also a constant balance of business and science, and the business of science. It is my responsibility to ensure that I contribute to the success of our business. Our technologies and our business have to be sustainable to allow for further evolution and advancement to ultimately enable our customers to succeed. Twenty years ago, I chose this path ... and 20 years ago, I chose to keep going no matter how exhausted I am. I am passionate. I am dedicated. #IAmAScientist."

– Mary Ann Santos, Senior Technical Specialist for Synthetic Biology, Thermo Fisher Scientific

Why did you choose to pursue science?

"I always had an interest in science, but the enthusiasm of my instructors and undergraduate research adviser really solidified my path. Ultimately I hope to be able to give younger students the same opportunities and investment I was given." – Krista Armbruster

If you could change society's perception about science, what would you say?

"The scientific community is working to advance science! Scientists do not go to work every day to prove conspiracy theories, harm the general public, or destroy the world. Most things that scientists do are overwhelmingly helpful to the general public. Look around at everything you interact with in your daily life and there is a scientist behind it! That food you ate? Tested for safety by a scientist. That car you drive? The fuel was tested to be the least harmful to the environment by a scientist. That toy that your child is playing with? The material was developed to be nontoxic and safe to play with by a scientist! We are good people and want to see the world evolve to become a better, safer place!" – Heather Hanson

Do you have any advice for the next generation of scientists?

"Quit focusing on the negatives. If you want to do it, start the journey. A journey can take as long as it needs to. So what if yours is longer than someone else's. It's the destination that counts." – Lisa Ramey

Read featured stories and submit your own at thermofisher.com/tellusyourstory



Connect your lab

Securely access your lab anytime, anywhere

It's your life's work — we'll help you protect it

Simplify your lab operations and data management with Thermo Fisher Cloud to help boost your lab's productivity.



COLLECT



STORE



MONITOR



SHARE



ANALYZE



Get connected at:
thermofisher.com/connect
 Powered by Thermo Fisher Cloud



For Research Use Only. Not for use in diagnostic procedures.

© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. B-27 is a trademark of Southern Illinois University. Luminex is a registered trademark of Luminex Corp. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. KYMRIA is a trademark of Novartis. C0240857 IN 0118