Molecular biology

PCR Enzyme 솔루션

Invitrogen™ SuperScript™ Reverse Transcriptases
Invitrogen™ Platinum™ DNA Polymerases

invitrogen

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
Traditional cDNA synthesis workflow with DNase I


Figure 1. High efficiency with degraded RNA. RT-qPCR of degraded RNA (RNA integrity number (RIN) 1–3) from human cells and plant tissues with different RTs and Applied Biosystems™ TaqMan® Assays. Delta C_t values ($\Delta C_t = C_t - C_{t \text{ SuperScript IV}}$) show that SuperScript IV RT, with its standard 10 min protocol, delivered higher cDNA yields and lower C_t values than the recommended protocols for SuperScript III reagent and other suppliers' RTs.


Figure 2. Tolerance to inhibitors. This data shows first-strand cDNA synthesis of an RNA ladder in the presence of reaction inhibitors. cDNA fragments were resolved by alkaline gel electrophoresis. All RTs except for SuperScript IV RT were severely affected with trace amounts of inhibitors.

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SuperScript IV First-Strand Synthesis System

- 합성에 필요한 모든 Component가 포함
- 추가적으로 Hela RNA, Control gene primer set 제공
- 구성품 : SuperScript® IV RT, 5X RT buffer, 0.1 M DTT, 10 mM dNTP mix, Oligo(dT)20, Random hexamers, Ribonuclease Inhibitor, E. coli RNase H, DEPC-treated water, Total HeLa RNA, Control Primer set

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SuperScript IV One-Step RT-PCR System

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

Figure 3. High sensitivity and reliable target detection from low amounts of input RNA. A 0.43 kb fragment was successfully amplified using serial dilution from 1,000 ng to 0.01 pg of Invitrogen™ Universal Human Reference RNA (UHRR) and the SuperScript IV UniPrime One-Step RT-PCR System. The molecular weight marker is the Thermo Scientific™ GeneRuler™ 100 bp Plus DNA Ladder, ready-to-use. NTC: no-template control.

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
Figure 4. RT-PCR cycling under two annealing conditions. Nine targets with varying calculated annealing temperatures were amplified from 10 ng UHRR using a universal annealing temperature of 60°C (left), or the annealing temperatures calculated with the Tm calculator for Invitrogen™ Platinum™ SuperFi™ DNA Polymerase (right). The molecular weight marker is the GeneRuler 1 kb Plus DNA Ladder, ready-to-use.

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



Figure 5. Fidelity comparison of commercially available enzymes and Taq enzyme. PCR amplicons (3.9 kb) obtained using different DNA polymerases were fragmented with a MuA transposase. Unique molecular identifiers (UMIs) of 12 random nucleotides were introduced during fragmentation to tag each product individually. After next-generation sequencing, reads from the same UMI family were aligned to call errors. Errors were identified only when present in all reads in the UMI family; otherwise they were discarded as sequencing errors. The fidelity values were normalized to Taq polymerase fidelity.

Benchtop stability for high-throughput applications


No need to calculate primer annealing temperature

High success with amplification of 20–40 kb sequences

Figure 6. Amplification of long fragments. A 20 kb target from human gDNA and 30–40 kb targets from *E. coli* gDNA were successfully amplified using Platinum SuperFi II DNA Polymerase.

No additives needed for GC-rich amplification

Figure 7. Robust amplification of GC-rich targets. Fifteen targets with a range of GC content were amplified from human gDNA without any supplementary buffer additives.

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Fast DNA synthesis at 15 sec/kb


Figure 8. Fast cycling reduces PCR run time. A 0.5 kb fragment was amplified for 35 cycles, using Platinum II Taq Hot-Start DNA Polymerase and hot-start DNA polymerases from other suppliers (A–F).

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Two protocols to fit your needs


Figure 10. Direct and lysis protocols. Two protocols are available to amplify target DNA directly from crude samples. The direct protocol offers a shorter workflow, whereas the lysis protocol allows flexibility and long-term storage.

High sensitivity and specificity


Figure 9. High sensitivity and reliable amplification from low amounts of input DNA. A 0.5 kb fragment was amplified from a range of human gDNA input, using Platinum II Taq Hot-Start DNA Polymerase.

Multiple PCR targets in the same reaction


Figure 11. Efficient multiplexing in direct PCR. Five fragments of 0.1–1.1 kb were amplified from mouse tail in singleplex to 5-plex reactions, using the lysis protocol.

Ordering information

SuperScript Family

Product Label	Size	SKU#
SuperScript™ II Reverse Transcriptase	2,000 units	18064022
	10,000 units	18064014
	4 x 10,000 units	18064071
SuperScript™ First-Strand Synthesis System for RT-PCR	50 reactions	11904018
SuperScript™ One-Step RT-PCR System with Platinum™ Taq DNA Polymerase	100 reactions	10928042
SuperScript™ VIL0™ cDNA Synthesis Kit	50 reactions	11754050
	50 reactions	11755050
SuperScript™ VIL0™ Master Mix	250 reactions	11755250
	2,000 units	18080093
SuperScript™ III Reverse Transcriptase	10,000 units	18080044
	4 x 10,000 units	18080085
	50 reactions	18080051
SuperScript™ III First-Strand Synthesis System	50 reactions	18080400
SuperScript™ III First-Strand Synthesis SuperMix	50 reactions	11752050
SuperScript™ III First-Strand Synthesis SuperMix for qRT-PCR	250 reactions	11752250
	25 reactions	12574018
SuperScript™ III One-Step RT-PCR System with Platinum™ Taq DNA Polymerase	100 reactions	12574026
	25 reactions	12574030
SuperScript™ III One-Step RT-PCR System with Platinum™ Taq High Fidelity DNA Polymerase	100 reactions	12574035
	25 reactions	18080200
SuperScript™ III CellsDirect™ cDNA Synthesis Kit	100 reactions	18080300
	2,000 units	18090010
SuperScript™ IV Reverse Transcriptase	10,000 units	18090050
	4 x 10,000 units	18090200
	50 reactions	18091050
SuperScript™ IV First-Strand Synthesis System	200 reactions	18091200
	50 reactions	18091150
SuperScript™ IV First-Strand Synthesis System with ezDNase™ Enzyme	200 reactions	18091300
	25 reactions	12594025
SuperScript™ IV One-Step RT-PCR System	100 reactions	12594100
	25 reactions	12595025
SuperScript™ IV One-Step RT-PCR System with ezDNase™	100 reactions	12595100
	25 reactions	12596025
SuperScript™ IV UniPrime One-Step RT-PCR System	100 reactions	12596100
	500 reactions	12596500
	25 reactions	12597025
SuperScript™ IV UniPrime One-Step RT-PCR System, Colored	100 reactions	12597100
	500 reactions	12597500
	50 reactions	11756050
SuperScript™ IV VIL0™ Master Mix	500 reactions	11756500
	50 reactions	11766050
SuperScript™ IV VIL0™ Master Mix with ezDNase™ Enzyme	500 reactions	11766500
	50 reactions	11750150
SuperScript™ IV CellsDirect™ cDNA Synthesis Kit	500 reactions	11750350
	500 reactions	11750550

Platinum Family

Product Label	Size	SKU#
Platinum™ Taq DNA Polymerase	120 reactions	10966018
	300 reactions	10966026
	600 reactions	10966034
Platinum™ Taq DNA Polymerase High Fidelity	100 reactions	11304011
	500 reactions	11304029
Platinum™ Taq Green Hot Start DNA Polymerase	600 reactions	11966034
Platinum™ PCR SuperMix High Fidelity	100 reactions	12532016
Platinum™ Hot Start PCR Master Mix (2X)	50 reactions	13000012
	200 reactions	13000013
Platinum™ Green Hot Start PCR Master Mix (2X)	50 reactions	13001012
	200 reactions	13001013
	1,000 reactions	13000014
Platinum™ SuperFi™ DNA Polymerase	100 units	12351010
	500 units	12351050
Platinum™ SuperFi™ PCR Master Mix	100 reactions	12358010
	500 reactions	12358050
	5 x 500 reactions	12358250
Platinum™ Multiplex PCR Master Mix	50 reactions	4464268
	250 reactions	4464269
Platinum™ II Taq Hot-Start DNA Polymerase	100 reaction	14966001
	500 reactions	14966005
Platinum™ II Hot-Start PCR Master Mix (2X)	50 reactions	14000012
	200 reactions	14000013
	1,000 reactions	14000014
Platinum™ II Hot-Start Green PCR Master Mix (2X)	50 reactions	14001012
	200 reactions	14001013
	1,000 reactions	14001014
Platinum™ SuperFi II DNA Polymerase	100 reactions	12361010
	500 reactions	12361050
Platinum™ SuperFi II PCR Master Mix	100 reactions	12368010
	500 reactions	12368050
	5 x 500 reactions	12368250
Platinum™ SuperFi II Green PCR Master Mix	100 reactions	12369010
	500 reactions	12369050
	5 x 500 reactions	12369250
Platinum™ Direct PCR Universal Master Mix	100 reactions	A44647100
	500 reactions	A44647500
	4 x 500 reactions	A44647200

Find out more at [thermofisher.com](https://www.thermofisher.com)

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For Research Use Only. Not for use in diagnostic procedures. © 2024 Thermo Fisher Scientific Inc.
All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

invitrogen