DATA SHEET

# Harvestainer BPC Microcarrier Separation Systems

Small-scale, self-contained, single-use solution for harvesting and separating microcarrier cell cultures

Thermo Scientific<sup>™</sup> Harvestainer<sup>™</sup> BioProcess Containers (BPCs) are closed, single-use microcarrier separation systems that help increase product yields over traditional methods, while reducing clean-in-place and steam-in-place requirements and water-for-injection usage.

Using a Harvestainer BPC helps ensure full containment of the cell culture supernatant in a closed, single-use system for user safety. Additionally, the sophisticated design of these BPCs reduces hold-up volumes to increase recovery efficiencies, yielding more product with less effort. Simple post-use disposal reduces cleaning costs and process cycle times.

For smaller-scale applications, we developed the Thermo Scientific<sup>™</sup> Harvestainer<sup>™</sup> Microcarrier Separation System, an integrated microcarrier separation system designed for when <12 L of microcarrier beads require separation.

The Harvestainer BPCs are available with the Thermo Scientific<sup>™</sup> BioTitan<sup>™</sup> Retention Device. This universal tubing-retention solution was designed to provide the best method for retaining flexible tubing on a barbed fitting, and it helps eliminate the risk of leaks and failures at the tubing connection point.

### **Applications**

The Harvestainer system is designed to retain ≥90 µm of microcarriers while allowing the supernatant to pass through to the next process stage. The Harvestainer system passes USP Class VI Biological Tests for Plastics.



### Key features and benefits

- In-line microcarrier separations
- Preassembled
- Scalable from 2D (3 L and 12 L) to 3D (25 L and 50 L) systems
- Secondary containment
- User safety with closed, single-use system
- High yield (cell recovery ≥85%)
- Reduction of cycle time and the manual process of liquid decanting into separate microcarriers



#### **Design features**

For small-scale applications, when <12 L of microcarrier beads require separation, the 3 L or 12 L Harvestainer system is ideal for in-line microcarrier separation. Each system consists of a preassembled 2D BPC and tray in a complete single-use unit, where the tray acts as the secondary containment device. The Harvestainer system consists of four parts:

- BioProcess Container (BPC)—composed of three layers; the outer two layers are constructed of Thermo Scientific<sup>™</sup> CX5-14 film and the inner layer is constructed of polyester mesh
- Inlet and drain lines—for easy system connection, these lines are made of weldable 3/8 x 5/8 in. C-Flex<sup>™</sup> tubing with a 3/8 in. quick-connect body on the inlet line and a 3/8 in. quick-connect insert on the drain line
- **Support rods**—the Harvestainer system support rods are designed to angle the Harvestainer BPC chamber for optimal drainage and improve recovery rates
- Top support and bottom containment trays—the support and containment trays are made of polyethylene terephthalate glycol (PETG) material. The top tray is designed to support the Harvestainer BPC chamber, while the bottom tray acts as a storage tray, bottom support, and secondary container





Now with ultimate
connection protection
thermofisher.com/biotitan

Description	3 L	12 L	
Tray dimensions (H x W x D)	7.11 x 36.61 x 55.25 cm (2.80 x 15.20 x 21.75 in.)	7.59 x 57.09 x 81.28 cm (2.99 x 22.48 x 32.0 in.)	
Chamber dimensions (H x W)	46.94 x 28.96 cm (18.5 x 11.4 in.)	70.36 x 50.04 cm (27.7 x 19.7 in.)	
Chamber weight	0.11 kg (0.25 lb)	0.25 kg (0.56 lb)	
Chamber surface area	2,303 cm <sup>2</sup> (357 in <sup>2</sup> )	5,909 cm <sup>2</sup> (916 in <sup>2</sup> )	
Mesh surface area	1,000 cm² (155 in²)	2,710 cm² (420 in²)	
Tray material thickness	0.18 cm (0.050 in.)	0.18 cm (0.050 in.)	
Tray material type	PETG		
BPC inlet lineC-Flex tubing; ID x OD: 9.53 x 16.0 mm (3/8 x 5/8 in.)Polycarbonate quick connect: 9.53 mm (3/8 in.) MPC body and MPC cap			
BPC drain line	C-Flex tubing; ID x OD: 9.53 x 16.0 mm (3/8 x 5/8 in.) Polycarbonate quick connect: 9.53 mm (3/8 in.) MPC insert and MPC plug		

## thermo scientific

### Physical and mechanical testing

Testing was conducted to validate the functionality and performance of Harvestainer systems.

- Time and volume of filtration through mesh
- Ability of mesh to retain microcarriers  $\ge$ 90  $\mu$ m
- Cell yield in terms of percent recovery; specification to pass is a recovery of ≥85%
- Cell population viability before and after separation from beads; specification to pass is a viability loss of <3%
- Cell density and cell viability comparable in 2D and 3D systems
- Secondary containment
- Class VI testing data for all subcomponents

The validation report is available upon request.

Table 2. Physical and mechanical test results byBPC size.

Test	Results	
	3 L	12 L
Microcarrier bead retention	≥90 µm	≥90 µm
Flow rates	≤2.5 L/min	≤3.3 L/min
Shipping (ISTA 2A)	Pass	Pass
System integrity	Pass	Pass
Microcarrier bead capacity*	3 L	12 L
Secondary containment	Yes	Yes
Cell yield	≥85%	≥85%
Cell viability	>97%	>97%
Sterility assurance level (25–40 kGy)	10-6	10 <sup>-6</sup>
Maximum internal BPC pressure rating	0.5 psi	0.5 psi

 $^{\ast}$  Use the following calculation to determine the appropriately sized Harvestainer system: Grams of carriers x swell factor of the carriers = total mL/1,000 = total L

### **Ordering information**

Product	Cat. No.
3 L Harvestainer Microcarrier Separation System	SH31078.01
3 L Harvestainer Microcarrier Separation System with BioTitan Retention Device	SH31268.01
12 L Harvestainer Microcarrier Separation System	SH31078.02
12 L Harvestainer Microcarrier Separation System with BioTitan Retention Device	SH31268.02

### Find out more at thermofisher.com/harvestainer



For Research Use or Further Manufacturing. Not for diagnostic use or direct administration into humans or animals. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. C-Flex is a trademark of Saint-GobainPerformance Plastics Corporation. **EXT1011 0721**