

## Handheld Raman spectroscopy for qualitative and quantitative analysis

As the need to quickly identify and quantify more materials increases, pharmaceutical and biotechnology manufacturers require instruments that decrease laboratory sample testing and enable more at-line decisions. The Thermo Scientific™ TruScan™ RM analyzer with TruTools™ embedded chemometrics package provides users with the flexibility to build customized qualitative and quantitative methods for complex material analysis problems.

Applications include:

### QA/QC

- Multiple component discrimination
- Coverage for additional raw materials
- Dosage form identification

### PAT

- Solvent distillation
- Blend analysis and endpoint determination

### Falsified/Substandard Medicines

- Better discrimination of authentic vs. falsified medicine
- API content quantification
- Identification of substandard medicines

### TruScan RM Key Features

- Meets cGMP and 21 CFR Part 11 requirements
- Rugged design; chemical and drop resistant
- Weighs less than 2 lbs (0.9 kg)

The TruScan RM analyzer is the foundation of the system, providing rapid and reliable material identity verification at the point of need. With its state-of-the-art optics and patented, multivariate residual analysis, the analyzer measures a broad range of solid and liquid materials. Method development is fast and simple, requiring minimal samples for creation of a robust qualitative model.

The TruScan RM analyzer with TruTools becomes a more powerful system, enabling users to create customized, advanced qualitative and quantitative methods for



deployment on the analyzer. No longer confined to the lab or benchtop spectrometer, users can conduct advanced chemometric analyses anywhere in the plant.

TruTools leverages Solo, a chemometrics software package from Eigenvector Research Inc. that allows users to develop models that can be deployed onto the TruScan RM analyzer as follows:

- Using TruScan RM with TruTools, the user customizes acquisition parameters including laser power, exposure time, and number of co-adds, then exports SPC files
- Once SPC files are imported into Solo, the user can define spectral range, preprocessing and algorithm type, then process data to generate custom models
- Using Eigenvector's Model\_Exporter, Solo models are exported and deployed on the TruScan RM analyzer via WebAdmin
- Once deployed, TruTools methods are selected through standard menus on the TruScan RM analyzer

The analyzer's embedded math engine runs customized models including:


Model	Type	Description
PLS	Quantitative	Supports quantification of up to 10 chemicals
PCA	Qualitative	Designed for pass/fail results based on a single class
PLSSQ	Qualitative	Suitable for putting limits on quantitative results; produces a pass or fail screen
PLSDA	Qualitative	Used to classify a group of chemicals or identify a chemical from a group of up to 10 chemicals; produces an identification screen

## Training and Support


Our subject matter experts provide training and technical consultation from method development and validation to general operator usage. Once up and running, we provide support anywhere in the world.

Prepared templates and documentation include:

- IQ/OQ/PQ
- SOP Templates
- Statements of Compliance



A PCA TruTools method screen result



A PLS TruTools method screen result.

## Specifications

Raman Spectrum Range	250 to 2875 cm <sup>-1</sup>
Spectral Resolution	8 to 10.5 cm <sup>-1</sup> (FWHM) across range
Laser (excitation wavelength)	785 nm +/-0.5 nm, 2 cm <sup>-1</sup> line width, stability <0.1 cm <sup>-1</sup>
Laser Output Power	250 mW +/-25 mW or customizable (low/med/high) with TruTools
Collection Optics	NA=0.33, 18mm working distance; 0.2 to 2.5 mm spot size
Exposure	Automatic modes (12 ms minimum) or customizable up to 10,000 ms with TruTools
Battery	Rechargeable internal lithium ion battery > 3 hours operation
External Power Supply	DC Wall Adapter, 100-240 V AC 50/60 Hz
Weight	2 lb (0.9 kg)
Size	8.2 in x 4.2 in x 1.7 in (20.8 cm x 10.7 cm x 4.3 cm)
Operating Temperature	-20°C to +40°C (continuous)
Connectivity	Ethernet
Ports	Up to 10 simultaneous ports
Software Version	Requires TruScan RM v.2.6 or later
Chemometrics Package	Works with Eigenvector Solo + Model_Exporter v8.1
Operating Systems and Browsers	Microsoft® Windows® 7,8,10, Internet Explorer® 11, Microsoft Edge™ 25; Google Chrome™ 51
Barcode Supported Symbologies	Most linear and 2D standards
Biometrics	Fingerprint reader for easy login
Measurement Accessories	Vial holder, universal tablet holder, cuvette holder
Compliance	FDA 1040, 21 CFR Part 11, CE certification, Ph. Eur. 8.7

**No. America**  
Boston, USA

**Central and So. America**  
Sao Paulo, Brazil

**Europe, Middle East, Africa**  
Munich, Germany

**Asia Pacific**  
Mumbai, India  
Shanghai, China  
Tokyo, Japan

sales.chemid@thermofisher.com

Find out more at [thermofisher.com/quality](http://thermofisher.com/quality)

**ThermoFisher**  
SCIENTIFIC

**For Research Use Only. Not for use in diagnostic procedures.** © 2016 Thermo Fisher Scientific Inc. All rights reserved. Solo and Model\_Exporter are trademarks of Eigenvector Research Inc. Microsoft, Windows, Internet Explorer and Microsoft Edge are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Chrome is a trademark of Google Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. **2101.1016**