



TDE Series ultra-low temperature chest freezers

TDE Series ultra-low temperature chest freezers

Thermo Scientific™ TDE Series –80°C ultra-low temperature (ULT) freezers, powered by the H-Drive system, are designed for efficiency. As a result, the series delivers a foundational ULT experience with the reliability and performance needed to protect your work—and the sustainability to meet global compliance.

Designed to keep your samples safe

- Worry-free storage; with low peak variation, samples are protected throughout the entire cabinet¹
- Fast door opening recovery; as little as 1 minute for chest freezers
- Secure freezer settings with sample security
- General-purpose models for nondiagnostic general laboratory use, and Class II Medical Device models for storage and preservation of blood and blood products at freezing temperatures²
- Backed by our Global Service and Support team

Limit environmental impact and lower HVAC costs

- Energy use as low as 6.6 kWh/day, depending on size and temperature
- Less heat emission lowers HVAC costs
- SNAP and EU F-gas regulation compliant, natural hydrocarbon refrigerants
- Environmentally friendly, polyurethane foam with non-HFC blowing agent is compliant with SNAP and EU F-gas regulations



User-friendly and quiet

- H-Drive Information Center (HIC) touchscreen user interface
- Quiet enough to reside directly in the laboratory; at just 50 dBA, the sound is comparable to a home refrigerator³
- 5-year parts and labor standard warranty*

* Restrictions and exclusions apply. Warranty is void if the unit is not maintained as set forth in the operation and service manuals.

1. Typical freezer data based on internal testing with freezer setpoint at –80°C and ambient temperature at 20°C. Performance data vary across freezer size and electrical specification.

2. Verification of medical devices in the United States can be found on: <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfrl/rl.cfm>.

3. All data based on internal testing unless otherwise noted. Data on file.


A rich heritage of innovation

From our first ULT freezer in 1940 to today, the combined experience of REVCO™, Forma Scientific™, Jewett™, Harris™, and other families of products continues to drive our ability to deliver industry-leading innovations that help accelerate science and health care.

Easy-to-access controls with setpoint security

- Easy to use: plug in and go, no setup required
- HIC combines easy-to-read temperature with simple capacitive touch buttons
- Advanced security features: alarm settings and setpoint security passcode

Heart—the heart is the health status for the freezer



The diagram shows a dark blue control panel with a digital display showing '-80 °C'. Above the display is a green heart icon. To the right of the display are four touch buttons: a plus sign (+), a checkmark (✓), a minus sign (-), and a gear icon (⚙️). Below the minus sign is a lock icon. To the left of the display is a white arrow pointing right. Lines connect these elements to their respective descriptions.

Easy-to-read alerts

- Alarm bell**—indicates visual and audible alarms
- Door**—illuminates during a door ajar alarm
- Wrench**—generic service warning
- Thermometer**—indicates when the cabinet temperature exceeds either warm alarm or cold alarm setpoints

Easy-to-navigate controls

- Plus**—increases the value of the selected setting
- Check mark**—saves a change to the selected value
- Minus**—decreases the value of the selected setting
- Settings**—represents the various settings

Setpoint security—protect your settings with a 3-digit numeric code

Ready for connectivity

TDE Series ULT freezer models are compatible with the Thermo Scientific™ DeviceLink™ Connect monitoring system. This enables you to have a connected experience with your unit, so you always know the safety of your samples and the performance of your refrigerator or freezer—anytime, anywhere. The units are also compatible with Thermo Scientific™ Smart-Vue™ Pro wireless monitoring systems.

Get your ULT freezers connected

The DeviceLink Connect system connects to your **instrument**. The InstrumentConnect app connects to **you**. Now you can receive real-time information on your equipment even when you're not in the lab.



User-friendly cabinet design

- Ergonomic design: single-hand operation with lockable door handle for easy and secure access to samples
- Flat top workspace for quick access to hands-free operation
- Casters for easy mobility; all are lockable to safely secure unit into place
- Easy-to-access USB for simplified service

Help improve connectivity with automated data collection, storage, and mobile access

- The DeviceLink Connect monitoring solution allows 24/7 access to data and alarms from monitored lab equipment
- Data are securely stored in the cloud and sent to users via email and mobile phone notifications
- Streamlined dashboard-style reports reduce time spent interpreting data


DeviceLink Connect kit ordering information:

| Description | Cat. No. |
|--|--------------|
| DeviceLink Connect Remote Monitoring System (Kits*) for ULT freezers | DLCKTULTS101 |

* Kit includes a DeviceLink Connect system, 1-year basic monitoring subscription, power cable, commissioner USB, monitoring hardware, and an applicable serial or USB cable.

For more information on DeviceLink kit, please visit thermofisher.com/dlc.

For more information on Smart-Vue Pro systems, please visit thermofisher.com/smartvuepro.



TDE Series –80°C chest ULT freezers

Thermo Scientific™ TDE Series –80°C chest ULT freezers, powered by the H-Drive system, enable the same excellent experience as our upright ULT freezers, with an additional benefit of enhanced performance. The combination of our chest ULT freezers' rapid door open recovery and tight peak variation make them perfect for frequent entry or long-term storage applications.



Performance data

| Model series | Energy consumption* (kWh/day) | Peak variation from –80°C setpoint (°C)** | 60-second door opening recovery time (minutes) | Warm-up time (–80°C to –50°C) (minutes)† | Noise level (dBA) |
|--------------|-------------------------------|---|--|--|-------------------|
| TDEC06386 | 6.6 | +5.1/–3.6 | 1 | 141 | 50.9 |
| TDEC25286 | 9.4 | 5.2/–3.4 | 3 | 134 | 55.5 |
| TDEC39686 | 11 | 5.2/–2.8 | 1 | 162 | 52.9 |

Specifications for chest freezers (–50°C to –86°C)

| General-purpose model | Capacity | | Electrical | Amps/breaker (plug type) | Interior dimensions H x W x D in. (cm) | Exterior dimensions H x W x D in. (cm) | Shipping weight lbs. (kg) |
|-----------------------|-------------|-------------------|---------------|--------------------------|--|---|---------------------------|
| | Cu. ft. (L) | Cryoboxes (2 in.) | | | | | |
| TDEC06386FA | 3 (85) | 63 | 115V/60Hz | 10.5/20 (NEMA 5–20) | 16.4 X 18.4 x 18.4 (41.7 x 46.7 x 46.7) | 43.7 x 28.8 x 32.7 (111.0 x 73.2 x 83.1) | 380 (172) |
| TDEC25286FA | 13 (360) | 252 | 115V/60Hz | 16/20 (NEMA 5–20) | 27.8 x 42.4 x 18.5 (70.6 x 107.7 x 47.0) | 40.5 x 72.3 x 33.1 (102.9 x 183.6 x 84.1) | 710 (322) |
| TDEC25286FD | | | 203–230V/60Hz | 12/15 (NEMA 6–15) | | | |
| TDEC39686FA | 20 (566) | 396 | 115V/60Hz | 16/20 (NEMA 5–20) | 27.8 x 66.4 x 18.5 (70.6 x 168.7 x 47.0) | 40.5 x 96.3 x 33.1 (102.9 x 244.6 x 84.1) | 860 (390) |
| TDEC39686FD | | | 203–230V/60Hz | 12/15 (NEMA 6–15) | | | |

Specifications for medical device chest freezers (–50°C to –86°C)

| Medical device model | Capacity | | Electrical | Amps/breaker (plug type) | Interior dimensions H x W x D in. (cm) | Exterior dimensions H x W x D in. (cm) | Shipping weight lbs. (kg) |
|----------------------|-------------|-------------------|---------------|--------------------------|--|---|---------------------------|
| | Cu. ft. (L) | Cryoboxes (2 in.) | | | | | |
| TDEC06386LA | 3 (85) | 63 | 115V/60Hz | 10.5/20 (NEMA 5–20) | 16.4 X 18.4 x 18.4 (41.7 x 46.7 x 46.7) | 43.7 x 28.8 x 32.7 (111.0 x 73.2 x 83.1) | 380 (172) |
| TDEC25286LA | 13 (360) | 252 | 115V/60Hz | 16/20 (NEMA 5–20) | 27.8 x 42.4 x 18.5 (70.6 x 107.7 x 47.0) | 40.5 x 72.3 x 33.1 (102.9 x 183.6 x 84.1) | 710 (322) |
| TDEC25286LD | | | 203–230V/60Hz | 12/15 (NEMA 6–15) | | | |
| TDEC39686LA | 20 (566) | 396 | 115V/60Hz | 16/20 (NEMA 5–20) | 27.8 x 66.4 x 18.5 (70.6 x 168.7 x 47.0) | 40.5 x 96.3 x 33.1 (102.9 x 244.6 x 84.1) | 860 (390) |
| TDEC39686LD | | | 203–230V/60Hz | 12/15 (NEMA 6–15) | | | |

* Energy consumption: Typical data based on internal testing with freezer setpoint at –80°C and 20°C ambient temperature.

** Peak variation from setpoint: Typical data based on internal testing with freezer setpoint at –80°C and 20°C ambient temperature.

† Warm-up time: Typical freezer data based on internal testing with empty freezer setpoint at –80°C and ambient temperature at 20°C.

TDE Series chest ULT freezers— options and accessories







Accessories

| Description | Model series | | | |
|---|--------------|-----------|--------------------|--------------------|
| | Installation | TDEC06386 | TDEC25286 | TDEC39686 |
| LN ₂ Backup System, maintains temperature down to -80°C with liquid nitrogen | Factory | 332353G01 | 332353G05 | 332353G05 |
| | Field | FLN743A* | FLN744A*/FLN744D** | FLN744A*/FLN744D** |
| CO ₂ Backup System, maintains temperature down to -70°C with CO ₂ | Factory | 332353G02 | 332353G06 | 332353G06 |
| | Field | FCO743A* | FCO744A*/FCO744D** | FCO744A*/FCO744D** |
| Chart Recorder (Inkless), 7-day inkless recorder Range: -115°C to +50°C (+5°C resolution) Diameter: 6 in. (15.2 cm) | Factory | 332791G01 | 332166G01 | 332166G01 |
| | Field | FCR741 | FCR742 | FCR742 |
| Chart Recorder (Ink), 7-day ink recorder Range: -115°C to +50°C (+5°C resolution) Diameter: 6 in. (15.2 cm) | Factory | 332791G02 | 332166G02 | 332166G02 |
| | Field | FCR741INK | FCR742INK | FCR742INK |

* A = 115V/60Hz (NEMA 5-20 plug).

** D = 208-230V/60Hz (NEMA 6-15 plug).

| Description | Cat. No. |
|---|-----------|
| Replacement Air Filter Chest Freezers (quantity of 5) | 332297H02 |
| Replacement Backup Battery | 400159 |
| Chart Paper, for Inkless Recorder, 50/pack | 90520H07 |
| Chart Paper, for Ink Recorder, 50/pack | 17020 |
| Cryo Gloves, medium | 4425 |
| Cryo Gloves, large | 4426 |

| | Description/dimensions H x W x D in. (cm) | Freezer capacity cu. ft. (liters) | Boxes/racks (plates) | Racks/freezer (plates) | Boxes/freezer (plates) | Cat. No. |
|---|--|--------------------------------------|-------------------------|---------------------------|---------------------------|----------|
|  | 2 in. Box Rack - Side Access Rack (Tower), 15.2 x 5.5 x 5.5 (38.6 x 14.0 x 14.0) | 3 (85) | 7 | 9 | 63 | 398184 |
|  | 3 in. Box Rack - Side Access Rack (Tower), 16.0 x 5.5 x 5.6 (40.6 x 14.0 x 14.2) | 3 (85) | 5 | 9 | 45 | 398185 |
|  | 2 in. Box Rack - Side Access Rack (Tower), 26.5 x 5.5 x 5.5 (67.3 x 14.0 x 14.0) | 13 (360) | 12 | 21 | 252 | 398186 |
| | | 20 (566) | 12 | 33 | 396 | |
|  | 3 in. Box Rack - Side Access Rack (Tower), 26.3 x 5.5 x 5.6 (66.8 x 14.0 x 14.2) | 13 (360) | 8 | 21 | 168 | 398187 |
| | | 20 (566) | 8 | 33 | 264 | |
|  | Microplate Rack, 26.5 x 5.5 x 5.5 (67.3 x 14.0 x 14.0) | 13 (360) | 27 | 21 | 567 | 398188 |
| | | 20 (566) | 27 | 33 | 891 | |
|  | Wire Basket with Rollers* 12.3 x 16.5 x 11.8 (31.2 x 41.9 x 29.9) | 13 (360) | NA | 5 | NA | 114048 |
| | | 20 (566) | NA | 9 | NA | |

* Basket rail set is required: Cat. No. 189001 for 360 L chest and Cat. No. 189000 for 566 L chest.

Services

Certificate of calibration

Control, alarm, and temperature recording systems are systematically calibrated at the factory for accurate performance. An optional certificate of calibration to a NIST (National Institute of Standards and Technology, formerly National Bureau of Standards)-traceable instrument includes a model number and serial number certificate (single-point calibration). Multiple-point calibration is also available. Calibration certificates must be requested when ordering. For details, contact your sales representative.

Validation testing program

Validation testing services ordered in advance of shipment may accelerate in-house installation qualification (IQ) and operation qualification (OQ) programs. Validation services include documentation and performance reports, compressor and subsystem component specifications, microprocessor flow diagrams, software specifications, and more. For more details, contact your sales representative.

| Services description | Cat. No. |
|---------------------------------------|----------|
| Calibration certificate | 260029 |
| Certificate of conformance | 260043 |
| Unit-specific calibration certificate | 260049 |

Note: Selected models require calibration certificate for NSF 456 vaccine standard labeling.



Also available: Thermo Scientific™ TSX™ Universal Series ULT freezers

Designed to meet the highest protection and sustainability standards for storage at ultra-low temperatures

Find out more at thermofisher.com/chestULT

thermo scientific

For Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific uses or applications. © 2022, 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ENERGY STAR and the ENERGY STAR mark are trademarks owned by the U.S. Environmental Protection Agency. ACT is a trademark of My Green Lab. COL35994 0524