



The pen-like test
for presumptive
detection:

3M™ Clean-Trace™ Surface Protein Swabs

3M™ Clean-Trace™ Surface Protein Plus

is designed to detect residues of protein and reducing agents on surfaces using the biuret reaction.

- Semi-quantitative protein test with swab-click-read simplicity
- Easily accesses hard-to-reach areas
- Results available in 10 minutes or less if the surface contamination is high
- No instrument required
- Based on color change for ease of interpretation



3M™ Clean-Trace™ Surface Protein (Allergen) test


is designed to detect residues of protein and reducing agents on surfaces and in process water using the biuret reaction to indicate presents of potential allergens.

- Highly sensitive: 3 µg of proteins
- Fast: testing takes only 15 minutes
- Swab, click and incubate (55 °C, 15 minutes), read
- Storage at ambient temperature
- Color change makes it easy to read result

Tested allergen	Minimal detection
Egg	10–20 µg
Peanut	25–50 µg
Milk	10–20 µg
Gluten	40 µg
Soy	20–40 µg

3M™ Clean-Trace™ Surface Protein Instant

- Economical alternative for customers who predominantly test flat surfaces
- Instant results



Thermo Fisher Scientific Australia and
New Zealand is an authorised distributor of 3M.