


Industry Information

Pharma

Solutions for the
Pharmaceutical
Industry



Solutions for the Pharma

The most common critical parameters of pharmaceutical processes are pH, oxygen and conductivity. Measurement deviations can have a large impact on the final production result.

The pharmaceutical industry has therefore adapted tools for continuously improving measurement reliability, knowledge and data acquisition. Summarized under the term "Process Analytical Technology (PAT)".



Sensors

Sensors


As a minimum, pharmaceutical sensors must be sterilizable (SIP) and cleanable in place (CIP). To get worthwhile information on their everyday performance, sensors are ideally equipped with a chip that stores relevant quality parameters. Knick Memosens sensors offer this information in an easy to read sensor radar chart (graphic)

The sensor radar chart – available in numerous Knick analyzers – is an ideal tool to monitor and control the sensor performance in real time. The radar chart supports the operator in making maintenance decisions and helps to improve the lifetime of the sensor and to reduce operational costs. An additional "load graph" displays the exposure of the sensor depending on pH and temperature.

Last but not least, Knick offers pharmaceutical certifications with its sensors. These certifications include FDA, USP88, ATEX and others.



Industry




Typical process definition in pharmaceutical industry for R & D to industrial production process.

	Digital Sensor	Measurement	Material	Design	Certifications
	SE 555	pH pH/ORP	Glass Glass, Pt	CIP, SIP, Autoclavable	Cytotoxicity, ADI/TSE, FDA-approved O-ring, ATEX
	SE 565	ORP	Glass, Pt	CIP, SIP, Autoclavable	FDA-approved O-ring, ATEX
	SE 546	pH	PEEK, ISFET	SIP, Autoclavable	Cytotoxicity, FDA-approved O-ring, ATEX
	SE 605H	Conductivity	Stainless Steel, PEEK, EPDM	CIP, SIP, Autoclavable	FDA, USP Class VI, R _a < 0.4 μm, ATEX
	SE 680	Toroidal Conductivity	PEEK	CIP, SIP	FDA, USP<87>Grade 1, USP<88>classVI, EHEDG, ATEX
	SE 706 SE 707	Oxygen Oxygen (Traces)	Stainless Steel, PTFE	CIP, SIP, Autoclavable	ADI/TSE, FDA + USP, ATEX
	SE 740	Optical Oxygen	Stainless Steel, Silicone	CIP, SIP	FDA, R _a < 0.4 μm, ADI/TSE

Solutions for the Pharma



Management



Sensor Management

A crucial requirement in any pharmaceutical environment is the ability to track any changes in the process or measurement and to make sure that data cannot be manipulated. To help customers achieve this, Memosens sensors can be easily connected to an asset management tool developed by Knick: MemoSuite Advanced.

MemoSuite Advanced offers the option to automatically calibrate up to 10 sensors, following GMP standards. This feature allows setting limits that must be met by the sensor during calibration. Failure to meet these limits results in the rejection of the sensor.


Industry



MemoSuite Advanced generates GMP calibration reports that can be automatically uploaded to any higher-level system.

Furthermore MemoSuite Advanced contains a database with all relevant sensor data. All relevant data of the sensor life cycle can be graphically displayed. MemoSuite Advanced is designed according to FDA 21 CFR Part 11, which assures that electronic data cannot be manipulated and electronic signatures are secure.

Solutions for the Pharma



Fittings

Static and Retractable Fittings

Static and retractable fittings provide safety and convenience when the task is to measure pH, ORP, dissolved oxygen and conductivity. Especially for pharmaceutical applications, the ability to clean and sterilize are of utmost importance, as well as the availability of a wide range of materials.

The optimal choice for operation in the pharmaceutical industry without requirement of automation is the ARI106 H, an inexpensive and robust inline fitting for 12-mm sensors.

For extra safety and optimal – automatic – operation, Knick recommends SensoGate H retractable fittings. They are steam-sterilizable, validated for all wetted parts including the rinsing/calibration chamber.

Many “hygienic” fittings on the market have no lock-gate function with only one sealing to the process. This means that when the immersion tube is moved (in no matter what direction), the lower rinsing chamber and “inlet” and “outlet” are flooded with process medium until the next shut-off. The user is forced to use valves and apply elaborate rinsing and steam sterilization cycles to re-establish sterile conditions. SensoGate provides a patented 2-chamber system which safely seals off the process from “inlet” and “outlet”. The customer only needs connections for steam and rinsing agents. The green gasket safely seals off the chamber when the probe is moved – in both directions. Afterwards, the process medium is safely removed from the rinsing chamber using rinse liquid and steam. This is supported by the patented flow body which increases the velocity of the rinsing agents at the inner surfaces.

Industry


Sensogate WA 130 H

Sensogate WA 131 H

Sensogate WA 131 MH

Recommended static and retractable fittings

Fitting	Type	Wetted materials	Design	Sensor type	Certifications
ARI 106 H	Inline	PEEK / 1.4404 1.4435 / 1.4404 PVDF / 1.4404 PVDF / PVDF 1.4539 / 1.4404	CIP, SIP, Autoclavable	12 mm with PG 13.5 or 12 mm pressurizable, length: 120 ... 250 mm	Material trace- ability 3.1, sur- face roughness < 0.8/0.4 µm, FDA, USP<88> class VI
Sensogate WA 130 H Sensogate WA 131 H Sensogate WA 131 MH	Inline retractable fitting (automatic or manual)	1.4404 / 1.4435	CIP, SIP, dairy pipe DN 50 ... 100, Ingold socket 25 mm, clamp 1.5 ... 3.5, Varivent DN 50 ... 80, BioControl DN 50, DN 65	12 mm with PG 13.5 or 12 mm pressurizable, length: 250 mm	ATEX, material trace- ability 3.1, sur- face roughness < 0.8/0.4 µm, inline steam sterilizability, FDA, USP<88> class VI

Solutions for the Pharma



MemoRail



Stratos



Protos



Portavo

Analyzers

Process Analyzers

A wide range of analyzers ensures that each application in the pharmaceutical process can be covered: Starting from the small and easy to install MemoRail for fermenters and cleaning machines to the robust, full-featured modular Protos for batch control applications. Additionally Knick offers Portavo, the only portable analyzers with Memosens technology.

Again, data integrity is a main critical element in any pharmaceutical environment. This is why Knick has implemented the Audit Trail function in many of its analyzers. Audit trail is designed according to FDA 21 CFR Part 11 and guarantees that any change to the analyzer or its data is monitored, stored and traceable.

The analyzers, such as Protos, feature reliable backups and complete data recording in accordance with FDA regulation 21 CFR Part 11. This regulation is basically composed of "Electronic Records" and "Electronic Signature"; both parts interlock directly.

For every change in the menu of the analyzer, there is an exact recording of what was modified and who carried out this modification. For safety reasons, the Audit Trail is stored solely in a specially coded data memory by Knick. Readability from PC is not affected. This Audit Trail memory is ready equipped with a PC program that makes the logged Audit Trail data readable and can also decode them again in "log coding" mode. The data can easily be exported to other programs like for example Microsoft Excel.


As an extraordinary function, Knick analyzers offer a sensor radar chart for pH measurement.

Industry

```

1 Audit Trail-Log
2 Measuring Point BOILER UF0815
3 Front 3400S-011 1527181
4 Base 3400S-021 0006152
5 PHU 3400-110 0007525
6 OXY 3400-061 0000005
7
8 01.12.17 10:43:39 Login: Admin (ADMIN) OK
9 01.12.17 10:43:41 > Configuration (Specialist) OK
10 01.12.17 10:43:57 > Module PHU 3400-110 (Specialist) OK
11 01.12.17 10:43:59 > Configuration (Specialist) OK
12 01.12.17 10:44:00 < Menue Selection OK
13 01.12.17 10:44:03 Login: Admin (ADMIN) OK
14 01.12.17 10:44:05 > Calibration OK
15 01.12.17 10:44:12 > Module PHU 3400-110 (Specialist) OK
16 01.12.17 10:44:20 > Calimatic OK
17 01.12.17 10:44:57 Warn same Buffer OK
18 01.12.17 10:45:06 > Calimatic OK
19 01.12.17 10:45:12 Cal-Protocol OK
20 Last Calibration 01.12.17 10:44 OK
21 User: ADMIN OK
22 Sensor Type: SE555 OK
23 Serial Number: 0815 OK
24 Cal-Mode: Calimatic OK
25 Zero: 7.02 pH OK
26 Slope: 58.9 mV/pH OK
27 Impedance Glass (25.0 °C): 825.9 M0hm OK
28 Impedance Ref (25.0 °C): 119.4 k0hm OK
29 1st Buffer Value: 7.00 pH OK
30 Electrode Potential: 1 mV OK
31 Cal Temperature: 25 °C OK
32 Response Time: 19 s OK
33 01.12.17 10:45:12 < Module PHU 3400-110 OK
34 01.12.17 10:45:16 < Calibration OK
35
36 07.15.17 10:42:15 < Module PHU 3400-110 OK
37 662b026 17#6: 17 2 OK
38 09J 16mb619f06: 52 oC OK
39 ET6C1006 606en179J: 1 #A OK
40 12f 80164 17#6: 1.00 pH OK
41 1mb619f06 661 (52.0 °C): 171.4 k0hm OK
42 1mb619f06 6192 (52.0 °C): 852.1 M0hm OK



```



Recommended analyzers for pharmaceutical applications:

Device	Size	Output	Mounting	Memosens	Certifications	Audit Trail
MemoRail	12.5 mm	4 ... 20 mA	Cabinet (DIN rail)	Yes	ATEX (2), UL	No
MemoRail Modbus	17.5 mm	Modbus RTU	Cabinet (DIN rail)	Yes		No
Stratos (various versions)	148x148x117 mm	4 ... 20 mA, HART, FF, DP, PA	Wall, Pipe	Yes	ATEX (1), FM, CSA, NEPSI	Yes
Protos	213x163x160 mm	4 ... 20 mA, FF, PA	Wall, Pipe	Yes	ATEX (1), FM, CSA, NEPSI	Yes

Solutions for the Pharma



Analyzers

GMP Calibration

Manufacturing pharmaceutical products requires maximum accuracy for the acquisition of measured values. Failure effects must be eliminated as far as possible by calibrations / adjustments.


Limits (tolerances) can be specified for the parameters determined during a calibration. If after several calibrations the zero point or slope values lie outside the specified tolerances, this indicates that the sensor drift is inadmissibly high. A message will be generated and the sensor should be replaced. If the parameters lie within the specified tolerances, no adjustment will be made. This allows tracking the sensor parameters before and after a batch and thus ensuring consistent QM documentation of the production process.

Knick analyzers offer a sensor radar chart for pH measurement as a graphic representation of current sensor parameters – with slope, zero point, reference impedance, glass impedance, response time, calibration timer, and sensor wear.

Cal SOP Calibration Procedure with Portavo

The Cal SOP calibration procedure sets new standards for use in the pharmaceutical sector. pH buffers can be individually selected for up to 3 calibration points. The third buffer acts as a verification buffer and allows users to define the maximum allowed deviation (delta pH). The system indicates deviations from the SOP tolerance and automatically suggests an “adjustment”. This increases the safety of the SOP process. Also new are a multilevel user management with four user profiles and different access rights as well as a sensor verification feature, allowing sensors and devices to be directly associated. This guards against mistakes when operating the analyzer.

Industry



Schematic overview of a typical production process in pharma industry

Configuration

- Sensor verification
 - ↳ Check model Info
 - ↳ Sensor type
 - ↳ Check TAG Info
 - ↳ TAG
 - ↳ Check group Off

Back

Cal data record

14.09.2015 13:29

pH 0.00 7.00 14.00

Zero point pH 6.915 +5.1 mV

Slope 58.6 mV 99.0 %

Response time 11 s

Delta pH pH 0.008

Discard Apply

PIN-Code

Enter your log-in data.

User ADMIN

Pin code ****

Back Continue

Process Analytics

- Industrial transmitters
- Fittings
- Automatic cleaning and calibration systems
- Sensors
- Portables
- Laboratory meters

Knick The Art of Measuring

Knick has been among the leading manufacturers of electronic measurement devices for process analysis for more than 70 years. Today, the company alignment is still focused on a high technical level and an intensely innovative orientation.

The current product range of Knick includes unique device series also for extremely difficult applications. As a system provider, Knick can guarantee expert analysis of your measurement requirements and provide qualified advice for complete measurement loops.

**Knick
Elektronische Messgeräte
GmbH & Co. KG**

Beuckestraße 22, 14163 Berlin
Telefon: +49 30 80191-0
Telefax: +49 30 80191-200
info@knick.de · www.knick.de