

Sani-Cloth® Chlor: Instructions for Use

Two in one sporicidal cleaning & disinfectant wipes

High-Level disinfection. For non-invasive medical devices and non-porous hard surfaces. It can be used to pre-clean or decontaminate medical devices prior to sterilisation. Where surfaces are soiled, pre-cleaning is required. For professional use only.



Preparation



Dispense




- + Dispense the wipe at a 45° angle, pulling away from you.
- + Close the lid when not in use.

Clean



- + If the surface is visibly soiled, clean the surface following an S-shaped motion from top to bottom.
- + Use additional wipes as needed.
- + Use a fresh wipe for each new surface.

Disinfect



- + Using a fresh wipe, disinfect the surface in an S-shaped motion from top to bottom, avoiding cross-contamination.
- + Use a fresh wipe for each new surface.

Contact Time



- + Ensure the surface remains visibly wet for 2 minutes for complete disinfection.
- + Use an additional wipe if necessary.
- + Allow to air dry.

Disposal



- + Dispose of wipe after a single use.
- + Packaging is recyclable



Sani-Cloth® Chlor: Instructions for Use

Two in one sporicidal cleaning & disinfectant wipes

Caution & other information:

- + Store upright in a cool, dry place.
- + Wear protective gloves. Avoid contact with eyes, skin & clothing.
- + Refer to the expiry date on the packaging.
- + Approx.1 wipe per 1m².
- + IF ON SKIN: Wash with plenty of water. IF IN EYES: Rinse cautiously for several minutes. If irritation occurs, discontinue use.
- + Dispose of contents and container in accordance with local regulations

Contra-indications:

- + Do not use for terminal disinfection of invasive medical devices.
- + Any serious incident that occurs in relation to the device should be reported to the ANZ Distributor
- + Compatibility with medical devices should be checked before use



Product Efficacy

Sani-Cloth Chlor	Test Method	Organism	Contact time	
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c)	2 min (d)
	EN 16615	Staphylococcus aureus	60 sec (c)	2 min (d)
	EN 16615	Enterococcus hirae	60 sec (c)	2 min (d)
	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
Yeast	EN 16615	Candida albicans	60 sec (c)	2 min (d)
	EN 13624	Candida albicans	30 sec (c)	
	EN 13624	Candida auris	30 sec (c)	
Enveloped Virus	EN 14476	Vaccinia Virus	30 sec (c)	60 sec (d)
	EN 14476	Human Corona Virus		60 sec (d)
Non-enveloped Virus	EN 14476	Norovirus	30 sec (c)	60 sec (d)
	EN 14476	Adenovirus	60 sec (c)	60 sec (d)
	EN 14476	Poliovirus	60 sec (c)	60 sec (d)
Spores	ASTM E 2197-02	C.Difficile	4 min	

*Clean Conditions (c) Dirty Conditions (d)

Authorised distributor

In Australia:

For customer service, call 1300-735-292
To email an order, ordersanz@thermofisher.com
To order online: thermofisher.com

In New Zealand:

For customer service, call 0800-933-966
To email an order, ordersanz@thermofisher.com
To order online: thermofisher.com

ThermoFisher
SCIENTIFIC

Find out more about how our Layered Approach to Infection Prevention can help you.

PDI Ltd, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK

11.09.2023 V1 THERMO FISHER SCIENTIFIC

+ PDI
BE THE DIFFERENCE™

135 rev. 300623