

IGF-1 R (β -Subunit) Ab-5 (Clone 1-2)

Mouse Monoclonal Antibody

Cat. #MS-645-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200 μ g/ml) (Purified Ab with BSA and Azide)

Cat. #MS-645-P1ABX or -PABX (0.1ml or 0.2ml at 1.0mg/ml) (Purified Ab without BSA and Azide)

Description: IGF-1R consists of α - and β -subunits, which are disulfide-linked in a β - α - α - β configuration in the mature receptor. The α -subunit is completely extracellular, while the β -subunit spans the membrane and the intracellular portion has intrinsic tyrosine kinase activity.

Comments: Ab-5 is the antibody-of-choice for Western blotting of β -subunit of IGF-1R. Ab-3 reacts only weakly with hybrids of Insulin Receptor and IGF-1R while Ab-5 reacts very well.³ So, Ab-3, in conjunction with Ab-5, has proved useful in the isolation of hybrid receptors.³

Mol. Wt. of Antigen: 95kDa

Epitope: aa 1323-1337

Species Reactivity: Human^{1,2} and Rat.^{2,3} Others-not known.

Clone Designation: 1-2

Ig Isotype / Light Chain: IgG_{2b} / κ

Immunogen: A 15-mer peptide corresponding to aa 1323-1337 (Tyr-RKNERALPLPQSSTC) from the C-terminal of β -subunit of human IGF-1R.^{2,3}

Applications and Suggested Dilutions:

- Affinity Purification of HYBRID and IGF-1 R^{2,3} (For conjugation, order Ab without BSA)
- Immunoprecipitation (Not optimum)¹
- Western Blotting (Well-suited for β -subunit)² (Ab 1-2 μ g/ml for 2hrs at RT)

The optimal dilution for a specific application should be determined by the investigator.

Positive Control: Placenta or breast carcinoma

Cellular Localization: Cell membrane

Supplied As:

200 μ g/ml of antibody purified from ascites fluid by Protein A chromatography. Prepared in 10mM PBS,

pH 7.4, with 0.2% BSA and 0.09% sodium azide. Also available without BSA and azide at 1mg/ml.

Storage and Stability:

Ab with sodium azide is stable for 24 months when stored at 2-8°C. Antibody WITHOUT sodium azide is stable for 36 months when stored at below 0°C.

Key References:

1. Soos MA, et. al. J Biol Chem 1992, 267(18):12955-63.
2. Soos MA, et. al. Biochem J, 1993, 290:419-26.
3. Soos MA, et. al. Adv Expl Med Biol, 1993, 343:145-57.

Limitations and Warranty:

Our products are intended FOR RESEARCH USE ONLY and are not approved for clinical diagnosis, drug use or therapeutic procedures. No products are to be construed as a recommendation for use in violation of any patents. We make no representations, warranties or assurances as to the accuracy or completeness of information provided on our data sheets and website. Our warranty is limited to the actual price paid for the product. NeoMarkers is not liable for any property damage, personal injury, time or effort or economic loss caused by our products.

Material Safety Data:

This product is not licensed or approved for administration to humans or to animals other than the experimental animals. Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. Appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. The material contains 0.09% sodium azide as a preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material as indicated above. The National Institute of Occupational Safety and Health has issued a bulletin citing the potential explosion hazard due to the reaction of sodium azide with copper, lead, brass, or solder in the plumbing systems. Sodium azide forms hydrazoic acid in acidic conditions and should be discarded in a large volume of running water to avoid deposits forming in metal drainage pipes.

For Research Use Only

