

Progesterone Receptor (SP2)

Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9102-R7 (7.0ml)

Please note this data sheet has been changed effective November 17, 2016

FOR IN VITRO DIAGNOSTIC USE**INDICATIONS AND USE****Intended Use**

NeoMarkers Rabbit Monoclonal anti-Human Progesterone Receptor Antibody (Clone SP2) is an immunohistochemical (IHC) assay intended for laboratory use for the qualitative detection of Progesterone receptor (PR) antigen by light microscopy in sections of formalin fixed, paraffin embedded normal and neoplastic tissues on a Lab Vision automated slide stainer. It is indicated as an aid in assessing the likelihood of response to therapy as well as in the prognosis and management of breast cancer patients.

Summary and Explanation

Immunohistochemistry has been used to detect specific antigens in cells or tissue since 1950⁶. The use of enzymes and peroxidase as markers for immunohistochemistry was reported by Nakane and Pierce in 1967²⁶. The increased sensitivity of the avidin-biotin-peroxidase detection system over the enzyme labeled antibody method was documented by Hsu *et.al.* in 1981¹⁷.

Determination of PR status for all primary breast carcinomas was recommended by the NIH in 1979, in order to better determine appropriate therapy. In 1985 both the NIH and the American Cancer Society independently published reports in support of determining hormone receptor status as an aid in the management of breast cancer. A number of methodologies to assess PR status have been in use. FDA – cleared therapies include cytosol receptor assay (SBA/DCC) analyzed by Scatchard plot (1981), histochemical analysis of tissue using fluorescent microscopy, histochemical analysis of frozen tissue using anti-PR rat monoclonal antibody conjugate (1988), and enzyme immuno-assay (EIA) also using anti-PR rat monoclonal antibody conjugate (1988)²⁹. The immunohistochemical detection of PR has been described in cultured human breast cells¹⁸, some human breast cancer tissues¹⁹, human endometrium³⁵, some endometrial cancers³, some low-grade endometrial stromal sarcomas²⁸, some sweat gland tumors⁴², some thyroid cancers⁸, some gastric cancers³², some prostatic carcinomas²⁴ and some female human bladders³³.

Principle of Procedure

Lab Vision's NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) binds to PR in the paraffin embedded tissue section. The specific antibody is localized by a biotin conjugated secondary antibody formulation that recognizes rabbit immunoglobulins. This step is followed by the addition of an avidin/streptavidin enzyme conjugate that binds to the biotin present on the secondary antibody. The specific antibody secondary antibody avidin/streptavidin enzyme complex is then visualized with a precipitating enzyme reaction product, which is readily detected by light microscopy. Each step is incubated for a precise time and at room temperature. At the end of each incubation step, the Lab Vision automated slide stainer (Lab Vision Autostainer) washes the sections to stop the reaction and remove unbound material that would interfere the desired reaction in subsequent steps.

Clinical cases should be evaluated within the context of the performance of appropriate controls. Lab Vision recommends the inclusion of a positive tissue control fixed and processed in the same manner as the patient specimen placed on every slide run in addition to the case tissue. In addition to staining with NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2), a second slide should be stained with Negative Control for Rabbit IgG. For the test to be considered valid, the positive control tissue should exhibit nuclear staining of the tumor cells, normal breast glands or endometrial glands and stroma. These components should be negative when stained with Negative Control for Rabbit IgG. In addition, it is recommended that a negative tissue control slide be included for every batch of samples processed and run on Lab Vision Autostainer. This negative tissue control should be stained with NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) to ensure that the antigen enhancement and other treatment procedures did not create false positive staining.

MATERIALS AND METHODS**Reagents**

1. #RM-9102-S, (or -S0, -S1) is 1ml (or 0.1ml, 0.5ml) concentrated tissue culture supernatant containing rabbit anti-human monoclonal antibody directed against PR antigen, with 0.05% Sodium Azide.

2. #RM-9102-R7 (Ready to Use) is provided as a rabbit anti-human monoclonal antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.


A prokaryotic recombinant protein corresponding to 412-562 amino acid of human progesterone receptor molecule was used as the immunogen. New Zealand white rabbits were immunized with the recombinant antigen. The sera were tested by immunoassay and immunohistochemical staining. The rabbit with the best titer in the immunoassay and IHC was selected for a final intravenous boost four days before removal of the spleen. Fusions were performed using lymphocytes from an immunized rabbit and the fusion partner (240E-w). Supernatants were tested for the presence of antibody, specific for the immunogen, by ELISA. Immunohistochemistry and western blotting was also used as screening assays. The hybridomas were sub-cloned by limit dilution. The final antibody is produced from hybridoma culture supernatants using serum free media with no further purification. The antibody is diluted in 10mM phosphate buffered saline with 0.3% carrier protein and 0.05% sodium azide. The total protein concentration of the reagent is approximately 0.4 mg/ml. Specific antibody concentration is approximately 160µg/ml in the concentrate format and 1.6µg/ml in the ready-to-use format. There are no known irrelevant antibodies in the preparation.

Reconstitution, Mixing, Dilution, Titration

NeoMarkers Ready-To-Use antibodies (#RM-9102-R7) have been optimized for use on a Lab Vision Autostainer in combination with Lab Vision UltraVision Detection Systems (see Instructions for Use), and should not require further dilution. No reconstitution or mixing is required. Further dilution may result in loss of sensitivity. The user must validate any such change.



Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision



Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcorn, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049

Progesterone Receptor (SP2)

Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9102-R7 (7.0ml)

Please note this data sheet has been changed effective November 17, 2016

NeoMarkers concentrated antibodies (#RM-9102-S, -S1, -S0) must be diluted in accordance with the staining procedure when used with Lab Vision UltraVision Detection Systems (see Instructions for Use for recommended dilution). Use of non-Lab Vision systems other than recommended systems and protocols require validation by the user.

Differences in tissue processing and technical procedures in the laboratory may produce significant variability in results and require regular use of controls, see Quality Control Procedure section.

Materials and Reagents Needed But Not Provided

The following reagents and materials may be required but are not provided:

1. NeoMarkers Negative Control for Rabbit IgG
2. Lab Vision Antibody Diluent (Catalog # TA-125-UD)
3. Lab Vision Autostainer Systems
4. Microtome
5. Microscope slides, treated to enhance tissue adherence
6. Positive and negative tissue controls
7. Drying oven capable of maintaining a temperature of 70°C +/- 5°C
8. Bar code labels (appropriate for negative control and primary antibody being tested)
9. 10% neutral buffered formalin
10. Staining jars or baths
11. Staining Dishes
12. Timer
13. Xylene
14. Ethanol or reagent alcohol
15. Deionized or distilled water
16. Epitope recovery/tissue pretreatment reagents
17. Lab Vision UltraVision Immunohistochemistry Detection Systems
18. Blocking Reagents
19. Chromogens
20. Counter stains and mounting media
21. Light microscope

Storage and Handling

Store NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) at 2-8°C. This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user. Lab Vision's NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) can be used immediately after removal from storage.

Indications of Instability

When properly stored, the reagent should be stable to the dating indicated on the label. NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) has designed the PR primary antibody (Clone SP2) to have a 24 month stability from the date of manufacture. The user must honor the expiration date on the label. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with unknown specimens. Positive controls assure that the specimen staining was carried out correctly. Negative controls are used to assess non-

specific staining, which must be taken into consideration when interpreting results. Whenever positive control material shows a decrease in staining, it is a possible indication of reagent instability and Lab Vision customer Service (1-510-991-2800 / 1-800-828-1628) should be contacted immediately.

Specimen Collection and Preparation for Analysis

Formalin-fixed, paraffin-embedded tissues are suitable for use with Lab Vision's NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) when used with Lab Vision's UltraVision Detection Systems and Lab Vision Autostainer.

For Formalin-fixed, Paraffin-Embedded Tissues:

Tissue specimens should be preserved by 10% buffered formalin (pH 7.4) as the tissue fixative followed by paraffin embedding. Paraffin-embedded tissue sections should be cut at 3-5 micrometer thick and mounted on treated glass slides to enhance tissue adherence. Tissue slides should be dried at 60-70°C for 1-2 hours. Alternatively, tissue slides can be dried overnight at 37°C. Cool to room temperature if slides are to be stored.

Manual Deparaffinization Procedure:

1. De-wax slides by immersing in xylene for 3 x 5 minutes.
2. Hydrate slides in 100%, 100%, 95 %, 80% ethanol for 3 minutes each, then immerse slides in tap water for 5 minutes. Do not allow slides to dry.

Endogenous Peroxidase Quenching (for horseradish peroxidase detection method):

Immerse slides in 3% hydrogen peroxide solution for 10 minutes, then wash slides in PBS for 2 x 3 minutes. Do not allow slides to dry.

Pretreatment/Antigen Retrieval:

Heat Induced Epitope Retrieval (HIER) is recommended:

- a. Place slides in Citrate Buffer, pH 6.0 (Catalog # AP-9003-500).
- b. Heat samples to near boiling (95°C-98°C) for 10-20 minutes. Some may require longer heating times and/or higher temperatures. Cool slides in buffer at room temperature for at least 20 minutes before proceeding.
- c. Rinse in PBS at least 3 x 1 minute before proceeding. Do not allow slides to dry.


INSTRUCTIONS FOR USE

NeoMarkers primary antibodies have been developed for use on a Lab Vision Autostainer in combination with Lab Vision UltraVision Detection Systems and accessories. The recommended procedures for NeoMarkers Anti-PR antibody (Clone SP2) are as follows:

- Tissue Section Pretreatment: Citrate Buffer, pH 6.0, 10 minutes at 95°C-98°C, cool 20 minutes at room temperature.
- Apply Primary Antibody:
 - Dilution of Concentrated Antibody (#RM-9102-S, -S1, -S0): 1:100 to 1:400 in antibody diluent (we



Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision



Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcorn, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049

Progesterone Receptor (SP2)

Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9102-R7 (7.0ml)

Please note this data sheet has been changed effective November 17, 2016

- recommend users to apply the diluted antibody to the slide as soon as possible)
- or
- Ready-To-Use Antibody (#RM-9102-R7): No further dilution, reconstitution, or mixing is required.
 - Primary Antibody Incubation Time: 30 minutes at room temperature
 - Visualization: To detect antibody, follow the instructions provided with the visualization system. The following are instructions specific to Lab Vision's UltraVision Detection Kit #TP-015-HA:
 - Cover tissue with 4-5 drops of UltraVision biotinylated goat anti-polyvalent secondary antibody, incubate at room temperature for 10 minutes, then wash in PBS for 3 x 3 minutes.
 - Add drops of streptavidin peroxidase to cover section, incubate at room temperature for 10 minutes, then wash in PBS for 3 x 3 minutes.
 - Add 20-50µl (1 drop) of AEC Chromogen to 1 ml of AEC Substrate and mix well. Apply mixture to tissue section within 20 minutes of mixing. Incubate tissue section for 5-20 minutes, depending on the desired stain intensity.
 - Counterstain with hematoxylin according to your laboratory's procedures.
 - Positive Control Tissue: Breast carcinoma known to be positive for progesterone receptor.

Using Lab Vision Autostainer:

- Apply slide barcode label, which corresponds to antibody protocol to be performed. Or, program Lab Vision Autostainer for the protocol you wish to use (See Autostainer instruction manual).
- Load slides and scan barcodes. Lab Vision Autostainer will indicate which reagents are necessary.
- Load reagents as indicated by Lab Vision Autostainer reagent map.
- Start the staining run.
- At the completion of the run, remove slides from the automated slide stainer.
- Counterstain if necessary.
- For DAB, dehydrate, clear and Coverslip with permanent mounting media.
- For AEC or Fast Red, do not dehydrate and clear. Mount with aqueous mounting medium.
- The stained slides should be read within two to three days of staining and are stable for at least two years if properly stored at room temperature (15 to 25°C).

2. Take reasonable precautions when handling reagents. Use disposable gloves when handling suspected carcinogens or toxic materials, for example xylene or formaldehyde.
3. Avoid contact of eyes, skin and mucous membranes with reagents. If reagent comes in contact with sensitive areas, wash with copious amounts of water.
4. Patient specimens and all materials should be handled as if capable of transmitting infection and disposed of with proper precautions.
5. Never pipette by mouth.
6. Avoid microbial contamination of reagents as this could produce incorrect results.
7. Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.
8. The reagents have been optimally diluted and further dilution may result in loss of antigen staining. Any such change must be validated.

Warnings

The reagents contain sodium azide. Avoid contact with skin and mucous membranes. Flush affected areas with copious amounts of water. Seek immediate medical attention if reagents are splashed in eye or ingested.

Sodium azide may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up. Clean exposed metal surfaces with 10% sodium hydroxide.

QUALITY CONTROL PROCEDURE**Positive Tissue Control**

A positive tissue control must be run with every staining procedure performed. The positive staining cells or tissue components (nuclear staining of tumor cells and normal gland cells) are used to confirm that the antibody was applied and the instrument functioned properly. This tissue may contain both positive and negative staining cells or tissue components and serve as both the positive and negative control tissue. Control tissues should be fresh autopsy, biopsy or surgical specimens prepared or fixed as soon as possible in a manner identical to the test sections. Such tissues may monitor all steps of the procedure, from tissue preparation through staining. Use of a tissue section fixed or processed differently from the test specimen will provide control for all reagents and method steps except fixation and tissue processing. This positive tissue control should only be utilized for monitoring the correct performance of processed tissues and test reagents, not as an aid in determining a specific diagnosis of patient samples. If the positive tissue controls fail to demonstrate positive staining, results with the test specimens should be considered invalid.

Negative Tissue Control


Use a tissue control known to be negative for progesterone receptor and to be fixed, processed and embedded in a manner identical to the patient sample(s) with each staining run to verify the specificity of NeoMarkers anti-PR antibody (Clone SP2) for demonstration of PR, and to provide an indication of specific background staining (false positive staining). Also a variety of different cell types in most tissue

SAFETY ISSUES**Precautions**

1. This antibody is intended for *in vitro* diagnostic use.



Richard Allan Scientific Co.
 Subsidiary of Thermo Fisher Scientific
 4481 Campus Drive
 Kalamazoo, MI 49008, USA
 US Toll Free: 1 (800) 522-7270
 Phone: 1 (269) 544-5600
 www.thermofisher.com/labvision



Thermo Shandon Limited
 Subsidiary of Thermo Fisher Scientific
 Tudor Road, Manor Park
 Runcorn, Cheshire, WA7 1TA, UK
 Tel: +44 (0) 1928 534 050
 Fax: +44 (0) 1928 534 049

Progesterone Receptor (SP2)

Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9102-R7 (7.0ml)

Please note this data sheet has been changed effective November 17, 2016

sections can be used by the laboratorian as internal negative control to verify NeoMarkers anti-PR antibody (Clone SP2) performance specifications. There should be no positive staining of the negative tissue control. If specific staining occurs in the negative tissue control sites, results with the patient specimens should be considered invalid.

Negative Reagent Control

A negative reagent control should be run for every specimen to aid in the interpretation of results. A negative reagent control is one that is not specific for progesterone receptor and is used in place of the primary antibody to check nonspecific staining. This provides an indication of nonspecific background staining for each slide. Any staining observed on the specimen is probably due to non-specific protein binding or non-specific binding of other reagents.

Assay Verification

Prior to initial use of this antibody in the user's laboratory or if there is a change of lot number, the specificity of the antibody should be verified by staining a number of positive and negative tissues with known performance characteristics. These quality control tests should be repeated for each new lot or whenever there is a change of a lot number of one of the reagents in a matched set or a change in assay parameters. Refer to the quality control procedure previously outlined in this section of the product insert and to the quality control recommendations of the CAP certification program for immunohistochemistry⁵ as well as the NCCLS IHC guideline³⁰. Assay verification on a daily basis may be accomplished through the proper use of the above-mentioned positive and negative controls. In addition, it is recommended that on a monthly basis, the PR positive tissue control be stained and compared to the same tissue control stained in previous month. Comparison of controls stained at monthly intervals serves to monitor the assay stability, sensitivity, specificity, and reproducibility. All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

Interpretation of Results

The Lab Vision automated immunostaining procedure results in a colored reaction product to precipitate at the antigen sites localized by NeoMarkers anti-PR antibody (Clone SP2). The cellular staining pattern for NeoMarkers anti-PR antibody (Clone SP2) is nuclear. Cytoplasmic staining is considered to be background, non-specific staining, if seen. Nuclei of positive tissue control should be positively stained, whereas nuclei of negative tissue control should not be stained. If the positive tissue control fails to give appropriate positive staining, or if specific staining is found in the negative tissue control, results with the test specimens should be considered invalid. Lab Vision recommends using a semi-quantitative scoring system for staining intensity, e.g., 0-3 for negative, weak, moderate and strong staining.

LIMITATIONS

General Limitations


1. Immunohistochemistry (IHC) is a multistep diagnostic process that requires specialized training in the selection of the appropriate reagents, tissue selection, fixation,

processing, preparation of the IHC slide, and interpretation of the staining results. Excessive or incomplete counterstaining may compromise proper interpretation of results.

2. Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities with the tissue.
3. Excessive or incomplete counterstaining may compromise proper interpretation of results.
4. The clinical interpretation of any positive staining, or its absence must be evaluated within the context of clinical history, morphology and other histopathological criteria. The clinical interpretation of any staining, or its absence must be complemented by morphological studies and proper controls as well as other diagnostic tests. This antibody is intended to be used in a panel of antibodies. It is the responsibility of a qualified pathologist to be familiar with the antibodies, reagents and methods used to produce the stained preparation. Staining must be performed in a certified licensed laboratory under the supervision of a pathologist who is responsible for reviewing the stained slides and assuring the adequacy of the positive and negative controls.
5. Lab Vision provides antibodies and reagents at optimal dilution for use when the provided instructions are followed. Any deviation from recommended test procedures may invalidate expected results. Appropriate controls must be employed and documented. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results.
6. The product is not intended for use in flow cytometry, where performance characteristics have not been determined.
7. Reagents may demonstrate unexpected reactions in previously untested tissues. The possibility of unexpected reactions even in tested tissue groups cannot be completely eliminated because of biological variability of antigen expression in neoplasms, or other pathological tissues¹⁶. Contact Lab Vision with documented unexpected reactions.
8. Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen (HbsAg) may exhibit nonspecific staining with horseradish peroxidase³².
9. When used in blocking steps, normal sera from the same animal source as the secondary antisera may cause false negative or false positive results because of autoantibodies or natural antibodies.
10. False positive results may be seen due to non immunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (e.g. erythrocytes), or endogenous peroxidase activity (Cytochrome C) or endogenous biotin (liver, breast,



Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision



Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcorn, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049

Progesterone Receptor (SP2)

Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9102-R7 (7.0ml)

Please note this data sheet has been changed effective November 17, 2016

brain kidney) depending on the type of immunostain used²⁶.

11. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissues assayed.

Specific Limitations

1. Lab Vision's NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) has been optimized for a 30 minute incubation time, with antigen enhanced tissue. Due to variation in tissue processing, however, it may be necessary to increase or decrease the PR Primary Antibody incubation time on individual specimens. Users who deviate from the recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.
2. Primary antibody incubation time depends on the degree of tissue fixation and may range from 4 to 32 minutes. Lab Vision recommends 30 minutes for use with its detection kits. For further information on fixation variables refer to Immunomicroscopy: A Diagnostic Tool for the Surgical Pathologist⁴³.
3. The following normal tissue was not tested: parathyroid. The user should determine appropriate staining in the above tissue prior to interpretation of staining information.
4. NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) negative result does not exclude the presence of PR. Negative reactions in breast carcinomas may be due to loss or marked decrease of expression of antigen. Therefore, it is recommended that this antibody be used in a panel of antibodies including progesterone receptor.

SUMMARY OF EXPECTED RESULTS

Immunoreactivity of NeoMarkers rabbit monoclonal anti-human PR antibody (Clone SP2) was demonstrated by a study using clinical specimens that showed appropriate staining in formalin fixed, paraffin embedded breast carcinoma tissue.

Performance Characteristics:

A Positive staining result is defined as more than 10% of tumor cells with stained nuclei of any intensity.

Method Comparison:

When compared to a predicate device, NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) showed agreement of 95.20% (238/250), with a 95% confidence interval of 91.77-97.5%⁴⁵. Positive percent agreement is 98.83% (169/171), with a 95% confidence interval of 95.84-99.86%⁴⁵. Negative percent agreement is 87.34% (69/79), with a 95% confidence interval of 77.95-93.76%⁴⁵.

In addition, Cano et. al⁴⁶ investigated the PR status on 40 paraffin sections from breast cancer patients (39 women and 1 man, ages 41–82 yr) using NeoMarkers' rabbit monoclonal antibody anti-PR (Clone SP2) without antigen retrieval. The results were compared with assessment by the classic method (with antigen retrieval) with a mouse monoclonal anti-PR antibody (Clone 1A6). Total agreement

was observed: The PR detection in paraffin sections using SP2 antibody showed a sensitivity of 100% (18/18), a specificity of 100% (22/22), and an accuracy of 100% (40/40) when compared to antibody 1A6.

Specificity:

Specificity of NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) was determined by a study that tested formalin fixed, paraffin embedded normal tissues. The 30 normal tissues examined included adrenal, bone marrow, breast, cerebrum, cervix, colon, endometrium, esophagus, heart, kidney, liver, lung, mesothelium, ovary, pancreas, pituitary, testis, thyroid, prostate, stomach, small intestine, salivary gland, skeletal muscle, skin, placenta, cerebellum, spleen, tonsil, thymus, peripheral nerve. NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) showed specific staining of nucleus in normal breast, uterus, and cervix and pituitary tissues. This is expected, as it is known progesterone receptors are present in these tissues. Lab Vision also tested 40 neoplastic tissues using the same method as for normal tissue testing. The tissues examined included lung, esophagus, stomach, small intestine, colon, rectum, liver, pancreas, kidney, prostate, thyroid, smooth muscle, soft tissue, bone, leg, lip, skin, nerve, brain, lymphatic system, testis, placenta, bladder, ovary, uterine, and breast. All cases tested negative for PR with the exception of ovary, uterine and breast. NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) showed no specific positive staining in tumor tissue that is not expected to have PR. Sensitivity is dependent upon the preservation of the antigen. Any improper tissue handling during fixation, sectioning, embedding or storage, which alters antigenicity, weakens PR detection by NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2).

Reproducibility:


Reproducibility of staining with NeoMarkers rabbit monoclonal anti-PR antibody (Clone SP2) was determined by staining 10 slides containing the same tissue for Intra-run and 10 slides each day for ten days for Inter-run reproducibility. Results between slides showed no variation in staining intensity. Users should verify within run reproducibility results by staining several sets of serial sections with low, medium, and high antigen density in a single run (Intra-run) or on different days (Inter-run).

REFERENCES

1. Barnes, DM. Et.al. Immunohistochemical determination of oestrogen receptor: comparison of different methods of assessment of staining and correlation of clinical outcome of breast cancer patients. British J. Cancer, 74: 1445-1451, 1996.
2. Bevit DJ, et. Al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. J. Pathol. 183: 228-232m 1997.
3. Chambers, J.t., et.al. Immunohistochemical evaluation of estrogen and progesterone receptor content in 183 patients with endometrial carcinoma. Part II: Correlation between biochemical and immunohistochemical methods and survival. Am. J. Clin. Pathol. 94: 255-260, 1990.
4. Clinical Laboratory Improvement Amendments of 1988: Final Rule. 57 FR 7163, February 28, 1992.



Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision



Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcom, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049

Progesterone Receptor (SP2)

Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9102-R7 (7.0ml)

Please note this data sheet has been changed effective November 17, 2016

5. College of American Pathologists (CAP) Certification Program for immunohistochemistry. Northfield, IL. <http://www.cap.org>. (800) 323-40410.
6. Coons, A. H., and M. H. Kaplan. Localization of antigen in tissue cells. II. Improvements in method of detection of antigen by means of fluorescent antibody. J Exp. Med. 91: 1-13, 1950.
7. Deamant, F.D. et. Al. Estrogen receptor immunochemistry as a predictor of site of origin in metastatic breast cancer. Applied Immunohistochemistry, 1: 188-192, 1993.
8. Bur, M. et. Al. Estrogen and progesterone receptor detection in neoplastic and non-neoplastic thyroid tissues. Mod. Pathol. 6(4): 469-472, 1993.
9. Elias, J.M. Immunohistochemistry, a practical approach to diagnosis. ASCP Press, Chicago, 1990.
10. Elias, J.M. et. Al. Special Report: Quality control in immunohistochemistry. Am. J. Clin. Path. 92: 36, 1989.
11. Esteban, J.M. et. Al. Prognostic relevance of carcinoembryonic antigen and estrogen receptor status in breast cancer patients. Cancer 74: 1575-1583, 1994.
12. Green, S. et. Al. Human oestrogen receptor cDNA: sequence, expression and homology to v-erb-A. Nature, 320: 134-139, 1986.
13. Greene, G. L. and Press, M. F. Immunochemical evaluation of estrogen receptor and progesterone receptor in breast cancer. In "Immunological Approaches to the Diagnosis and Therapy of Breast Cancer" Roberto L. Ceriani, editor. Plenum Press, New York, 119-135, 1987.
14. Greene, G. L. et. Al. Sequence and expression of human estrogen receptor complementary DNA. Science. 231: 115-1154, 1986.
15. Hata, H. and Kuramoto, H. Immunocytochemical determination of estrogen and progesterone receptors in human endometrial adenocarcinoma cells (Ishikawa cells). J. Steroid Biochem. Molec. Biol. 42: 201-210, 1992.
16. Herman GE, Elfont EA. The taming of immunohistochemistry: the new era of quality control. Biotech Histochem. 66: 194-199, 1991.
17. Hsu, S. -M., et. al. Use of avidin-biotin-peroxidase complex (ABC) in immunoperoxidase techniques: A comparison between ABC and unlabeled antibody (PAP) procedures. J. Histochem. Cytochem. 29: 577-580, 1981.
18. Malet, C. et. al. Estradiol and progesterone receptors in cultured normal human breast epithelial cells and fibroblasts. J. Clin. Endocrinol Metab. 73(1): 8-17, 1991.
19. Zidan, A. et. al. Oestrogen and progesterone receptor assessment in core biopsy specimens of breast carcinoma. J. Clin. Pathol. 50: 27-29, 1997.
20. Koretz, K. et. al. Metachromasis of 3-amino-9ethylcarbazole (AEC) and its prevention in immunoperoxidase techniques. Histochemistry. 86: 471-478, 1987.
21. Kuukasjarvi, T. et. al. Loss of estrogen receptor in recurrent breast cancer is associated with poor response to endocrine therapy. J. Clin. Oncol. 14: 2584-2589, 1996.
22. Molino, A. et. al. Estrogen receptors in 699 primary breast cancers: A comparison of immunohistochemical and biochemical methods. Breast Cancer Research and Treatment, 34: 221-228, 1995.
23. Money, S.R. , et. al. Immunocytochemical localization of estrogen and progesterone receptors in human thyroid Surger. 106: 975-979, 1989.
24. Mobbs, B.G. et. al. Immunohistochemical localization of progesterone receptor in benign and malignant human prostate. Prostate. 16(3): 245-251, 1990.
25. Nadjji M, and Morales AR. Immunoperoxidase: part 1. The technique and its pitfalls. Lab. Med. 14: 767, 1983.
26. Nakane, P. K. and Pierce, G. B. Enzyme labeled antibodies: preparation and applications for the localization of antigens. J. Histochem. 14: 929-931, 1967.
27. Nedergaard, L. et. al. Immunohistochemical study of estrogen receptors in primary breast carcinomas and their lymph node metastases including comparison of two monoclonal antibodies> APMIS, 103: 20-24, 1995.
28. Nyholm, H.C.J. et. al. Biochemical and immunohistochemical estrogen and progesterone receptors in adenomatous hyperplasia and endometrial carcinoma: Correlations with stage and other clinopathologic features. Am. J. Obstet. Gynecol. 167: 1334-1342, 1992.
29. ODE Guidance documents, Center for Devices and Radiological Health. Draft: Premarketing approval review criteria for premarket approval of estrogen (ER) or progesterone (PGR) receptors in vitro diagnostic devices using steroid hormone binding (SBA) with dextran coated-charcoal (DCC) separation, histochemical receptor binding assays, or solid phase enzyme immunoassay (EIA) methodologies. Updated 7-29-1996.
30. O'Leary, T. J. et. al. Quality assurance for immunohistochemistry: proposed guideline. MM4-P. National committee for clinical laboratory standards (NCCLS). Wayne, PA. 1-46, 1997.
31. Omata, M. et. al. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immunohistochemistry. Am J. Clin. Pathol. 73: 626-632, 1980.
32. van Velthuysen, M-L.F. et. al. Expression of oestrogen receptor and loss of E-cadherin are diagnostic for gastric metastasis and breast carcinoma. Histopathology. 46(2): 153, 2005.
33. Rizk, D.E. et. al. Evidence of progesterone receptors in the mucosa of the urinary bladder. Scand. J. Urol. Nephrol. 35(4): 305-309, 2001.
34. PDQ. Treatment. Health Professionals. Breast Cancer. Downloaded from http://icisun.nci.nih...t_cancer_Physician.html on 4-10-1997.
35. Bergqvist, A. et. al. Immunohistochemical analysis of oestrogen and progesterone receptors in endometriotic tissue and endometrium. Human Reproduction. 8: 1915-1922, 1993.
36. Roche, P.C. Immunohistochemical stains for breast cancer. Mayo Clin. Proc. 69: 57-58, 1994.
37. Sabini, G. et. al. Steroid hormone receptors in endometrial stromal sarcomas. A biochemical and immunohistochemical study. Am. J. Clin. Pathol. 97: 381-386, 1992.
38. Sheehan, D. C. and B. B. Hrapachak. Theory and Practice of Histotechnology. 1980. The C.V. Mosby Company, St. Louis.



Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision



Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcorn, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049

Progesterone Receptor (SP2)**Catalog # RM-9102-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)****Catalog # RM-9102-R7 (7.0ml)****Please note this data sheet has been changed effective November 17, 2016**

39. Shi, S.R. et. al. Antigen retrieval immunohistochemistry under the influence of pH using monoclonal antibodies. J. Histochem. Cytochem. 43: 193-201, 1995.
40. Shi, S.R. et. al. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based on microwave oven heating of tissue sections. J. Histochem. Cytochem. 39: 741, 1991
41. Stierer, M. et. al. A prospective analysis of immunohistochemically determined hormone receptors and nuclear features as predictors of early recurrence in primary breast cancer. Breast cancer Research and treatment, 36: 11-21, 1995.
42. Busam, K.J. et. al. Epidermal growth factor, estrogen, and progesterone receptor expression in primary sweat gland carcinomas and primary and metastatic mammary carcinomas. Mod. Pathol. 12(8): 786-93, 1999.
43. Taylor RT and Cote RJ. Immunomicroscopy: A Diagnostic Tool for the Surgical Pathologist. W.B. Saunders company, Philadelphia, 1994.
44. Yokazaki, H. et. al. Estrogen receptors in gastric adenocarcinoma: a retrospective immunohistochemical analysis. Virchows Archiv. A Pathol. Anat. 413: 297-302, 1988
45. Clopper, CJ & Pearson, ES. The use of confidence or fiducial limits illustrated in the case of the binomial. Biometrika, 26: 404-413, 1934
46. Cano, G. et. al. Estimation of hormone receptor status in fine-needle aspirates and paraffin-embedded sections from breast cancer using the novel rabbit monoclonal Antibodies SP1 and SP2. Diagn. Cytopathol. 29(4): 207-211, 2003.

Manufacturer

Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision

**Authorized
Representative in E.U.:**


Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcorn, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049
sales.ap.uk@thermofisher.com

Other Countries:

Please contact your local Representative



Richard Allan Scientific Co.
Subsidiary of Thermo Fisher Scientific
4481 Campus Drive
Kalamazoo, MI 49008, USA
US Toll Free: 1 (800) 522-7270
Phone: 1 (269) 544-5600
www.thermofisher.com/labvision



Thermo Shandon Limited
Subsidiary of Thermo Fisher Scientific
Tudor Road, Manor Park
Runcorn, Cheshire, WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049