

Compensation and instrument beads

Compensation beads for flow cytometry

Emission profiles of fluorophores are broad, which can result in overlapping profiles that require compensation for signal correction. Compensation can be set using beads, particularly when cell samples are limited or when a positive population is needed.

The latest generation of compensation beads

Build flow cytometry panels with more accurate compensation using new Invitrogen™ UltraComp eBeads™ Plus Compensation Beads. When a fluorophore-conjugated antibody is added to the beads, both positive and negative populations result. This bimodal distribution can be used for single-color compensation controls in multicolor flow cytometry experiments. In addition to all the features of the first-generation Invitrogen™ UltraComp eBeads™ product, UltraComp eBeads Plus Compensation Beads now offer:

- Increased species reactivity including rabbit- and human-origin antibodies (Figure 1)
- Compatibility with fluorophores excited by ultraviolet (355 nm), violet (405 nm), blue (488 nm), green (532 nm), yellow-green (561 nm), and red (633–640 nm) lasers
- Better compensation resolution for antibodies conjugated with Invitrogen™ eBioscience™ Super Bright 780, Brilliant Violet 711, or Brilliant Violet 786 dyes




Figure 1. Staining of UltraComp eBeads Plus Compensation Beads with 14 different Invitrogen™ PE dye-conjugated monoclonal antibodies, including one of each subclass commonly used in flow cytometry. Beads were stained with 0.25 µg of each antibody and analyzed by flow cytometry. Each histogram represents one staining antibody.

Invitrogen™ antibody compensation beads

	UltraComp eBeads™ Plus beads	UltraComp eBeads™ beads	OneComp eBeads™ beads	AbC™ Total Antibody Compensation Bead Kit*	ArC™ Amine Reactive Compensation Bead Kit	GFP BrightComp eBeads™ beads
Application	Immunophenotyping				Cell viability assay	GFP expression; beads are present at 3 levels of GFP-like intensity
Reactivity	Human, rabbit, hamster, mouse, and rat antibodies	Hamster, mouse, and rat antibodies with recognition of the kappa and lambda chains		Hamster, mouse, rabbit, and rat antibodies	Invitrogen™ LIVE/DEAD™ fixable dead cell stains*	GFP isoforms
Format	One vial: dispense as a single drop			1 vial positive beads, 1 vial negative beads		One vial: dispense as a single drop
Laser compatibility	Compatible with most standard lasers, UV to 633 nm; improved for polymer dye use from the violet laser	Compatible with most standard lasers, UV to 633 nm	Compatible with most standard lasers, but not with UV or violet lasers	Compatible with most standard lasers, UV to 633 nm		488 nm
Quantity	25 tests or 100 tests					25 tests
Cat. No.	01-3333-41 01-3333-42	01-2222-41 01-2222-42	01-1111-41 01-1111-42	A10513 A10497	A10628 A10346	A10514

* Also applicable to similar amine-reactive dyes.

Counting beads

Absolute cell counting is a method for quantifying cell concentration or the absolute count of cells in a sample. Benefits of our absolute counting beads include:

- Wide range of fluorophores to fit a broad spectrum (Figure 2)
- Simple protocols that work with multiple cell types
- Increased consistency and reliable results
- Accommodates most cell sizes with increased percentage of singlets




Figure 2. CountBright Plus beads can be used with a broader range of fluorophores. CountBright Plus beads (red) can be detected simultaneously with cells stained with Invitrogen™ CD19 APC-eFluor™ 780 antibody (pink) in lysed whole blood when excited with an IR laser (808 nm) with an 840/20 nm emission filter.

Invitrogen™ absolute counting beads

	CountBright™ Plus beads*	AccuCheck™ beads		LIVE/DEAD™ BacLight™ Bacterial Viability and Counting Kit**
Parameters measured	Cell concentration in sample	<ul style="list-style-type: none"> • Cell concentration in sample • Pipetting accuracy 		<ul style="list-style-type: none"> • Viability • Bacterial concentration in sample
Sample type	Any type	Whole blood		Bacteria
Bead size	4 µm	Bead A 6.40 µm	Bead B 6.36 µm	6 µm
Range	Ex: UV–800 nm Em: 385–860 nm	Bead A Ex: 488 nm Em: 575–585 nm	Bead B Ex: 635 nm Em: 660–680 nm	Ex: 488 nm Em: 617 nm, 498 nm
Cat. No.	C36995	PCB100		L34856

* The original CountBright Absolute Counting Beads are still available, but not compatible with IR-excitable fluorophores.

** Stains all cells, so a pure bacterial sample is required for accurate results.

Calibration and size beads

Instrument calibration is critical to collecting and analyzing accurate experimental data. Our beads are designed to help ensure robust flow cytometer performance.

Invitrogen™ calibration beads

	Size calibration		Instrument control	Alignment control
Product	Flow Cytometry Size Calibration Kit	Flow Cytometry Sub-micron Particle Size Reference Kit	Rainbow Calibration Particles	Alignflow™ Flow Cytometry Alignment Beads
Use	Size reference	Size reference	Routine calibration of flow cytometers	Calibrate laser alignment
Emission	No fluorescence	Green fluorescence	400–680 nm	3 types: 400–470 nm (for UV lasers), 515–660 nm (for blue lasers), or 645–680 nm (for red lasers)
Bead size	6 sizes, 1.0–15 µm range	6 sizes, 0.02–2.0 µm range	3.0–3.4 µm	2 sizes, 2.5 µm or 6.0 µm diameter
Cat. No.	F13838	F13839	A34305	2.5 µm: A16502, A16500, A16501 6.0 µm: A16505, A16503, A16504

Contact your sales representative:

Find out more about flow cytometry beads and controls at [thermofisher.com/flow-controls](https://www.thermofisher.com/flow-controls)