



ThermoFisher
SCIENTIFIC

Invitrogen GeneArt Gene Synthesis API


The world leader in serving science

Contents

- 1 What is an API
- 2 Who can benefit
- 3 What we offer
- 4 How to get started
- 5 Appendix



What is an API?



- *An API (Application Programming Interface) is a set of rules that let programs talk to each other, exposing data and functionality across the internet in a consistent format.*
- *With GeneArt APIs we enable programs on the customer's side to allow access to services and algorithms through thermofisher.com to automate design and ordering for clones and fragments.*

Who can benefit from GeneArt APIs?

GeneArt APIs

Will Bring Benefit To:



Customers ordering large amounts of sequences frequently.



Customers having an automated pipeline to design genes before ordering.



Third-party developers offering commercial platforms for gene design or open source libraries for bioinformaticians.



Design & Diagnostics APIs

- Check sequences against production acceptance criteria and get a detailed description of the found problems
- Optimize sequences for expression in a selected host organism and/or producibility



Upload & Ordering APIs

- Validate a project.
Get a brief problem summary for the rejected constructs and prices for the accepted ones
- Upload a project.
Validation is done in the background, only fully accepted projects can be uploaded
- Submit an uploaded project to the cart.
Checkout is to be done manually on the thermofisher.com website
- Review the status of the uploaded project

Supported products: Strings, High-Q Strings, Gene Syntheses into pMX vectors and direct Gene Synthesis into one of the 9 standard expression vectors.

How to get started

- Each user need be enabled manually. Send a request to geneartapi@thermofisher.com to get instructions.
- To try the API in action use the provided Postman collection with examples for requests.
 - ✓ Install Postman tool (postman.com)
 - ✓ Import “GeneArtApis April2020.postman_collection.json” as a collection.
 - ✓ Select individual requests from the tree on the left.
- To see detailed specifications and to start implementation
 - ✓ Go to Swagger editor (editor.swagger.io)
 - ✓ Upload one of the documentation files: UploadAndOrdering.yaml or DesignAndDiagnosticsBulk.yaml
 - ✓ Automatically generate client code in the language of your choice



ThermoFisher
S C I E N T I F I C

Appendix

Upload & Ordering API example: construct upload

- Request body consists of the two main parts: Authentication and Project.
- Authentication information consists of the thermofisher.com username and access token provided when access is requested.
 - For the purpose of testing one can use example access data, but no confidential sequences should be uploaded with it.
- Project contains list of constructs, each with name, sequence, product type and (optionally) details.
- I.e. for “geneSynthesisPrime” product details can include vector name, delivery quantity and some other options.
- Most important value in Response is project id, which can be later used to add project to cart or review status.

More information can be found in [UploadAndOrdering.yaml](#) file.

The screenshot displays a REST client interface for a POST request to the URL `https://www.thermofisher.com/order/gene-design-ordering/api/upload/v1/`. The request body is shown in JSON format, containing authentication details and a project with two constructs. The response body is also shown in JSON format, indicating the project ID and name.

```
POST https://www.thermofisher.com/order/gene-design-ordering/api/upload/v1/ URL

Params Authorization Headers (11) Body ● Pre-request Script Tests Settings

● none ● form-data ● x-www-form-urlencoded ● raw ● binary ● GraphQL JSON ▼

1 {
2   "authentication": {
3     "username": "synbio.customer@gmail.com",
4     "token": "vc0BsQAZLS4LJnInD5bmj6LvIT1mBvrF6ZB8G05B"
5   },
6   "project": {
7     "name": "my-new-project",
8     "constructs": [
9       {
10        "name": "somePrimeConstruct",
11        "sequence": "AAGAGATTGATGATAGTCATGGGTGCTGAGCTGAGCGCGCTGATCATAGCGGACTTT",
12        "product": "geneSynthesisPrime"
13      },
14      {
15        "name": "exprVectorAndPPrep",
16        "sequence": "CCGGCGTCGCATTTGTCCGCATTGCGGTGGCCCTCTATACCGTAATCGCTCGGCC",
17        "product": "geneSynthesisPrime",
18        "details": {
19          "geneSynthesisPrime": {
20            "vector": {
21              "name": "pcDNA3.4-TOPO"
22            },
23            "deliveryQuantity": "H2O 250ml"
24          }
25        }
26      }
27    ]
28   }
29 }
```

```
Body Cookies (5) Headers (17) Test Results

Pretty Raw Preview Visualize JSON ▼

1 {
2   "project": {
3     "projectId": "2020AAA0GX",
4     "name": "my-new-project",
5     "constructs": [
```



ThermoFisher
S C I E N T I F I C

For Research Use or Further Manufacturing. Not for diagnostic use or direct administration into humans or animals. © 2020 Thermo Fisher Scientific Inc.
All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. 0520

The world leader in serving science