



Find high-quality OEM enzymes optimized for molecular assay development

Select the category of enzymes you need for different applications.
Then answer a series of simple questions to get the best OEM
enzymes for your research and experiments.

For Research Use Only. Not for use in diagnostic procedures.

© 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. 1024

PCR DNA polymerases

What is your application?

PCR DNA polymerases

What is your application?

Need hot start PCR?

PCR DNA polymerases

What is your application?

Need hot start PCR?

Brand	Enzymes		
	Ampli <i>Taq</i> ™ DNA Polymerase	Ampli <i>Taq</i> ™ 360 DNA Polymerase	<i>Taq</i> DNA Polymerase, recombinant
	Applied Biosystems™	Applied Biosystems™	Thermo Scientific™
	TaqMan™ probe compatible	Yes	Yes
	Proofreading activity (3'-5' exonuclease)	No	No
	Lyo-ready/glycerol free	On request	On request
	Inhibitor tolerance	•	•
	Extention rate	• •	• •
	Benchtop stability	No	No
	GC-rich target amplification	Standard	★ High
Fidelity vs Taq Pol	Standard (1x)	Standard (1x)	Standard (1x)

★ Uses GC enhancer

PCR DNA polymerases

What is your application?

Need hot start PCR?

Need lyo-ready/glycerol free?

Brand	Enzymes				
	Platinum™ II <i>Taq</i> Hot-Start DNA polymerase	Platinum™ <i>Taq</i> DNA Polymerase, DNA-free	Platinum™ <i>Taq</i> DNA Polymerase	Ampli <i>Taq</i> Gold™ 360 DNA Polymerase	Ampli <i>Taq</i> Gold™ DNA Polymerase
	Invitrogen™	Invitrogen™	Invitrogen™	Applied Biosystems™	Applied Biosystems™
	TaqMan™ probe compatible	Yes	Yes	Yes	Yes
	Proofreading activity (3'-5' exonuclease)	No	No	No	No
	Lyo-ready/glycerol free	Yes	On request	Yes	On request
	Inhibitor tolerance	●●●	●	●	●
	Extention rate	***	**	**	**
	Benchtop stability	Yes	Yes	Yes	Yes
	GC-rich target amplification	★ High	Standard	★ High	Standard
	Fidelity vs Taq Pol	Standard (1x)	Standard (1x)	Standard (1x)	Standard (1x)

★ Uses GC enhancer

PCR DNA polymerases

What is your application?

Need hot start PCR?

Need lyo-ready/glycerol free?

	Enzymes		
	Platinum™ II <i>Taq</i> Hot-Start DNA Polymerase	Platinum™ <i>Taq</i> DNA Polymerase	Ampli <i>Taq</i> Gold™ DNA Polymerase
	Invitrogen™	Invitrogen™	Applied Biosystems™
	TaqMan™ probe compatible	Yes	Yes
	Proofreading activity (3'-5' exonuclease)	No	No
	Lyo-ready/glycerol free	Yes	Yes
	Inhibitor tolerance	●●●	●
	Extention rate	***	**
	Benchtop stability	Yes	Yes
	GC-rich target amplification	★ High	Standard
Fidelity vs Taq Pol	Standard (1x)	Standard (1x)	Standard (1x)

★ Uses GC enhancer

PCR DNA polymerases

What is your application?

Need hot start PCR?

Need lyo-ready/glycerol free?

	Enzymes	
	Platinum™ <i>Taq</i> DNA Polymerase, DNA-free	Ampli <i>Taq</i> Gold™ 360 DNA Polymerase
Enzymes		
Brand	Invitrogen™	Applied Biosystems™
TaqMan™ probe compatible	Yes	Yes
Proofreading activity (3'-5' exonuclease)	No	No
Lyo-ready/glycerol free	On request	On request
Inhibitor tolerance	●	●
Extention rate	**	**
Benchtop stability	Yes	Yes
GC-rich target amplification	Standard	High
Fidelity vs Taq Pol	Standard (1x)	Standard (1x)

PCR DNA polymerases

What is your application?

Relative fidelity compared to Standard *Taq* Polymerase

	Enzymes							
	Platinum™ SuperFi II DNA Polymerase	Platinum™ SuperFi™ DNA Polymerase	Phusion™ Plus DNA Polymerase	Phusion™ Hot start II DNA Polymerase	Phusion™ High-Fidelity DNA Polymerase	Phusion™ U Hot start DNA Polymerase	Platinum™ <i>Taq</i> DNA Polymerase, High Fidelity	AccuPrime™ <i>Taq</i> DNA Polymerase, High Fidelity
Brand	Invitrogen™	Invitrogen™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Invitrogen™	Invitrogen™
Applications	End-point PCR	End-point PCR	End-point PCR	End-point PCR	End-point PCR	End-point PCR	End-point PCR	End-point PCR
Hot start	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
TaqMan™ probe compatible	No	No	No	No	No	No	No	No
Proofreading activity (3’-5’ exonuclease)	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Lyo-ready/ glycerol free	On request	Yes	On request	Yes	On request	On request	On request	On request
Inhibitor tolerance	●●●	●●●	●●●	●●	●●	●●	●	●
Extention rate	●●●	●●●	●●●	●●●	●●●	●●●	●●	●●
Benchtop stability	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
GC-rich target amplification	High	High	High	High	High	High	Standard	Standard
Fidelity vs Taq Pol	High (300x)	High (300x)	Medium (100x)	Medium (52x)	Medium (52x)	Low (25x)	Low (6x)	Low (9x)

PCR DNA polymerases

What is your application?

Relative fidelity compared to Standard *Taq* Polymerase

	Enzymes		
	Platinum™ <i>Taq</i> DNA Polymerase, High Fidelity	Phusion™ U Hot start DNA Polymerase	AccuPrime™ <i>Taq</i> DNA Polymerase, High Fidelity
	Invitrogen™	Thermo Scientific™	Invitrogen™
	Yes	Yes	Yes
	No	No	No
	Yes	Yes	Yes
	On request	On request	On request
	•	••	•
	Standard > 30-60 sec/kb	Fast < 10-30 sec/kb	Standard > 30-60 sec/kb
	Yes	Yes	Yes
	Standard < 65% GC	High > 65% GC	Standard < 65% GC
	Low (6x)	Low (25x)	Low (9x)

PCR DNA polymerases

What is your application?

Relative fidelity compared to Standard *Taq* Polymerase

	Enzymes		
	Phusion™ Plus DNA Polymerase	Phusion™ High-Fidelity DNA Polymerase	Phusion™ Hot start II DNA Polymerase
Brand	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Hot start	Yes	No	Yes
TaqMan™ probe compatible	No	No	No
Proofreading activity (3'-5' exonuclease)	Yes	Yes	Yes
Lyo-ready/glycerol free	On request	On request	Yes
Inhibitor tolerance	●●●	●●	●●
Extention rate	Fast < 10-30 sec/kb	Fast < 10-30 sec/kb	Fast < 10-30 sec/kb
Benchtop stability	Yes	No	Yes
GC-rich target amplification	High > 65% GC	High > 65% GC	High > 65% GC
Fidelity vs Taq Pol	Medium (100x)	Medium (52x)	Medium (52x)

PCR DNA polymerases

What is your application?

Relative fidelity compared to Standard *Taq* Polymerase

	Enzymes	
	Platinum™ SuperFi II DNA Polymerase	Platinum™ SuperFi™ DNA Polymerase
Brand	Invitrogen™	Invitrogen™
Hot start	Yes	Yes
TaqMan™ probe compatible	No	No
Proofreading activity (3'-5' exonuclease)	Yes	Yes
Lyo-ready/glycerol free	On request	Yes
Inhibitor tolerance	•••	•••
Extention rate	Fast < 10-30 sec/kb	Fast < 10-30 sec/kb
Benchtop stability	Yes	Yes
GC-rich target amplification	High > 65% GC	High > 65% GC
Fidelity vs Taq Pol	High (300x)	High (300x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

	Enzymes				
	DreamTaq DNA Polymerase	AmpliTaq™ DNA Polymerase	AmpliTaq™ 360 DNA Polymerase	Taq DNA Polymerase, recombinant	Phusion™ High-Fidelity DNA Polymerase
	Thermo Scientific™	Applied Biosystems™	Applied Biosystems™	Thermo Scientific™	Thermo Scientific™
	TaqMan™ probe compatible	Yes	Yes	Yes	No
	Proofreading activity (3'-5' exonuclease)	No	No	No	Yes
	Lyo-ready/glycerol free	On request	On request	On request	On request
	Inhibitor tolerance	•	•	•	••
	Extention rate	••	••	••	•••
	Benchtop stability	No	No	No	No
GC-rich target amplification	Standard	Standard	High	Standard	High
Fidelity vs Taq Pol	Standard (1x)	Standard (1x)	Standard (1x)	Standard (1x)	Medium (52x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

Working with 65+% GC-rich targets?

Brand	Enzymes						
	Platinum™ SuperFi II DNA Polymerase	Platinum™ SuperFi™ DNA Polymerase	Phusion™ Plus DNA Polymerase	Phusion™ Hot start II DNA Polymerase	Phusion™ U Hot start DNA Polymerase	Platinum™ <i>Taq</i> DNA Polymerase, High Fidelity	AccuPrime™ <i>Taq</i> DNA Polymerase, High Fidelity
	Invitrogen™	Invitrogen™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Invitrogen™	Invitrogen™
	TaqMan™ probe compatible	No	No	No	No	No	No
	Proofreading activity (3'-5' exonuclease)	Yes	Yes	Yes	Yes	Yes	Yes
	Lyo-ready/glycerol free	On request	Yes	On request	Yes	On request	On request
	Inhibitor tolerance	●●●	●●●	●●●	●●	●●	●
	Extention rate	●●●	●●●	●●●	●●●	●●●	●●
	Benchtop stability	Yes	Yes	Yes	Yes	Yes	Yes
	GC-rich target amplification	High	High	High	High	High	Standard
	Fidelity vs Taq Pol	High (300x)	High (300x)	Medium (100x)	Medium (52x)	Low (25x)	Low (6x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

Working with 65+% GC-rich targets?

	Enzymes				
	Platinum™ SuperFi II DNA Polymerase	Platinum™ SuperFi™ DNA Polymerase	Phusion™ Plus DNA Polymerase	Phusion™ Hot start II DNA Polymerase	Phusion™ U Hot start DNA Polymerase
	Invitrogen™	Invitrogen™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
	TaqMan™ probe compatible	No	No	No	No
	Proofreading activity (3'-5' exonuclease)	Yes	Yes	Yes	Yes
	Lyo-ready/glycerol free	On request	Yes	On request	Yes
	Inhibitor tolerance	●●●	●●●	●●	●●
	Extention rate	●●●	●●●	●●●	●●●
	Benchtop stability	Yes	Yes	Yes	Yes
	GC-rich target amplification	High	High	High	High
	Fidelity vs Taq Pol	High (300x)	High (300x)	Medium (100x)	Medium (52x)
					Low (25x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

Working with 65+% GC-rich targets?

	Enzymes	
	Platinum™ <i>Taq</i> DNA Polymerase, High Fidelity	AccuPrime™ <i>Taq</i> DNA Polymerase, High fidelity
	Invitrogen™	Invitrogen™
	No	No
	Yes	Yes
	On request	On request
	•	•
	••	••
	Yes	Yes
	Standard	Standard
	Low (6x)	Low (9x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

Working with 65+% GC-rich targets?

	Enzymes						
	Platinum™ II <i>Taq</i> Hot-Start DNA polymerase	Platinum™ <i>Taq</i> DNA Polymerase, DNA-free	Platinum™ <i>Taq</i> DNA Polymerase	Ampli <i>Taq</i> Gold™ 360 DNA Polymerase	Ampli <i>Taq</i> Gold™ DNA Polymerase	Dream <i>Taq</i> ™ Hot start DNA Polymerase	Phire Hot start II DNA Polymerase
	Invitrogen™	Invitrogen™	Invitrogen™	Applied Biosystems™	Applied Biosystems™	Thermo Scientific™	Invitrogen™
	TaqMan™ probe compatible	Yes	Yes	Yes	Yes	Yes	No
	Proofreading activity (3'-5' exonuclease)	No	No	No	No	No	Yes
	Lyo-ready/glycerol free	Yes	On request	Yes	On request	Yes	On request
	Inhibitor tolerance	***	*	*	*	*	**
	Extention rate	***	**	**	**	**	***
	Benchtop stability	Yes	Yes	Yes	Yes	Yes	Yes
	GC-rich target amplification	High	Standard	High	Standard	Standard	Standard
	Fidelity vs Taq Pol	Standard (1x)	Standard (1x)	Standard (1x)	Standard (1x)	Standard (1x)	Low (2x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

Working with 65+% GC-rich targets?

	Enzymes		
	Platinum™ II <i>Taq</i> Hot-Start DNA polymerase	Platinum™ <i>Taq</i> DNA Polymerase	Ampli <i>Taq</i> Gold™ 360 DNA Polymerase
	Invitrogen™	Invitrogen™	Applied Biosystems™
	TaqMan™ probe compatible	Yes	Yes
	Proofreading activity (3'-5' exonuclease)	No	No
	Lyo-ready/glycerol free	Yes	On request
	Inhibitor tolerance	***	*
	Extention rate	***	**
	Benchtop stability	Yes	Yes
	GC-rich target amplification	High	High
	Fidelity vs Taq Pol	Standard (1x)	Standard (1x)

PCR DNA polymerases

What is your application?

Need Hot-Start PCR?

Need high fidelity PCR?

Working with 65+% GC-rich targets?

	Enzymes			
	Platinum™ <i>Taq</i> DNA Polymerase, DNA-free	Ampli <i>Taq</i> Gold™ DNA Polymerase	Dream <i>Taq</i> ™ Hot start DNA Polymerase	Phire Hot start II DNA Polymerase
	Invitrogen™	Applied Biosystems™	Thermo Scientific™	Invitrogen™
	TaqMan™ probe compatible	Yes	Yes	No
	Proofreading activity (3'-5' exonuclease)	No	No	Yes
	Lyo-ready/glycerol free	On request	On request	Yes
	Inhibitor tolerance	*	*	**
	Extention rate	**	**	***
	Benchtop stability	Yes	Yes	Yes
	GC-rich target amplification	Standard	Standard	Standard
	Fidelity vs Taq Pol	Standard (1x)	Standard (1x)	Low (2x)

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

	Enzymes	
	M-MLV Reverse Transcriptase	RevertAid Reverse Transcriptase
Brands	Invitrogen™	Thermo Scientific™
Lyo-ready/glycerol free	Yes	Yes
Sensitivity	Low	Low
Inhibitor resistance	Low	Low
Reaction speed	Low	Low
Optimal reaction temp	37°–42°C	37°–42°C
Trancript length	•	•
GC-rich target amplification	Standard	Standard
RNase H activity	Yes	Yes
Template switch	No	No

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

	Enzymes	
	SuperScript™ III Reverse Transcriptase	SuperScript™ II Reverse Transcriptase
Brands	Invitrogen™	Invitrogen™
Lyo-ready/glycerol free	Yes	On request
Sensitivity	Medium	Medium
Inhibitor resistance	Medium	Medium
Reaction speed	Medium	Medium
Optimal reaction temp	50°C	50°C
Trancript length	●●	●●
GC-rich target amplification	High	High
RNase H activity	No	No
Template switch	No	Yes

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

	Enzymes
	Lyo-ready Super-Script™ III Flash Reverse Transcriptase
Brands	Invitrogen™
Lyo-ready/glycerol free	Yes
Sensitivity	High
Inhibitor resistance	High
Reaction speed	High
Optimal reaction temp	60°–70°C
Trancript length	••
GC-rich target amplification	High
RNase H activity	No
Template switch	No

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

	Enzymes		
	M-MLV Reverse Transcriptase	RevertAid Reverse Transcriptase	RevertAid H Minus Reverse Transcriptase
Brands	Invitrogen™	Thermo Scientific™	Thermo Scientific™
Lyo-ready/glycerol free	Yes	Yes	Yes
Sensitivity	Low	Low	Low
Inhibitor resistance	Low	Low	Low
Reaction speed	Low	Low	Low
Optimal reaction temp	37°–42°C	37°–42°C	37°–42°C
Trancript length	•	•	•
GC-rich target amplification	Standard	Standard	Standard
RNase H activity	Yes	Yes	No
Template switch	No	No	No

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

	Enzymes				
	SuperScript™ III Reverse Transcriptase	SuperScript™ II Reverse Transcriptase	Maxima Reverse Transcriptase	Maxima H Minus Reverse Transcriptase	
	Brands	Invitrogen™	Invitrogen™	Thermo Scientific™	Thermo Scientific™
	Lyo-ready/glycerol free	Yes	On request	Yes	Yes
	Sensitivity	Medium	Medium	Medium	Medium
	Inhibitor resistance	Medium	Medium	Medium	Medium
	Reaction speed	Medium	Medium	Medium	Medium
	Optimal reaction temp	50°C	50°C	50°C	50°C
	Trancript length	●●	●●	●●●	●●●
	GC-rich target amplification	High	High	Standard	Standard
RNase H activity	No	No	Yes	No	
Template switch	No	Yes	No	Yes	

Reverse transcriptases

What is your RT-PCR / RT-qPCR application?

Inhibitor tolerance?

	Enzymes
	SuperScript™ IV Reverse Transcriptase
Brands	Invitrogen™
Lyo-ready/glycerol free	Yes
Sensitivity	High
Inhibitor resistance	High
Reaction speed	High
Optimal reaction temp	50°C
Trancript length	●●●
GC-rich target amplification	High
RNase H activity	No
Template switch	Yes

Modifying enzymes

What is your analyte?

Modifying enzymes

What is your analyte?

Choose your enzyme.

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	Poly(A) Polymerase, Yeast	IVT, Cloning, Labeling
Thermo Scientific™	T7 RNA Polymerase	IVT
Thermo Scientific™	Terminal Deoxynucleotidyl Transferease	Cloning, NGS library preparation
Thermo Scientific™	SP6 RNA Polymerase	IVT, NGS library preparation
Thermo Scientific™	T3 RNA Polymerase	IVT, NGS library preparation

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	RNase A/T1 Mix	NA preparation
Thermo Scientific™	RNase A, DNase and protease-free	NA preparation, NGS library preparation
Thermo Scientific™	RNase I	NA preparation
Thermo Scientific™	Ribonuclease H (RNase H)	PCR, Cloning, NGS library preparation
Thermo Scientific™	RNase T1	NA preparation, NGS library preparation

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	T4 Polynucleotide Kinase	Cloning, NGS library preparation, Labeling
Thermo Scientific™	T4 RNA Ligase	Cloning, NGS library preparation
Thermo Scientific™	S1 Nuclease	Cloning, NGS library preparation
Thermo Scientific™	FastAP Thermosensitive Alkaline Phosphatase	Cloning, PCR

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	RiboLock RNase Inhibitor	PCR, isothermal amplification, IVT
Invitrogen™	RNaseOUT™ Recombinant Ribonuclease Inhibitor	PCR, isothermal amplification, IVT
Applied Biosystem™	RNase Inhibitor	PCR, isothermal amplification, IVT
Invitrogen™	Suprase-In™ RNase Inhibitor	PCR, isothermal amplification, IVT

Modifying enzymes

What is your analyte?

Choose your enzyme.

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	DNA Polymerase I	Cloning, NGS library preparation, Labeling
Thermo Scientific™	Klenow Fragment	Cloning, NGS library preparation, Labeling
Thermo Scientific™	Klenow Fragment, exo-	NGS library preparation, Labeling
Thermo Scientific™	T4 DNA Polymerase	IVT, Cloning
Thermo Scientific™	T7 DNA Polymerase	Cloning, NGS library preparation, Labeling
Thermo Scientific™	Terminal Deoxynucleotidyl Transferease	Cloning, NGS library preparation
Invitrogen™	Lyo-ready Bst DNA Polymerase	DNA preamplification
Thermo Scientific™	phi29 DNA Polymerase	DNA preamplification
Thermo Scientific™	EquiPhi29™ DNA Polymerase	DNA preamplification

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	Endonuclease V, T.maritima	Cloning, NGS library preparation
Thermo Scientific™	Endonuclease IV, E.coli	PCR, NGS library preparation
Thermo Scientific™	Exonuclease I	PCR, NGS library preparation
Thermo Scientific™	Exonuclease VII	PCR, NGS library preparation
Thermo Scientific™	Exonuclease III	PCR, NGS library preparation
Thermo Scientific™	T7 Gene 6 Exonuclease	PCR, Cloning
Thermo Scientific™	Lambda Exonuclease	PCR, Cloning
Thermo Scientific™	S1 Nuclease	Cloning, NGS library preparation
Thermo Scientific™	DNase I, RNase-free	PCR, IVT, NGS library preparation
Thermo Scientific™	Uracil-DNA Glycosylase, heat-labile	PCR, NGS library preparation
Thermo Scientific™	Uracil DNA Glycosylase	PCR, NGS library preparation
Thermo Scientific™	Micrococcal Nuclease	Cloning, PCR

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Application
Thermo Scientific™	FastDigest Restriction Enzymes	Cloning, PCR, IVT
Thermo Scientific™	Conventional Restriction Enzymes	Cloning, PCR, IVT
Invitrogen™	Anza Restriction Enzymes	Cloning, PCR, IVT

Modifying enzymes

What is your analyte?

Choose your enzyme.

Brand	Enzyme	Enzyme Categories	Application
Thermo Scientific™	T4 Polynucleotide Kinase	Phosphatases & Kinases	Cloning, NGS library preparation, Labeling
Thermo Scientific™	Pyrophosphatase, inorganic	Phosphatases & Kinases	IVT, NGS library preparation
Thermo Scientific™	T4 beta-glucosyltransferase	Transferase	PCR, NGS library preparation, Methylation analysis
Thermo Scientific™	T4 DNA Ligase	Ligase	Cloning, NGS library preparation
Thermo Scientific™	CpG Methyltransferase (M.SssI)	Transferase	PCR
Thermo Scientific™	FastAP Thermosensitive Alkaline Phosphatase	Phosphatases & Kinases	Cloning, PCR
Thermo Scientific™	Single-Strand Binding Protein (SSB)	Binding Proteins	PCR, NGS library preparation
Invitrogen™	<i>E. coli</i> DNA Ligase	Ligase	Cloning, NGS library preparation
Thermo Scientific™	T4 DNA Ligase	Ligase	Cloning, NGS library preparation
Thermo Scientific™	<i>Thermus thermophilus</i> DNA Ligase	Ligase	NGS, DNA Repair

Modifying enzymes

What is your analyte?

Brand	Enzyme	Enzyme Categories	Application
Thermo Scientific™	T4 Polynucleotide Kinase	Phosphatases & Kinases	Cloning, NGS library preparation, Labeling
Thermo Scientific™	Terminal Deoxynucleotidyl Transferease	Polymerases	Cloning, NGS library preparation
Thermo Scientific™	S1 Nuclease	DNA Repair Enzymes, Exo- & Endonucleases	Cloning, NGS library preparation
Thermo Scientific™	FastAP Thermosensitive Alkaline Phosphatase	Phosphatases & Kinases	Cloning, PCR

Isothermal amplification

What is your analyte?

Isothermal amplification

What is your analyte?

What is your application?

Isothermal amplification

What is your analyte?

What is your application?

Enzymes						
	SuperScript™ IV RT-LAMP Master Mix	RiboLock RNase Inhibitor	RnaseOUT™ Recombinant Ribonuclease Inhibitor	RNase Inhibitor	SUPERase-In™ RNase Inhibitor	Lyo-ready Bst DNA Polymerase
Brand	Invitrogen™	Thermo Scientific™	Invitrogen™	Applied Biosystems™	Invitrogen™	Invitrogen™
Format	Master mix	Standalone enzyme	Standalone enzyme	Standalone enzyme	Standalone enzyme	Standalone enzyme, kit
Reaction temperature	65°C	<60°C	<65°C	<60°C	<60°C	37°–65°C
Lyo-ready/glycerol free	On request	Yes	Yes	On request	No	Yes

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes		
	RiboLock RNase Inhibitor	RNaseOUT™ Recombinant Ribonuclease Inhibitor	SUPERase-In™ RNase Inhibitor
	Thermo Scientific™	Invitrogen™	Invitrogen™
	Standalone enzyme	Standalone enzyme	Standalone enzyme
	<60°C	<65°C	<60°C
Brand			
Format			
Reaction temperature			
Lyo-ready/glycerol free	Yes	Yes	No

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes	
	RNase Inhibitor	T7 RNA Polymerase
	Applied Biosystems™	Thermo Scientific™
	Standalone enzyme	Standalone enzyme
	<60°C	37°C
Brand		
Format		
Reaction temperature		
Lyo-ready/glycerol free	On request	High concentration

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes		
	RNase Inhibitor	T7 RNA Polymerase	Ribonuclease H (RNase H)
	Applied Biosystems™	Thermo Scientific™	Thermo Scientific™
	Standalone enzyme	Standalone enzyme	Standalone enzyme
	<60°C	37°C	37°C
Brand			
Format			
Reaction temperature			
Lyo-ready/glycerol free	On request	High concentration	High concentration

Isothermal amplification

What is your analyte?

What is your application?

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes			
	Lyo-ready Bst DNA Polymerase	SuperScript™ IV RT-LAMP Master Mix	Bsm DNA Polymerase, large fragment	
	Brand	Invitrogen™	Invitrogen™	Thermo Scientific™
	Format	Standalone enzyme, kit	Master mix	Standalone enzyme
	Reaction temperature	37–65°C	65°C	37°–60°C
Lyo-ready/glycerol free	Yes	On request	Yes	

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes			
	phi29 DNA Polymerase	EquiPhi29™ DNA Polymerase	EquiPhi29™ DNA Amplification Kit	Lyo-ready Bst DNA Polymerase
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Invitrogen™
	Standalone enzyme	Standalone enzyme	Kit	Standalone enzyme
	37°C	37°–45°C	42°C	37°–65°C
Lyo-ready/glycerol free	High concentration	Yes	On request	Yes

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes	
	phi29 DNA Polymerase	EquiPhi29™ DNA Polymerase
Brand	Thermo Scientific™	Thermo Scientific™
Format	Standalone enzyme	Standalone enzyme
Reaction temperature	37°C	37°–45°C
Lyo-ready/glycerol free	High concentration	Yes

Isothermal amplification

What is your analyte?

What is your application?

		Enzymes				
		Lyo-ready Bst DNA Polymerase	T4 UvsX	T4 Gene 32 Protein	T4 UvsY	Lyo-ready RPA Kit
Brand		Invitrogen™	Invitrogen™	Invitrogen™	Invitrogen™	Invitrogen™
Format		Standalone enzyme	Standalone enzyme	Standalone enzyme	Standalone enzyme	Kit
Reaction temperature		37°–65°C	42°C	42°C	42°C	42°C
Lyo-ready/glycerol free		Yes	Yes	Yes	Yes	Yes

Isothermal amplification

What is your analyte?

What is your application?

 HDA

Enzymes	
Lyo-ready Bst DNA Polymerase	
Brand	Invitrogen™
Format	Standalone enzyme
Reaction temperature	37°–65°C
Lyo-ready/glycerol free	Yes

Isothermal amplification

What is your analyte?

What is your application?

	Enzymes
	Klenow Fragment, exo–
Brand	Thermo Scientific™
Format	Standalone enzyme
Reaction temperature	37°C
Lyo-ready/glycerol free	High concentration

Isothermal amplification

What is your analyte?

What is your application?

Enzymes	
Lyo-ready Bst DNA Polymerase	
Brand	Invitrogen™
Format	Standalone enzyme
Reaction temperature	37°–65°C
Lyo-ready/glycerol free	Yes

NGS library prep

What is your analyte?

NGS library prep

What is your analyte?

Choose your sample type.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes
	T4 Polynucleotide Kinase
Brand	Thermo Scientific™
Description	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes
	Ribonuclease H (RNase H)
Brand	Thermo Scientific™
Description	Ribonuclease H (RNase H) specifically degrades the RNA strand in RNA-DNA hybrids.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes						
	SuperScript™ IV Reverse Transcriptase	RiboLock RNase Inhibitor	DNA Polymerase I	Ribonuclease H (RNase H)	dNTP Mix	Random Hexamer Primer	Oligo(dT) ₁₈ Primer
	Brand	Invitrogen™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	SuperScript IV Reverse Transcriptase (RT) is a proprietary MMLV mutant with superior robustness and reliability in RT reactions. SuperScript IV RT is designed to provide reliable, consistent, and fast cDNA synthesis in the presence of inhibitors found in a wide variety of samples that cause other currently available RTs to perform inefficiently.	RiboLock RNase Inhibitor inhibits the activity of RNases A,B and C by binding them in a noncompetitive mode. The enzyme is used to prevent RNA from degradation by RNases.	DNA Polymerase I, a template-dependent DNA polymerase, catalyzes 5'→3' synthesis of DNA. The enzyme also exhibits 3'→5' exonuclease (proofreading) activity, 5'→3' exonuclease activity, and ribonuclease H activity. The enzyme is used for second cDNA strand synthesis.	Ribonuclease H (RNase H) specifically degrades the RNA strand in RNA-DNA hybrids.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	Random Hexamer Primer is a mixture of single-stranded random hexanucleotides with 5'- and 3'-hydroxyl ends. Random Hexamer primers are used for cDNA synthesis.	Oligo(dT) ₁₈ Primer is a synthetic single-stranded 18-mer oligonucleotide with 5'- and 3'-hydroxyl ends. Oligo(dT) primers are used for cDNA synthesis starting from the poly(A) tails of mRNAs.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes			
	T4 Polynucleotide Kinase	T4 DNA Polymerase	Klenow Fragment	Klenow Fragment, exo-
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.	T4 DNA Polymerase is a template-dependent DNA polymerase with strong 3'-5' exonuclease activity. T4 DNA Polymerase is used for blunting of DNA ends: removal of 3'-overhangs.	Klenow Fragment is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity and 3'→5' exonuclease activity, but lacks 5'→3' exonuclease activity of DNA polymerase I. The enzyme is used for DNA blunting by fill-in 5'-overhangs.	Klenow Fragment, exo- , is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity, but lacks the 3'→5' and 5'→3' exonuclease activities of DNA Polymerase I. The enzyme is used for filling 5'-overhangs of dsDNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes	
	T4 DNA Ligase	T4 RNA Ligase
	Thermo Scientific™	Thermo Scientific™
Description	T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA. The enzyme is used for addition of adapters to dsDNA.	T4 RNA Ligase catalyzes the ATP-dependent intra- and intermolecular formation of phosphodiester bonds between 5'-phosphate and 3'-hydroxyl termini of oligonucleotides, single-stranded RNA and DNA. The enzyme is used for addition of adapters to RNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes			
	Phusion™ Plus DNA Polymerase	dNTP Mix	Exonuclease I	Exonuclease VII
	Brand	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	Phusion Plus DNA Polymerase is a hot-start, high-fidelity DNA polymerase that brings together protein fusion technology and universal primer annealing. The enzyme is used for NGS library amplification. The enzyme should be used with GC enhancer.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	Exonuclease I (ExoI) degrades single-stranded DNA in a 3'→5' direction. The enzyme is used for primer removal after NGS library amplification.	Exonuclease VII is a strict single-strand directed enzyme with 5'→3' and 3'→5' exonuclease activities, making it the only bi-directional exonuclease with single-stranded specificity. The enzyme is used for oligonucleotides/primers removal after NGS library amplification.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes
	T4 Polynucleotide Kinase
	Thermo Scientific™
Brand	
Description	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes
	Ribonuclease H (RNase H)
Brand	Thermo Scientific™
Description	Ribonuclease H (RNase H) specifically degrades the RNA strand in RNA-DNA hybrids.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

		Enzymes					
Brand	SuperScript™ IV Reverse Transcriptase	RiboLock RNase Inhibitor	DNA Polymerase I	Ribonuclease H (RNase H)	dNTP Mix	Random Hexamer Primer	Oligo(dT) ₁₈ Primer
	Invitrogen™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	SuperScript IV Reverse Transcriptase (RT) is a proprietary MMLV mutant with superior robustness and reliability in RT reactions. SuperScript IV RT is designed to provide reliable, consistent, and fast cDNA synthesis in the presence of inhibitors found in a wide variety of samples that cause other currently available RTs to perform inefficiently.	RiboLock RNase Inhibitor inhibits the activity of RNases A,B and C by binding them in a noncompetitive mode. The enzyme is used to prevent RNA from degradation by RNases.	DNA Polymerase I, a template-dependent DNA polymerase, catalyzes 5'→3' synthesis of DNA. The enzyme also exhibits 3'→5' exonuclease (proofreading) activity, 5'→3' exonuclease activity, and ribonuclease H activity. The enzyme is used for second cDNA strand synthesis.	Ribonuclease H (RNase H) specifically degrades the RNA strand in RNA-DNA hybrids.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	Random Hexamer Primer is a mixture of single-stranded random hexanucleotides with 5'- and 3'-hydroxyl ends. Random Hexamer primers are used for cDNA synthesis.	Oligo(dT)18 Primer is a synthetic single-stranded 18-mer oligonucleotide with 5'- and 3'-hydroxyl ends. Oligo(dT) primers are used for cDNA synthesis starting from the poly(A) tails of mRNAs.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

		Enzymes				
Brand	DNase I, RNase-free	Single-Stranded DNA Binding Protein (SSB)	T4 Polynucleotide Kinase	T4 DNA Polymerase	Klenow Fragment	Klenow Fragment, exo-
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	DNase I (RNase-free) is an endonuclease that digests single- and double-stranded DNA.	Single-Stranded DNA Binding Protein (SSB) binds with high affinity in a cooperative manner to single-stranded DNA. After binding single-stranded DNA, SSB destabilizes helical duplexes.	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.	T4 DNA Polymerase is a template-dependent DNA polymerase with strong 3'-5' exonuclease activity. T4 DNA Polymerase is used for blunting of DNA ends: removal of 3'-overhangs.	Klenow Fragment is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity and 3'→5' exonuclease activity, but lacks 5'→3' exonuclease activity of DNA polymerase I. The enzyme is used for DNA blunting by fill-in 5'-overhangs.	Klenow Fragment, exo-, is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity, but lacks the 3'→5' and 5'→3' exonuclease activities of DNA Polymerase I. The enzyme is used for filling 5'-overhangs of dsDNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes	
	T4 DNA Ligase	T4 RNA Ligase
	Thermo Scientific™	Thermo Scientific™
Description	T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA. The enzyme is used for addition of adapters to dsDNA.	T4 RNA Ligase catalyzes the ATP-dependent intra- and intermolecular formation of phosphodiester bonds between 5'-phosphate and 3'-hydroxyl termini of oligonucleotides, single-stranded RNA and DNA. The enzyme is used for addition of adapters to RNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes			
	Phusion™ Plus DNA Polymerase	dNTP Mix	Exonuclease I	Exonuclease VII
	Brand	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	Phusion Plus DNA Polymerase is a hot-start, high-fidelity DNA polymerase that brings together protein fusion technology and universal primer annealing. The enzyme is used for NGS library amplification. The enzyme should be used with GC enhancer.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	Exonuclease I (ExoI) degrades single-stranded DNA in a 3'→5' direction. The enzyme is used for primer removal after NGS library amplification.	Exonuclease VII is a strict single-strand directed enzyme with 5'→3' and 3'→5' exonuclease activities, making it the only bi-directional exonuclease with single-stranded specificity. The enzyme is used for oligonucleotides/primers removal after NGS library amplification.

NGS library prep

What is your analyte?

Choose your sample type.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

		Enzymes				
Brand	Uracil-DNA Glycosylase	Endonuclease IV, E.coli	Endonuclease V, T.maritima	Bsm DNA Polymerase, large fragment	T4 Polynucleotide Kinase	T4 DNA Polymerase
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	Uracil-DNA Glycosylase (UDG, UNG) catalyzes the hydrolysis of the N-glycosylic bond between uracil and sugar, leaving an apyrimidinic site in uracil-containing single or double-stranded DNA. The enzyme is used for damaged DNA repair (removes uracil).	Endonuclease IV recognizes apurinic/apyrimidinic (AP) sites of dsDNA and cleaves the phosphodiester bond 5' to the lesion generating a hydroxyl group at the 3'-terminus. The enzyme is used to repair DNA ends.	Endonuclease V is a 3'-endonuclease involved in DNA repair, which initiates removal of deaminated bases from damaged DNA, including uracil, hypoxanthine, and xanthine.	Bsm DNA Polymerase, Large Fragment, is an equivalent to Bst DNA polymerase, which catalyzes 5'→3' synthesis of DNA and lacks 5'→3' and 3'→5' exonuclease activities. The enzyme is used for DNA end repair.	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.	T4 DNA Polymerase is a template-dependent DNA polymerase with strong 3'-5' exonuclease activity. T4 DNA Polymerase is used for blunting of DNA ends: removal of 3'-overhangs.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes	
	DNase I, RNase-free	Single-Stranded DNA Binding Protein (SSB)
	Thermo Scientific™	Thermo Scientific™
Brand		
Description	DNase I (RNase-free) is an endo-nuclease that digests single- and double-stranded DNA.	Single-Stranded DNA Binding Protein (SSB) binds with high affinity in a cooperative manner to single-stranded DNA. After binding single-stranded DNA, SSB destabilizes helical duplexes.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes			
	T4 Polynucleotide Kinase	T4 DNA Polymerase	Klenow Fragment	Klenow Fragment, exo-
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.	T4 DNA Polymerase is a template-dependent DNA polymerase with strong 3'-5' exonuclease activity. T4 DNA Polymerase is used for blunting of DNA ends: removal of 3'-overhangs.	Klenow Fragment is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity and 3'→5' exonuclease activity, but lacks 5'→3' exonuclease activity of DNA polymerase I. The enzyme is used for DNA blunting by fill-in 5'-overhangs.	Klenow Fragment, exo-, is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity, but lacks the 3'→5' and 5'→3' exonuclease activities of DNA Polymerase I. The enzyme is used for filling 5'-overhangs of dsDNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

Brand	Enzymes
	T4 DNA Ligase
	Thermo Scientific™
Description	T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA. The enzyme is used for addition of adapters to dsDNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes		
	Phusion™ Plus DNA Polymerase	Phusion U Hot start DNA Polymerase	dNTP Mix
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	Phusion Plus DNA Polymerase is a hot-start, high-fidelity DNA polymerase that brings together protein fusion technology and universal primer annealing. The enzyme is used for target enrichment.	Phusion U DNA polymerase can incorporate dUTP and read through uracil present in DNA templates. The enzyme is used to amplify DNA containing uracil.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes		
	Phusion™ Plus DNA Polymerase	Exonuclease I	Exonuclease VII
Brand	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	Phusion Plus DNA Polymerase is a hot-start, high-fidelity DNA polymerase that brings together protein fusion technology and universal primer annealing. The enzyme is used for target enrichment.	Exonuclease I (ExoI) degrades single-stranded DNA in a 3'→5' direction. The enzyme is used for primer removal after NGS library amplification.	Exonuclease VII is a strict single-strand directed enzyme with 5'→3' and 3'→5' exonuclease activities, making it the only bi-directional exonuclease with single-stranded specificity. The enzyme is used for oligonucleotides/primers removal after NGS library amplification.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes		
	DNase I, RNase-free	Single-Stranded DNA Binding Protein (SSB)	MuA Transposase
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	DNase I (RNase-free) is an endo-nuclease that digests single- and double-stranded DNA.	Single-Stranded DNA Binding Protein (SSB) binds with high affinity in a cooperative manner to single-stranded DNA. After binding single-stranded DNA, SSB destabilizes helical duplexes.	MuA Transposase catalyzes trans-position reaction in vitro.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes			
	T4 Polynucleotide Kinase	T4 DNA Polymerase	Klenow Fragment	Klenow Fragment, exo-
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the gamma-phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides, or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. The enzyme is used for 5'-phosphorylation and removal of 3'-phosphate groups.	T4 DNA Polymerase is a template-dependent DNA polymerase with strong 3'-5' exonuclease activity. T4 DNA Polymerase is used for blunting of DNA ends: removal of 3'-overhangs.	Klenow Fragment is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity and 3'→5' exonuclease activity, but lacks 5'→3' exonuclease activity of DNA polymerase I. The enzyme is used for DNA blunting by fill-in 5'-overhangs.	Klenow Fragment, exo-, is the large fragment of DNA polymerase I. It exhibits 5'→3' polymerase activity, but lacks the 3'→5' and 5'→3' exonuclease activities of DNA Polymerase I. The enzyme is used for filling 5'-overhangs of dsDNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes
	T4 DNA Ligase
	Thermo Scientific™
Brand	
Description	T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA. The enzyme is used for addition of adapters to dsDNA.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes	
	Phusion™ Plus DNA Polymerase	dNTP Mix
	Thermo Scientific™	Thermo Scientific™
Description	Phusion Plus DNA Polymerase is a hot-start, high-fidelity DNA polymerase that brings together protein fusion technology and universal primer annealing. The enzyme is used for target enrichment.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.

NGS library prep

What is your analyte?

Choose your sample type.

Which step?

	Enzymes			
	Phusion™ Plus DNA Polymerase	dNTP Mix	Exonuclease I	Exonuclease VII
	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Description	Phusion Plus DNA Polymerase is a hot-start, high-fidelity DNA polymerase that brings together protein fusion technology and universal primer annealing. The enzyme is used for target enrichment.	dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	Exonuclease I (ExoI) degrades single-stranded DNA in a 3'→5' direction. The enzyme is used for primer removal after NGS library amplification.	Exonuclease VII is a strict single-strand directed enzyme with 5'→3' and 3'→5' exonuclease activities, making it the only bi-directional exonuclease with single-stranded specificity. The enzyme is used for oligonucleotides/primers removal after NGS library amplification.

Nucleotides & reagents

What is your format?

Nucleotides & reagents

What is your format?

Reagents					
dNTP Mix (25 mM each)	dNTP Mix (10 mM each)	dNTP/dUTP Mix	dNTP Mix (10 mM ea)	GeneAmp™ dNTP Blend with dUTP (12.5 mM)	dNTP Mix (2.5 mM)
Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Invitrogen™	Applied Biosystems™	Invitrogen™
Nucleotide	Nucleotide	Nucleotide	Nucleotide	Nucleotide	Nucleotide
“dNTP Mixes contain dATP, dCTP, dGTP and dTTP, each at different final concentrations (from 2.5 mM to 25 mM). The Mix offers the possibility to reduce the number of pipetting steps and the risk of reaction set-up errors.”					
Thermo Scientific dNTP Mix (25 mM) contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 25 mM.	Thermo Scientific dNTP Mix (10 mM) contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	Thermo Scientific dNTP/dUTP Mix is a solution containing dATP, dCTP and dGTP, each at a final concentration of 2 mM and dUTP of 4 mM.	Thermo Scientific dNTP Mix (10 mM) contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.	GeneAmp dNTP Blend, 12.5 mM with dUTP contains 2.5 mM each of dATP, dCTP, dGTP and 5.0 mM of dUTP	Invitrogen 2.5 mM dNTP Mix consists of a solution of all four nucleotides, dATP, dCTP, dGTP, and dTTP, each at a concentration of 2.5 mM.

Nucleotides & reagents

What is your format?

	Reagents						
	dNTP Set, 100 mM Solutions	dNTP Set (100 mM)	dUTP Solution (100 mM)	dGTP Solution (100 mM)	dTTP Solution (100 mM)	dATP Solution (100 mM)	dCTP Solution(100mM)
Brand	Thermo Scientific™	Invitrogen™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™	Thermo Scientific™
Type	Nucleotide	Nucleotide	Nucleotide	Nucleotide	Nucleotide	Nucleotide	Nucleotide
Description	Nucleotides, molecular biology grade. highly purified dNTPs for direct use in enzymatic reactions. Thermo Fisher Scientific is one of the few primary manufacturers of nucleotides in the industry. All Thermo Scientific™ nucleotides are supplied in aqueous solutions titrated to pH 7.5 with NaOH.						
Applications	dNTPs can be used in PCR, qPCR, cDNA synthesis, high-fidelity and long PCR, isothermal amplification, next-generation sequencing, DNA labeling, and other molecular biology applications.						