

Phusion High-Fidelity DNA Polymerases

Thermo Scientific™ Phusion™ High-Fidelity DNA Polymerases offer very high fidelity, speed, and yield for all PCR applications.

General instructions

- Due to the unique nature of Phusion DNA polymerases, always use the T_m calculator on our website to determine optimal annealing temperature (thermofisher.com/tmcalculator).
- Use 98°C for denaturation.
- Use 15–30 sec/kb for extension. Do not exceed 1 min/kb.
- Use Phusion DNA polymerases at 0.5–1.0 U per 50 µL reaction volume. Do not exceed 2 U per 50 µL reaction volume.
- Use 200 µM of each dNTP.
- If uracil is present in the dNTP mix or DNA template, use Thermo Scientific™ Phusion™ U Hot Start DNA Polymerase.

Note: Phusion DNA polymerases produce blunt-end DNA products.

Choosing the right Phusion product

		Phusion High-Fidelity DNA Polymerase (Cat. No. F530S)	Phusion Hot Start II High-Fidelity DNA Polymerase (Cat. No. F549S)	Phusion Flash High-Fidelity DNA Polymerase (Cat. No. F548S)	Phusion U Hot Start DNA Polymerase (Cat. No. F555S)	Phusion U Multiplex PCR Master Mix (Cat. No. F562S)
Characteristics	Blunt or 3'-A end	Blunt	Blunt	Blunt	Blunt	Blunt
	Target length, genomic/phage DNA	≤16/20 kb	≤16/20 kb	≤16/20 kb	≤7.5/20 kb	≤2.5/2.5 kb
	Hot start	No	Yes	Yes	Yes	Yes
	Recommended extension time	15–30 sec/kb	15–30 sec/kb	15 sec/kb	15–30 sec/kb	15–30 sec/kb
	Fidelity vs. <i>Taq</i>	52x	52x	25x	25x	NA
	dUTP tolerance	No	No	No	Yes	Yes
Formats	Enzyme*	✓	✓		✓	
	Green buffer**	✓	✓			
	Master mix†	✓	✓	✓	✓	✓
	Complete kit‡	✓				

* DNA polymerase, buffer, DMSO, and MgCl₂.

** DNA polymerase supplied with Phusion Green Buffer, which includes a density reagent and two tracking dyes for direct loading on gel.

† 2X master mix format.

‡ All the necessary PCR components, including control template and primers.

Reaction setup

Component	50 μ L reaction	20 μ L reaction	Final concentration
5X Phusion buffer*	10 μ L	4 μ L	1X
10 mM dNTPs*	1 μ L	0.4 μ L	200 μ M each
Primer A	x μ L	x μ L	0.5 μ M
Primer B	y μ L	y μ L	0.5 μ M
Template DNA	z μ L	z μ L	–
DMSO (optional)	(1.5 μ L)	(0.6 μ L)	(3%)
Phusion DNA polymerase	0.5 μ L	0.2 μ L	0.02 U/ μ L
Water	To 50 μ L total	To 20 μ L total	–

* If you are using any of the Phusion PCR master mix products, add 25 or 10 μ L of the 2X master mix (depending on the final reaction volume). Do not add dNTPs.

Cycling instructions for Phusion and Phusion Hot Start II High-Fidelity DNA Polymerases

Cycle step	2-step protocol		3-step protocol		Cycles
	Temperature	Time	Temperature	Time	
Initial denaturation	98°C	30 sec	98°C	30 sec	1
Denaturation	98°C	5–10 sec	98°C	5–10 sec	25–35
Annealing*	–	–	X°C*	10–30 sec	
Extension	72°C	15–30 sec/kb	72°C	15–30 sec/kb	
Final extension	72°C	5–10 min	72°C	5–10 min	1
	4°C	Hold	4°C	Hold	

* Depends on the primer T_m values. Use the T_m calculator at thermofisher.com/tmcalculator

Cycling instructions for Phusion Flash High-Fidelity PCR Master Mix

Cycle step	2-step protocol		3-step protocol		Cycles
	Temperature	Time	Temperature	Time	
Initial denaturation	98°C	10 sec	98°C	10 sec	1
Denaturation*	98°C	0 or 1 sec	98°C	0 or 1 sec	30
Annealing**	–	–	50–72°C	5 sec	
Extension	72°C	15 sec/kb	72°C	15 sec/kb	
Final extension	72°C	1 min	72°C	1 min	1
	4°C	Hold	4°C	Hold	

* A very short denaturation step is recommended. If the PCR instrument used does not accept 0 sec as a value, then a 1 sec value can be programmed.

** Depends on the primer T_m values. Use the T_m calculator at thermofisher.com/tmcalculator

Find out more at thermofisher.com/phusion