



# Performance comparison of the updated and original formulations of Maxima H Minus Reverse Transcriptase



# First-strand cDNA synthesis using Maxima H Minus Reverse Transcriptase (RT)



## Ability to synthesize cDNAs of different lengths at a range of temperatures.

cDNA was synthesized using 1 µg/µL of Invitrogen™ RNA Millennium™ marker at different temperatures using both the updated and original formulations of Thermo Scientific™ Maxima™ H Minus Reverse Transcriptase. Reaction products were resolved on alkaline agarose gels.

### Original formulation—


Maxima H Minus Reverse Transcriptase formulated with Triton™ X-100 detergent

### Updated formulation—

Maxima H Minus Reverse Transcriptase reformulated with an alternative detergent

There is no difference in efficiency between the updated and original formulations of Maxima H Minus Reverse Transcriptase.

# Expression analysis with targets of different lengths from human RNA, using Maxima H Minus RT




**Versatility of Maxima H Minus Reverse Transcriptase across a broad range of target lengths.**

cDNA was synthesized from targets ranging from 0.5 to 9.4 kb, using 100 ng of human total RNA and the updated and original formulations of Maxima H Minus Reverse Transcriptase. Synthesized cDNA was PCR-amplified with Thermo Scientific™ Phusion™ Hot Start II High-Fidelity PCR Master Mix.

There is no difference in amplification across a broad range of targets of human RNA, between the updated and original formulations of Maxima H Minus Reverse Transcriptase.

# Expression analysis with targets of different lengths from plant RNA, using Maxima H Minus RT



**Original formulation—**  
Maxima H Minus Reverse Transcriptase  
formulated with Triton X-100 detergent


**Updated formulation—**  
Maxima H Minus Reverse Transcriptase  
reformulated with an alternative detergent

## Versatility of Maxima H Minus Reverse Transcriptase across a broad range of target lengths.

cDNA was synthesized from targets ranging from 0.382 to 6.6 kb, using 100 ng of plant RNA (**A**, wheat; **B**, tobacco) and the updated and original formulations of Maxima H Minus Reverse Transcriptase. Synthesized cDNA was PCR-amplified with Phusion Hot Start II High-Fidelity PCR Master Mix.

There is no difference in amplification across a broad range of targets of plant RNA between the updated and original formulations of Maxima H Minus Reverse Transcriptase.

# Detection sensitivity of Maxima H Minus Reverse Transcriptase with human RNA



**Original formulation (O)—**  
Maxima H Minus Reverse Transcriptase formulated with Triton X-100 detergent

**Updated formulation (U)—**  
Maxima H Minus Reverse Transcriptase reformulated with an alternative detergent


## Detection sensitivity with different amounts of input human RNA.

cDNA was synthesized from various targets (1 kb, 3 kb, 7.8 kb) with inputs of 10 pg to 1 µg of Invitrogen™ HeLa-S3 Total RNA, using the updated and original formulations of Maxima H Minus Reverse Transcriptase. Synthesized cDNA was PCR-amplified with Phusion Hot Start II High-Fidelity PCR Master Mix.

The updated formulation of Maxima H Minus Reverse Transcriptase retains detection sensitivity equivalent to that of the original formulation of Maxima H Minus Reverse Transcriptase.



# Maxima H Minus Reverse Transcriptase sensitivity testing with plant RNA



**Original formulation—**  
Maxima H Minus Reverse Transcriptase  
formulated with Triton X-100 detergent

**Updated formulation—**  
Maxima H Minus Reverse Transcriptase  
reformulated with an alternative detergent


## Detection sensitivity with different amounts of input plant RNA.

cDNA was synthesized from various targets (0.382–6.66 kb) with inputs of 10 pg to 1 µg of total plant RNA (**A**, wheat; **B**, tobacco), using the updated and original formulations of Maxima H Minus Reverse Transcriptase. Synthesized cDNA was PCR-amplified with Phusion Hot Start II High-Fidelity PCR Master Mix.

The updated formulation of Maxima H Minus Reverse Transcriptase retains detection sensitivity equivalent to that of the original formulation of Maxima H Minus Reverse Transcriptase.

# Performance of Maxima H Minus cDNA Synthesis Master Mix with different targets in RT-qPCR


## Applied Biosystems™ TaqMan® Assay



## Performance of Maxima H Minus cDNA Synthesis Master Mix across different TaqMan® Assay targets.

cDNA was synthesized from 100 ng of HeLa-S3 total RNA input and analyzed with 9 different gene-specific primer sets in TaqMan Assays, using Applied Biosystems™ TaqMan® Fast Advanced Master Mix. Higher normalized fold changes indicate higher cDNA yields.

## Applied Biosystems™ PowerTrack™ SYBR™ Green Assay




## Performance of Maxima H Minus cDNA Synthesis Master Mix across different PowerTrack SYBR Green assay targets.

cDNA was synthesized from 100 ng of HeLa-S3 total RNA input and analyzed with 5 different gene-specific primer sets in SYBR Green assays, using Applied Biosystems™ PowerTrack™ SYBR™ Green Master Mix. Higher normalized fold changes indicate higher cDNA yields.


The updated formulation of Maxima H Minus cDNA Synthesis Master Mix retains performance equivalent to that of the original formulation of Maxima H Minus cDNA Synthesis Master Mix.

# Dynamic range of RT-qPCR using Maxima H Minus cDNA Synthesis Master Mix

**A** TaqMan Assay (*UBE2B* gene)



**B** SYBR Green Assay (*GAPDH* gene)



## Dynamic range of Maxima H Minus cDNA Synthesis Master Mix.

cDNA was synthesized from gene targets **(A)** *UBE2B* and **(B)** *GAPDH* with 100 ng to 100 pg of HeLa-S3 total RNA, using the updated (green) and original (gray) formulations of Maxima H Minus cDNA Synthesis Master Mix. cDNA was analyzed with **(A)** TaqMan Fast Advanced Master Mix and **(B)** SYBR Green Master Mix.

*UBE2B*: E = 102.9%,  $R^2 = 0.998$  (updated formulation); E = 112.9%,  $R^2 = 0.995$  (original formulation).

*GAPDH*: E = 109.9%,  $R^2 = 0.994$  (updated formulation); E = 110.3%,  $R^2 = 0.994$  (original formulation).

The updated formulation of Maxima H Minus cDNA Synthesis Master Mix retains dynamic range equivalent to that of the original formulation of Maxima H Minus cDNA Synthesis Master Mix.



# Thank you

