



Thermo Scientific Sera-Mag Speedbeads

Carboxylate-Modified Magnetic Particles

Thermo Scientific Sera-Mag Speedbeads Carboxylate-Modified Magnetic Particles provide the speed, precision, and increased binding capacity required for sample preparation, proteomics, nucleic acid isolation and clinical immunoassay applications.



This SEM image shows the unique cauliflower-like surface of the Sera-Mag Speedbeads which dramatically increases the overall surface area available for binding. Sera-Mag SpeedBeads also have a second layer of magnetite within the particle, resulting in a 2x faster increase in speed in response to magnetic field.



- Fast reaction kinetics increases throughput and precision
- Low non-specific binding improves assay accuracy
- Unique cauliflower-like surface increases overall binding capacity
- Carboxylic groups on the surface permit easy covalent coupling using simple carbodiimide chemistry
- Salt-tolerance and slow settling rate provides excellent colloidal stability in the absence of a magnetic field
- Uniform diameter provides high surface area and excellent lot-to-lot reproducibility
- Surfactant-free particles require no washing
- Stability in buffer systems and detergents allows versatility in reagent and sample preparation

Thermo Scientific Sera-Mag SpeedBeads Carboxylate-Modified Particles are nominal 1 μ m magnetic particles of uniform size and feature a second layer of magnetite.

As a result, Sera-Mag SpeedBeads respond much faster to a magnetic field to separate quickly and completely from suspensions. It also ensures shorter assay times in clinical diagnostic tests as well as faster particle movement through viscous solutions.

Available in Sera-Mag SpeedBeads and original Sera-Mag versions, these particles provide exquisite sensitivity and low non-specific binding for greater accuracy.

These carboxylate-modified magnetic particles also have a distinct cauliflower-like surface that adds to the overall surface area available for binding and without any reduction in particle size.

Product Specifications


Thermo Scientific Sera-Mag SpeedBeads Carboxylate-Modified Magnetic Particles

Specifications

	Sera-Mag SpeedBeads	Sera-Mag
Particle Composition	Double layer of magnetite	Single layer of magnetite
Particle Density	~2.0 g/cm ³	~1.7 g/cm ³
Content		
Magnetite	~60%	~40%
Acid	High acid content with parking areas ranging from ~2 - 5 Å ² per carboxyl group	
Product Attributes		
Nominal diameter	1 µm	1 µm
Concentration	Supplied at approximately 5% solids (50 mg/mL)	
Fill volume	15 mL, 100 mL, and 1,000 mL bottles	
pH stability	Stable in detergents, biological buffer systems and lysis or elution buffer systems (pH 1 to 13), including high salt solutions and DMF or DMSO	
Additives	0.05% sodium azide	0.05% sodium azide
Package includes	Certificate of Analysis and Packaging Insert	
Storage and handling	Unless otherwise stated, refrigerate (2-8 °C) product when not in use but do not freeze. Store upright and keep bottle tightly sealed. Mix product with gentle inversion by hand, roller or vortex mixer.	

Product Specifications	Type	Parking Area (Å ² per group)
4515-2105-050250 (15 mL) ¹	Sera-Mag SpeedBead	2 - 5
4515-2105-050350 (100 mL) ¹	Sera-Mag SpeedBead	2 - 5
6515-2105-050250 (15 mL) ¹	Sera-Mag SpeedBead	2 - 5
6515-2105-050350 (100 mL) ¹	Sera-Mag SpeedBead	2 - 5
2415-2105-050250 (15 mL) ¹	Sera-Mag	2 - 5
2415-2105-050350 (100 mL) ¹	Sera-Mag	2 - 5
4415-2105-050250 (15 mL) ¹	Sera-Mag	2 - 5
4415-2105-050350 (100mL) ¹	Sera-Mag	2 - 5

¹The 6515 and 4415 particles have surfaces that are more hydrophobic while the 4515 and 2415 have surfaces that are more hydrophilic. Please ask about our sample evaluation packs. Some performance feedback has indicated that the 6515 and 4415 particles, which are manufactured with a lower level of EDAC than the 4515 and 2415 particles, work better with nucleic acid isolation.



Applications

Sera-Mag SpeedBeads Carboxylate-Modified Particles are useful as a solid phase support for many applications including sample preparation and nucleic acid isolation. Prior to downstream applications, samples can be mixed with our particles, whereupon biomolecules of interest are covalently attached to carboxyl groups on the particles. Isolation or purification of the biomolecules occurs through magnetic separation. The carboxyl groups on these particles are activated by the water soluble carbodiimide 1-ethyl-3-(3 dimethylaminopropyl)-carbodiimide (EDAC). Once activated, the groups react with the free amino groups of the adsorbed protein to form amide bonds. If exposure to EDAC is found to be potentially harmful to the protein, then we recommend using our two-step covalent coupling procedure which prevents this. To learn more about covalently coupling proteins and other ligands to Sera-Mag SpeedBeads and Sera-Mag particles, review TN-027 "Recommended Adsorption and Covalent Coupling Procedures". This and other relevant documents can be obtained at www.thermoscientific.com/particletechnology

© 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

PS2000.1_06/10

Clinical Diagnostics

46360 Fremont Blvd.
Fremont, CA
94538

USA: 1-800-232-3342
International: 1-510-979-5488

www.thermoscientific.com/particletechnology
info.microparticles@thermofisher.com

Thermo
SCIENTIFIC