



Reducing Antibiotic Exposure

B·R·A·H·M·S PCT (Procalcitonin) sensitive KRYPTOR
An Effective Tool to Aid in Antibiotic Stewardship and
Assessing the Risk of Bacterial Infection

B·R·A·H·M·S PCT (Procalcitonin)

Expanded Insight for Reducing Antibiotic Use

The Centers for Disease Control and Prevention (CDC) estimates that at least 30% of antibiotics prescribed in the United States are unnecessary.¹ The use of antibiotics is believed to be a primary cause of the spread of antibiotic-resistant bacteria.²

Thermo Scientific™ B·R·A·H·M·S PCT™ sensitive KRYPTOR™ provides clinicians in the emergency department, intensive care unit and hospital wards with a sensitive, specific STAT biomarker that provides newly expanded insight to aid clinicians in making decisions regarding antibiotic therapy, as well as differential diagnosis of bacterial infection and sepsis.



1/3 of Antibiotic Prescriptions are Unnecessary¹

The Surviving Sepsis Campaign International Guidelines for Management of Severe Sepsis and Septic Shock (2016) suggests that PCT measurements can be used to support shortening the duration of antimicrobial therapy in sepsis patients.³



B·R·A·H·M·S PCT: The Quality Standard

When using PCT assays to support clinical decisions, quality and experience count. Clinicians worldwide rely on B·R·A·H·M·S PCT quality to make confident patient care decisions due to:


- **Evidence:** 4,500+ publications in both the U.S. and Europe demonstrate the clinical utility of PCT.
- **Experience:** B·R·A·H·M·S PCT has been in clinical use since 1996.
- **Precision:** Demonstrated high concordance at clinically relevant cut-offs.
- **Standardization:** Standardized, comparable clinical cut-offs, independent of platform.
- **Adoption:** PCT is included in antibiotic stewardship guidelines issued by IDSA (Infectious Disease Society of America) and the Surviving Sepsis Campaign.^{3,4}

Insight for Initiating Antibiotics

As many as 75% of patients with acute respiratory-tract infections are treated with antibiotics, despite a mainly viral cause for these infections.⁵ B·R·A·H·M·S PCT aids clinicians in the emergency room or inpatient hospital settings in determining whether antibiotics are appropriate for patients with suspected or confirmed lower respiratory tract infections (LRTI), defined as community-acquired pneumonia (CAP), acute bronchitis, and acute exacerbation of chronic obstructive pulmonary disease (AECOPD).


B·R·A·H·M·S PCT has been shown to reduce antibiotic prescription rate and exposure duration in LRTI.

Duration of antibiotic exposure and antibiotic prescription rates were significantly reduced in the PCT group in comparison to the standard of care group for community-acquired pneumonia (CAP) (n=925), acute exacerbations of COPD (n=228) and bronchitis (n=151) in the ProHosp trial.⁵




Guiding Antibiotic Therapy Decisions for LRTI

PCT Cut-offs



Aiding Differential Diagnosis

As a sensitive and specific biomarker of the inflammatory response to bacterial infection⁶ B·R·A·H·M·S PCT aids clinicians in determining a critically ill patient's risk of progression to severe sepsis and septic shock.



7,8,9,10




The short turn around time of B·R·A·H·M·S PCT sensitive KRYPTOR assay supports clinicians' decisions.

Decisions regarding antibiotic therapy should NOT be based solely on B·R·A·H·M·S PCT concentrations, but only in conjunction with clinical signs and symptoms and other diagnostic evidence.

Insight for Safely Discontinuing Antibiotics

Paired with clinical assessment, B·R·A·H·M·S PCT also aids decisions about whether to discontinue antibiotic therapy for patients with LRTI, or with suspected or confirmed sepsis.



PCT-supported therapy has been shown to reduce inpatient antibiotic exposure by


↓ 35% for LRTI patients⁵ ↓ 23% for critically ill ICU patients¹¹

without negative effects for mortality or length of stay.^{5,11-21}

Aiding Assessment of Mortality Risk

Following ICU admission, evaluating serial B·R·A·H·M·S PCT measurements over consecutive days aids in assessing the host response to antibiotic therapy and the risk of all-cause mortality. When the infection is controlled, PCT will decline daily.²² If the PCT level has not declined, the patient and therapy should be reassessed.

A baseline PCT measurement greater than 2.0 ng/mL on Day 0 is an additional factor to consider when evaluating PCT measurements on subsequent days.²³



Decisions regarding antibiotic therapy should NOT be based solely on B·R·A·H·M·S PCT concentrations, but only in conjunction with clinical signs and symptoms and other diagnostic evidence.

Understanding PCT Kinetics


The utility of B·R·A·H·M·S PCT as a tool for assessing the risk of bacterial infection stems from its unique kinetics in response to severe bacterial infection.

PCT levels increase 3-6 hours after bacterial insult and return to normal as the infection is resolved.^{7,8,9}

- **Approximate half-life of 24 hours**
- **High specificity and sensitivity for bacterial infection**
- **Indicator for disease severity and treatment response**

PCT Kinetics

PCT values rise in relation to sepsis severity, providing clinicians with a valuable tool for assessing patients suspected of sepsis.



The U.S. Multicenter MOSES Study²³

The ability of B·R·A·H·M·S PCT to support mortality risk assessment over multiple days was demonstrated in a major multi-site U.S. study. The Procalcitonin Monitoring Sepsis Study (MOSES) included 858 adult patients with sepsis recruited across 13 investigational sites in the U.S. Key findings of the study included:

- The change of PCT over time aids in assessing the cumulative 28-day risk of all-cause mortality for patients diagnosed with severe sepsis or septic shock.
- A 2-fold increased risk of death is seen for patients showing a decrease in PCT less than or equal to 80 percent during the first four days following diagnosis of severe sepsis or septic shock compared to those who experienced a decrease in PCT greater than 80 percent. Mortality was the same for men and women.
- The initial PCT level (≤ 2.0 ng/mL or > 2.0 ng/mL) provided important additional information about the mortality risk when evaluating the patient's clinical course with PCT measurements on subsequent days.

Comparing PCT to Other Biomarkers

The sensitivity and specificity of B·R·A·H·M·S PCT to the host response to severe bacterial infection, together with its rapid rise after an infectious insult, offer clinical advantages that complement existing biomarkers.

Lactate

Lactate (lactic acid) is produced due to inadequate tissue perfusion, a defining parameter of late sepsis. However, lactate is not specific for bacterial infection. Clinical conditions including microcirculatory dysfunction, shunting, regional blood flow maldistribution, exaggerated aerobic or anaerobic glycolysis, hypovolemia or arterial hypotension can increase lactate levels.²⁴ In addition, lactate does not rise until late in the course of sepsis.²⁵

For patients evaluated in the ED for a suspected infection, the combination of lactate and B·R·A·H·M·S PCT measurements, together with clinical data and vital signs, provide complementary information for risk stratification.²⁶


C-Reactive Protein (CRP)

CRP secretion is triggered by cytokines (IL-6, IL-1, TNF- α) in response to acute or chronic inflammation associated with bacterial, viral, or fungal infection, and conditions such as obesity and tissue injury. It has no correlation to Sepsis-related Organ Failure Assessment (SOFA) score and its kinetics are slow, peaking 36 to 50 hours after causal challenge.^{27,28,29} In recent years, CRP has not been recommended because of its lack of specificity for systemic bacterial infection and its suppression when corticosteroids are used.¹⁰


Blood Cultures

Only 30–50% of patients with a clinical diagnosis of severe sepsis or septic shock have positive blood cultures.³⁰ One U.S. study found that more than 50% of patients diagnosed with sepsis had normal white blood cell count (WBC).³¹ Challenges with relying on blood cultures for assessing systemic infection include the delay for response, decreased sensitivity in patients already on antibiotics, and potential for false-positive results due to sample contamination.

The unique kinetics of PCT make it a valuable complement to other biomarkers of sepsis.³²



PCT's sensitivity and specificity as a marker of systemic bacterial infection make it a valuable complement to traditional biomarkers used in sepsis risk assessment.⁸



Important Considerations When Interpreting PCT Results

Increased B•R•A•H•M•S PCT levels may not always be related to systemic bacterial infection.^{9,33,34,35} They may also be associated with:

- Injuries including major trauma, burns and heat stroke

- Acute medical conditions such as biliary pancreatitis, chemical pneumonitis, viral hepatitis and/or decompensated severe cirrhosis (Child-Pugh Class 3), prolonged or severe cardiogenic shock, prolonged severe organ perfusion anomalies, and post-cardiac arrest

- Unusual infectious diseases including invasive fungal infections and acute plasmodium falciparum malaria

- Active medullary C-cell carcinoma, small cell lung carcinoma, and bronchial carcinoid

- Following interventions such as surgery with extracorporeal circulation, treatment with drugs stimulating release of pro-inflammatory cytokines or resulting in anaphylaxis, peritoneal or hemodialysis

- Neonates during the first three days of life. B•R•A•H•M•S PCT values should be interpreted using a specific nomogram during the first 72 hours following birth

B•R•A•H•M•S PCT results should be evaluated in context of all laboratory findings and the total clinical status of the patient. In cases where the laboratory results do not agree with the clinical picture or history, additional tests should be performed.



References

1. Fleming-Dutra KE, et al. Prevalence of inappropriate antibiotic prescriptions among us ambulatory care visits, 2010-2011, JAMA. 2016;315(17):1864-1873. doi:10.1001/jama.2016.4151.
2. Centers for disease control and prevention. Antibiotic resistance threats in the United States, 2013. <http://www.cdc.gov/drugresistance/threat-report-2013/>.
3. Surviving sepsis campaign: International guidelines for management of severe sepsis and septic shock: 2016.
4. Barlam, TF, et al. Implementing an antibiotic stewardship program: Guidelines by the Infectious Diseases Society of America and the Society for Healthcare Epidemiology of America. *Clinical infectious diseases* 62.10 (2016): e51-e77.
5. Schuetz P, et al., Effect of procalcitonin-based guidelines vs standard guidelines on antibiotic use in lower respiratory tract infections: The ProHOSP randomized controlled trial. *Jama* 2009; 302(10): 1059-1066.
6. Brunkhorst FM, et al. Kinetics of procalcitonin in iatrogenic sepsis. *Intensive Care Med* 1998;24(8):888-889.
7. Harbarth S, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med* 2001;164(3):396-402.
8. Müller B, et al. Calcitonin precursors are reliable markers of sepsis in a medical intensive care unit. *Critical Care Medicine* 2000;28(4):977-983.
9. Meisner M. Procalcitonin: Biochemistry and clinical diagnosis. UNI-MED Verlag AG; 2010.
10. Morgenthaler NG, et al. Detection of procalcitonin (PCT) in healthy controls and patients with local infection by a sensitive ILMA. *Clin Lab* 2002;48(5-6):263-270.
11. Bouadma L, et al., Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial): A multicentre randomised controlled trial. *Lancet* 2010; 375(9713): 463-474.
12. Briel M, et al. Procalcitonin-guided antibiotic use vs a standard approach for acute respiratory tract infections in primary care. *Arch Intern Med*. 2008; 168(18): 2000-7.
13. Burkhardt O, et al. Procalcitonin guidance and reduction of antibiotic use in acute respiratory tract infection. *Eur. Resp. J.* 2010; 36(3): 601-7.
14. Christ-Crain M, et al. Procalcitonin guidance of antibiotic therapy in community-acquired pneumonia. *Am. J. Resp. Crit. Care Med.* 2006; 174: 84-93.
15. Hochreiter M, et al. Procalcitonin to guide duration of antibiotic therapy in surgical intensive care patients: a randomized prospective controlled trial. *Crit Care* 2009; 13(3), R83.
16. Kristoffersen KB, et al. Antibiotic treatment interruption of suspected lower respiratory tract infections based on a single procalcitonin measurement at hospital admission—a randomized trial. *Clin. Microbiol. Infect.* 2009; 15(5): 481-7.
17. Long W, et al. Procalcitonin guidance for reduction of antibiotic use in low-risk outpatients with community-acquired pneumonia. *Respirology*. 2011; 16(5): 819-24.
18. Nobre V, et al. Use of procalcitonin to shorten antibiotic treatment duration in septic patients: a randomized trial. *Am. J. Resp. Crit. Care Med.* 2008; 177: 498-505.
19. Schroeder S, et al., Procalcitonin (PCT)-guided algorithm reduces length of antibiotic treatment in surgical intensive care patients with severe sepsis: Results of a prospective randomized study. *Langenbecks Arch Surg* 2009;394(2): 221-6.
20. Stolz D, et al. Antibiotic treatment of exacerbations of copd: a randomized, controlled trial comparing procalcitonin-guidance with standard therapy. *Chest* 2007; 131(1): 9-19.
21. Stolz D, et al. Procalcitonin for reduced antibiotic exposure in ventilator-associated pneumonia: a randomised study. *Eur. Respir. J.* 2009; 34(6): 1364-75.
22. Soni NJ, et al. Procalcitonin-guided antibiotic therapy: a systematic review and meta-analysis. *J Hosp Med* 2013;8(9):530-540.
23. Moses clinical trial data. On file Thermo Fisher Scientific.
24. Thomas-Rueddel DO, et al. Hyperlactatemia is an independent predictor of mortality and denotes distinct subtypes of severe sepsis and septic shock. *J Crit Care* 2015;30(2):439.e1-439.e6.
25. Freund Y. Serum lactate and procalcitonin measurements in emergency room for the diagnosis and risk-stratification of patients with suspected infection. *Biomarkers* 2012;17(7):590-596.
26. Blomkalns AL. Sick or not sick? Evolving biomarkers for severe bacterial infection. 2007;7(November). Emergency Medicine Cardiac Research and Education Group.
27. Vigushin DM. Metabolic and scintigraphic studies of radioiodinated human C-reactive protein in health and disease. *J Clin Invest* 1993;91(4):1351-1357.
28. Pepys MB, et al. C-reactive protein: A critical update. *J Clin Invest* 2003;111(12):1805-1812.
29. Standage SW, et al. Biomarkers for pediatric sepsis and septic shock. *Expert Rev Anti Infect Ther.* 2011;9(1):71-79.
30. The ProCESS Investigators. A randomized trial of protocol-based care for early septic shock. *N Engl J Med.* 2014;370(18):1683-1693.
31. Seigel TA, et al. Inadequacy of temperature and white blood cell count in predicting bacteremia in patients with suspected infection. *J Emerg Med.* 2012;42(3):254-259.
32. Meisner M. Procalcitonin: Experience with a new diagnostic tool for bacterial infection and systemic inflammation. *J Lab Med* 1999;23:263-72.
33. Meisner M, et al. Postoperative plasma concentrations of procalcitonin after different types of surgery. *Intensive Care Med.* 1998;24(7):680-684.
34. Stocker M, et al. Neonatal Procalcitonin Intervention Study (NeoPInS): Effect of Procalcitonin-guided decision making on duration of antibiotic therapy in suspected neonatal early-onset sepsis: a multi-centre randomized superiority and non-inferiority intervention study. *BMC Pediatrics* 2010 10:89.
35. Reith HB, et al. Procalcitonin in early detection of postoperative complications. *Dig Surg.* 1998;15(3):260-265.

Clinical Diagnostics

Germany:

Thermo Fisher Scientific +49 (0)3302 883 0
B·R·A·H·M·S GmbH +49 (0)3302 883 100 fax
Neuendorfstr. 25 info.brahms@thermofisher.com
16761 Hennigsdorf www.procalcitonin.com
Germany

USA:

Thermo Fisher Scientific (800) 232 3342
8365 Valley Pike customerservice.diagnostics.mtn@thermofisher.com
Middletown, VA 22645 www.thermofisher.com/aboutsepsis

Find out more at thermofisher.com/aboutsepsis

Availability of product in each country depends on local regulatory marketing authorization status.

© 2019 Thermo Fisher Scientific Inc. All rights reserved. B·R·A·H·M·S PCT and all other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. All data regarding specifications, terms and pricing correspond to the existing knowledge at the time of the printing. We are not responsible for any errors, misprints or changes. Reprint, also in parts, solely with prior written consent of B·R·A·H·M·S GmbH.

ThermoFisher
SCIENTIFIC