Pacific Hemostasis®

Abnormal Coagulation Reference Plasma (ACRP)

- I. **Intended Use**

Pacific Hemostasis® Abnormal Coagulation Reference Plasma (ACRP) is an assayed reference plasma intended for use in the quantitative determination of coagulation proteins.
- II. **Summary and Principles**

Assayed reference plasma is needed to construct and validate standard curves that are used to measure the activity of patient plasmas. To provide maximum accuracy, ACRP is referenced against a freeze dried assayed standard. Some analytes are also referenced to World Health Organization (WHO) or International Society on Thrombosis and Haemostasis/Scientific and Standardization Committe (ISTH/SSC) standards. The assigned value of the standard directly influences the final results, so a high quality standard helps ensure accurate patient results.
- III. **Reagent**

For *in vitro* diagnostic use.

Composition: ACRP is manufactured from human plasma collected with <0.1% sodium citrate anticoagulant. ACRP has been adjusted to yield coagulation factors in the abnormal range. The plasmas are rapidly processed, buffered, and freeze dried to ensure stability.

Store unopened vials at -25°C or below. Reconstitute with 1.0 mL of distilled water. Swirl gently until no material adheres to the glass and let stand undisturbed for 15 minutes at room temperature. Do not invert vial or mix vigorously. After proper reconstitution, ACRP is stable for 2 hours at 2 to 8°C. Gently mix contents prior to each use.

Erratic values, product color variations, or lack of vacuum in the vials could indicate product deterioration. However, poor performance could also be due to other factors within the test system.

Caution: Each unit of source material used in the preparation of this product has been tested by an FDA licensed method and found non-reactive for HBsAg and negative for antibodies to HIV and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS, or other infectious diseases. This product, like all materials of human origin, should be handled as potentially infectious biological material.
- IV. **Test Procedure**

Materials Provided:

Abnormal Coagulation Reference Plasma (ACRP), 10x1 mL.
- Materials Required, But Not Provided:
- Distilled water
- Pipette capable of accurately delivering 1.0 mL
- Use the reconstituted ACRP to perform quantitative coagulation assays according to routine test procedures. Follow instructions provided by instrument and reagent manufacturers.
- VI. **Results**

Assay values are provided for each lot. However, actual results depend largely on the standard, reagents, instrument, and individual laboratory protocols in use. Each laboratory should evaluate these products with their reagent system and methodology to establish acceptable limits.
- VI. **Limitations**

All plasmas are subject to the limitations of the test system. Variables such as temperature, reagent stability, instrument performance, and individual technique can influence final results. Always follow instrument and reagent manufacturers guidelines.
- VII. **Expected Values**

Actual values recovered depend on individual technique, instrument, standard, and reagent used.
- VIII.Performance Characteristics
- Abnormal Coagulation Reference Plasma is formulated to have values approximately 50% of normal. ³
- IX. **References**
 1. *Design History File*
 2. Websites
www.isth.org
www.nibsc.ac.uk
www.who.org

Ordering Information			
Cat. No.	Description	Contents	
100603	Abnormal Coagulation Reference Plasma (ACRP)	10 x 1 mL	

FISHER DIAGNOSTICS ® LIMITED WARRANTY

Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD's sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Pacific Hemostasis

Abnormal Coagulation Reference Plasma (ACRP)

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Abnormal Coagulation Reference Plasma (ACRP), 10 x 1 mL.

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All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Pacific Hemostasis Abnormal Coagulation Reference Plasma (ACRP)

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