




PepMap Neo Columns Enabling brilliant performance

Is your peptide mapping analysis delivering the consistent performance you need?

Developing methods that are robust and reproducible so they can handle analytical variability can be time-consuming and frustrating. The solution? Thermo Scientific™ EASY-Spray™ and double nanoViper™ PepMap Neo™ Columns, which have been designed and manufactured for consistent column-to-column performance to remove these obstacles. They also connect seamlessly with our Thermo Scientific™ Vanquish™ Neo UHPLC System, EASY-Spray and NanoSpray Flex™ Ion Sources—eliminating the risk of leaks, failures, and poor connections.

The consistency of these columns has been reported in [Thermo Scientific Technical Note 74152](#). Reproducible results were obtained with multiple columns over a large sample cohort. Four 75 µm I.D. × 750 mm columns were used for HeLa proteome profiling with a 2000 ng loading and a 240 min gradient. The column-to-column variation for peptide and protein identifications was below 5% and 1%, respectively. Similar results were obtained with the number of quantifiable peptides and proteins in each sample, as shown on the next page.

Excellent column performance



Reproducible identification and quantification using PepMap Neo columns (ES75750PN, 75 µm I.D. × 750 mm): HeLa peptides and proteins over four EASY-Spray PepMap Neo columns while using the Vanquish Neo UHPLC system (direct injection workflow) coupled with the Thermo Scientific™ Orbitrap™ Exploris 480 Mass Spectrometer.

Learn more at thermofisher.com/lowflowHPLCcolumns
and thermofisher.com/VanquishNeo