

AI-driven TMT Data Analysis using CHIMERY5 and INFERY5 in Proteome Discoverer

Mark Sanders¹, Mick Greer¹, David M. Horn¹, Kai Fritzsche²

¹Thermo Fisher Scientific, San Jose, CA; ²Thermo Fisher Scientific, Bremen, Germany

ABSTRACT

The most popular acquisition strategies for multiplexed quantitation using tandem mass tags (TMT) are broadly grouped into MS2 and MS3 based approaches. Here we demonstrate the benefit of adding machine learning (ML) algorithms to TMT data analysis tools. Addition of INFERY5 and CHIMERY5¹ increase identified and quantified proteins. Furthermore, we recommend data analysis workflows that best suit MS2 and MS3 based methods.

INTRODUCTION

Multiplexed proteomics quantitation using Thermo Scientific™ Tandem Mass Tags™ (TMT™) is a powerful tool to measure differences in cellular states. Much effort has been applied to instrument development and development of smart acquisition strategies such as real time search (RTS) and synchronous precursor selection MS3 (SPS-MS3) to improve the quantitative performance of TMT workflows. At the same time, there has been significant progress leveraging deep-learning and AI to enhance the depth of proteomics data analysis. For instance, Here we demonstrate the application of deep learning algorithms (INFERY5 and CHIMERY5) to increase identified and quantified TMT labeled peptides in MS2 and MS3 based acquisition methods.

MATERIALS AND METHODS

Sample preparation and data acquisition were carried out as described in Fürtwangler *et al*, MCP, 2022 (PXD029320)², Koenig *et al*, Proteomics, 2022 (PXD031277)³, and Paulo *et al*, JASMS, 2021 (PXD020815)⁴. Data from these studies were used here. CHIMERY5 and INFERY5 algorithms were implemented as nodes in the Proteome Discoverer 3.0 software, except for data from PXD031277 which was processed in Thermo Scientific™ Proteome Discoverer™ 3.1 using the updated INFERY5 partition model that now supports phosphorylation modifications. CHIMERY5 was deployed in the cloud and received prepared data from a local instance of Proteome Discoverer 3.0. Once processed by the CHIMERY5 service, results were transmitted to the local instance for final processing and results viewing.

CHIMERY5 AND INFERY5 WORKFLOWS

Figure 1. Representative TMT Processing Workflows in Proteome Discoverer





Figure 1 shows the processing workflows used for TMT data analysis. Note: ptmRS is included here because it was used for phosphorylation modification localization confidence (75% threshold), however this can be removed if modification localization is not required. CHIMERY5 is a complete search engine which includes FDR validation, whereas INFERY5 rescoring receives PSMs from SEQUEST and adds additional figures of merit (i.e. spectral angle) to each PSM which can be used by Percolator for the FDR calculation. An overview of both approaches can be seen in Figure 2.

Figure 2. Features of CHIMERY5 and INFERY5



- Search Engine using ML based spectral prediction.
- Mokapot FDR validation introduced in PD3.1

- Re-scores PSMs from SEQUEST using ML spectral prediction.
- Provides novel metrics to percolator (spectral angle etc.)

IMPROVEMENTS for MS2 and MS3 TMT DATA ANALYSIS

While CHIMERY5 and INFERY5 show improvements over SEQUEST in most cases Figures 3&4 demonstrate that CHIMERY5 offers the greatest improvement for MS2 based TMT data while INFERY5 and CHIMERY5 show comparable performance for MS3 based approaches. These results show the impact of data complexity on the results. For more complex data such as the MS2 data CHIMERY5 can accurately deconvolute and predict complex chimeric spectra (Figure 5) while INFERY5 is limited to the PSMs assigned by SEQUEST (via Precursor Detector node), which assign fewer chimeric spectra. However, for MS3 based approaches (Figure 4) INFERY5 and CHIMERY5 perform similarly.

Figure 3 results show the Hypro16 standard described by Paulo *et al.*⁴ which consists of increasing amounts of yeast spiked into a constant background of human peptides. It is important to note in these data closeout was not used during real time search. Addition of closeout would restrict the number of peptides per protein selected for activation, therefore increasing the breadth of coverage for the targeted proteome. hrMS2 shows similar performance to RTS-MS3 however RTS-MS3 generates more accurate ratios (data not shown) and with closeout could yield more than hrMS2 in some cases.

Figure 3. INFERY5 and CHIMERY5 increase results vs. SEQUEST HT alone





Figure 4. INFERY5 and CHIMERY5 increase results vs. SEQUEST HT for single cell data



PHOSPHOPEPTIDE PREDICTION

CHIMERY5 and INFERY5 models in PD 3.1 now support prediction of phosphorylation modifications and neutral losses. Figure 5 demonstrates the capability of CHIMERY5 to deconvolute and identify several PSMs from chimeric spectra. Figure 6 shows a complex chimeric spectra which contains 7 TMTpro-labeled phosphopeptide PSMs. The upper half of the mirror plot shows the fragment ion matches assigned by SEQUEST for the 4 most abundant PSMs. The lower shows the matching INFERY5 predicted fragment spectra. All PSMs which are assigned and predicted in Figure 6 are shown highlighted in the table below. For clarity only the first 4 assigned/predicted PSMs are displayed in the mirror plot.

Figure 5. CHIMERY5 ID's up to 7 PSMs for a single MS2 scan (high-load phospho sample)




Figure 6. Predicted and Observed TMT labeled phosphopeptide PSMs




Figure 7. CHIMERY5 improves analysis of phospho TMTpro data vs SEQUEST HT




Figure 7 shows the increase in identification performance of TMT labeled phosphopeptides using CHIMERY5 compared to SEQUEST HT. Larger gains in PSMs and IDs can be attributed to thorough assignment of chimeric spectra (Figure 5) from these complex data (200 µg starting material, before enrichment, 45 min gradient and 1.3 Th isolation.)³

CONCLUSIONS

INFERY5 and CHIMERY5 can be used to increase identified and quantified TMT peptides.

- INFERY5 and CHIMERY5 show greater performance than SEQUEST HT for complex data
- CHIMERY5 should be included for MS2 data analysis
- For MS3 based TMT acquisition, INFERY5 Rescoring should be included
- Enhanced CHIMERY5 model supports TMTpro and phosphopeptide prediction and showed a substantial increase in the number of identified TMTpro modified phosphopeptides.

REFERENCES

- Zolg, DP, Gessulat, S, Paschke, C, et al. (2021) INFERY5 rescoring: Boosting peptide identifications and scoring confidence of database search results. Rapid Commun Mass Spectrom.
- Furtwangler, Benjamin et al. (2022) Real-Time Search-Assisted Acquisition on a Tribrid Mass Spectrometer Improves Coverage in Multiplexed Single-Cell Proteomics MCP, Volume 21, Issue 4
- Koenig, C., et al.(2022). Optimal analytical strategies for sensitive and quantitative phosphoproteomics using TMT-based multiplexing. *Proteomics*, 22, e2100245.
- Paulo JA., et al.(2021) Hypro16: A Two-Proteome Mixture to Assess Interference in Isobaric Tag-Based Sample Multiplexing Experiments. *JASMS* 2021 6;32(1):247-254.

TRADEMARKS/LICENSING

© 2023 Thermo Fisher Scientific Inc. All rights reserved. SEQUEST is a trademark of the University of Washington. CHIMERY5 and INFERY5 are trademarks of MSAID. TMTpro is a trademark of Proteome Sciences. Tandem Mass Tag is a trademark of Electrophoretics Limited. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

PO2023-58EN