

Real time library search for the confident annotation of metabolites on Thermo Ascend Tribrid instrument

Rahul Deshpande¹, Bashar Amer¹, Jingjing Huang¹, David Bergen¹, Susan Bird¹, and Thomas Moehring², ¹Thermo Fisher Scientific, San Jose, USA, ²Thermo Fisher Scientific, Bremen, Germany

ABSTRACT

Purpose: Development of a compound class specific deep scan for confident annotation of relevant analytes using a high quality curated High Resolution Accurate Mass (HRAM) MSⁿ spectral library in conjunction with Real Time Library Search (RTLS).

Methods: A Thermo Scientific™ Vanquish™ Horizon UHPLC system coupled to a Thermo Scientific™ Orbitrap Ascend™ Tribrid™ mass spectrometer is used for collecting all MS and MSⁿ data. For data acquisition, MS/MS is always collected with precursor ions detected in the survey MS scan within 1.2 second cycle time. High order MSⁿ (3-5) is collected using the built-in spectral library creation templates using LC-MS. The generated MSⁿ tree data are processed using Thermo Scientific™ Mass Frontier™ 8.0 and Thermo Scientific™ Compound Discoverer™ 3.3 software. Flavonoid standards were used for the library creation and putative flavonoids were annotated in tea extract. Thermo Scientific™ mzVault™ 2.3 software was used for creating the spectral library to be used for RTLS.

Results: Use of standard spectral library and RTLS allows more unknown flavonoid compounds being identified from the natural products. Nearly 100 more putative flavonoids were annotated using this workflow from tea extracts.

INTRODUCTION

Confident compound annotation is essential for translating untargeted metabolomics data into meaningful biological information. Mining of relevant information from untargeted metabolomic data remains a challenge. Local mass spectral libraries consisting of MSⁿ data from authentic standards provide the opportunity to confidently annotate known and related unknown compounds in biological samples. Here we describe the process of creating and curating a high-quality spectral library of high-resolution accurate mass data using authentic flavonoid standards and the built-in library builder method template on the a Thermo Scientific™ Orbitrap Ascend™ Tribrid™ Mass Spectrometer (MS) and its utility for the confident annotation of flavonoids in biological samples using Real Time Library Search (RTLS) to automatically guide MSⁿ data acquisition on the fly.

MATERIALS AND METHODS

Sample Preparation

The Flavonoid Library consisting of 40 standards was obtained from MetaSci (Toronto, ON, CA). The library consists of 4 vials each containing 10 authentic standards which could be directly injected for LC-MS analysis. Three varieties of green tea, black tea and herbal tea were obtained from the local market. Three replicates of each were extracted using 240 mL of hot water in commercially available disposable coffee cup. For blank, hot water without tea in similar coffee cup was used. Pooled samples were prepared by adding 100 µL of each extract. This pooled sample was used for QC.

HPLC Conditions

A Thermo Scientific™ Vanquish™ UHPLC system performed separations. The column was a Thermo Scientific™ Accucore™ Vanquish™ C18+ UHPLC column (2.1 x 150mm, 1.5µm).




Figure 1. Thermo Scientific™ Orbitrap Ascend™ Tribrid™ mass spectrometer. Dual Ion Routing Multipoles improves ion transfer, capture and dissociation leading to higher number of MSⁿ scans.




Figure 2. Example flavonoid structures showing shared structural similarity

Mass Spectrometer

Thermo Scientific™ Orbitrap Ascend™ Tribrid™ Mass Spectrometer was used.

MS Conditions-Spectral Library Building

The inbuilt template for library building on the instrument control software using LC-MS was used for the flavonoid standards. The template allows for the collection of MSⁿ spectral data. Data was also collected with a modified template allowing for collection of MS² data utilizing different fragmentation strategies (HCD, CID and UVPD).

MS Conditions-RTLS

The instrument control software template for using RTLS on the tribrid (named Met-IQ) was used for data acquisition on the pooled sample. The template allows for similarity search of compounds present in the library and a trigger for MS³ for compounds similar to those in the library.

MS Conditions-Data Acquisition

Deep Scan Acquire-X acquisition workflow was used for analysis. This workflow automatically creates an exclusion list from background ions (Blank) and an inclusion list of metabolites of interest from the reference sample (Pooled) with iterative updating in between each injection after performing MS². Full scan MS¹ (200-1200 m/z) was performed on all the tea samples.




Figure 3. Thermo Scientific™ AcquireX Deep Scan mode for intelligent data acquisition to maximize the number of relevant compounds interrogated by MS/MS.

RESULTS

Spectral Library

Thermo Scientific™ Compound Discoverer™ 3.3 software was used to process the LC-MS data from the standards. The spectral data with the various fragmentation strategies from the standards was transferred to a mzVault™ library. The mzVault™ library containing the spectral fragmentation information of the flavonoid standards was populated with the structural and chemical information of individual compounds. This library was used for the RTLS on the instrument.

Data Processing

Thermo Scientific™ Compound Discoverer™ 3.3 software was used to process the LC-MS data from the samples. Differential analysis and compound annotation revealed relative differences between the various tea samples. Several annotation tools are available in the Compound Discoverer 3.3 software including database searching and spectral library matching against the mzCloud™ spectral library at the MSⁿ level. In parallel it has the ability to search custom mass lists with corresponding chemical structures. Here the data was processed using the Arita Flavonoid Database which has the mass and structural information of around 6500 Flavonoids.

Data Filtering

The data obtained from Compound Discoverer™ 3.3 was filtered for showing the utility of RTLS to trigger more relevant compounds (structurally similar to flavonoids) for MSⁿ fragmentation. The filter applied extracted only compounds in the Arita Flavonoid Database which had a subsequent MS³ fragmentation.




Figure 4. Characterization of structurally similar unknown compounds using Real Time Library Search (RTLS)




Figure 5. A) Using RTLS reduces the number of MS³ scans allowing more MS² scans for non-targeted metabolite annotation. B) Use of RTLS increases the number of putative flavonoid identified in the tea samples

RTLS for compound specific deep scan

Two major challenges in untargeted metabolomics are

1. Extracting biologically relevant information
2. Confident compound annotation with structural information

The workflow described here addresses both challenges by utilizing RTLS. RTLS can be used to identify and confidently annotate analytes belonging to certain biological classes. The annotation is done by getting structural information utilizing MSⁿ and multiple fragmentation methods in combination with using Compound Discoverer™ 3.3. The use of RTLS allows more MS² scans for annotation of additional non-flavonoid metabolites.

CONCLUSION

A workflow for utilization of RTLS is for compound class (flavonoids) specific identification of analytes in biological matrices (Tea) is described. It involves

1. Creation of spectral library using authentic flavonoid standards.
2. Utilizing the created spectral library for RTLS of related compounds in biological matrices (tea).

More putative flavonoids were identified in tea using RTLS in comparison with samples without acquired without RTLS.

The challenge of obtaining deep and relevant metabolome coverage can be addressed using intelligent data acquisition strategies such as AcquireX in combination with RTLS.

TRADEMARKS/LICENSING

© 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

For Research Use/Purposes Only – Not for Diagnostic Procedures.
PO2023-125EN.

