
Thermo Fisher S C I E N T I F I C

Best practices to simplify environmental sample analysis by ICP-MS

Jeff Bown

Sr. Application Scientist Trace Elemental Analysis





Agenda



Environmental analysis Overview and challenges



Best practice for the elemental analysis workflow Sample/standard preparation and sample introduction



ICP-MS instrument innovations that simplify analysis Hardware and software features, interference correction **ThermoFisher** SCIENTIFIC

Instrument optimization routines Instrument checks and routine maintenance

Troubleshooting tips and tricks Sensitivity, accuracy, precision, and carryover



 Environmental analysis overview and challenges



Why is elemental analysis required for the environment?

Federal, state, and local laws require the analysis of trace metals to ensure the safety of drinking water, surface and ground waters, and soils for consumption, recreation, farming, building, and urban development. Elemental analysis is also required for industrial wastes for proper disposal.

• Federal laws and regulations

- Safe Drinking Water Act (SDWA)
- National Primary and Secondary Drinking Water Regulations (NPDWR & NSDWR)
- Clean Water Act (CWA)
- Resource Conservation and Recovery Act (RCRA)
- Comprehensive Environmental Response, Compensation, and Liability Act (CERCLA)
- Lead Copper Rule (LCR)
- Unregulated Contaminant Monitoring Rule (UCMR)







U.S. EPA drinking water regulations

National Primary Drinking Water Regulations

Max Contaminant Levels (MCL)

Contaminant	MCL (mg/L)	MCLG (mg/L)
Antimony	0.006	0.006
Arsenic	0.010	0
Barium	2	2
Beryllium	0.004	0.004
Cadmium	0.005	0.005
Chromium (total)	0.1	0.1
Copper	1.3	1.3
Lead	0.015	0
Mercury (inorganic)	0.002	0.002
Selenium	0.05	0.05
Thallium	0.002	0.0005

National Secondary Drinking Water Regulations Secondary Max Contaminant Levels (SMCL)

Contaminant	SMCL (mg/L)	Effects Above SMCL	
Aluminium	0.05 – 0.20	Coloring of water	
Chloride	250	Salty taste	
Copper	1.0	Metallic taste; blue-green staining	
Fluoride	2.0	Tooth discoloration	
Iron	0.3	Rusty color, metallic taste, reddish/orange staining	
Manganese	0.05	Black to brown color and staining; bitter metallic taste	
Silver	0.10	Skin discoloration	
Sulfate	250	Salty taste	
Zinc	5	Metallic taste	

U.S. EPA regulations

Unregulated Contaminant Monitoring Rules (UCMR) Inorganic Contaminant Levels

Rule	Contaminant	Minimum Reporting Level (µg/L)
	Vanadium	0.2
	Molybdenum	1.0
UCMR 3	Cobalt	1.0
(2012 - 2016)	Strontium	0.3
	Chromium (total)	0.2
	Chromium - 6	0.03
UCMR 4	Germanium	0.3
(2017 - 2021)	Manganese	0.4
UCMR 5 (2022 - 2026)	Lithium	9.0

Lead Copper Rule

Contaminant	MCL (mg/L)	MCLG (mg/L)
Copper	1.3	1.3
Lead	0.015	0

Environmental analysis

Environmental analysis encompasses a broad range of applications throughout industry

Thermo Fisher

SCIEN



Challenges when analyzing environmental samples

- Complex sample matrices
- Variety of samples
- Regulation and compliance
- Trace and ultra-trace detection requirements
- Comprehensive QC standards and protocols
- Failed sample analyses and re-reruns
- Sample throughput
- Quick turnaround of results
- Data transfer and management
- Reporting requirements
- Traceability and documentation
- Maintenance and troubleshooting
- Training new analysts



Thermo Fi








- Why is environmental analysis challenging?
- Where do we begin to address these challenges?
- ❑ How can we prevent these challenges?
- When do we call service or applications support?

Let's start with the sample matrix...





Thermo Fisher

SCIENTIFIC

Regulation and compliance add another layer of challenges

- Detection limit requirements
 - National Primary Drinking Water Regulations
 - National Secondary Drinking Water Regulations
 - Unregulated Contaminant Monitoring Rule (UCMR)
 - Different state/municipal regulations

- Analysis according to EPA or other industry methods
 - Specific quality control protocols
 - Method validation
 - Numerous QC standards and samples analyzed prior to samples
 - Control limit criteria

- Audit, data management, and reporting requirements
- Documentation and traceability
- Data reporting
- Data package audit
- Onsite audit
- Data transfer to LIMS
- Data security



Addressing the challenges in environmental analysis

General best practices to streamline workflow

- Sample and standard preparation
- Sample handling
- Contamination prevention
- Sample introduction automation

Instrument innovations

- Hardware design and features
- Software features
- Interference correction
- Instrument optimization routines
- Troubleshooting and maintenance tips
 - Troubleshoot failures due to sensitivity, accuracy, precision, and carryover














 Best practices for the elemental analysis workflow



Elemental analysis workflow

Processes in the elemental analysis workflow



Thermo Fisher

How can these processes be done more efficiently?

What are some best practices to ensure accuracy and quality data?

General best practices for the elemental analysis workflow

Key for achieving trace and ultra-trace detection required in environmental analysis



Be aware of all contamination sources.



Use clean and compatible apparatus.



Minimize handling and transfer steps.



Measure weights and volumes with accuracy.



Use high purity reagents.



Use high purity reagent water.



Have separate sample and standard preparation areas.



Apply proper skill, be consistent, and pay attention to detail.

Tips and tricks for sample and standard preparations

Along with the general best practices, apply the following for standard and sample preparations

Apparatus

- Use plastic, avoid glass especially for ultra-trace detection limits
 - E.g., PTFE, PFA, PMP, FEP, HDPE, LDPE, PP
- Use mechanical pipettes with disposable plastic tips
- Use Class A volumetric flasks



Standards

- Use high purity, NIST traceable, ISO certified stock standards designated for ICP-MS or ICP-OES analyses
- Use multi-elements stock standards for most preparations
- Use custom standards, if possible
- Use single element standards to prepare the Internal Standard

solution



Reagents

- Use ASTM Type I water (resistivity 18.2 MΩ-cm) or ultrapure water
- Use high purity (e.g., Optima[™]) grade concentrated acids
- Ensure water purification system delivers ASTM Type 1 water
- Ultrapure water is not Deionized water!



Tips and tricks for sample and standard preparation

Tips and tricks to ensure accurate weights and measurements

Analytical Balances

- Calibrate yearly or as needed
- Check at least weekly using Class 1 standard weights and spot check daily, document all checks
- Store balances on a heavy table away from windows, heat, high traffic areas, doorways, vibration



Mechanical Pipettes

- Calibrate yearly or as needed
- Check weekly by measuring increasing volumes of water on an analytical balance
- Spot check daily, document checks
- Use colorless pipette tips
- Always hold pipette upright when drawing up and dispensing liquid
- Pull up and dispense liquid slowly to avoid air bubbles and liquid from going up to pipette causing damage



Tips and tricks for sample and standard preparations

Tips and tricks to streamline sample handling and prevent contamination

Handling & Transfers

- Never place pipette tip directly into container of stock standard or high purity concentrated acid as this will cause contamination.
- Pour an aliquot of stock standard and concentrated acids into disposable plastic beakers (e.g., 5 mL PP) to pipette from when preparing standard solutions.
- Use plastic (e.g., Teflon) wash bottles to store and dispense dilute acid solutions (e.g., 1% HNO₃) for preparations.
- Avoid multiple transfers by preparing calibration standard solutions in autosampler tubes. Ensure autosampler tubes are Class A, metal free, and thoroughly cleansed prior to use.



Labcon MetalFree™ centrifuge tube, Class A, made from ultra clean resins, with additive free cap

Teflon wash bottle, best for ppt level preparations


Sample preparation

Goals of an optimized sample preparation process

- > To convert sample into a solution suitable for introduction to an ICP-OES or ICP-MS
- Decompose the sample matrix, completely or partially
- Complete solution and retention of analytes at measurable concentrations
- Prevent loss of analytes
- Minimize sample contamination
- Reduce digestion time to satisfy laboratory throughput and turnaround requirements



Environmental Sample (e.g., ground/surface water, wastewater, soils, sludge)







Digestion Method (hot plate, hot block, and microwave)

Digestate (clean, colorless solution)

Digestion methods



Hot Block Acid Digestion

Advantages

- Reduced sample handling
 and transfers
- Reduced exposure to contamination
- All plastic parts, no metal
- Elimination of issues associated with use of glassware

Disadvantages

- Long digestion time (hours)
- Incomplete digestion
- Loss of analyte
- Exposure to atmosphere
- High reagent consumption
- Constant monitoring

Hot Plate Acid Digestion

Disadvantages

Exposure to contamination

High reagent consumption

Numerous handling steps

Incomplete digestion

Constant monitoring

Loss of analyte

Inefficient

Long digestion time (hours)

Advantages

- Simple set-up, needs minimal and common apparatus
- Uncomplicated procedures
- Higher sample weights
- High number of samples
- Low initial investment


Thermo Fisher S C I E N T I F I C

Digestion method


Microwave assisted acid digestion

- Advantages
 - Faster digestion (e.g., 20 minutes)
 - Quality digestion
 - Reduced exposure to contamination
 - Reduced reagent consumption
 - Reduced loss of analyte
 - Overall efficiency
- Disadvantages
 - Higher initial investment
 - Limited number of samples
 - Ease of set-up





Microwave digestion – US EPA Method 3051a



73

Digest up to 24 solid matrix samples simultaneously

ThermoFis

- Digestion time10 minutes
- Samples ready for analysis in under 30 minutes compared to over 2 hours with traditional hot plate/block digestion
- Use disposable liners to eliminate vessel washing

Sample introduction automation

Autosamplers and autodilution systems

Teledyne CETAC



Teledyne CETAC SDX_{HPLD} Liquid Dilution System with ASX-560 Autosampler

Elemental Scientific (ESI)




Autosamplers and autodilution systems

Discrete sampling valves for optimized uptake and washout

- FAST is a reliable, high-throughput, automated sample introduction system using a discrete sampling valve
- Speeds up total analysis time by reducing sample uptake time and rinse time between samples







ESI prepFAST autodilution system

- PrepFAST is an inline dilution system that fully automates laboratory dilutions while providing high sample throughput
- High precision autodilution up to 400x
- Autocalibration from one or multiple stock standards
- Syringe-driven internal standard addition
- Vacuum or syringe sample loading on prep*FAST*
- Autodilution performed when Qtegra records an out-of-range analyte concentration or internal standard failure






prep*FAST* valve module on an iCAP RQ and SC-2DX autosampler cart with S400V syringe pump


Autodilution system for automatic standard preparation

Automatic standard preparation reduces preparation time and systematic error



Streamlined elemental analysis workflow







ICP-MS instrument innovations that simplify analysis



Instrument solutions for environmental analysis

A portfolio of innovative instruments for simplified environmental analysis

FA **Atomic Absorption Spectrometry Inductively Coupled Plasma Optical Inductively Coupled Plasma Mass** (AAS) **Emission Spectroscopy (ICP-OES) Spectrometry (ICP-MS)** Flame & Graphite **Graphite Furnace** Flame Furnace Thermo Scientific[™] Thermo Scientific™ iCAP™ RQ ICP-MS iCAP™ TQ ICP-MS Thermo Scientific[™] iCE[™] 3000 Series AAS Thermo Scientific[™] PRO[™] Series ICP-OES Improved detection capability Lower investment and complexity ✓ Fast, multi-element analysis High sensitivity and wide dynamic range ✓ Robustness for high matrix samples ✓ GFAA offers high sensitivity for key elements

Thermo Fisher

Removal of advanced spectral interferences


Hyphenated techniques

✓ Flexibility, performance, and ease of use

29 jeff.bown@thermofisher.com | 22-August-2022


✓ Flame offers fast, single element analysis

Thermo Scientific iCAP RQ ICP-MS



Thermo Fisher S C I E N T I F I C

Key components for performance, robustness, and ease of use




Sample introduction system and plasma interface

Quick connect sample introduction components








Plasma interface unique drop-down door





Unique interface designed for maximum coverage of sensitivity and dynamic range:

Thermo Fisher


.....

High Sensitivity 2.8 mm	High Matrix 3.5 mm	Robust 4.5 mm	Argon Gas Dilution
Below 0.1% TDS	Up to 0.2% TDS	Up to 0.5% TDS	Above 0.5% and up to 3-4% up to 25% possible
Pure solutions	Universal, e.g. digested samples, drinking waters etc.	Waste waters, hard drinking waters	Sea water, highly saline sample types

Easily implemented through dedicated grouping and autotunes in the software


High matrix samples

- Use of Argon Gas Dilution (AGD) for robust, long-term analysis of high matrix samples (e.g., sea water, brackish water, wastewater) in one analytical run sequence
 - Consistent internal standard recovery despite the changing sample matrices





Spectral Interferences in ICP-MS



Thermo Fis

Polyatomic interferences are formed when 2 or more isotopes combine to form species with the same m/z as the analyte ion, e.g., ⁴⁰Ar¹⁶O — ⁵⁶Fe; ⁴⁰Ar³⁵CI — ⁷⁵As; ⁴⁰Ar¹²C — ⁵²Cr

Spectral interference removal with QCell[™] Collision/Reaction Cell (CRC) Technology


Kinetic Energy Discrimination (KED)

- Single mode interference removal with He KED
- Method development is simplified as He KED
 eliminates interferences for most applications



Quadrupole set to filter out exact mass of target analyte

QCell in collision mode with pure He uses energy discrimination




With QCell, KED is complemented by a second active mechanism...

Addressing challenges through instrument innovations

Spectral interference removal with QCell CRC technology

Kinetic Energy Discrimination (KED) plus Low Mass Cutoff (LMCO)

• A unique characteristic of flatapoles, used in QCell, to remove lower mass precursor ions, that can form new interferences, from the ion beam



Target

Analyte

⁷⁵As⁺

Thermo Fisher

Comprehensive

interference

removal is

achieved

Addressing interferences with Low Mass Cutoff

Low Mass Cutoff eliminates lower mass ions from traveling through the Qcell and contributing to interferences.....

Mass	Interferences	Precursors	Mass	Interferences	Precursors	
⁵¹ V	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, N, O, S, Cl	⁵¹ V	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, N, O, S, CI	⁴⁰ Ar+, ⁴⁰ Ca+, ³⁵ Cl+, ³¹ P+, ³² S+, ¹⁴ N+, ¹
⁵⁶ Fe	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca	⁵⁶ Fe	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca	г, О , N,
⁶³ Cu	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ³¹ P ³² S	C, N, O, Na, P, S, Cl, Ar	⁶³ Cu	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ³¹ P ³² S	C, N, O, Na, P, S, Cl, Ar	
⁷⁵ As	⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ⁴⁰ Ar ³⁴ SH, ³⁷ Cl ² H	H, S, Cl, Ca, Ar	⁷⁵ As	⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ⁴⁰ Ar ³⁴ SH, ³⁷ Cl ² H	H, S, Cl, Ca, Ar	


....and reduces BECs further than KED alone

Target Analyte ⁷⁵As⁺ Thermo

Thermo Fisher

Handling interferences – KED with LMCO

Calibration curve for ⁷⁵As in a solution containing 1.2% v/v HCI



STD mode: Polyatomic interference leads to poor IDL and elevated BEC

KED mode: Polyatomic interference removed - IDL below 5 ppt


Interferences in ICP-MS

Advanced spectral interferences

- Isobaric Interference
 - · occurs when two elements have isotopes with the same nominal mass
 - + E.g., ^{58}Fe and $^{58}\text{Ni};\,^{204}\text{Pb}$ and $^{204}\text{Hg},\,^{40}\text{Ca}$ and ^{40}Ar
- Doubly charged ion interferences
 - formed from elements having a 2nd ionization potential lower than the ionization potential of argon (15.8 eV)


Thermo

- Appear at half the parent isotope mass
- 150Nd2+ on 75As+, 156Gd2+ on 78Se+



Advanced interference removal

Thermo Scientific[™] iCAP[™] Qnova Series ICP-MS



iCAP RQ ICP-MS

iCAP TQ ICP-MS

Built-in safety for handling reactive gases

Thermo Fisher

4 mass flow controllers with optimized flow rates

Additional quadrupole for superior interference removal


Single Quadrupole (SQ) ICP-MS

 Effective removal of polyatomic interferences with QCell technology using KED with LMCO

Triple Quadrupole (TQ) ICP-MS

- Effective removal of polyatomic interferences with Qcell technology using KED with LMCO
- Advaned interference removal of isobaric and doubly charged ion interferences

Triple quadrupole ICP-MS



First mass filtering quadrupole situated axially in front of a second quadrupole

Second quadrupole acting as a collision/reaction cell (CRC)

Q3 Third off-axis mass filtering quadrupole for mass analysis.

$$R \propto N = \frac{fl}{v_z}$$


R = Mass resolving power N = Number of RF cycles experienced by the ion f = RF frequency

- = Quadrupole length
- v_z = Initial ion velocity along quadrupole axis

Voo et al., J. Vac. Sci. Technol. A, 1997, 15 (4)


- ✓ Q1 isolates analyte and interferences that affect accurate measurement
- High Sensitivity (iMS) or High Resolution (1amu) as required by the application

How does triple quadrupole ICP-MS work?




Mass shift mode with oxygen

- Use reaction gas to selectively shift the analyte to a higher mass
- Measure the analyte isotope at its new, higher mass




On mass mode with oxygen

- Use reaction gas to selectively shift the interference to a higher mass
- Measure the analyte isotope at its original mass




Triple quadrupole reactions – interference removal


Calibration curve for ⁷⁵As in a solution containing 1.2% v/v HCI



STD mode: IDL = 227 ppt BEC = 3,718 ppt



KED mode: IDL = 4.1 ppt (55x better) BEC = 3.0 ppt (1240x better)



TQ-O₂ mode: IDL = 0.17 ppt (1340x better) BEC = 0.90 ppt (4130x better)

Simplified method development with a streamline workflow software

Thermo Scientific Qtegra[™] Intelligent Scientific Data Solution[™] (ISDS) Software



- Intuitive, streamline workflow platform
- A range of new features added for ease of use
- Built-in QC features for Method 200.8 and 6020B
- Same leading 21 CFR Part 11 compliance capabilities
- Common across Thermo Scientific ICP-OES and ICP-MS instruments




Thermo Fig

Streamline workflow software


Qtegra ISDS Software

- From Simple Workflows to Quality Results
- Easy to use, yet contains features for flexibility in method development



Streamline workflow software

Intuitive method development from starting from the top moving down




ThermoFisher

SCIENTIFIC

Simplified method development for triple quadrupole ICP-MS

- Triple quadrupole ICP-MS offers multiple interference modes for accurate analysis of your sample
- Problem: When faced with the measurement of a sample where interferences are expected, which is the best measurement mode???

- Which analyte isotope?
- Which gas? None, He, reactive gas?
- Which product ion?




Thermo Fi

Reaction Finder method development assistant


How does it work?

Select the element/isotope of interest with a single click



How Does Reaction Finder Method Development Assistant Work?

• Reaction Finder loads a list of corresponding default optimum settings from a stored database



Thermo Fisher


iMS ensures *complete interference removal* while *maximizing sensitivity*

Reaction Finder method development assistant

Flexibility and easy method set up

Identifier	Q3 Analyte	SQ/TQ	CR Gas Flow	CR Gas	Dwell time (s)	Channels	Spacing (u)	Q1 resolution	Q3 resolution
9Be (S-SQ-KED)		SQ	Normal	KED	0.1	1	0.1	High	Normal
23Na (S-SQ-KED)		SQ	Normal	KED	0.1	1	0.1	Normal	High
39K (S-SO-KED)		SQ	Normal	KED	0.1	1	0.1	Normal	Normal
48Ti 48Ti. 14N4. 1H10 (S-TQ-NH3)	48Ti.14N4.1H10	TQ	Normal	NHa	0.1	1	0.1	Normal	Normal
500. 500 100 (0. TO. 00)	500, 100	TO	1	0.	0.1	1	0.1	Normal	Normal
59Co (S-SQ-KED)		SQ	Normal	KED	0.1	1	0.1	Normal	Normal
75As 75As.160 (S-TQ-02)	75As.160	TQ	Normal	01	0.1	1	0.1	Normal	Normal
80Se 80Se. 16O (S-TQ-O2)	80Se.16O	TQ	Normal	01	0.1	1	0.1	Normal	Normal
111Cd 111Cd (S-TQ-iO2)	111Cd (110.904u) (de	TQ	High	01	0.1	1	0.1	Normal	Normal
138Ba 138Ba. 16O (S-TQ-O2)	138Ba.16O	TQ	Normal	01	0.1	1	0.1	Normal	Normal
205TI (S-SQ-KED)		SQ	Normal	KED	0.1	1	0.1	Normal	Normal
208Pb (S-SQ-KED)		SQ	Normal	KED	0.1	1	0.1	Normal	Normal
111Cd (S-SQ-KED)		SQ	Normal	KED 👻	0.1	1	0.1	Normal	Normal

- Flexibility to modify settings or add modes for key elements
 - Add other modes for method development
 - Choose another gas or product ion (list automatically displayed)



Intelligent Mode switching chooses the most optimal analysis order for each run

20 🌲


Advanced Parameters


Number of sweeps:

Instrument optimization routines and troubleshooting tips and tricks



Typical ICP-MS analysis workflow





Torch assembly

- Inspect injector/center tube and torch for matrix build up, blockages, cracks, fractures, and devitrification
 - these defects will affect sensitivity, precision, accuracy, and cause drift throughout the analysis
- **TIP:** Always have a spare torch assembly, clean and ready to use

iCAP RQ/TQ ICP-MS torch assembly



Thermol

Sample introduction system – torch assembly

Torch assembly - 2 types of torches available for ICP-MS analysis

Quartz torch

- Comes standard with the iCAP RQ ICP-MS
- Good for most aqueous applications consisting of dilute acid solutions
- Quartz has a high coefficient of linear expansion
- Disadvantages
 - Devitrification
 - Poor tolerance to high matrix
 - Not compatible with HF
 - Maintenance and replacement

Ceramic PLUS torch

- PLUS Performance, Lifetime,
 Ultraclean Spectrum
- Made from high purity and highperformance ceramic material
- Identical geometry as the standard quartz torch
- Benefits
 - Decrease in background for Si
 - Resistant to HF
 - Improved robustness for high matrix samples
 - Less maintenance and less frequent replacement



Quartz	Ceramic			
torch	torch			

Thermo Fisher S C I E N T I F I C

Sample introduction system

Check the spray chamber


Instrument checks

- There should be no droplets inside the spray chamber
 - > Droplets and condensation along the walls of the spray chamber cause signal instability



Baffled cyclonic spray chamber







Nebulizer


checks

- Inspect the nebulizer for deposits at the tip or any damage
- Ensure that the nebulizer is clean and that there are no blockages
 - Deposits and blockages restrict aerosol formation decreasing sensitivity, causing signal drift, and affecting accuracy and precision





Nebulizer - ensure the appropriate type of nebulizer is used for the application



Sample introduction system

Autosamplers and autodilution systems

- Inspect the autosampler lines and autosampler probe
 - Obstructions will cause longer uptake times and poor stability resulting to precision issues
- Ensure that the samples are loaded according to the method
 - Samples must be in the correct location and on the correct rack
- Remove autosampler caps, tops, and any covering from the samples
- Remove any items that will interfere with the movement of the sample probe
- Check the sample probe depth and ensure it is above any precipitate/solids that have settled
- Inspect autosampler wash station pump tubing for wear and tear
- Tip: Use an autosampler cover to prevent dust or dirt from depositing onto samples



Instrument checks

Instrument checks – interface and cones

Instrument • Interface and sample and skimmer cones


- The ICP-MS interface is the point where sample ions are transferred into the mass spectrometer
- Cones can be prone to build up of sample matrix
- Inspect sample and skimmer cones prior to analysis for blockage and wear around the orifices
- Ensure that the appropriate skimmer cone insert is placed at the back of the skimmer cone



Robust 4.5 mm	High matrix 3.5 mm	High sensitivity 2.8 mm

checks

Maintenance, tips, and tricks for sample and skimmer cones



- Sonicating with reagent water for 5 -10 minutes. This should be adequate to clean and restore performance. Conditioning is not necessary as coating of oxides should still be intact.
- If performance issues persist or for tough deposits, sonicate in 2% Citranox or 2% nitric acid for 5 10 minutes. Rinse cones and allow to air dry. Condition cones prior to analysis.
- Handle both cones with care, especially the skimmer cone as the tip is more delicate.
- If tips are chipped or the orifices are enlarged, replace cones as soon as possible.



Thermol

Sample cone




Skimmer cone

Routine maintenance

When does the peristaltic pump tubing need to be replaced?

nte	ensities							
	2	No	Date / Time	Label 🗸	7Li +¤	59Co +¤	115In +⊐	238U +¤
].[38	9/3/2019 4:06:25 PM	<ldentifier></ldentifier>	353,995	376,923	319,111	384,184
]		39	9/3/2019 4:07:20 PM	<ldentifier></ldentifier>	400,960	439,965	364,775	454,618
]	6	40	9/3/2019 4:08:40 PM	<ldentifier></ldentifier>	386,138	412,953	342,489	417,645
	•			1	410,788.1	440,256.5	368,845.4	444,588.5
	•			2	389,443.0	413,674.8	337,935.3	416,452.5
	•			3	358,183.4	384,926.6	320,684.9	391,895.2
				Mean:	386,138.1	412,952.6	342,488.5	417,645.4
				RSD [%]:	6.9	6.7	7.1	6.3
				SD:	26,457.6	27,672.0	24,400.9	26,366.9
					° .	•	0	0
					o	0	0	0
1	2	No	Date / Time	Label 🗸	7Li +¤	59Co +⊐	115In +¤	238U +¤
-		41	9/3/2019 4:09:37 PM	<ldentifier></ldentifier>	421,115	444,523	379,697	454,945
		42	9/3/2019 4:13:45 PM	<ldentifier></ldentifier>	7,165	1,104	371	115
-	10	43	9/3/2019 4:14:41 PM	<ldentifier></ldentifier>	6,993	906	270	3
	30	44	9/3/2019 4:15:59 PM	<ldentifier></ldentifier>	6,973	880	273	1
	10	46	9/3/2019 4:27:54 PM	<ldentifier></ldentifier>	321,208	337,602	280,372	343,619
3		47	9/3/2019 4:28:51 PM	d d an tiff and	316,370	330,756	274,959	338,767

Maintenance, tips, and tricks for the torch assembly



• Clean the quartz torch by:

- Soaking the ends of the torch up to where matrix has deposited in acid solutions (e.g., 5% nitric acid and 2% HCl) for at least 30 minutes or a few hours for persistent deposits.
- Rinse thoroughly with reagent water and allow to air dry completely.
- Do not sonicate the torch and injector or use a wire brush or scraping tools to remove deposits.
- Do not touch torch and injector with bare hands. Always wear gloves when handling torch and injector to prevent oil and moisture from contaminating/damaging the surface.
- Ensure that the argon flow rates are optimized for the application and set prior to plasma ignition. Incorrect settings may cause damage, such as melting of the torch.
- After analysis, rinse the torch by running the blank solution followed by reagent water for a few minutes to prevent formation of matrix/salts inside the injector.



Thermo

Maintenance, tips, and tricks for concentric glass nebulizers


Proactively prevent nebulizer blockage by:

- filtering particulates/suspended solids in the samples prior to analysis
- covering samples using an autosampler enclosure especially for long runs
- Clean the nebulizer by:
 - Soaking in acid solutions (e.g.,10% nitric acid or aqua regia)
 - For heavy deposits, soak the nebulizer for several hours in more concentrated acid (e.g., 20% HNO₃) solution and rinse thoroughly with reagent water
- Do not sonicate or insert a wire through the tip of the nebulizer to remove blockage!
- Do not touch the delicate tip of the nebulizer and do not handle aggressively, store in its original packaging when not in use.
- Monitor the nebulizer back pressure to detect blockages. Record back pressure daily to track upward or downward trends.
- After analysis, rinse the nebulizer by running the blank solution followed by reagent water for a few minutes to prevent sales, sample matrix, etc., from forming inside the capillary. Allow the nebulizer to run dry.
- Disconnect the sample line to prevent liquid from being drawn up to the nebulizer when not in operation.



Concentric glass nebulizer

Maintenance, tips, and tricks for the cyclonic spray chamber



- Soaking in acid solutions (e.g., 5% HNO₃ and 2% HCI) for a few hours or overnight for persistent contamination.
- Rinse with reagent water and allow to air dry completely.
- Do not touch spray chamber with bare hands and do not use a wire brush for cleaning.
- Clean new spray chambers following the same procedure. Although the spray chamber is new, there may be dust or dirt settled inside.
- For samples containing HF, always use a PFA spray chamber.
- After analysis, rinse the spray chamber by running the blank solution followed by reagent water for several minutes to prevent sample deposits from forming inside the spray chamber when the solvent dries out.



Thermo

Dirty spray chamber



Clean spray chamber

Troubleshooting tips and tricks

Troubleshoot issues with sensitivity, precision, accuracy, and contamination/carry over



Sensitivity

Sensitivity issues are typically characterized by decrease or increase of signal and failure of continuing calibration standard (CCV) recoveries.

To Troubleshoot

Check the following:

- Nebulizer or injector blockage
- Sample and skimmer cone orifices for blockage/damage
- Use of nebulizer appropriate for sample matrix
- Dirty spray chamber
- Operating parameters, nebulizer and gas flows, power setting and pump speed
- Interferences and appropriate correction applied
- Old/expired calibration standards
- Analysis of second source standard for reference



Precision issues are typically characterized by high % RSD between sample replicates.

To Troubleshoot

Check the following:

- Worn peristaltic pump tubing
- Nebulizer or injector blockage
- Use of nebulizer appropriate for the sample matrix
- Dirty spray chamber
- Sufficient sample uptake time
- Sufficient rinse time between samples
- Operating parameters, gas flows, pump speed
- Use of the appropriate rinse solution for sample matrix

Troubleshooting Tips and Tricks

Troubleshoot issues with sensitivity, precision, accuracy, carryover and contamination



Accuracy

Accuracy issues are typically characterized by poor sample recoveries, failures in the analysis of CRMs and second source standards.

To Troubleshoot

Check the following:

- Nebulizer or injector blockage
- Use of nebulizer appropriate for sample matrix
- Dirty spray chamber
- Operating parameters, nebulizer and gas flows, power setting and pump speed
- Sufficient uptake time for sample matrix
- Interferences and appropriate correction applied
- Use of appropriate Internal Standard
- Old/expired calibration standards



Contamination and Carryover

Contamination issues are shown by high blanks and sample or standard recoveries. Carryover is characterized by high standard blanks (CCB) and decreasing sample replicates resulting to high % RSD.


To Troubleshoot

Check the following:

- Sufficient rinse time for sample matrix
- Appropriate rinse solution for sample matrix
- Dirty spray chamber
- Contaminated DI water supply and acids, use trace metal or higher-grade acid if possible
- For "sticky" elements (e.g., Hg, Mo, Sb), use longer rinse times. For Hg, use Au to help rinse out Hg.
- Clean work bench/environment free of dust and dirt

Best practices to simplify environmental sample analysis

Streamlining workflow helps to obtain fast, accurate results and quality data



Thermo Fisher

SCLENTLELC

New resource

Guide for environmental sample analysis by ICP-MS

If your laboratory is

- Experiencing analytical challenges, inaccurate results, & sample reruns
- · Seeking to streamline current methodologies and workflows
- Starting up or preparing for environmental sample analysis by ICP-MS

our eBook, "Guide for Environmental Sample Analysis by ICP-MS: Recommendations for Getting Started and Best Practices to Streamline Workflow," serves as a helpful resource

Topics include:

- Considerations and tips for selecting laboratory apparatus, equipment, reagents, and standard solutions
- Best practices for the entire elemental analysis workflow
- Recommended pre-calibration routines and instrument inspections
- General instrument maintenance and troubleshooting tips and tricks

https://www.thermofisher.com/us/en/home/global/forms/ind ustrial/environmental-sample-analysis-by-icp-ms-ebook.html



Guide for environmental sample analysis by ICP-MS:

Recommendations for getting started and best practices to streamline workflow

by Sabrina Antonio

Visit this page or scan the QR code to download the eBook





Thermo Fisher S C I E N T I F I C

72 jeff.bown@thermofisher.com | 22-August-2022