

Solvents

Acids

Purity Grade Chart

Safety Packaging

Chemical Storage Codes



Chemicals for Analytical Techniques

A reference guide to key solvents and acids for every application

"Chemicals for Analytical Techniques" is Fisher Chemical's latest reference guide providing extensive information on chemicals used in wide range of analytical techniques, from LC/MS to HPLC to gas chromatography to spectrophotometry. Use the safety reference guide to quickly access chemical storage codes or spill cleanup symbols. Read more about Fisher Chemical's commitment to not just safely pack each chemical, but ensuring the packaging is environmental friendly as well. Finally, don't waste time searching for the right solvent or acid, use this easy-to-identify product list. Fisher Chemical delivers quality, purity and certainty

Fisher Chemical offers more than 5,000 high-quality chemicals: dry reagents, acids and solutions and high-purity solvents. All chemicals are manufactured in ISO 9001:2008-certified facilities, and undergo rigorous quality assurance and testing procedures, ensuring excellent lot-to-lot and bottle-to-bottle consistency.

For more information about any of our products, or to learn about special offers, please contact your local Fisher Scientific representative or log on to www.fishersci.com/chemicals.

Table of Contents

Safety Quick Reference Guide	pages 1-2
Fisher Chemical Bottle Label	page 1
Catalog Bottle Label	page 2
ChemAlert* Storage Codes	page 2
Spill Cleanup Kit Symbols	page 2
Safety Packaging	pages 3-4
Solvents.....	pages 5-11
LC/MS.....	page 7
HPLC.....	page 8
GC	page 9
ACS	pages 10-11
Certified.....	page 11
Acids.....	page 12
High-Purity.....	page 13
TraceMetal* Analysis	page 14
ACS	pages 14-15
Certified.....	page 16
Periodic Table	page 17

Need more information about chemicals and bioreagents?

Our team of manufacturing chemists, chemical specialists, account representatives and customer service personnel are committed to bringing you the highest level of customer service, technical assistance and fast, accurate delivery.

Customer Service: **1-800-766-7000**

Chemical Technical Support: **1-800-227-6701**
Available in the United States, Puerto Rico and Canada.

Monday – Friday, 8:00 a.m. to 6:30 p.m. ET

Your Partners in Smart Chemistry Solutions

Our essential reagents and chemicals enable analysis, research and development, scale-up and manufacturing through:

Supply Chain Flexibility

Fisher Chemical product line stands for high-quality products and customer-focused, optimized, global solutions. We are committed to improving your transaction efficiency through the broadest and most flexible purchasing options including direct sales, catalogs and e-commerce.

Manufacturing Capabilities

Our blending, custom synthesis, distillation, crystallization, testing, packaging, filling, USP grade, multicompendial and GMP capabilities allow us to manufacture a broad range of products to assist on the various stages from drug discovery to production. Each year we produce and package millions of liters of high-purity liquid reagents and semi-bulk quantities of high-purity dry reagents. We specialize in high-volume packaging, as well as custom micro-scale weighing, while maintaining chemical purity. Through our extensive supply chain network and sourcing excellence we deliver customized chemical solutions to a wide range of customers in pharmaceutical, biotechnology and industrial markets.

Quality and Reliability

Our world-class manufacturing and unmatched expertise enable us to provide customers with the best quality, and deliver products when and where they are needed. Our reagents are manufactured, tested and packaged under ISO 9001:2008 certification.

Innovative Packaging

Our products come in a wide variety of innovative packaging designed for safety, environmental protection, convenient handling and storage, and preservation of product integrity. Our packaging is compliant with all government regulations.

Specialized Chemical Services

All Fisher Chemical products are available in the following custom services and solutions:

- Semi-bulk and bulk chemicals
- Special Solutions
- Tailored Solvents and Solvent Blends
- Testing Services
- Customized Packaging and Labeling

Please contact your Fisher Scientific sales representative for a custom quote.



Safely Serving Science

Safety and customer convenience are at the foundation of our business. From product selection to packaging innovations, we continue to strive to meet the changing needs of our customers, faster and more effectively.

We ship our products in the environmentally friendly, fully recyclable corrugated boxes produced by a Sustainable Forestry Initiative* (SFI)-certified manufacturer. (Learn more on p. 4.)

Our FisherLOCK* cap has been designed to reduce the potential for leakage during transportation and incorporates the easy-to-see color coding that follows National Fire Protection Association (NFPA) industry standards. (Learn more on p. 3.) This same color coding is also reflected on our pressure-sensitive labels, which are printed with chemical resistant ink to protect the important product safety information. (See below.)

Safety Quick Reference Guide

Fisher Chemical* Bottle Labels



- 1 **Bar Code:** an alphanumeric code widely used by industry for inventory and tracking control; this “three-of-nine” code meets the Health Industry Bar Code (HIBC) Supplier Labeling Standard as prepared by the Health Industry Business Communication Council (HIBCC)
- 2 Appropriate **Precautionary Signal Word** (i.e., DANGER, WARNING, CAUTION), which highlights the Potential Health Hazards, both acute and chronic, of the product
- 3 **Molecular Formula, Formula Weight and Product Chemical Abstract Service (CAS) Registry Number; Date Received** is left blank for the customer to write the date the bottle was received to assist in inventory control
- 4 **Maximum Limits of Impurities**
- 5 **Product Name and Purity Grade**
- 6 **Fisher Chemical ChemAlert* Guide:** A three-part guide that provides the Safety Code, the NFPA Code and a Storage Code
- 7 **NFPA Code:** Provides information relating to the severity of health, flammability, reactivity and related hazards that may be presented by short-term, acute exposure to material during handling under emergency conditions such as spills, leaks and fires. Each section displays a numerical severity rating from four (most severe) to zero (least severe). The fourth section, as stipulated by NFPA, is left blank and is reserved for indicating any unusual reactivity with water or oxidizing properties
- 8 **Storage Code:** The storage code is spelled out for clarification and the color-coded bar provides an instant reference guide; refer to p. 2 for the five ChemAlert colors and their descriptions
- 9 **U.S. Department of Transportation (DOT) Shipping Information:** Includes the DOT Proper Shipping Name and U.N./N.A. Number, which identifies the United Nations/North American numerical designation for transportation hazards
- 10 **Catalog Number:** For identification reference
- 11 **Manufacturer’s Address and Telephone Number** (Fisher Chemical is listed)
- 12 **Package Size:** The amount of material per unit or case
- 13 **Lot Number**

Note: ChemAlert is an instant guide only. It should be supplemented by reading the rest of the label (provides detailed instructions in the event of accidental exposure, spill or fire, and applicable OSHA, DOT and ANSI* data), and the appropriate Material Safety Data Sheet and standard references.

Catalog Product Information Labels

Acetonitrile		HPLC	
Also meet 7 CS specifications			
Quantity	Packaging	Cat. No.	Each Case of
1L	Amber Glass	A998-1	93.24 6397.15
1L	Amber Glass/Safe-Cols*	A998SK-1	106.47 6453.48
2.5L	Amber Glass	A998-212	191.24 4580.73
4L	Amber Glass	A998-4	228.23 4548.19
4L	Amber Glass/Safe-Cols*	A998SK-4	240.74 4583.67
19L	NOWPak* I (Plastic Exterior)	A998NI-19	773.72

Also available in recycled FisherPak* and NOWPak* containers.

Methyl Cyanide, Ethanenitrile
Melting point: -34.3°C
 CH_3CN ; $\text{C}_2\text{H}_3\text{N}$; F.W. 41.05
 Submicron Filtered,
 [75-05-8]

Product Specifications
 Actual Lot Analysis is reported on label.

Assay (CH_3CN)	≥99.9%
Color (APHA)	≤10
Optical Absorbance:	
at 190nm	≤1.00
at 195nm	≤0.15
at 200nm	≤0.07
at 205nm	≤0.05
at 210nm	≤0.04
at 220nm	≤0.02
at 254nm	≤0.01
Fluorescence Background (as Quinine Sulfate)	To pass test
LC Gradient Suitability	To pass test
Refractive Index at 25°C	1.3405-1.3425
Residue after Evaporation	≤1ppm
Water	≤0.01%
Titratable Acid	≤0.008mEq/g
Titratable Base	≤0.0005mEq/g
Density at 25°C	0.775-0.780g/mL

- Product Name
- Company Logo
- Container Size
- Synonym(s), Formula, Formula Weight, CAS Registry Number
Also: Physical Data and Special Application Notes
- Product Specifications
- Type of Primary Container
- Catalog Number
- Purity Grade (Fisher Chemical* grades are described on p. 5)
- Price
- Storage Code Color (ChemAlert* Storage Codes are described on the right.)
- Spill Cleanup Kit Symbols (See right for more information.)

Accepted Nomenclature

- Chemical salts and compounds are listed under element name "Sodium Benzoate," not "Benzoate of Sodium"
- Prefixes for organic compounds – alpha, iso, meta, secondary and abbreviations – do not affect alphabetical order (Exceptions: cases where prefix is considered part of the name)
- Numerical prefixes, such as tri and tetra, are treated as an integral part of the name; the prefix mono is not used

Packaging for Safety, Convenience and Product Quality

The Fisher Chemical portfolio comes in a wide variety of innovative packaging designed for safety, environmental protection, convenient handling and storage, and preservation of product integrity. Our packaging is compliant with all government regulations. Shipping cartons covered by the U.S. Department of Transportation regulations also meet international hazardous materials regulations and standards.

Primary containers include:

- Plastic and glass bottles/jars
- Tin cans
- SafeTin* cans
- Aluminum cans
- Square poly bottles
- Steel drum/pails
- PolyPac* containers
- FisherPak*

Learn more about our packaging on pp. 3-4.

ChemAlert Storage Codes



RED (R): Flammable; *Store in an area segregated for flammable reagents*



BLUE (B): Health and hazard. Toxic if inhaled, ingested or absorbed through skin; *Store in secure area*



YELLOW (Y): Reactive and oxidizing reagents. May react violently with air, water or other substance; *Store away from flammable and combustible materials*



WHITE (W): Corrosive. May harm skin, eyes, mucus membranes; *Store away from red, yellow and blue coded reagents*



GRAY (G): Presents no more than moderate hazard in any of the categories above; *For general chemical storage*

Exception: *Reagent incompatible with other reagents of the same color bar; store separately*

Spill Cleanup Kit Symbols

To simplify ordering the appropriate kit, the following symbols are used to identify chemicals that require special cleanup materials and the type of cleanup kit required. These are one-time-use kits containing everything necessary for spill cleanup.



Mercury Warning; this product contains Mercury.



Acid Spill Emergency Cleanup Kit for neutralizing and absorbing up to 1L of acid. (**Cat. No. 18-061A**)




Caustic Spill Emergency Cleanup Kit for neutralizing and absorbing up to 1L of caustic material. (**Cat. No. 18-061C**)



Solvent Spill Emergency Cleanup Kit for neutralizing and absorbing up to 1L of solvent. (**Cat. No. 18-061B**)

Excellence in Packaging through Innovation



Excellence in Packaging through Innovation

The Fisher Chemical* portfolio comes in a wide variety of innovative packaging designed for safety, environmental protection, convenient handling and storage, and preservation of product integrity. Our packaging is compliant with all government regulations.

Plastic Bottles – Safety and Convenience

Available when chemical properties are compatible to minimize the risk of breakage, provide lighter weight packaging and enable more economical shipping.

Amber Glass Bottles – Quality and Reliability

Used to package photosensitive chemicals, to protect them from light.

Safe-Cote* PVC-Coated Bottles – Safety

These tough bottles provide the purity of glass and most of the benefits of plastic for storing and dispensing of solvents and acids. If the bottle breaks, glass fragments and liquids are more likely to remain trapped.

Aluminum Cans – Quality and Reliability

Designed to contain ethers, Fisher Chemical 1L and 4L aluminum cans feature a round, seamless construction, which reduces the possibility of leakage through a seam.

Lock in Quality and Safety with Fisher Chemical Products

Fisher Chemical products are now packaged with the innovative, new FisherLOCK* Cap, designed to help maintain the highest standards in product quality and chemical safety practices. With an exclusive color-coded design (patent pending), the FisherLOCK Cap is engineered to lock to the transfer bead of the bottle when applied during manufacturing to a filled bottle of laboratory chemicals. The cap design includes an interior ring, visible from various angles, which offers resistance until separated from the cap when opened.

The presence and location of that ring imparts tamper evidence to the configuration — without the drawbacks and complications of a plastic seal over the cap. The FisherLOCK Cap is easy-to-use, and supports product quality, reliability of use and safety-in-use of the chemical products you purchase.



The FisherLOCK Cap LOCKS in quality, safety, reliability and convenience:

QUALITY

- Provides a tight, tamper-evident, secure seal to ensure chemical contents arrive fresh and unopened
- Eliminates polyethylene glycol contamination that is possible with a plastic overseal

SAFETY

- Caps are designed to resist back-off during transport, reducing risk of leakage
- Color-coded rings indicate storage requirements and hazard categories, and enhance proper recognition, handling and storage, even before the bottle is removed from the case

RELIABILITY

- Rigorously tested for chemical compatibility
- Bottle threads are unchanged, allowing you to continue to attach the opened bottle to standard equipment

CONVENIENCE

- Cap design facilitates correct initial torque application during manufacturing, thus eliminating caps that may be hard to twist open
- Larger ridges on the exterior of the cap make it easier to open
- Caps readily reseal after initial opening

Going Green Has Never Been So Easy With...



Styrofoam*-free EcoSafPak* packaging. We now deliver all 4 x 4L, 1 x 4L, 6 x 1L, 6 x 500mL, 1 x 2.5L, 1 x 1L, 1 x 500mL, 2 x 4L, 2 x 1L and 2 x 500mL glass bottles (including the Safe-Cote* glass bottles) of your favorite products in eco-friendly EcoSafPak packaging manufactured by an SFI-certified manufacturer.

EcoSafPak packaging minimizes the use of non-recyclable material through the use of the revolutionary Hexacomb* design inserts.

The Hexacomb is:

- Fully recyclable
- Sturdy and durable
- Versatile

Improved Safety and Handling

- EcoSafPak is the only package of its kind to pass the requirements of the International Safe Transport Association "3A" test:
 - 17 Drops to simulate real-world parcel shipment handling
 - Shock testing
 - Vibration testing
 - Corrugated material is shock absorbing
 - Greater stability in the frame of the box
 - Staggered handholds for ease of transport

SFI's Nine Objectives:

- Sustainable forestry
- Responsible practices
- Reforestation and productive capacity
- Forest health and productivity
- Long-term forest and soil productivity
- Protection of water resources
- Protection of special sites and biological diversity
- Legal compliance
- Continual improvement



FIBER USED IN THIS PRODUCT LINE
MEETS THE SOURCING REQUIREMENTS
OF THE SFI PROGRAM
WWW.SFIPROGRAM.ORG

The SFI Program:

- is based on the premise that responsible environmental behavior and sound business decisions can co-exist
- is a comprehensive system of principles, objectives and performance measures developed by professional foresters, conservationists and scientists
- promotes the perpetual growing and harvesting of trees with the long-term protection of wildlife, plants, soil and water



Solvents

The Fisher Chemical* portfolio of high-purity solvents comprises more than 600 solvents suited to a wide range of laboratory applications, from LC/MS to HPLC to gas chromatography to spectrophotometry.

Manufactured in ISO 9001:2008-certified facilities, each Fisher Chemical solvent undergoes rigorous quality assurance and testing measures that ensure excellent lot-to-lot and bottle-to-bottle consistency.

Fisher Chemical Solvents Purity Grades

Grade	Definition	Application	Certificate Of Analysis
Optima* LC/MS	Solvents and additives of exceptionally high purity specially designed and use-tested on LC/MS instruments. Solvents impurity levels in ppb.	LC/MS, HPLC, GC and spectrophotometry	Provided with each shipment
Optima	Acids and solvents of extremely high purity. Acids are analyzed for 65 metals by ICP/MS; impurity levels in ppt. Solvent impurity levels in ppm. UV absorbance curves and sample chromatograms available on request. For Optima acids, a typical lot analysis is given in the catalog. Impurity levels in an actual lot received may vary from the amount listed.	HPLC, GC, plasma/ICP, spectrophotometry, and pesticide residue analysis	Provided with each shipment
GC <i>Resolv</i> *	Solvents with the highest purity and lot-to-lot consistency. Free of contaminants to the ppb level, including those listed in Contract Laboratory Program Target Compound List. Meet ACS specifications. Chromatogram available on request.	Gas Chromatography (GC)	Provided with each shipment
HPLC	Solvents manufactured specifically for use with HPLC instruments. Meet all ACS specifications. Submicron filtered.	HPLC and spectrophotometry procedures	Available on request
Certified ACS Plus	Acids that, in addition to meeting or exceeding the latest specifications of the ACS, are analyzed for more than 16 metals. Actual lot analysis on label.	Analytical applications with tighter metal specifications	Available on request
Certified ACS	Reagent chemicals that meet or exceed the latest ACS Specifications. Actual lot analysis on label.	Analytical applications requiring tight specifications	Available on request
Certified	Reagent chemicals for which the purity standard is established by Fisher Scientific. Purity is guaranteed to meet published maximum limits of impurities.	General analytical procedures	Available on request
Ultra Trace Elemental Analysis	Solvents manufactured for use with plasma/ICP instruments. Impurity levels in ppt. Packaged in acid-cleaned polyethylene bottles.	Plasma/ICP, environmental testing, trace-metal analysis	Available on request
Environmental Grade	Solvents for use in HPLC analysis, trace-organic analysis, and environmental testing. Packaged in precleaned glass bottles in M3.5 (Class 100) clean-room. Shipped with Certificate of Analysis.	HPLC, trace-organic analysis, environmental testing	Provided with each shipment
Pesticide	Solvents for use in analysis of pesticide residue. Meet or exceed ACS standards of purity for pesticide residue analysis.	GC with electron capture detector (ECD), pesticide residue analysis	Available on request
USP/NF/FCC/EP/BP/JP	Reagent chemicals that meet or surpass specifications of the United States Pharmacopeia (USP), the National Formulary (NF), the Food Chemicals Codex (FCC), the European Pharmacopeia (EP), the British Pharmacopeia (BP) and/or the Japanese Pharmacopeia (JP).	Food and drug laboratories, biological testing	20L or greater; 10kg or greater
Spectranalyzed*	Solvents for use in spectrophotometry. Also meet ACS specifications. Actual lot analysis on label.	Ultraviolet and visible wavelength detectors (UV-Vis)	Available on request
Biotechnology	Solvents and reagents that have been specially purified and assayed for biotechnology applications.	Electrophoresis, molecular biology, sequencing, and peptide and oligonucleotide synthesis	Available on request
Scintanalyzed*	Solvents, fluors and prepared cocktails for liquid scintillation counting. Includes nonflammable, nontoxic, biodegradable ScintiSafe* cocktails.	Liquid scintillation counting	Available on request
Electronic	Solvents manufactured to ensure low levels of metal contamination. Meet Semiconductor Equipment and Materials Institute (SEMI) requirements. Actual lot analysis on label.	Electronics and circuit board manufacturing	Available on request
TraceMetal*	Acids manufactured to achieve low metal contamination measurable in ppm to ppb range. Each lot is analyzed for more than 55 metals by ICP/MS. For TraceMetal acids, a typical lot analysis is given in the catalog. Impurity levels in an actual lot received may vary from amounts listed	Primarily used in digestion of samples prior to instrument (ICP) analysis	Provided with each shipment
Histology	Solvents and products that are specially prepared for use in the histology laboratory setting. Solvents are filtered for tissue processing applications.	Tissue processing, clinical or histology procedures	Available on request
Laboratory and Technical	Chemicals of reasonable quality and purity for use in procedures where no official standards are required.	Manufacturing and general laboratory use	Available on request

Fisher Chemical* High Purity Solvents Selection Guide

Solvent	Available Grades [†]	UV Cutoff (nm)	Boiling Point (°C)	Density (g/mL, 25°C)	Refractive Index (25°C)	Melting Point (°C)	Polarity Index (P')	Eluotropic Value on Silica (ε°)	Viscosity (cP, 20°C)	Flash Point (°C)	Mol. Wt.
Acetone	GOHPSEF	330	56.1	0.7857	1.3568	-94.3	5.1	0.53	0.36	20	58.08
Acetonitrile	OHP	190	81.6	0.7780	1.3415	-50.0	5.8	0.52	0.36	2	41.05
1-Butanol	HS	215	117.7	0.8098	1.3972	-88.6	3.9	----	2.98	35	74.12
Chloroform	OMHPS	245	61.7	1.4840	1.4445	-63.3	4.1	0.26	0.58	none	119.38
Cyclohexane	HPS	202	80.7	0.7740	1.4247	-6.5	0.2	0.03	0.90	-20	84.16
N,N-Dimethylformamide	S	268	153.0	0.9440	1.4280	-61.0	6.4	----	0.92	58	73.09
Dimethyl Sulfoxide	SH	262	189.0	1.1014	1.4783	18.5	7.2	----	2.24	87.8	78.13
Ethyl Acetate	OHPSPH	255	77.1	0.8940	1.3695	-83.9	4.4	0.38	0.45	-4	88.11
Ethyl Ether	HPS	218	34.6	0.7134	1.3500	-116.3	2.8	0.43	0.24	-45	74.12
Glycerol	MSPH	205	290.0	1.2613	1.4746	18.2	----	----	----	193	92.09
Heptane	HSO	197	98.4	0.6838	1.3855	-90.6	0.2	0.01	0.40	-4	100.20
Hexanes	GOHPSPH	195	69.0	0.6630	1.3759	-95.3	0.1	0.01	0.31	-23	86.18
Isooctane	OHPS	205	99.2	0.6919	1.3895	109.5	0.1	0.01	0.50	28	114.23
Methanol	GPtOHPSLw	205	64.7	0.7915	1.3288	-97.8	5.1	0.73	0.55	12	32.04
Methylene Chloride	GOHPS	233	39.5	1.3180	1.4215	-96.7	3.1	0.32	0.30	N/A	84.93
N-Methylpyrrolidinone	M	275	202.2	1.03	1.469	-24.4	----	----	1.67	95	99.13
Pentane	HPS	190	36.1	0.6264	1.3555	-129.7	0.0	0.00	0.22	-49	72.15
Petroleum Ether	OP	----	35-60	0.6400	1.3610	----	0.1	----	----	-18	----
2-Propanol	OHPSLwPh	205	82.3	0.7855	1.3772	-90.0	3.9	0.63	2.40	-12	60.10
Tetrahydrofuran	OHS	210	66.1	0.8892	1.4060	-108.3	4.0	0.35	0.55	-14	72.11
Toluene	OHSSc	285	110.6	0.8660	1.4940	-95.0	2.4	0.22	0.59	-4	92.14
Water	OH	----	100.0	0.9982	1.3330	0.0	10.2	----	1.00	N/A	18.02

[†]G = GC *Resolv** Pt = Purge & Trap O = Optima* M = Molecular Biology Grade H = HPLC Grade P = Pesticide Grade S = Spectranalyzed* E = Electronic Grade
Sc = Scintanalyzed* Lw = Low water

Ph = Pharmacopeia, one or more of the following: USP, NF, FCC, EP, BP, JP

Liquid Chromatography/ Mass Spectrometry (LC/MS)

The joining of liquid chromatography (LC) with mass spectrometry (MS) has become an indispensable tool for various fields of research. The value of LC/MS derives from its ability to combine separation chemistry with selective mass ion detection. As instrumentation advances lead to ever-lower analyte detection limits, it is crucial for the chromatographer to consider the level of purity when selecting appropriate solvents for use in the LC/MS mobile phase. The Fisher Chemical* product line offers **Optima* LC/MS Grade** Solvents designed to meet the required purity level of advanced LC/MS systems.

Key Features:

- Lower contamination from plasticizer peaks
- Exceptionally low metal ion content
- Fewer background peaks
- Higher signal intensity low LC/ D UV response



Thermo Scientific* Q Exactive* Benchtop Orbitrap* LC-MS/MS

offers "Quanfirmation" – identification, quantitation, confirmation.

For more product information, go to:
www.thermoscientific.com/ms

Product Description	Grade	Size	Packaging	Cat. No.
OPTIMA LC/MS				
SOLVENTS				
Acetonitrile	Optima LC/MS	500mL	Amber Glass	A955-500
Acetonitrile	Optima LC/MS	1L	Amber Glass	A955-1
Acetonitrile	Optima LC/MS	2.5L	Amber Glass	A955-212
Acetonitrile	Optima LC/MS	4L	Amber Glass	A955-4
Methanol	Optima LC/MS	500mL	Amber Glass	A456-500
Methanol	Optima LC/MS	1L	Amber Glass	A456-1
Methanol	Optima LC/MS	2.5L	Amber Glass	A456-212
Methanol	Optima LC/MS	4L	Amber Glass	A456-4
2-Propanol	Optima LC/MS	500mL	Amber Glass	A461-500
2-Propanol	Optima LC/MS	1L	Amber Glass	A461-1
2-Propanol	Optima LC/MS	2.5L	Amber Glass	A461-212
2-Propanol	Optima LC/MS	4L	Amber Glass	A461-4
Water	Optima LC/MS	500mL	Amber Glass	W6-500
Water	Optima LC/MS	1L	Amber Glass	W6-1
Water	Optima LC/MS	2.5L	Amber Glass	W6-212
Water	Optima LC/MS	4L	Amber Glass	W6-4
BLENDS				
0.1% FA in Water	Optima LC/MS	500mL	Amber Glass	LS118-500
0.1% FA in Water	Optima LC/MS	1L	Amber Glass	LS118-1
0.1% FA in Water	Optima LC/MS	2.5L	Amber Glass	LS118-212
0.1% FA in Water	Optima LC/MS	4L	Amber Glass	LS118-4
0.1% TFA in Water	Optima LC/MS	500mL	Amber Glass	LS119-500
0.1% TFA in Water	Optima LC/MS	1L	Amber Glass	LS119-1
0.1% TFA in Water	Optima LC/MS	2.5L	Amber Glass	LS119-212
0.1% TFA in Water	Optima LC/MS	4L	Amber Glass	LS119-4
0.1% FA in Acetonitrile	Optima LC/MS	500mL	Amber Glass	LS120-500
0.1% FA in Acetonitrile	Optima LC/MS	1L	Amber Glass	LS120-1
0.1% FA in Acetonitrile	Optima LC/MS	2.5L	Amber Glass	LS120-212
0.1% FA in Acetonitrile	Optima LC/MS	4L	Amber Glass	LS120-4
0.1% TFA in Acetonitrile	Optima LC/MS	500mL	Amber Glass	LS121-500
0.1% TFA in Acetonitrile	Optima LC/MS	1L	Amber Glass	LS121-1
0.1% TFA in Acetonitrile	Optima LC/MS	2.5L	Amber Glass	LS121-212
0.1% TFA in Acetonitrile	Optima LC/MS	4L	Amber Glass	LS121-4
ADDITIVES				
Formic Acid	Optima LC/MS	50mL	Poly Bottle	A117-50
Formic Acid	Optima LC/MS	10x1mL	Ampules	A117-10X1AMP
Formic Acid	Optima LC/MS	1mL	Ampule	A117-1AMP
Formic Acid	Optima LC/MS	0.5mL	Ampule	A117-05AMP
Formic Acid	Optima LC/MS	2mL	Ampule	A117-2AMP
Trifluoroacetic Acid	Optima LC/MS	50mL	Amber Glass	A116-50
Trifluoroacetic Acid	Optima LC/MS	10x1mL	Ampules	A116-10X1AMP
Trifluoroacetic Acid	Optima LC/MS	1mL	Ampule	A116-1AMP
Trifluoroacetic Acid	Optima LC/MS	0.5mL	Ampule	A116-05AMP
Trifluoroacetic Acid	Optima LC/MS	2mL	Ampule	A116-2AMP
Acetic Acid	Optima LC/MS	50mL	Poly Bottle	A113-50
Acetic Acid	Optima LC/MS	1mL	Ampule	A113-1AMP
Acetic Acid	Optima LC/MS	10 x 1mL	Ampules	A113-10X1AMP
Ammonium Acetate	Optima LC/MS	50g	Amber Glass	A114-50
Ammonium Formate	Optima LC/MS	50g	Amber Glass	A115-50

High-Performance Liquid Chromatography (HPLC)

HPLC is a technique that can separate a mixture of compounds and is used in biochemistry and analytical chemistry to identify, quantify and purify the individual components of the mixture. Fisher Chemical* Optima* Grade Solvents meet extremely high purity (ppm) levels. They are manufactured for use when contaminant-free performance is essential — HPLC, GC, plasma/ICP, spectrophotometry, environmental testing and other analytical applications.

Key Features:

- Contaminant-free to ppb and ppm levels
- Supplied in specially cleaned bottles
- Certificate of Analysis provided with each shipment



UltiMate 3000 LC System

Provides excellent chromatographic performance while maintaining easy and reliable operation. Select from a wide variety of systems for applications at all scales.

For more product information, go to:

www.thermoscientific.com/chromatography

Product Description	Grade	Size	Packaging	Cat. No.
OPTIMA				
Acetone	Optima	1L	Amber Glass	A929-1
Acetone	Optima	4L	Amber Glass	A929-4
Acetone	Optima	4L	Glass/Safe-Cote*	A929SK-4
Acetonitrile	Optima	1L	Amber Glass	A996-1
Acetonitrile	Optima	4L	Amber Glass	A996-4
Acetonitrile	Optima	4L	Amber Glass/Safe-Cote	A996SK4
Chloroform (Approx. 50ppm Amylene as Preservative)	Optima	4L	Amber Glass	C297-4
Ethyl Acetate	Optima	4L	Amber Glass	E196-4
Ethyl Acetate	Optima	4L	Safe-Cote	E196SK-4
Hexanes	Optima	1L	Amber Glass	H303-1
Hexanes	Optima	4L	Amber Glass	H303-4
Hexanes	Optima	4L	Amber Glass/Safe-Cote	H303SK-4
n-Hexane,95%,	Optima	1L	Amber Glass	H306-1
n-Hexane,95%,	Optima	4L	Amber Glass	H306-4
n-Hexane,95%,	Optima	4L	Amber Glass/Safe-Cote	H306SK-4
Methanol	Optima	1L	Amber Glass	A454-1
Methanol	Optima	1L	Amber Glass	A454-4
Methanol	Optima	4L	Amber Glass/Safe-Cote	A454SK-4
Methylene Chloride	Optima	2.5L	Amber Glass	D151-1
Methylene Chloride	Optima	4L	Amber Glass	D151-4
Methylene Chloride	Optima	4L	Amber Glass/Safe-Cote	D151SK-4
Petroleum Ether	Optima	4L	Amber Glass	E120-4
Petroleum Ether	Optima	4L	Amber Glass/Safe-Cote	E120SK-4
2-Propanol	Optima	4L	Amber Glass	A464-4
2-Propanol	Optima	4L	Amber Glass/Safe-Cote	A464SK-4
Tetrahydrofuran	Optima	1L	Amber Glass	T427-1
Tetrahydrofuran	Optima	4L	Amber Glass/Safe-Cote	T427SK-4
Tetrahydrofuran	Optima	4L	Amber Glass/Safe-Cote	T427SK-4
Toluene	Optima	4L	Amber Glass	T291-4
Toluene	Optima	4L	Amber Glass/Safe-Cote	T291SK-4
Water	Optima	4L	Amber Glass	W7-4
Water	Optima	4L	Amber Glass/Safe-Cote	W7SK-4
BLENDS				
Acetonitrile with 0.01% Trifluoroacetic Acid (v/v)	HPLC	4L	Amber Glass	HB98104
Acetonitrile with 0.035% Formic Acid (v/v)	HPLC	4L	Amber Glass	HB98214
Acetonitrile with 0.035% Trifluoroacetic Acid (v/v)	HPLC	4L	Amber Glass	HB98114
Acetonitrile with 0.05% Formic Acid (v/v)	HPLC	4L	Amber Glass	HB98224
Acetonitrile with 0.05% Trifluoroacetic Acid (v/v)	HPLC	4L	Amber Glass	HB98124
Acetonitrile with 0.1% Formic Acid (v/v)	HPLC	4L	Amber Glass	HB98234
Acetonitrile with 0.1% Formic Acid and 0.01% Trifluoroacetic Acid (v/v)	HPLC	4L	Amber Glass	HB98344
Acetonitrile with 0.1% Trifluoroacetic Acid (v/v)	HPLC	4L	Amber Glass	HB98134

Gas Chromatography (GC)

Gas Chromatography is a common type of technique used to separate and analyze compounds that can be vaporized without decomposition. GC *Resolv** Grade Solvents are manufactured for GC applications, and are free of contaminants to the ppb level. Chromatogram on GC *Resolv* Grade Solvents is available on request. Pesticide Grade solvents are for use in the analysis of pesticide residue on GC with electron capture detector (ECD).

Product Description	Grade	Size	Packaging	Cat. No.
GC RESOLV/ PESTICIDE				
Acetone	GC <i>Resolv</i>	4L	Amber Glass	A928-4
n-Hexane	GC <i>Resolv</i>	4L	Amber Glass	H307-4
Methanol	GC <i>Resolv</i>	4L	Amber Glass	A457-4
Methylene Chloride	GC <i>Resolv</i>	4L	Amber Glass	D154-4
Acetone	Pesticide	4L	Amber Glass	A40-4
Chloroform (With Amylene Preservative)	Pesticide	4L	Amber Glass	C603-4
Cyclohexane	Pesticide	4L	Amber Glass	C553-4
Ethyl Acetate	Pesticide	4L	Amber Glass	E191-4
Ethyl Ether	Pesticide	4L	Amber Glass	E199-4
Hexanes	Pesticide	4L	Amber Glass	H300-4
Isooctane	Pesticide	4L	Amber Glass	O297-4
Methanol	Pesticide	4L	Amber Glass	A450-4
Methylene Chloride	Pesticide	4L	Amber Glass	D142-4
Pentane	Pesticide	4L	Amber Glass	P400-4
Petroleum Ether	Pesticide	4L	Amber Glass	P480-4
2-Propanol	Pesticide	4L	Amber Glass	A519-4



TRACE GC Ultra Gas Chromatographs

With the most complete range of proprietary inlets, sensitive detection systems, smart accessories, and ancillary devices, the Ultra platform features unique technologies that raise the standards of speed, sensitivity and separation in gas chromatography.

For more product information, go to:

www.thermoscientific.com/chromatography

ACS

ACS Grade Solvents meet or exceed the specifications set by the American Chemical Society. ACS Grade Solvents are suitable for analytical applications requiring tight specifications.

There are two grades among this classification:

1. **Certified ACS Plus Grade** Solvents meet or exceed the latest specifications of the ACS and are suitable for applications with tighter metal specifications
2. **Certified ACS Grade** Solvents meet or exceed the latest ACS Specifications

Product Description	Grade	Size	Packaging	Cat. No.
ACS				
1,2-Dichloroethane	Certified ACS	500mL	Amber Glass	E175-500
1,2-Dichloroethane	Certified ACS	4L	Amber Glass	E175-4
1,2-Dichloroethane	Certified ACS	20L	Steel Pail	E175-20
1-Butanol	Certified ACS	500mL	Amber Glass	A399-500
1-Butanol	Certified ACS	1L	Amber Glass	A399-1
1-Butanol	Certified ACS	4L	Amber Glass	A399-4
1-Butanol	Certified ACS	4L	SafeTin	A399S-4
2-Propanol	Certified ACS Plus	500mL	Amber Glass	A416-500
2-Propanol	Certified ACS Plus	1L	Amber Glass	A416-1
2-Propanol	Certified ACS Plus	4L	Amber Glass	A416-4
2-Propanol	Certified ACS Plus	4L	Poly Bottle	A416P-4
2-Propanol	Certified ACS Plus	4L	SafeTin	A416S-4
2-Propanol	Certified ACS Plus	4L	Amber Glass/Safe-Cote*	A416SK-4
Acetone	Certified ACS	500mL	Amber Glass	A18-500
Acetone	Certified ACS	1L	Amber Glass	A18-1
Acetone	Certified ACS	4L	Amber Glass	A18-4
Acetone	Certified ACS	4L	Poly Bottle	A18P-4
Acetone	Certified ACS	4L	SafeTin	A18S-4
Acetone	Certified ACS	4L	Amber Glass/Safe-Cote	A18SK-4
Acetonitrile	Certified ACS	1L	Amber Glass	A21-1
Acetonitrile	Certified ACS	4L	Amber Glass	A21-4
Chloroform (Approx. 0.75% Ethanol as Preservative)	Certified ACS	500mL	Amber Glass	C298-500
Chloroform (Approx. 0.75% Ethanol as Preservative)	Certified ACS	1L	Amber Glass	C298-1
Chloroform (Approx. 0.75% Ethanol as Preservative)	Certified ACS	4L	Amber Glass	C298-4
Chloroform (Approx. 0.75% Ethanol as Preservative)	Certified ACS	4L	SafeTin	C298S-4
Chloroform (Approx. 0.75% Ethanol as Preservative)	Certified ACS	4L	Amber Glass/Safe-Cote	C298SK-4
Ethyl Acetate	Certified ACS	500mL	Amber Glass	E145-500
Ethyl Acetate	Certified ACS	1L	Amber Glass	E145-1
Ethyl Acetate	Certified ACS	4L	Amber Glass	E145-4
Ethyl Acetate	Certified ACS	4L	SafeTin	E145S-4
Ethyl Acetate	Certified ACS	4L	Amber Glass/Safe-Cote	E145SK-4
Hexanes	Certified ACS	500mL	Amber Glass	H292-500
Hexanes	Certified ACS	1L	Amber Glass	H292-1
Hexanes	Certified ACS	4L	Amber Glass	H292-4
Hexanes	Certified ACS	4L	Amber Glass/Safe-Cote	H292SK-4

ACS (continued)

Product Description	Grade	Size	Packaging	Cat. No.
ACS				
Methanol	Certified ACS	500mL	Amber Glass	A412-500
Methanol	Certified ACS	1L	Amber Glass	A412-1
Methanol	Certified ACS	4L	Amber Glass	A412-4
Methanol	Certified ACS	4L	Poly Bottle	A412P-4
Methanol	Certified ACS	4L	Amber Glass/Safe-Cote*	A412SK-4
Methylene Chloride (Stabilized)	Certified ACS	500mL	Amber Glass	D37-500
Methylene Chloride (Stabilized)	Certified ACS	1L	Amber Glass	D37-1
Methylene Chloride (Stabilized)	Certified ACS	4L	Amber Glass	D37-4
Methylene Chloride (Stabilized)	Certified ACS	4L	Amber Glass/Safe-Cote	D37SK-4
N,N-Dimethylformamide	Certified ACS	500mL	Amber Glass	D119-500
N,N-Dimethylformamide	Certified ACS	1L	Amber Glass	D119-1
N,N-Dimethylformamide	Certified ACS	4L	Amber Glass	D119-4
N,N-Dimethylformamide	Certified ACS	4L	SafeTin	D119S-4
Toluene	Certified ACS	500mL	Amber Glass	T324-500
Toluene	Certified ACS	1L	Amber Glass	T324-1
Toluene	Certified ACS	4L	Amber Glass	T324-4
Toluene	Certified ACS	4L	SafeTin	T324S-4
Toluene	Certified ACS	4L	Amber Glass/Safe-Cote	T324SK-4
Xylene	Certified ACS	500mL	Amber Glass	X5-500
Xylene	Certified ACS	1L	Amber Glass	X5-1
Xylene	Certified ACS	1 gal.	Poly Bottle	X5P-1GAL
Xylene	Certified ACS	4L	Amber Glass	X5-4
Xylene	Certified ACS	4L	SafeTin	X5S-4
Xylene	Certified ACS	4L	Amber Glass/Safe-Cote	X5SK-4

Certified

Certified Grade Solvents are suitable for general analytical procedures, for which the purity standard is established by Fisher Scientific. Purity is guaranteed to meet published maximum limits of impurities.

Product Description	Grade	Size	Packaging	Cat. No.
CERTIFIED				
Tetrahydrofuran	Certified	500mL	Amber Glass	T397-500
Tetrahydrofuran	Certified	1L	Amber Glass	T397-1
Tetrahydrofuran	Certified	4L	Amber Glass	T397-4
Tetrahydrofuran	Certified	4L	Amber Glass/Safe-Cote	T397SK-4

Acids

The Fisher Chemical* portfolio of high-purity acids and bases comprises more than 400 products suited to a wide range of laboratory applications, including trace element analysis in sample preparation, HPLC and spectrophotometry.

Manufactured in ISO 9001:2008-certified facilities, all Fisher Chemical acids and bases undergo rigorous quality assurance and testing measures that ensure excellent lot-to-lot and bottle-to-bottle consistency.

Select the grade that is best suited to your application:

Fisher Chemical Acid Purity Grades

Grade	Application	Features
Optima*	HPLC, GC, AA, plasma/ICP/ICP-MS, spectrophotometry, pesticide residue analysis	Acids of extremely high purity. Lowest metal content of any commercially available acid. Acids are analyzed for 65 metals by ICP/MS; impurity levels in ppt. Packaged in specially manufactured and cleaned PFA and FEP fluoropolymer bottles (Optima Ammonium Hydroxide is packaged in specially leached HDPE bottles.)
TraceMetal*	Primarily used in digestion of samples prior to instrument (ICP) analysis	Acids manufactured under proprietary distillation to achieve low metal contamination measurable in ppm to ppb range. Each lot is analyzed for more than 55 metals by ICP/MS. Now packaged in poly bottles for safer and easier handling in the laboratory and limited or no breakage during transportation.
HPLC	HPLC and spectrophotometry	Specifically formulated for use in liquid chromatography, Fisher Chemical HPLC Grade Acids ensure minimum baseline drift and maximum retention reproducibility. Meet all ACS specifications.
Certified ACS Plus	Analytical applications with tighter metal specifications	Acids that, in addition to meeting or exceeding the latest specifications of the ACS, are analyzed for more than 16 metals.
Certified ACS	Analytical applications requiring tight specifications	Reagent chemicals that meet or exceed the latest ACS specifications. Actual lot analysis on label.
USP/NF/FCC/EP/BP/JP	Food and drug laboratories, biological testing	Solvents meet or surpass specifications of the United States Pharmacopeia (USP), the National Formulary (NF), the Food Chemicals Codex (FCC), the European Pharmacopeia (EP), the British Pharmacopeia (BP) and/or the Japanese Pharmacopeia (JP).
Reagent, Laboratory and Technical	Manufacturing and general laboratory use	Chemicals of reasonable quality and purity for use in procedures where no official standards are required.

Evolution 600 UV-Vis Spectrophotometer

Thermo Scientific Evolution 600 incorporates a high-performance optical design, versatile VISION* software packages, and the highest quality accessories for the most demanding materials science and research applications.

For more product information, go to:
www.thermoscientific.com/uv-vis



High-Purity

Optima* Grade Acids offer the lowest metal content and the highest purity. They are the perfect choice for environmental testing, plasma analysis, electronic research and other exacting procedures requiring the utmost purity, precision and accuracy.

Key Features:

- Analyzed for up to 65 metals to ppt level by ICP-MS to provide unequaled purity (Certificate of Analysis supplied with each bottle)
- Packaged in specially manufactured and cleaned PFA and FEP fluoropolymer bottles (Optima Ammonium Hydroxide is packaged in specially leached HDPE bottles)
- Produced and packaged in Class 100 cleanroom environment to ensure maximum purity
- Established product lines provide excellent lot-to-lot consistency for reproducible results

Applications: HPLC, GC, plasma/ICP, spectrophotometry, pesticide residue analysis

Product Description	Grade	Size	Packaging	Cat. No.
OPTIMA ACIDS				
Acetic Acid, Glacial	Optima	250mL	PFA bottle	A465250
Acetic Acid, Glacial	Optima	500mL	PFA bottle	A465500
Acetic Acid, Glacial	Optima	1L	PFA bottle	A4651
Ammonium Hydroxide	Optima	250mL	PFA bottle	A470250
Ammonium Hydroxide	Optima	500mL	PFA bottle	A470500
Ammonium Hydroxide	Optima	1L	PFA bottle	A4701
Hydrobromic Acid	Optima	1L	PFA bottle	A4711
Hydrobromic Acid	Optima	2L	FEP bottle	A4712
Hydrobromic Acid	Optima	250mL	PFA bottle	A471250
Hydrobromic Acid	Optima	500mL	PFA bottle	A471500
Hydrochloric Acid	Optima	1L	PFA bottle	A4661
Hydrochloric Acid	Optima	2L	FEP bottle	A4662
Hydrochloric Acid	Optima	250mL	PFA bottle	A466250
Hydrochloric Acid	Optima	500mL	PFA bottle	A466500
Hydrofluoric Acid	Optima	1L	PFA bottle	A4631
Hydrofluoric Acid	Optima	2L	FEP bottle	A4632
Hydrofluoric Acid	Optima	250mL	PFA bottle	A463250
Hydrofluoric Acid	Optima	500mL	PFA bottle	A463500
Hydrogen Peroxide	Optima	500mL	FEP bottle	P170500
Nitric Acid	Optima	1L	PFA bottle	A4671
Nitric Acid	Optima	2L	FEP bottle	A4672
Nitric Acid	Optima	250mL	PFA bottle	A467250 [†]
Nitric Acid	Optima	500mL	PFA bottle	A467500
Perchloric Acid	Optima	1L	PFA bottle	A4691
Perchloric Acid	Optima	250mL	PFA bottle	A469250
Perchloric Acid	Optima	500mL	PFA bottle	A469500
Sulfuric Acid	Optima	2L	FEP bottle	A4682
Sulfuric Acid	Optima	250mL	PFA bottle	A468250
Sulfuric Acid	Optima	500mL	PFA bottle	A468500
Sulfuric Acid	Optima	1L	PFA bottle	A4681

[†]Products might not be available in all regions. Contact your local sales representative for details.



XSERIES 2 ICP-MS

The Thermo Scientific* XSERIES 2 ICP-MS offers outstanding productivity in a quadrupole ICP-MS for both routine and high performance analytical work.

For more product information, go to:

www.thermoscientific.com/ms



iCAP 6500 ICP Spectrometer

Thermo Scientific* iCAP 6500 ICP Spectrometer is small in size, yet big on performance, offering the best detection capability of any optical ICP, with enhanced productivity and the optimum in application flexibility.

For more product information, go to: www.thermoscientific.com/ms

TraceMetal* Analysis

TraceMetal Grade Acids are used primarily in digestion of samples prior to instrument (ICP) analysis. Each lot is analyzed for more than 55 metals by ECP-MS. Proprietary distillation achieves metal concentrations in the sub-ppb range.

Selected Fisher Chemical* TraceMetal Grade Acids and Bases are now available in a *poly* bottle. These bottles consist of a new high-density polyethylene made with proprietary resin that provides up to 80% less metallic extractables compared to glass. This new packaging material for TraceMetal Acids and Bases provides more convenience and stability during transport, plus increased resistance to breakage compared to conventional glass packaging.

Packaging Advantages

- Poly bottles provide reduced weight for safer and easier handling in the laboratory and limited or no breakage during transportation
- FisherLOCK* tamper evidence cap for chemical safety and drip lip pour feature for safer pouring in the laboratory
- Fully recyclable package and poly bottle for reduced waste disposable
- Poly bottles take up less space than glass bottles, providing better chemical storage space utilization in the laboratory
- With the poly bottle, there is no longer a need to remove the PVC coating on glass to allow for disposal in landfills or for recycling glass



Product Description	Grade	Size	Packaging	Cat. No.
TraceMetal				
Acetic Acid, Glacial	TraceMetal	2.5L	Poly Bottle	A507P212
Acetic Acid, Glacial	TraceMetal	500mL	Poly Bottle	A507P500
Ammonium Hydroxide	TraceMetal	500mL	Poly Bottle	A512P500
Hydrochloric Acid	TraceMetal	2.5L	Poly Bottle	A508P212
Hydrochloric Acid	TraceMetal	500mL	Poly Bottle	A508500
Hydrofluoric Acid	TraceMetal	4L	HDPE Bottle	A5134
Hydrofluoric Acid	TraceMetal	500mL	LDPE Bottle	A513500
Nitric Acid	TraceMetal	2.5L	Poly Bottle	A509P212
Nitric Acid	TraceMetal	500mL	Poly Bottle	A509P500
Perchloric Acid	TraceMetal	2.5L	Poly Bottle	A511P212
Perchloric Acid	TraceMetal	500mL	Poly Bottle	A511P500
Sulfuric Acid	TraceMetal	2.5L	Poly Bottle	A510P212
Sulfuric Acid	TraceMetal	500mL	Poly Bottle	A510P500
Water	Ultra Trace Elemental Analysis	500mL	LDPE Bottle	W9500
Water	Ultra Trace Elemental Analysis	1L	LDPE Bottle	W91
Water	Ultra Trace Elemental Analysis	2L	LDPE Bottle	W92

ACS

ACS Grade Acids meet or exceed the specifications set by the American Chemical Society. ACS Grade Acids are suitable for analytical applications requiring tight specifications.

There are two grades among this classification:

1. **Certified ACS Plus Grade:** Acids that, in addition to meeting or exceeding the latest specifications of the ACS, are analyzed for more than 16 metals; actual lot analysis on label
2. **Certified ACS Grade:** Acids meet or exceed the latest ACS Specifications

Product Description	Grade	Size	Packaging	Cat. No.
ACS ACIDS				
Acetic Acid, Glacial	Certified ACS	6 x 2.5L	Glass	A38C212 [†]
Acetic Acid, Glacial	Certified ACS	2.5L	Safe-Cote*/Glass	A38S212 [†]
Acetic Acid, Glacial	Certified ACS	2.5L	Glass	A38212
Acetic Acid, Glacial	Certified ACS	500mL	Glass	A38500
Acetic Acid, Glacial	Certified ACS	500mL	Safe-Cote/Glass	A38S500
Acetic Acid, Glacial	Certified ACS	2.5L	Safe-Cote/Glass	A38SI212
Ammonium Hydroxide	Certified ACS Plus	6 x 2.5L	Glass	A669C212 [†]
Ammonium Hydroxide	Certified ACS Plus	2.5L	Glass	A669212
Ammonium Hydroxide	Certified ACS Plus	500mL	Safe-Cote/Glass	A669S500
Ammonium Hydroxide	Certified ACS Plus	2.5L	Safe-Cote/Glass	A669S212
Ammonium Hydroxide	Certified ACS Plus	500mL	Poly Bottle	A669500
Formic Acid	Certified ACS	4L	SafeTin	A118P4
Formic Acid	Certified ACS	500mL	Poly Bottle	A118P500
Formic Acid	Certified ACS	100mL	Poly Bottle	A118P100
Hydrochloric Acid	Certified ACS Plus	6 x 2.5L	Glass	A144C212 [†]
Hydrochloric Acid	Certified ACS Plus	2.5L	Safe-Cote/Glass	A144S212
Hydrochloric Acid	Certified ACS Plus	500mL	Glass	A144500
Hydrochloric Acid	Certified ACS Plus	2.5L	Glass	A144212
Hydrochloric Acid	Certified ACS Plus	500mL	Safe-Cote/Glass	A144S500
Hydrochloric Acid	Certified ACS Plus	2.5L	Safe-Cote/Glass	A144SI212 [†]
Nitric Acid	Certified ACS Plus	500mL	Glass	A200500 [†]
Nitric Acid	Certified ACS Plus	500mL	Safe-Cote/Glass	A200S500 [†]
Nitric Acid	Certified ACS Plus	2.5L	Safe-Cote/Glass	A200SI212
Nitric Acid	Certified ACS Plus	2.5L	Safe-Cote/Glass	A200S212
Nitric Acid	Certified ACS Plus	2.5L	Glass	A200212
Nitric Acid	Certified ACS Plus	6 x2.5L	Glass	A200C212
o-Phosphoric Acid 85%	Certified ACS	2.5L	Glass	A242212
o-Phosphoric Acid 85%	Certified ACS	2.5L	Safe-Cote/Glass	A242SK212
o-Phosphoric Acid 85%	Certified ACS	500mL	Amber Glass	A242500
o-Phosphoric Acid 85%	Certified ACS	4L	Amber Glass	A2424
o-Phosphoric Acid 85%	Certified ACS	1L	Amber Glass	A2421
Potassium Hydroxide	Certified ACS	500g	Poly Bottle	P250500
Potassium Hydroxide	Certified ACS	10kg	Poly Pail	P25010
Potassium Hydroxide	Certified ACS	1kg	Poly Bottle	P2501
Potassium Hydroxide	Certified ACS	3kg	Poly Bottle	P2503
Sodium Hydroxide	Certified ACS	1kg	Poly Bottle	S3181
Sodium Hydroxide	Certified ACS	100g	Poly Bottle	S318100

[†]Products might not be available in all regions. Contact your local sales representative for details.

ACS (continued)

Product Description	Grade	Size	Packaging	Cat. No.
ACS ACIDS				
Sodium Hydroxide	Certified ACS	3kg	Poly Bottle	S3183
Sodium Hydroxide	Certified ACS	500g	Poly Bottle	S318500
Sodium Hydroxide	Certified ACS	5kg	Poly Pail	S3185
Sulfuric Acid	Certified ACS Plus	2.5L	Glass	A300212
Sulfuric Acid	Certified ACS Plus	500mL	Glass	A300500
Sulfuric Acid	Certified ACS Plus	2.5L	Safe-Cote*/Glass	A300SI212
Sulfuric Acid	Certified ACS Plus	500mL	Safe-Cote/Glass	A300S500
Sulfuric Acid	Certified ACS Plus	6 x 2.5L	Glass	A300C212†
Sulfuric Acid	Certified ACS Plus	2.5L	Safe-Cote/Glass	A300S212
Trichloroacetic Acid	Certified ACS	500g	Amber Glass	A322500
Trichloroacetic Acid	Certified ACS	100g	Amber Glass	A322100
Trichloroacetic Acid	ACS Certified	3kg	Glass	A3223

*Products might not be available in all regions. Contact your local sales representative for details.

Certified

Certified Grade Acids are suitable for general analytical procedures, for which the purity standard is established by Fisher Scientific. Purity is guaranteed to meet published maximum limits of impurities.

Product Description	Grade	Size	Packaging	Cat. No.
CERTIFIED				
Hydrochloric Acid Solution, 1N	Certified	4L	Poly Bottle	SA484
Hydrochloric Acid Solution, 1N	Certified	500mL	Poly Bottle	SA48500
Hydrochloric Acid Solution, 1N	Certified	1L	Poly Bottle	SA481
Hydrochloric Acid Solution, 2N	Certified	500mL	Poly Bottle	SA431500
Hydrochloric Acid Solution, 6N	Certified	1L	Amber Glass	SA561
Hydrochloric Acid Solution, 6N	Certified	4L	Amber Glass	SA564
Hydrochloric Acid Solution, 6N	Certified	500mL	Amber Glass	SA56500
Nitric Acid Concentrate	Certified	100mL	Amber Glass	SA95
Nitric Acid Solution, 0.1N (N/10)	Certified	1L	Amber Glass	SA941
Perchloric Acid Solution, N/10 (0.1N)	Certified	500mL	Amber Glass	SP339500
Potassium Hydroxide Solution, 0.1N (In Isopropanol)	Certified	500mL	Amber Glass	ST110500
Propionic Acid	Certified	500mL	Amber Glass	A258500
Sulfuric Acid Solution (0.02N (N/50))	Certified	4L	Poly Bottle	SA2264
Sulfuric Acid Solution (0.02N (N/50))Certified	Certified	1L	Poly Bottle	SA2261

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
1a	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
H 1 1.008 Hydrogen	He 2 4.003 Helium	Li 3 6.941 Lithium	Be 4 9.012 Beryllium	B 5 10.811 Boron	C 6 12.011 Carbon	N 7 14.007 Nitrogen	O 8 15.999 Oxygen	F 9 18.998 Fluorine	Ne 10 20.180 Neon	Na 11 22.990 Sodium	Mg 12 24.305 Magnesium	Al 13 26.982 Aluminum	Si 14 28.086 Silicon	P 15 30.974 Phosphorous	S 16 32.066 Sulfur	Cl 17 35.453 Chlorine	Ar 18 39.948 Argon
K 19 39.098 Potassium	Ca 20 40.078 Calcium	Sc 21 44.956 Scandium	Ti 22 47.88 Titanium	V 23 50.942 Vanadium	Cr 24 51.996 Chromium	Mn 25 54.938 Manganese	Fe 26 55.847 Iron	Co 27 58.933 Cobalt	Ni 28 58.69 Nickel	Cu 29 63.546 Copper	Zn 30 65.39 Zinc	Ga 31 69.723 Gallium	Ge 32 72.61 Germanium	As 33 74.922 Arsenic	Se 34 78.96 Selenium	Br 35 79.904 Bromine	Kr 36 83.80 Krypton
Rb 37 85.468 Rubidium	Sr 38 87.62 Strontium	Y 39 88.906 Yttrium	Zr 40 91.224 Zirconium	Nb 41 92.906 Niobium	Mo 42 95.94 Molybdenum	Tc 43 (98) Technetium	Ru 44 101.07 Ruthenium	Rh 45 102.906 Rhodium	Pd 46 106.42 Palladium	Ag 47 107.868 Silver	Cd 48 112.411 Cadmium	In 49 114.82 Indium	Sn 50 118.71 Tin	Sb 51 121.75 Antimony	Te 52 127.60 Tellurium	I 53 126.905 Iodine	Xe 54 131.29 Xenon
Cs 55 132.905 Cesium	Ba 56 137.327 Barium	La 57 138.906 Lanthanum	Hf 72 178.49 Hafnium	Ta 73 180.948 Tantalum	W 74 183.85 Tungsten	Re 75 186.207 Rhenium	Os 76 190.2 Osmium	Ir 77 192.22 Iridium	Pt 78 195.08 Platinum	Au 79 196.967 Gold	Hg 80 200.59 Mercury	Tl 81 204.383 Thallium	Pb 82 207.2 Lead	Bi 83 208.980 Bismuth	Po 84 (209) Polonium	At 85 (210) Astatine	Rn 86 (222) Radon
Fr 87 (223) Francium	Ra 88 226.025 Radium	Ac 89 227.028 Actinium	Rf 104 (261) Rutherfordium	Db 105 (261) Dubnium	Sg 106 (263) Seaborgium	Bh 107 (262) Bohrium	Hs 108 (265) Hassium	Mt 109 (266) Meitnerium	Ds 110 (271) Darmstadtium	Rg 111 (272) Roentgenium	Cn 112 (277) Copernicium	Uut 113 (281) Ununtrium	Uuq 114 (285) Ununquadium	Uup 115 (287) Ununpentium	Uuh 116 (289) Ununhexium	Uus 117 (291) Ununseptium	Uuo 118 (293) Ununoctium

La 57 138.906 Lanthanum	Ce 58 140.115 Cerium	Pr 59 140.908 Praseodymium	Nd 60 144.24 Neodymium	Sm 62 150.36 Samarium	Eu 63 151.965 Europium	Gd 64 157.25 Gadolinium	Tb 65 158.925 Terbium	Dy 66 162.50 Dysprosium	Ho 67 164.93 Holmium	Er 68 167.26 Erbium	Tm 69 168.934 Thulium	Yb 70 173.04 Ytterbium	Lu 71 174.967 Lutetium
Ac 89 227.028 Actinium	Th 90 232.038 Thorium	Pa 91 231.036 Protactinium	U 92 238.029 Uranium	Pu 94 244 Plutonium	Am 95 243 Americium	Cm 96 247 Curium	Bk 97 247 Berkelium	Cf 98 251 Californium	Es 99 252 Einsteinium	Fm 100 257 Fermium	Md 101 258 Mendelevium	No 102 259 Nobelium	Lr 103 260 Lawrencium

Lanthanide Series

Actinide Series

KEY TO TABLE

- Alkali Metals
- Alkaline Earth Metals
- Transition Elements
- Non-metals
- Other Metals
- Halogens
- Inert Gases
- Lanthanide Series
- Actinide Series
- Aa - Solid
- Aa - Gas
- Aa - Liquid
- Aa - Synthetically Prepared

Symbol: **H**
Atomic Weight: 1.008
Element Name: Hydrogen

AMERICAS
Canada

Fisher Scientific Ltd.
112 Colonnade Road
Nepean
Ontario K2E 7L6 Canada
Tel: 800-234-7437
Fax: 800-463-2996
www.fishersci.ca

United States

Fisher Scientific
300 Industry Drive
Pittsburgh, PA
Post Code: 15275
Tel: 800-766-7000
Fax: 800-926-1166
www.fishersci.com

Fisher HealthCare
9999 Veterans Memorial Drive
Houston, TX
Post Code: 77038
Tel: 800.640.0640
Fax: 800.290.0290
www.fisherhealthcare.com

Fisher Scientific Global Export
Atlanta Office
3970 Johns Creek Court
Suite 500
Suwanee, GA
Post Code: 30024
Tel: 770-871-4725
Fax: 770-871-4726
www.fishersci.com

ASIA
China

Fisher Scientific China
Toll-Free Number: 400 881 5117
sales.china@thermofisher.com
www.fishersci.com.cn

Shanghai Corporate Office
2F, Building #7, No 27 Xin Jinqiao Rd
Shanghai 201206
Tel: +86 21 6192 5688
Fax: +86 21 6100 2128

Beijing Office
Units 702-715, 7th Floor
Tower West, Yonghe Plaza
No. 28 Andingmen
East Street
Beijing, China, 100007
Tel: +86 10 8419 3588
Fax: +86 10 8419 3580

Guangzhou Office
Room 2405-2406
Times Property Center
No. 410-412 Middle DongFend Rd
Guangzhou, China, 510030
Tel: +86 20 8314 5288
Fax: +86 20 3877 1941

India

Thermo Fisher Scientific India Pvt Ltd.
403-404, Delphi "B" Wing
Hiranandani Business Park
Powai, Mumbai 400 076
Tel: +91 (0) 22 6680 3000
Fax: +91 (0) 22 6680 3001/02
qfc.customer@thermofisher.com
www.fishersci.in

Japan

Fisher Scientific Japan Ltd.
NBF Ueno-Bldg. 10F
4-24-11 Higashi-Ueno
Taito-ku Tokyo 110-0015
Tel: +81 (0) 3 5826 1600
support@fishersci.co.jp
www.fishersci.co.jp

Korea

Fisher Scientific Korea
Sambu Bldg. 13F
676 Yoksam-Dong
Kangnam-Ku
Seoul 140-210 South Korea
Tel: +82 2 527 0300
Fax: +82 2 527 5662
sales@fishersci.co.kr
www.fishersci.co.kr

Malaysia

Fisher Scientific Malaysia Sdn Bhd
No. 3 Jalan Sepadu 25/123
Taman Perindustrian Axis Seksyen 25
40400 Shah Alam
Selangor Darul Ehsan Malaysia
Tel: +60 3-5122-8888
Fax: +60 3-5121-8899
enquiry.my@thermofisher.com
www.fishersci.com.my

Singapore

Fisher Scientific Singapore
Fisher Scientific Pte Ltd.
8 Pandan Crescent
LL4, #05-05 UE Tech Park
128464 Singapore
Tel: +65 6873 6006
Fax: +65 6873 5005
www.fishersci.com.sg

EUROPE
Austria

Fisher Scientific (Austria) GmbH
Rudolf von Alt-Platz 1
A-1030 Wien
Tel: +49 (0) 800 20 88 40
Fax: +49 (0) 800 20 66 90
info.austria@thermofisher.com
www.at.fishersci.com

Belgium

Fisher Scientific
BP 567, B-7500 Tournai 1
Tel: +32 056 260 260
Fax: +32 056 260 270
be.fisher@thermofisher.com
www.be.fishersci.com

Czech Republic

Fisher Scientific, spol. s r.o.
Kosmonautu 324
530 09 Pardubice
Tel: +420 466 798 230
Fax: +420 466 435 008
info.cz@thermofisher.com
www.thermofisher.cz

Denmark

Fisher Scientific Biotech Line A/S
Postboks 60
Industrivej 3
3550 Slangerup
Denmark
Tel: +45 70 27 99 20
Fax: +45 70 27 99 29
kundeservice@thermofisher.com
www.fishersci.dk

Finland

Fisher Scientific Oy
Ratastie 2
FI-01620 Vantaa, Finland
Tel: +358 802 76 237
Fax: +358 802 76 235
fisher.fi@thermofisher.com
www.fishersci.fi

France

Fisher Scientific
Parc d'Innovation BP 50111
67403 Illkirch Cedex
Tel: +33 (0) 3 88 67 53 20
Fax: +33 (0) 3 88 67 11 68
fr.commande@thermofisher.com
www.fr.fishersci.com

Germany

Fisher Scientific GmbH
Im Heiligen Feld 17
Schwerte
D-58239 Germany
Tel: +49 (0) 2304 932 5
Fax: +49 (0) 2304 932 950
info.germany@thermofisher.com
www.de.fishersci.com

Ireland

Fisher Scientific Ireland
Suite 3, Plaza 212
Blanchardstown Corporate Park 2
Ballycoolin
Dublin 15
Tel: +353 01 885 5854
Fax: +353 01 899 1855
fsie.sales@thermofisher.com
www.ie.fishersci.com

Italy

Fisher Scientific
Parc D'Innovation, BI Sebastien Brand,
BP 50111 - F67403 Illkirch Cedex
Tel: +39 02 953 28 258
Fax: +39 02 953 27 374
it.fisher@thermofisher.com

Luxembourg

Fisher Scientific
BP 567
B-7500 Tournai
Tel: +352 056 260 260
Fax: +352 056 260 270
be.fisher@thermofisher.com
www.be.fishersci.com

Netherlands

Fisher Scientific
P.O. Box 4
1120 AA Landsmeer
Scheepsbouwersweg 1B
1121 PC Landsmeer
Tel: +31 20 487 7000
Fax: +31 20 487 7070
nl.info@thermofisher.com
www.nl.fishersci.com

Norway

Fisher Scientific
Frysjaeveien 33E
0884 Oslo
Tel: +47 22 95 59 59
Fax: +47 22 95 59 40
fisher.no@thermofisher.com
www.fishersci.no

Portugal

Fisher Scientific
Rua Pedro Alvares Cabral, n°24, 3°D
Edificio Euro
Infantado 2670-391 Loures
Portugal
Tel: +351 21 425 33 50/4
Fax: +351 21 425 33 51
pt.fi.sher@thermofisher.com
www.pt.fishersci.com

Spain

Fisher Scientific, S.L.
Luis I, 9
28031 Madrid
Tel: +34 902 239 303
Fax: +34 902 239 404
es.fisher@thermofisher.com
www.es.fishersci.com

Sweden

Fisher Scientific GTF AB
Box 9193
400 94 Göteborg
Sweden
Tel: +46 (0) 31 352 32 00
Fax: +46 (0) 31 352 32 50
fisher.se@thermofisher.com
www.fishersci.se

Switzerland

Fisher Scientific AG
Wilstrasse 57
CH-5610 Wohlen
Switzerland
Tel: +41 56 618 4111
Fax: +41 56 618 4141
info.ch@thermofisher.com
www.ch.fishersci.com

United Kingdom

Fisher Scientific UK Ltd.
Bishop Meadow Road
Loughborough
Leicestershire LE11 5RG UK
Tel: +44 150 923 1166
Fax: +44 150 923 1893
info@fisher.co.uk
www.fisher.co.uk

MIDDLE EAST AND AFRICA
Fisher Scientific Global Export

MEA Office
Fisher Scientific Global Export
DAFZA
Office Block 4A #504
P.O. Box 54710
Dubai, UAE
Tel: +971 4609 1594
Fax: +971 4204 5564
www.fishersci.com

OCEANIA
Australia

Thermo Fisher Scientific
5 Caribbean Drive
P.O. Box 9092
Scoresby, Vic 3179
Tel: +61 1300 735 292
Fax: +61 1800 067 639
auinfo@thermofisher.com
www.thermofisher.com.au

New Zealand

Thermo Fisher Scientific
244 Bush Road
Albany, North Shore City 0632
New Zealand
Tel: +64 9980 6700
Fax: +64 9980 6788
nzinfo@thermofisher.com
www.thermofisher.co.nz

©2011 Thermo Fisher Scientific Inc.
All rights reserved. Litho in U.S.A.
11_4604 MA-JJ/LY 10M-IW 12/2011
BN12151110