



Remco Swart

Pittcon

March 8-12, 2015


Outline

- Thermo Scientific[™] Vanquish[™] UHPLC System
- Thermo Scientific[™] TSQ Endura[™] Triple Quadrupole Mass Spectrometer (MS) and Thermo Scientific[™] Quantiva[™] Triple Quadrupole Mass Spectrometer
- Thermo Scientific[™] Dionex[™] Chromeleon[™] Chromatography
 Data System (CDS) Software 7.2 for LCMS

Why a New UHPLC?


Customers told us they need:

- More resolution
- More sensitivity
- Easier and less method development
- More productivity
- Easier interaction and operation
- More robust and more reliable
- More sample capacity
- More access to experts
- More and fast service
- New, fresh & exciting instruments
- Better results in shorter time





The New Thermo Scientific Vanquish UHPLC System



We heard you! And we:

- Started with a complete new design in line with MS portfolio
- Put in 30 years of experience and know-how
- Developed a new and unique technology
- Closed gaps and pushed beyond previous limitations
- Boosted robustness

- Developed a new column technology
- Built the best UPHLC around the column/user
- Further leveraged our gold standard CDS software
- Boosted serviceability beyond established concepts



Vanquish System...Provides More Solutions!

Vanquish UHPLC is not only a system – it is a problem-solving solution philosophy



Better integration with the world's best mass spectrometers

Maximum speed and resolution with new Thermo Scientific™ Accucore™ Vanquish™ columns








Revolutionary ease-ofuse experience by new Chromeleon CDS features and workflows More throughput capacity
with the new Thermo
Scientific Vanquish
Charger module



Vanquish System is...





Vanquish Binary Pump

Feature	Benefit
1500 bar / 22 000 psi Up to 5000 μL/min	Flexibility for all working conditions
New parallel dual piston principle	very small pressure ripple => column and sensitivity
Adaptive Thermal Effect Compensation (ATEC™)	best flow and gradient performance => RT stable
2 × 3 solvent channels	9 different solvent combinations
Fingertight check valve design	Easy to service





Vanquish Split Sampler

Feature	Benefit
New switching valve	no maintenance
Adjustable GDV	easy method transfer
Automatic sample pre-compression	 reduced physical stress on column = longer lifetime better retention time reproducibility
New air-to-air cooling concept	no condensation
Increased sample capacity	4 racks or well plates
Automated barcode reader	Rack tracking (esp. with Charger)






Vanquish Charger Module

- Fully integrated robotic unit with less than 1 min cycle time
- Hosts up to 9 additional racks/ deep well plates or 20 shallow well plates
- Up to 9000 samples temperature controlled
- Precise temperature control using new designed air stream cooling concept

$$> 4 - 40 \, {}^{\circ}\text{C}$$

Integrated barcode reader







The Vanquish Column Compartment

Feature	Benefit
Active Pre-Heater	Perfect temperature control (QC)
Two thermostatting modes	Easy method transfer
Three independent column compartments	Method scouting, different temperature zones





Vanquish Column Compartment for MS integration

Feature	Benefit
Column compartment is vertical	The column outlet is closer to the MS ion source
Position of the column compartment can be to the left or to the right	Better integration with all Thermo Scientific Mass Spectrometers







Vanquish Diode Array Detector

- LightPipe technology for an unmatched detection experience
 - > 10 mm standard flow cell
 - > 60 mm high sensitivity flow cell
- Best signal-to-noise performance
- Lowest baseline noise
- Ultra-wide dynamic range
 - ➤ Up to 3000 mAU
 - For simultaneous detection of highly concentrated compounds and trace impurities
- Variable slit width for optimum resolution
 1 8 nm
- Low peak dispersion
- Reduced RI and thermal effects





Vanquish System Features & Benefits for LC-MS Users

Improved separation performance

FEATURE	BENEFIT
Better resolutionMinimal peak dispersion	Less matrix effectsMore consistent quantificationMore consistent sensitivity
 Faster analysis 	Improved LCMS productivity
 Rock solid retention time stability 	More reliable identificationMore reliable quantification
 New core enhanced 1.5 µm columns 	Less method developmentLess peaks of interest in the void volume




Challenges in Quantitative Mass Spectrometry



Large Sample Numbers



Low Concentration






TSQ Endura MS and TSQ Quantiva MS



TSQ Endura MS

Extreme Quantitative Value

- Best-in-class sensitivity
- Unprecedented usability
- Ultimate robustness



TSQ Quantiva MS

Extreme Quantitative Performance

- Attogram sensitivity
- Unprecedented usability
- Exceptional robustness



Sensitivity, Speed & Robustness for High-Throughput Laboratories





	TSQ Quantiva MS	TSQ Endura MS
Mass Range	10-1850	10-3400
Max SRM Number	30,000 SRMs	30,000 SRMs
SRM/Sec	500 SRMs/sec	500 SRMs/sec
Ion Optics	Active Ion Management (AIM)	S- LENS with Beam Blocker Technology
Quadrupole Design	 Ion Max NG source Electrodynamic ion funnel ion beam guide with neutral blocker 6 mm HyperQuad quadrupoles with asymmetric RF drive 	4mm Quadrupoles with Asymmetric RF
Reserpine Specification	100,000 :1 S/N for 1 pg Reserpine	10,000 :1 S/N for 1 pg Reserpine

TSQ Quantiva MS

Extreme quantitative performance

- Designed for the most challenging assays.
- For scientists needing to stay at the forefront of analytical technology


TSQ Endura MS

Extreme quantitative value

- Designed for non-stop operation.
- For scientist who need to run routine samples day in and day out.



Low Sample Concentration



Demanding assays which require the absolute lowest limit of detection



TSQ Quantiva MS

Powered by AIM (Active Ion Management) technology, the TSQ Quantiva MS is the world's most sensitive triple quadrupole MS, detecting compounds at the ppt level.

Large Sample Numbers



Challenge

Need to analyze multi-residue analytes in very short run times, using UHPLC with robust performance 24/7




Vanquish System with Charger and TSQ Quantiva MS

The AIM technology includes an ultra fast collision cell designed for fast SRMS with no loss in sensitivity. Neutral clusters eliminated for robustness.



Routine Quantitation



Challenge

Robust and reliable quantitation



TSQ Endura MS

LODs and LOQs unrivaled in its class. It delivers this best-in-class quantitation run-after-run and dayafter-day regardless of sample type or matrix.



TSQ Quantiva MS: Ease-of-Use







Vanquish System and TSQ MS with Chromeleon CDS Control

Analysts don't have time to waste learning new software or optimizing methods.

The entire software workflow has been redesigned to remove un-needed parameters.



Introduction

- How can you include MS in your laboratory CDS?
 - Currently it's on a workstation and there is no way to work remotely
- Chromeleon[™] 7.2 Chromatography Data System is the first CDS to support MS instrument control and data processing with all main front-end separation techniques (GC, LC, IC) in an enterprise environment







Quantiva & Endura - Fully Integrated within Chromeleon CDS

Chromeleon 7.2 CDS supported workflows

- Compliance, scalability, suite of quantification features
- MS-specific data organization at component level
- Context menu to open data in MS applications (e.g. Thermo Scientific™ Xcalibur™ Software)



Key benefits

- Only need to install, validate and learn one software package
- Can use all compliance and processing features of Chromeleon CDS
- Enhanced data security
- Easily start analyses with eWorkflows[™]
- Additional MS-specific detection algorithms (Genesis, ICIS)







MS Data Visualization



Component and Channel Centric Data Viewing




MS Data Views – Total Ion Chromatogram (TIC)






MS Data Views – MS Quantitation Channel




MS Data Views – Mass Spectrum Plot






MS Data Views - Component Traces





Data Visualization – Multiple Components / Injections





MS Data Processing




MS Data Processing

- With multiple traces and possibly hundreds of components, how can you quickly and easily process your MS data?
- How can you rapidly and accurately setup your component table?
- Existing data processing tools in Chromeleon CDS are now complimented by a suite of MS-specific tools for fast and efficient MS data processing, including:
 - Easy extracted ion chromatogram (XIC) creation
 - Additional MS-specific detection algorithms
 - Full processing control on a component-by-component basis
 - MS Library searching



MS Data Processing


- Apply all existing calibration modes and options to MS
- Easy extracted ion chromatogram creation
- Additional MS-specific detection algorithms
 - Genesis, ICIS
- Full processing control on a component-by-component basis
 - Additional MS-specific parameters in component table quantitation and confirming ions, CAS numbers, molecular formulas, etc.





MS Data Processing – Extracted Ion Chromatograms


Manual extraction of XIC





MS Data Processing – Extracted Ion Chromatograms

Graphical extraction




MS Data Processing – Extracted Ion Chromatograms

Graphical extraction



Integrated NIST Library Search


 Use NIST libraries for spectral library screening and ad-hoc library search






Import from Compound Databases and NIST Libraries

- Quickly populate processing method component tables
- Use library searches
 - NIST
 - User libraries
- Use predefined compound data bases
 - Thermo Scientific[™]
 TraceFinder[™] software
- Use Xcalibur software raw files with component info





MS Data Reporting




MS Data Reporting

- Can be a very time-consuming activity!
- Multiple components and traces make generic reporting tricky, if not impossible
- How do you effectively report such large quantities of data?
- Several new additions to Chromeleon CDS reporting tools
 - New MS based reporting objects
 - Mass Spectral Plot
 - Component Traces Plot
 - MS related report variables and tables
 - New tools for reporting of large amounts of multi-channel data
 - Consolidated Report Tables
 - Autorepeat Area Enhancements



MS Data Reporting

MS peaks often only appear in certain detection channels (filters) – how do you create a report for your peaks of interest?





MS Data Reporting – Consolidated Report Tables

Chromeleon CDS Consolidated Report Tables


Automatically report for the Channel with the highest response

No.	Peak Name	Area	Channel	
		counts*min		
TIC	TIC	Highest Response	Highest Response	
3	PCB#18	15392	TIC_F01	
10	PCB#28	19480	TIC_F01	
11	PCB#33	17945	TIC_F01	
25	PCB#52	5455	TIC_F03	
34	PCB#95	3640	TIC_F06	
50	PCB#70	2364	TIC	
79	PCB#156	3650	TIC_F06	
145	PCB#126	2757	TIC_F06	
179	PCB#153	4781	TIC_F07	
219	PCB#188	1100	TIC_F10	
225	PCB#170	889	TIC_F10	

 Useful for reporting large quantities of MS data in a generic, compact table



Antibiotics 1 ppb








36 Antibiotics Targeted Screening

Accucore Vanquish C18 column, 1.5µm x 2.1mm x 100mm

Overlaid selected reaction monitoring chromatograms showing detection of 36 antibiotics within a 5 minute detection window, binary Vanquish system and TSQ Vantage MS



Mobile phase A:

0.1% formic acid / Water

Mobile phase B: 0.1% formic

 $\begin{array}{ccc} & \text{acid / MeOH} \\ \text{Flow rate:} & 400 \, \mu\text{L/min} \\ \text{Column temperature:} & 40 \, ^{\circ} \, \text{C, active} \end{array}$

preheating 2 µL


Injection volume

Table 1. LC gradient conditions

Time (min)	%В
0.000	10
4.375	90
5.000	90
5.125	10
8.750	10

47 Drugs in 4 Minutes

Overlaid selected reaction monitoring chromatograms showing detection of 47 drugs within a 4 minute detection window, binary Vanquish system and TSQ Vantage MS



Mobile phase A: 10 mM ammonium

acetate in water

Mobile phase B: 0.1% formic acid in

methanol

Flow rate: $500 \,\mu\text{L/min}$

Column temperature: 50 $^{\circ}$ C, with active

preheating

Injection volume: 2 μL

Table 1. LC gradient conditions

Time (min)	%B
0.00	10
0.16	10
2.88	90
3.20	90
3.28	10
5.60	10



Summary – Bringing It All Together







Chromeleon CDS has all the tools to help you process and report your MS data in the shortest time!





Thank You!

Any Questions?