**ThermoFisher**  
S C I E N T I F I C


# Chromeleon 7.2 CDS eWorkflows Simply Intelligent

Matthias Hrobarsch  
Chromeleon Support Specialist  
Thermo Fisher Scientific, Dreieich

The world leader in serving science

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- Kategorie eWorkflows


The screenshot displays the Chromeleon Console interface. The main window title is 'Chromeleon Console'. The menu bar includes 'Back', 'Create', 'File', 'Edit', 'View', 'Tools', and 'Help'. The left sidebar shows a tree view with categories: 'eWorkflows', 'HPLC', 'Labor 1', 'Labor 2', 'Instruments', 'Data', and 'eWorkflows' (highlighted). The main content area shows the configuration for the selected eWorkflow '02 geschlossener eWorkflow'. The description is 'UHPLC Methode Produkt\_A für die Schichtarbeiter keiner Eingaben mehr notwendig auuser der Probenanzahl alles ist vorgegeben'. The type is 'HPLC' and the status is 'Approved'. Below this, there is a 'Launch' button and a table with columns: '#', 'Instrument Name', 'Instrument Status', 'Sequence Status', and 'Queue Status'. The table contains two rows: Row 1: '1', 'UltiMateRSLC', 'Idle', 'Sequence Status' is empty, 'Queue Status' is 'Pending sequences: 1'. Row 2: '2', 'Vanquish', 'Idle', 'Sequence Status' is empty, 'Queue Status' is empty. The status bar at the bottom indicates 'eWorkflow '02 geschlossener eWorkflow' selected' and 'Matthias: Full Access'.

#	Instrument Name	Instrument Status	Sequence Status	Queue Status
1	UltiMateRSLC	Idle		Pending sequences: 1
2	Vanquish	Idle		

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
**Filter:** Product A Tests und Type HPLC

The screenshot displays the 'eWorkflows' interface. At the top, there is a 'Filter' button with a funnel icon and a dropdown menu. Below this, a 'More filters' section is expanded, showing four filter criteria: 'Data Vault' (set to 'All'), 'State' (set to 'All'), 'Tag' (set to 'Product A Tests'), and 'Type' (set to 'HPLC'). A red rectangular box highlights these four filter settings. Below the filter section, the 'ChromeleonLocal' section is visible, listing three workflow items: 'CUT Templates Stage 1 (10 units)', 'CUT Templates Stage 2 (+20 units)', and 'ICH Accuracy'.

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
The screenshot displays the eWorkflow Editor interface. The main window is titled "eWorkflow Editor - 02 geschlossener eWorkflow <Modified>". The "Sequence General" tab is active, and its settings are highlighted with a red circle. The settings include:

- Data Vault:** ChromeleonLocal
- Path:** 000\_Spring Tour 2014/SEQ
- Sequence Name:** Fruehschicht - {gen.currentTime; "yyyy-mm-dd"}
- Comment:** durch den eWorkflow Fruehschicht erstellt am {gen.currentTime; "yyyy-mm-dd"}

Below the main settings, there are checkboxes for "Append injection(s) if sequence already exists", "Include header block", "Include bracket block", and "Include footer block".

The "Insert Formula" dialog is open, showing the "General" category selected. The "Variables" list includes "Current Time", which is selected. The "Formula" field contains "gen.currentTime". The "Header" field contains "Current Time". The "Unit" field is empty. The "Format" dropdown is set to "yyyy-mm-dd". The "Preview" field shows the resulting formula: "{gen.currentTime; 'yyyy-mm-dd'}".

At the bottom of the main window, there is a table for "Custom Variables and Formulas" with columns for Name, Description, Type, and Sequence Property Value. The "Custom Variables" tab is selected.

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
The screenshot shows the 'Sequence General' configuration window with the following fields and callouts:

- Speicherort:** Points to the 'Data Vault' dropdown menu, which is set to 'ChromeleonLocal'.
- Aktive Ordnerstruktur:** Points to the 'Path' field, which contains '000\_Spring Tour 2014/SEQ/{gen.currentTime; "yyyy-mm"}'.
- Dynamischer Sequenzname:** Points to the 'Sequence Name' field, which contains 'Fruehschicht -{gen.currentTime; "yyyy-mm-dd"}'.
- Kommentar aus Variablen:** Points to the 'Comment' field, which contains 'durch den eWorkflow Frühschicht erstellt am {gen.currentTime; "yyyy-mm-dd"}'.

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eWorkflow Editor - A\_GC\_positions\_test <Modified>

File | eWorkflow General | Sequence General | Sequence Layout

Max. Samples per Bracket: 1 | Max. Brackets per Sequence: 5 | No. of alternate Brackets: 1 |  Use Bracket Block after Sequence Header

Chromatogram	No. of Inj.	Type	Name	Level	Spike Group	Position	Volume
Sequence Header - 1 item							
n.A.	1	Calibration Standard	calibration#n_#p_#i_#y_#r			R:A1	2.200
Sample Block - 1 item							
n.A.	2	Unknown	sample #n_#p_#i_#y_#r			<defined at launch>	1.400
Bracket - 1 item							
n.A.	1	Unknown	bracket#n_#p_#i_#y_#r			6	0.500
Alternating Bracket 1 - 1 item							
n.A.	1	Unknown	bracket2_#n_#p_#i_#y_#r			7	6.100
Sequence Footer - 1 item							
n.A.	1	Check Standard	footer#n_#p_#i_#y_#r			test2	3.100

Sequence Preview

Number of Samples: 1 | vanquish on DEGER-C-NSPRUNK | Sampler start position: R:A2 | Try Launch

#	Chromatogram	Type	Name	Level	Spike Group	Position	Volume	Instrument Method
1	None	Calibration Standard	calibration1_R:A1_2_2_1_1_1			R:A1	2.20	GC6890_InstrumentMethod
2	None	Calibration Standard	calibration2_R:A1_2_2_2_2_1			R:A1	2.20	GC6890_InstrumentMethod
3	None	Unknown	bracket3_R:A7_0.5_1_1_1			R:A7	0.50	GC6890_InstrumentMethod
4	None	Unknown	sample 4_R:A2_1.4_2_1_1			R:A2	1.40	GC6890_InstrumentMethod
5	None	Unknown	sample 5_R:A2_1.4_3_2_1			R:A2	1.40	GC6890_InstrumentMethod
6	None	Unknown	bracket2_6_912395_6.1_4_1_1			912395	6.10	GC6890_InstrumentMethod
7	None	Check Standard	footer7_test2_3.1_1_1_1			test2	3.10	GC6890_InstrumentMethod

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- Sequence Layout incl. eigener Variablen

eWorkflow Editor - 02 geschlossener eWorkflow <Modified>

File eWorkflow General Sequence General Sequence Layout

Max. Samples per Bracket: 1 Max. Brackets per Sequence: 99 No. of alternate Brackets: 0  Use Bracket Block after Sequence Header

Chromatogram	No. of Inj.	Name	Type	Position	Volume	*test	Instrument Method	Processing Method	Status	Inject Time	Weight	Dilution
Sequence Header - 1 item												
n.A.	1	Spülen	Blank	A01	5,000 [µl]		Produkt_A	Proben_Produkt_A	Idle		1,000	1,000
New eWorkflow Injection Template												
Sample Block - 2 items												
n.A.	1	<HIER PROBENNAME>	Unknown	<defined at launch>	5,000 [µl]	TEST1	Produkt_A	Proben_Produkt_A	Idle		1,000	1,000
n.A.	1	<HIER PROBENNAME>	Unknown	<defined at launch>	5,000 [µl]	TEST1	Produkt_A	Proben_Produkt_A	Idle		1,000	1,000
New eWorkflow Injection Template												
Bracket - 1 item												
n.A.	1	Spülen	Blank	A01	5,000 [µl]		Produkt_A	Proben_Produkt_A	Idle		1,000	1,000
New eWorkflow Injection Template												
Sequence Footer - 1 item												
n.A.	1	Spülen	Blank	A01	1,000 [µl]			Proben_Produkt_A	Idle		1,000	1,000
New eWorkflow Injection Template												


Sequence Preview

Number of Samples: 1  Sampler start position:  Try Launch

#	Chromatogram	Name	Type	Position	Volume	*test	Instrument Method	Processing Method	Status	Inject Time	Weight	Dilution	Comments
1	None	Spülen	Blank	A01	5,000		Produkt_A	Proben_Produkt_A	Idle		1,000	1,000	Blank.dad3000
2	None	Spülen	Blank	A01	5,000		Produkt_A	Proben_Produkt_A	Idle		1,000	1,000	Blank.dad3000
3	None	<HIER PROBENNAME>	Unknown	n.a.	5,000	TEST1	Produkt_A	Proben_Produkt_A	Idle		1,000	1,000	Linearity 2 run 1 da
4	None	<HIER PROBENNAME>	Unknown	n.a.	5,000	TEST1	Produkt_A	Proben_Produkt_A	Idle		1,000	1,000	Linearity 2 run 2 da
5	None	Spülen	Blank	A01	5,000		Produkt_A	Proben_Produkt_A	Idle		1,000	1,000	Blank.dad3000
6	None	Spülen	Blank	A01	1,000			Proben_Produkt_A	Idle		1,000	1,000	

# Sequenz erstellen via eWorkflow – Schritt 1/3

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
The screenshot displays the Chromeleon Console interface. The main window title is "Chromeleon Console". The menu bar includes "Back", "Create", "File", "Edit", "View", "Tools", and "Help". The left sidebar shows a tree view with "eWorkflows" selected, and a list of workflows including "02 geschlossener eWorkflow". The main area displays the details for the selected workflow, "02 geschlossener eWorkflow", with a description: "UHPLC Methode Produkt\_A für die Schichtarbeiter keiner Eingaben mehr notwendig auuser der Probenanzahl alles ist vorgegeben". The workflow type is "HPLC" and the status is "Approved". A "Launch" button is highlighted with a red box. Below the details, there is a table with columns: "#", "Instrument Name", "Instrument Status", "Sequence Status", and "Queue Status". The table contains two rows: Row 1: "1", "UltiMateRSLC", "Idle", "Sequence Status", "Pending sequences: 1"; Row 2: "2", "Vanquish", "Idle", "Sequence Status", "Queue Status". The status bar at the bottom indicates "eWorkflow '02 geschlossener eWorkflow' selected" and "Matthias: Full Access".

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- Anzahl Proben angeben - *Finish*

eWorkflow Wizard

**Sample Configuration**  
Define the desired number of samples and the corresponding start position

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The screenshot displays the Chromeleon Console interface. The main window shows a sequence titled "Fruehschicht -2017-06-12". The sequence is composed of 6 steps, each with a specific name, type, position, volume, and method. The steps are as follows:

#	UV_VIS_1	Name	Type	Position	Volume [µl]	*test	Instrument Method	Processing Method
1	None	Spülen	Blank	BA1	5,000		Produkt_A	Proben_Produkt_A
2	None	Spülen	Blank	BA1	5,000		Produkt_A	Proben_Produkt_A
3	None	<HIER PROBENNAME>	Unknown	RA1	5,000	TEST1	Produkt_A	Proben_Produkt_A
4	None	<HIER PROBENNAME>	Unknown	RA2	5,000	TEST1	Produkt_A	Proben_Produkt_A
5	None	Spülen	Blank	BA1	5,000		Produkt_A	Proben_Produkt_A
6	None	Spülen	Blank	BA1	1,000		Produkt_A	Proben_Produkt_A

Below the sequence table, there is a section for "Associated Items" which lists various files and templates associated with the sequence:

Name	Date Modified	Type	Comment
065565 Chromeleon 7.2 C...	18.07.2013 12:38:06 +02...	Attachment	
Ausdruck_Produkt_A	25.04.2014 13:25:26 +02...	Report Template	
Bildschirmansicht	25.04.2014 13:25:26 +02...	View Settings	
Proben_Produkt_A	25.04.2014 13:25:26 +02...	Processing Met...	
Produkt_A	25.04.2014 13:56:42 +02...	Instrument Met...	

The interface also shows a file tree on the left with the selected sequence "Fruehschicht -2017-06-12" highlighted. The bottom status bar indicates "Sequence 'Fruehschicht -2017-06-12' selected" and "Matthias: Full Access".

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