




# Advanced QA/QC characterization MS in QC : Multi Attribute Method

Global BioPharma Summit

# A Complex Problem: Drug Safety and Quality



## Safety

Is the product safe to use?  
(e.g. Immunogenic effects?)

## Potency

Does the drug have the expected effect?  
(e.g. CDR complementation)

## Knowledge

How do changes effect the therapeutic?  
(e.g. Oxidation)

## Quality

How do changes in process effect the product?  
(e.g. Glucose concentration on glycoforms)

# Drug Development Workflow: From Discovery to Production

## What?

## How much?



### MS in QC: Multi Attribute Method (MAM)

- Molecular Assessment
- Quality Control
- Attribute Science Groups
- Analytical and Automation Process Technologies
- Protein Attribute Chemistry
- Manufacturing

# Advantages of MAM



## Required Characterizations


## MAM Method

## Conventional Methods


	Pep Map-MS	SEC	CEX	rCE-SDS	nrCE-SDS	HILIC	ID ELISA	HCP ELISA
Aggregate Assessment	No	Yes	Indirect	Yes	Yes	No	No	No
Deamidation (Isomerization) Assessment	Yes	No	Indirect	No	No	No	No	No
Disulfide Isoform Assessment	Maybe	No	Indirect	No	Yes	No	No	No
Glycation Assessment	Yes	No	No	Yes	Yes	No	No	No
High Mannose Assessment	Yes	No	No	No	No	Yes	No	No
Methionine Oxidation Assessment	Yes	No	No	No	No	No	No	No
Signal Peptide Assessment	Yes	No	No	No	No	No	No	No
Unusual Glycosylation Assessment	Yes	No	Indirect	Maybe	Maybe	Yes	No	No
CDR Tryptophan Degradation Assessment	Yes	Indirect	No	No	No	No	No	No
Non-consensus Glycosylation Assessment	Yes	No	No	Maybe	Maybe	No	No	No
N-terminal pyroGlutamate Assessment	Yes	No	Indirect	No	No	No	No	No
C-terminal Lysine Assessment	Yes	No	Yes	No	No	No	No	No
Galactosylation Assessment	Yes	No	No	No	No	No	No	No
Dimer Assessment	No	Yes	No	No	No	No	No	No
Fragmentation (peptide bond) Assessment	Maybe	Maybe	No	Yes	Yes	No	No	No
Disulfide Reduction (DS Fragmentation) Assessment	Maybe	No	No	No	Yes	No	No	No
Host Cell Protein Assessment	Yes	No	No	No	No	No	No	Yes
Mutations/Misincorporations Assessment	Yes	No	No	No	No	No	No	No
Hydroxylysine Assessment	Yes	No	No	No	No	No	No	No
Thioether Assessment	Yes	No	No	No	No	No	No	No
Trisulfide Assessment	Maybe	No	No	No	No	No	No	No
Non-glycosylated Heavy Chain	Yes	No	No	No	No	No	No	No
DNA Assessment	No	No	No	No	No	No	No	No
Cysteine Adducts Assessment	Maybe	No	Maybe	No	No	No	No	No
C-terminal Amidation Assessment	Yes	No	Indirect	No	No	No	No	No
CDR Conformers (HIC Isoform) Assessment	No	No	Indirect	No	No	No	No	No
O-linked Glycans Assessment	Maybe	No	No	No	No	No	No	No
Fucosylation Assessment	Yes	No	No	No	No	No	No	No
Residual Protein A	Yes	No	No	No	No	No	No	No
Identity	Yes	No	Yes	No	No	No	Yes	No

# Multi-Attribute Method for QC


## Current Release Method



## Product Knowledge



## Multi Attribute Method



## Multi Attribute Method




### Confirmation

Confirm process by evaluating release product against a gold reference standard



### High Resolution

Using High Resolution/Accurate Mass MS instrumentation to directly measure CQA



### Quality by Design



Implement hybrid traditional (evaluation in QC) and QbD approach (monitoring of CQA and end-result of CPP space)

### Knowledge



Eliminate the need of traditional lot release methodologies while increasing product knowledge



# Materials

3 ug Trypsin Digested  
NIST mAb

Sample



**ThermoScientific™ Vanquish™ UHPLC**  
70 minute Gradient 250 uL/min  
0.1% FA in H<sub>2</sub>O and MeCN  
50° C Column Temp

Column



2.1x150mm 1.5 uM  
**ThermoScientific™ Accucore™  
Vanquish C18+**

uHPLC



**ThermoScientific™ QExactive™ HF**  
120k Resolution @ 200 m/z  
300 to 1800 m/z

MS



Discovery




**ThermoScientific™ Chromeleon™  
CDS 7.2 SR5**  
CFR 21 Part 11 Compliant Data  
Acquisition,  
Processing,  
and Reporting

Processing



**ThermoScientific™ BioPharma  
Finder™ 2.0**  
Component Detection, Peptide  
Mapping, and CQA Selection




# Methods

1. Product Denaturation with 7 M Guanidine HCl and reduction with 500 mM DTT (30 min)
2. Alkylation with 500 mM Iodoacetate (20 min)

5. LC-MS/MS Analysis – Top 5

6. Peptide Mapping Selection of Attributes to Monitor

R. Rogers, et al. MABs 2015, 7, 881-890






# Peptide Mapping

TIC

Intensity 3e9

Basepeak

Intensity 1.6e9



# Peptide Mapping

TIC

Intensity 3e9

Basepeak

Intensity 1.6e9

Light Chain

Heavy Chain

Unidentified






Sequence Coverage

Complete sequence coverage of heavy chain and light chain

## Feature Chromatogram

Detected features belonging to light chain, heavy chain, and unidentified species.




# Establishing CQAs



BioPharma Finder™

Predicted/Observed RT shift



Modification Summary  
Chart

SELECT Component Used for  
Quantitation

Predicted  
Component  
Extracted Ion  
Chromatogram

Observed Spectrum

<3 ppm mass error

Confidence

# MS in QC – Late Stage Discovery to Process Dev / QC

- Analytical group is primary driver and is responsible for
  - Identification and selection of target peptides / modifications
  - Generation of targeted component list (HR/AM MS1) for QC
  - Complete analytical method (processing and instrument methods, designing report templates, eWorkflows)
  - Development of peptide mapping assay SOPs (for QC)
  - MS list is transferred directly from the discovery experiment in BioPharma Finder 2.0 to Chromeleon
  - Checkboxes enable easy selection of components to transfer.
  - More can be added at anytime in the future









Results								
		Level	No.	Identification	Peptide Sequence	Modification	Site	Delta (ppm)
		Δa	Δa	Δa	Δa	Δa	Δa	Δa
				(NonBlanks)				
94		Component	1414	1:Q154-K168 = 1622.70199m[nonspecific]	QSGNSQESVTEQDSK	nonspecific		-1.18
95		Component	1416	1:Q154-K168 = 1622.70199m[nonspecific]	QSGNSQESVTEQDSK	nonspecific		-0.65
96		Component	1422	2:P294-R304 = 1441.65861m(N300+A2G2F)	PREEQYNSTYR	A2G2F	N300	-1.49
97		Component	1424	2:H60-R68 = 1128.56761m	HYNPSLKDR	None		-1.60
98		Component	1425	2:P294-R304 = 1441.65861m(N300+A2G1F)	PREEQYNSTYR	A2G1F	N300	-1.21
99		Component	1428	2:H60-R68 = 1128.56761m(K66+Glycation )	HYNPSLKDR	Glycation	K66	-1.38
100		Component	1430	2:P294-R304 = 1441.65861m(N300+A2G0F)	PREEQYNSTYR	A2G0F	N300	-1.24
101	<input checked="" type="checkbox"/>	Component	1431	2:H60-R68 = 1128.56761m(K66+Glycation )	HYNPSLKDR	Glycation	K66	-1.00
102		Component	1437	1:V149-K168 = 2134.96145m(D166+Isomeriz...	VDNALQSGN	D166		-0.46
103		Component	1457	1:V149-K168 = 2134.96145m(D166+Isomeriz...	VDNALQSGN			0.00
104		Component	1505	2:E296-R304 = 1188.50473m(N300+A3Ga3F)	EEQYNSTYR			-1.31
105		Component	1515	2:T259-K277 = 2080.99869m(C264+Carboxy...	TPEVTCVVVD			-0.58
106		Component	1517	2:E296-R304 = 1188.50473m(N300+A2Ga1G...	EEQYNSTYR	A2Ga1G1F	N300	-1.87
107		Component	1518	2:E296-R304 = 1188.50473m(N300+A2Ga2F)	EEQYNSTYR	A2Ga2F	N300	-1.59
108		Component	1519	2:E296-R304 = 1188.50473m(N300+A2Ga1G...	EEQYNSTYR	A2Ga1G1F	N300	-1.63
109		Component	1523	2:E296-R304 = 1188.50473m(N300+A3Ga1G...	EEQYNSTYR	A3Ga1G2F	N300	-1.08
110		Component	1524	2:E296-R304 = 1188.50473m(N300+A2Ga2F)	EEQYNSTYR	A2Ga2F	N300	-1.66
111		Component	1525	2:E296-R304 = 1188.50473m(N300+A2Ga2F)	EEQYNSTYR	A2Ga2F	N300	-1.74

- Export All components ▶
- Export Checked components ▶
- Create mgf File
- Run De Novo Processing
- Component Information
- Excel Workbook
- CSV
- Chromeleon

Compound Data Import										
Data Source										
PinPoint Data Path: C:\Users\VRKomatsu\k\Desktop\Step_c...										
	Name	RT	Charge	Precursor Mass	Isotope 1	Isotope 2	Isotope 3	Isotope 4	Isotope 5	
<input checked="" type="checkbox"/>	EEQFN[Hex9HexNAc2]STFR	11.617	2	1511.58154	1512.08337	1512.58484	1511.58179	1513.08618	1513.58752	
<input checked="" type="checkbox"/>	EEQFN[dHex1Hex5HexNAc4]STFR	11.722	2	1008.05670	1008.39136	1008.72565	1008.05695	1009.05988	1009.39417	
<input checked="" type="checkbox"/>	EEQFN[Hex8HexNAc2]STFR	11.722	2	1463.58411	1464.08606	1464.58752	1463.58447	1465.08887	1465.59021	
<input checked="" type="checkbox"/>	EEQFN[Hex8HexNAc2]STFR	11.722	2	1430.55505	1431.05688	1431.55835	1430.55542	1432.05681	1432.55816	
<input checked="" type="checkbox"/>	EEQFN[dHex1Hex6HexNAc3]STFR	11.730	2	954.03912	954.37372	954.70801	954.03937	955.04230	955.37653	
<input checked="" type="checkbox"/>	EEQFN[Hex7HexNAc2]STFR	11.800	2	1349.52865	1350.03052	1349.52893	1350.03198	1351.03333	1351.53467	
<input checked="" type="checkbox"/>	EEQFN[dHex1Hex5HexNAc3]STFR	11.870	3	900.0214	900.35614	900.02173	900.69043	901.02466	901.35889	
<input checked="" type="checkbox"/>	EEQFN[dHex1Hex4HexNAc4]STFR	11.888	2	1363.04480	1363.04473	1363.04473	1363.04473	1363.04473	1363.04473	
<input checked="" type="checkbox"/>	EEQFN[Hex6HexNAc2]STFR	11.905	2	1268.5022	1269.00415	1268.50256	1269.00549	1270.00696	1270.50830	
<input checked="" type="checkbox"/>	EEQFN[Hex4HexNAc4]STFR	11.975	2	1309.52881	1310.03064	1309.52917	1310.03210	1311.03345	1311.53479	
<input checked="" type="checkbox"/>	EEQFN[dHex1Hex3HexNAc4]STFR	12.001	2	873.35492	873.69851	873.35516	874.02380	874.35803	874.69226	
			3	1301.53125	1302.03320	1301.53162	1302.03467	1303.03601	1303.53735	
			3	868.02332	868.35791	868.02356	868.69220	869.02643	869.36066	



# Acquiring Data



# Data Analysis Processing Chromeleon

## Chromeleon 7.2



- Set up injection sequence
- Build LC and MS methods
- Target. Confirm. Integrate.
- Look for new peptides/impurities/features
- Automatically generate reports outputs/formulas

Simple. Logical.

**1. Create Injection List** (eWorkflow)

**2. Define Instrument Method** (LC and MS, eWorkflow)


**3. Targeted MS Processing**

Verify all expected peptides

**4. Non-Targeted MS Processing**


Identify new peptides/impurities/modifications  
Comparison to reference standard (A to B)

**5. Reporting**





# Building Targeted List of Critical Quality Attributes



# Adjust Processing Settings as Desired

Chromeleon 7.2



## MS Settings

high resolution enables accurate extraction and quantitation

*no offline recalibration necessary*

# Adjust Processing Settings as Desired

## Chromeleon 7.2



MS Settings

Composite Scoring

Detection MS Detection MS Component Table Calibration MS Settings MS Library Screen

Composite Scoring

Scoring result on checked criteria

Pass score if at least  criteria passed.

Fail score if less than  criteria passed.

Indeterminate score if neither passed nor failed.

2D Criteria

☐ Amount based peak identity verification

Reference Channel:  Confirmation Channel:  Tolerance(%): +/-

MS Criteria

☐ Confirming ion ratio passed

☒ Isotopic dot product  $\geq$

☒ Mass accuracy  $\leq$   PPM

☒ Peak apex alignment  $\leq$   min

No. of Conditions Satisfied for a Full Pass

No. of Conditions Satisfied for a Partial Pass

- Isotopic Ratios (Individual)
- Isotopic Correlation (Envelope)
- Mass Accuracy
- Retention Time Apex Alignment

# Adjust Processing Settings as Desired

## Chromeleon 7.2



MS Settings



Composite Scoring

Integration Settings


RT Extraction Settings

MS or MS/MS Filters

Ion Ratio Settings

MS Spectra

Smoothing



# Adjust Processing Settings as Desired

## Chromeleon 7.2



MS Settings

Composite Scoring

Integration Settings

Retention Time Settings

- Set absolute or relative retention time
- Set absolute or relative RT window
- Setting on specific to individual component
- First, greatest, nearest match within window

### Retention

Retention Time

Retention Time:  min

☐ Use this component as reference component

Interpretation

- ☒ Absolute Time  
☐ Time Distance  
☐ Time Ratio

Reference component:

☒ Requires Reference Peak

☐ Retention Time Standard Peak Area Ratio

Ratio Tolerance:

Window

Detect peak within retention time +/-  min


Interpretation

- ☒ Absolute  
☐ Relative

Component Match

- ☐ First  
☒ Greatest  
☐ Nearest  
☐ Spectrum only  
☐ Spectrum and Time

# Quantification of CQA: M255 Oxidation



## Components

Targeted components and their detection and quantitation status

## Experimental Stage

Current state of experiment from injection to reporting

## Ion Chromatogram

5 ppm extraction based on scan filters and optional smoothing

## Interactive Results

Continually updated results during processing and acquisition

## Mass Spectra

Visualize Full MS spectra for selected component

## Isotopic Correlation

Validation of peptide identify based on theoretical elemental composition




## Deamidation of N392

GFYPSDIAVEWESNGQPEN[Deamid]NYK



### Chromatographic Resolution

Temporal separation of many deamidated species



## Deamidation of N62

HYN[Deamid]PSLK


Main Peak



**Not Chromatographically Resolved**

No temporal separation of this deamidated species

Deamidation



# Separation by Resolution

## Deamidation of N62

HYN[Deamid]PSLK


Main Peak



**Not Chromatographically Resolved**

No temporal separation of this deamidated species

Deamidation



30k Resolution  
25 ppm Extraction

Completely incorrect quantitation

## Deamidation of N62

HYN[Deamid]PSLK



**Not Chromatographically Resolved**

No temporal separation of this deamidated species



30k Resolution  
25 ppm Extraction

Completely incorrect quantitation

Peak Name	Ret. Time min	Area counts*min
First Injection	NistAb02_MAM30k	NistAb02_MAM30k
MS Quantitation Peak	MS Quantitation Peak	MS_Quantitation
HYN[Deamid]PSLK	7.928	40595222.619
HYNPSLK	7.928	104371927.127
	TOTAL	144967149.7
	Rel Abundance	72.00%

## Deamidation of N62

HYN[Deamid]PSLK




### Not Chromatographically Resolved

No temporal separation of this deamidated species



30k Resolution  
25 ppm Extraction

Completely incorrect quantitation



## Deamidation of N62

HYN[Deamid]PSLK



**Not Chromatographically Resolved**

No temporal separation of this deamidated species

Peak Name	Ret.Time min	Area counts*min
First Injection	NistAb02_MAM30k	NistAb02_MAM30k
MS Quantitation Peak	MS Quantitation Peak	MS_Quantitation
HYNPSLK	7.928	96322458.682
HYN[Deamid]PSLK	8.014	82470.684
	TOTAL	96404929.37
	Rel Abundance	0.09%
HYNPSLK - 2		
✓ HYNPSLK		
✓ (+2) HYNPSLK		
HYN[Deamid]PSLK - 2		
✓ HYN[Deamid]PSLK		
? (+2) HYN[Deamid]PSLK		



30k Resolution  
5 ppm Extraction

Not enough true resolution, over  
extracts the data



## Deamidation of N62

HYN[Deamid]PSLK


Main Peak



**Not Chromatographically Resolved**

No temporal separation of this deamidated species

Deamidation



**30k Resolution**  
**5 ppm Extraction**

Not enough true resolution, over extracts the data

## Deamidation of N62

HYN[Deamid]PSLK


Main Peak



**Not Chromatographically Resolved**

No temporal separation of this deamidated species

Deamidation



**120k Resolution**  
**5 ppm Extraction**

Spectrometrically resolve deamidated peptides

## Deamidation of N62

HYN[Deamid]PSLK



Not Chromatographically Resolved

No temporal separation of this deamidated species

Peak Name	Ret.Time min	Peak Area counts*min
First Injection	NistAb02_MAM120k_1	NistAb02_MAM120k_1
MS Quantitation Peak	MS Quantitation Peak	MS Quantitation Peak
HYNPSLK	7.930	83977247.099
HYN[Deamid]PSLK	8.022	202432.341
	TOTAL	84179679.44
	Rel Abundance	0.24%
HYNPSLK - 2		
✓ HYNPSLK		
✓ (+2) HYNPSLK		
HYN[Deamid]PSLK - 2		
✓ HYN[Deamid]PSLK		
✓ (+2) HYN[Deamid]PSLK		



120k Resolution  
5 ppm Extraction

Spectrometrically resolve deamidated peptides

## Chromeleon 7.2



- Set up injection sequence
- Build LC and MS methods
- Target. Confirm. Integrate.
- Look for new peptides/impurities/features
- Automatically generate reports outputs/formulas

Simple. Logical.

1. Create Injection List

2. Define Instrument Method (LC and MS)


3. Targeted MS Processing

Verify all expected peptides

4. Non-Targeted MS Processing

Identify new peptides/impurities/modifications  
Comparison to reference standard (A to B)

5. Reporting



# Detection of New Features: Batch to Batch

- ✓ **Assign a reference injection (A) and compare as many new injection (B) as desired**

First Sample Prep (A)

Second Sample Prep (B)

- ✓ **Set Alignment and Framing Settings**

0.1% Base Peak      10 ppm   0.5 min

- ✓ **Automatic Alignment and Framing**

- ✓ **Reporting of Features**


- ✓ **Filtering of Detected Features**

Monoisotopic

Multiple Isotopes

Charge between 2 and 4

More than 10-fold Change



- ✓ **Validation of Results**

# Detection of New Features: Batch to Batch

- ✓ Assign a reference injection (A) and compare as many new injection (B) as desired

First Sample Prep (A)

Second Sample Prep (B)

- ✓ Set Alignment and Framing Settings

0.1% Base Peak      10 ppm   0.5 min

- ✓ Automatic Alignment and Framing

- ✓ Reporting of Features


- ✓ Filtering of Detected Features

Monoisotopic

Multiple Isotopes

Charge between 2 and 4

More than 10-fold Change



- ✓ Validation of Results



# Detection of New Features: Batch to Batch

- ✓ **Assign a reference injection (A) and compare as many new injection (B) as desired**

First Sample Prep (A)

Second Sample Prep (B)

- ✓ **Set Alignment and Framing Settings**

0.1% Base Peak      10 ppm 0.5 min

- ✓ **Automatic Alignment and Framing**

- ✓ **Reporting of Features**


- ✓ **Filtering of Detected Features**

Monoisotopic

Multiple Isotopes

Charge between 2 and 4

More than 10-fold Change



- ✓ **Validation of Results**

# Detection of New Features: Batch to Batch

- ✓ Assign a reference injection (A) and compare as many new injection (B) as desired

First Sample Prep (A)

Second Sample Prep (B)

- ✓ Set Alignment and Framing Settings

0.1% Base Peak      10 ppm   0.5 min

- ✓ Automatic Alignment and Framing

- ✓ Reporting of Features


- ✓ Filtering of Detected Features

Monoisotopic

Multiple Isotopes

Charge between 2 and 4

More than 10-fold Change



- ✓ Validation of Results

✓ In Spec: No New Features

# Detection of New Features: PRTC Syn Peptide Kit

New Peak Detection

ThermoFisher

SCIENTIFIC

Injection Details

Injection Name:

Acdaim15cm\_PRTC

Vial Number:

R-C4

Injection Type:

Unknown

Calibration Level:

Instrument Method:

QEHF\_MAM\_MkV

Processing Method:

SuperMAM

Injection Date/Time:

29/Mar/17 19:40

Run Time (min):

68.51

Injection Volume:

3.00

Dilution Factor:

1.0000

Sample Weight:

1.0000

Total Ion Chromatogram Table

No.

Injection Name

Maximum Signal Value counts  
TIC

1% of TIC counts  
TIC

0.1% of TIC counts  
TIC

1

Acdaim15cm

1580663200

15806632

1580663

2

Acdaim15cm

1490999700

14909997

1491000

3

Acdaim15cm\_PRTC

1438209500

14382095

1438210

4

Acdaim15cm\_PRTC

1459710600

14597106

1459711

Detection Settings

Status

Not In Spec

m/z min

300.0000

Algorithm

FRAME

m/z max

1800.0000

Alignment Bypass

False

Frame Width [ppm]

10.0000

Alignment Min Intensity

1000.00

Max Frames

5000

Alignment Bin Width

1.00

Frame Width [ppm]

10.0000

Frame Filter

FTMS + p ESI Full lock ms [300.00-1800.00]

Peak Intensity Threshold

1500000

Reference Injection No.

2

RT Start [min]

2.00

RT Stop [min]

70.52

Frame Width [min]

2.50

Frame Width [ppm]

10.0000

New Features

Frame ID

Component ID

Mass m/z

RT min

Ratio

Base Component

PR Element

PR Root

PR Size

Charge

Molecular Weight m/z

190

252

496.2873

8.79

99999.9000

252

0

190

3

2

990.55997

210

336

613.3156

8.14

99999.9000

336

0

210

3

2

1224.61655

286

299

422.7304

10.34

99999.9000

299

0

286

2

2

843.45817

316

357

493.7691

8.07

22834.7506

357

0

316

3

2

985.52372

331

219

658.3267

25.04

320923.5021

219

0

331

3

2

1114.63877

422

301

745.3943

26.10

99999.9000

301

0

422

3

2

1488.77402

444

330

773.8973

18.03

99999.9000

330

0

444

3

2

1545.78013

460

359

695.8323

12.55

5016.2316

359

0

460

3

2

1389.65012

548

317

301.1919

11.10

4675.0764

317

0

548

2

3

900.55376

551

334

801.4132

25.82

99999.9000

334

0

551

3

2

1600.91186

578

309

673.3038

34.85

111.9689

309

0

578

2

2

1144.59311

587

332

680.3748

38.89

99999.9000

332

0

587

3

2

1358.73496

655

358

498.8026

30.88

37933.1051

358

0



655

2

2

995.59073

## PRTC Peptides Detected



✓ Detection of new features in Nist Ab digest with 1% of Peptide RT Calibration peptides spiked in

Digest Only (A)  
Digest + 500 fmol/uL PRTC (B)


Alignment and Framing Settings

1.0% Base Peak 10 ppm

# Detection of New Features: Stress Study

## Two Weeks at pH 8


Forced stress study versus sample immediately after preparation



# Detection of New Features: Stress Study


## Two Weeks at pH 8

Forced stress study versus sample immediately after preparation





Sequence Details			Created On: 09/Nov/16 14:49:34	
Name:	HMSVictory	Instrument:	GEHF_TheGhost	Updated On: 02/Feb/17 11:08:45
Imported Data:	False	First Injection:	NISTmAb_Control1_DDA_1	
Processing Method:	SuperMAM	MS Acquisition Time [min]:	70.52	
Method Length [min]:	115.00	Total Time [hrs]:	63.25	
Data Vault:	ChromleonLocal	Created By:	Thermo	
No. of Injections:	33	Updated By:	Thermo	



Retention Time Stability									
No.	Name	Retention (DETECTED) min	Retention (DETECTED) min	Retention (DETECTED) min	Retention (DETECTED) min	Retention (DETECTED) min	Retention (DETECTED) min	Retention (DETECTED) min	Retention (DETECTED) min
15	Control_Digest2_Rep1	EEQYN[dHex1Hex3HexNAc4]STYR	DTLM[Oxid]SR	DTLMISR	DIQMTQSPSTLSASVGDR	FNWYVDGVEVHNAK	WQQGNVFSC[Carboxymethyl]SYMHEALHNHYTK	GFYPSDIAVEWESNGQPENNYK	GFYPSDIAVEWESN
16	Control_Digest2_Rep2	8.14	10.33	12.87	20.89	24.37	25.40	35.03	35.71
17	Control_Digest2_Rep3	8.14	10.35	12.88	20.87	24.38	25.43	35.05	35.72
		8.14	10.34	12.89	20.93	24.43	25.50	35.11	35.75
	coefficient of variation	0.0%	0.1%	0.1%	0.1%	0.1%	0.2%	0.1%	0.2%
		PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Peak Area Stability and Minimum Intensity									
No.	Name	Area counts*min	Area counts*min	Area counts*min	Area counts*min	Area counts*min	Area counts*min	Area counts*min	Area counts*min
15	Control_Digest2_Rep1	EEQYN[dHex1Hex3HexNAc4]STYR	DTLM[Oxid]SR	DTLMISR	DIQMTQSPSTLSASVGDR	FNWYVDGVEVHNAK	WQQGNVFSC[Carboxymethyl]SYMHEALHNHYTK	GFYPSDIAVEWESNGQPENNYK	GFYPSDIAVEWESN
16	Control_Digest2_Rep2	1.68E+07	1.80E+06	1.39E+08	2.58E+08	3.10E+08	2.78E+08	2.00E+08	3.12E+08
17	Control_Digest2_Rep3	1.64E+07	1.83E+06	1.38E+08	2.56E+08	3.08E+08	2.84E+08	2.01E+08	3.11E+08
		1.64E+07	1.84E+06	1.38E+08	2.55E+08	3.07E+08	2.85E+08	2.00E+08	3.01E+08
	coefficient of variation	1.3%	1.1%	0.7%	0.5%	0.4%	1.3%	0.4%	2.1%
		PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS
	specified min area	5.00E+06	5.00E+05	5.00E+07	1.00E+08	1.00E+08	1.00E+08	1.00E+08	7.00E+08
		PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Mass Accuracy									
No.	Name	Mass Accuracy ppm	Mass Accuracy ppm	Mass Accuracy ppm	Mass Accuracy ppm	Mass Accuracy ppm	Mass Accuracy ppm	Mass Accuracy ppm	Mass Accuracy ppm
15	Control_Digest2_Rep1	EEQYN[dHex1Hex3HexNAc4]STYR	DTLM[Oxid]SR	DTLMISR	DIQMTQSPSTLSASVGDR	FNWYVDGVEVHNAK	WQQGNVFSC[Carboxymethyl]SYMHEALHNHYTK	GFYPSDIAVEWESNGQPENNYK	GFYPSDIAVEWESN
16	Control_Digest2_Rep2	-2.32	0.72	-1.17	-0.71	-1.09	-1.31	-2.01	1.17
17	Control_Digest2_Rep3	-1.94	0.79	-0.51	-0.45	-0.80	-0.46	-2.11	1.17
		-1.94	0.79	-0.88	-0.45	-1.16	-0.59	-1.92	1.17
	maximum mass error	2.32	0.79	1.17	0.71	1.16	1.31	2.11	1.17
		PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Custom MS  
Report Template



## CQA Overview


**ThermoFisher**  
SCIENTIFIC

### Sequence Details

Name:	HMSVictory	Created On:	09/Nov/16 14:49:34
Instrument:	QEHF_TheGhost	Updated On:	20/Dec/16 17:07:26
Imported Data:	False		
First Injection:	NISTmAb_Control1_DDA_1		
Processing Method:	SuperMAM		
MS Acquisition Time [min]:	70.51		
Method Length [min]:	115.00		
Total Time [hrs]:	63.25		
Data Vault:	ChromleonLocal	Created By:	Thermo
No. of Injections:	33	Updated By:	Thermo


### Deamidations and Isomerizations

No.	Name	PENNKY Deam HC N391	FNWY IsoD HC D283	FNWY Deam HC N283	NQVVLK Deam HC N78	VVSV Deam HC N317	WQQ Deam HC N436	YCAR Deam HC N85	VDNAL Deam LC N151	SGTAS Deam LC N136
		%	%	%	%	%	%	%	%	%
11	Control_Digest1_Rep1	1.58	0.11	0.04	0.03	0.05	0.23	0.00	0.21	8.3
12	Control_Digest1_Rep2	1.59	0.11	0.04	0.03	0.05	0.24	0.00	0.30	7.9
13	Control_Digest1_Rep3	1.58	0.11	0.03	0.03	0.05	0.24	0.00	0.32	8.5
15	Control_Digest2_Rep1	1.54	0.11	0.04	0.04	0.05	0.23	0.00	0.40	7.9
18	Control_Digest3_Rep1	1.48	0.10	0.04	0.04	0.04	0.40	0.00	0.38	8.4
24	Stressed_Digest2_Rep1	10.89	22.34	6.66	0.04	11.78	6.80	0.12	1.94	31.6
25	Stressed_Digest2_Rep2	11.03	22.14	6.55	0.04	11.98	6.86	0.12	1.68	31.1
26	Stressed_Digest2_Rep3	11.25	22.41	6.31	0.04	11.45	6.84	0.12	1.65	30.7




### Oxidations

No.	Injection Name	DTLM Oxid HC M253	DIQM Oxid LC M4	DMIF Oxid HC M100	FNWY Oxid HC W279
		%	%	%	%
11	Control_Digest1_Rep1	1.4	0.3	1.0	0.01
12	Control_Digest1_Rep2	1.4	0.3	1.0	0.01
13	Control_Digest1_Rep3	1.5	0.3	1.1	0.01
15	Control_Digest2_Rep1	1.3	0.2	0.7	0.01
18	Control_Digest3_Rep1	1.2	0.2	0.7	0.01
24	Stressed_Digest2_Rep1	2.9	0.2	1.1	0.02
25	Stressed_Digest2_Rep2	2.9	0.2	1.2	0.02
26	Stressed_Digest2_Rep3	3.0	0.2	1.2	0.02



### Lysine Modifications

No.	Injection Name	SLSSPGK Lysine Trunc.	VSNK HC K329 Gly	APK LC K44 Gly
		%	%	%
11	Control_Digest1_Rep1	28.4	0.6	0.1
12	Control_Digest1_Rep2	28.2	0.5	0.1
13	Control_Digest1_Rep3	27.2	0.6	0.1
18	Control_Digest3_Rep1	26.6	0.6	0.2
24	Stressed_Digest2_Rep1	25.6	0.3	0.1
25	Stressed_Digest2_Rep2	25.9	0.3	0.1
26	Stressed_Digest2_Rep3	25.5	0.3	0.1



Sequence Overview / Glycan Profile / Sheet / Master CQA Table / Integration

Custom MS  
Report Template

# Glycan Profile

## Glycan Profile

**ThermoFisher**  
SCIENTIFIC


### Sequence Details

Name:	AppsNotail	Created On:	01/Feb/17 17:00:03
Instrument:	TheGhostandtheDarkness	Updated On:	03/Feb/17 15:17:33
Imported Data:	False		
First Injection:	Control_One_1		
Processing Method:	SuperMAM		
MS Acquisition Time (min):	70.52		
Method Length (min):	115.00		
Total Time (hrs):	15.33		
Data Vault:	ChromleonLocal	Created By:	Thermo
No. of Injections:	8	Updated By:	Thermo

### EEQYNSTYR Glycopeptides

No.	Name	G0F - GlcNAc	G0F	G1F	G2F	G2F + Hex	G0	G1	M5	G2	M6	M7	Degly	Total
1	Control_One_1	5.2	38.8	39.5	9.0	1.8	0.1	0.1	0.9	0.04	0.05	0.02	0.6	100.0
2	Control_One_2	5.2	39.0	39.1	9.1	1.8	0.1	0.1	0.9	0.03	0.05	0.02	0.6	100.0
3	Control_One_3	5.3	38.4	39.6	9.3	1.8	0.1	0.1	0.9	0.05	0.05	0.02	0.6	100.0
4	Control_Two_1	5.2	38.7	39.5	9.1	1.8	0.1	0.1	0.9	0.05	0.06	0.03	0.6	100.0
5	Control_Two_2	5.4	38.8	39.3	9.0	1.7	0.1	0.1	0.9	0.04	0.05	0.02	0.6	100.0
6	Control_Two_3	5.4	38.7	39.2	9.1	1.8	0.1	0.1	1.0	0.05	0.06	0.02	0.6	100.0
7	Control_Three_1	5.3	38.7	39.3	9.0	1.8	0.1	0.1	0.9	0.05	0.06	0.03	0.6	100.0
8	Control_Three_2	5.3	39.1	39.2	8.9	1.7	0.1	0.1	0.9	0.05	0.05	0.02	0.6	100.0


% Abundance



Custom MS  
Report Template

## Sequence Template

## Targeted / non-targeted Processing Method




- eWorkflow automates the acquisition, data processing, and reporting processes
- Ultimate goal of a QC implementation

# Pack Everything Up

## Chromeleon 7.2



**Contains:**  
Sequence  
Methods (Instrument  
and processing)  
Report Formats  
Studio Layouts  
Raw Data  
And more...




# Enterprise Based Solution

## Chromeleon 7.2



**Remote:**  
Data Vaults  
Data Processing  
Licensing  
Administration



## Multi-Attribute Method

Not just about reducing the number of tests, but a better way to validate purity



### Insightful

Use high quality deeper knowledge of product to improve from development to production through QbD



### Easy to Use

Simple LCMS methods without the need for advanced chromatography or gas phase separations



### High Resolution

Take advantage of another dimension of separation to see what is missed.



### Powerful

Replace numerous conventional lot release techniques while providing greater knowledge of product attributes



### Previously Filed

Software/hardware already used for several filings to begin clinical trials with the FDA



### Confident

Composite peptide scoring and separation of isotopes/charge states helps eliminate false quantitation