



Accelerating **Innovation** & Enhancing Productivity

Overcoming Challenges of Manually Tracking the Performance of Your IC Consumables

John Madden and David Moore

Thermo Scientific Consumables Tracking System – Introduction

- New feature available on the recently introduced Thermo Scientific™ Dionex™ Integrion™ HPIC™ system
- Benefits of consumables tracking will be presented
- With the consumables tracking system you can:
 - Automatically track consumables without human intervention
 - Confirm installation
 - Monitor usage
 - Store and retrieve performance metrics



Why Should We Care About Consumables Tracking?

- The heart of any chromatography system is the consumables
- Consumables affect overall system performance:
 - Eluent generation and purification is achieved in RFIC systems with an EGC and CR-TC
 - Separation is achieved with a guard and separator column set
 - Detection is enhanced with an eluent suppressor and carbonate removal device



What is Consumables Tracking?

- SW identifies consumable install and maintains record of:
 - Installation
 - Usage
 - Performance Metrics
- Data is stored directly on the consumable itself
 - Ensures data is transferred with the consumable item
 - Ensures traceability of the data
- Takes the guesswork out of consumables monitoring
 - Provides information on the specifications
 - Records installation, use, and performance
 - Alerts operator when consumables are mismatched



How Does it Work?

- Memory chip embedded on each consumable
 - Instrument controller detects presence of memory chip
 - Chip downloads stored data
- Current date and time is written to memory tag
 - Consumable age can be calculated
 - Data is not written until consumable is first used
- Instrument controller writes key performance metrics to memory tag
 - Data can be accessed by instrument and controlling SW to monitor performance



Consumable Specification Data

- Comprehensive set of specification data in each device
 - Product Identification
 - Serial and Lot Number
 - Best Use-by Date
 - Column Diameters
 - Other Device Specifications
- Quality Assurance Report (QAR) results also stored
 - Thermo Scientific™ Dionex™ Chromeleon™ CDS software directly compares actual results with reported QAR data
- Easily retrieve data from front panel tablet app or by using the Chromeleon CDS software



Usage Tracking

- On first usage, a date and time stamp is written to consumable
 - Age of the consumable automatically tracked
 - Used consumable can be identified prior to operation
- System records number of injections while consumable is in use
- Usage data follows the consumable if moved between systems
 - Keeps records accurate



- Consumables Tracking System monitors key performance metrics
 - Date of First Install
 - Number of Injections
 - Key Metric Maxima
 - Last Columns Paired With
 - Key Metric Tallies
 - Eluent Types
 - Ion Count
 - Key Metric Trend Data
- Data can be stored on up to nine different consumables
 - Separator and Guard Columns
 - Concentrator and Trap Columns
 - EGC, CR-TC, Electrolytic Suppressors, Chemical Suppressors and CRD Devices



Configuration Compatibility

- Consumables Tracking System checks compatibility
- Warning issued if incompatibility is detected
- Two configuration restrictions are monitored
 - Size mismatches, such as combining 4 mm and 2 mm consumables
 - Chemistry mismatches, such as anion and cation
- Chromeleon CDS software issues configuration mismatches at warning level
 - Warnings do not prevent operation of the instrument



Case Study: Isaac – Switching Between Cations and Anions

- Lab performs compliance testing for drinking water using EPA Method 300.1
- Lab also receives requests for cation analysis
- The IC system is converted from anion analysis to cation analysis about every two weeks
- Anion method uses a carbonate eluent
- Cation method uses a methanesulfonic acid eluent
- After switching, it usually takes a few hours to flush out the system before analysis can resume



Switching Between Cations and Anions – Challenges

- Switchover from cation to anion analysis can cause issues:
 - Suppressor could be left in place after system is converted
 - Poor peak shapes compromise results
 - Cation guard left in place during anion analysis



Switching Between Cations and Anions – Solution

- Dionex Integrion HPIC system
 - Initial setup configured for anion analysis
 - After two weeks of operation, system converted for cation analysis
 - System operated for 24 hours in cation mode before being converted back to anion mode
 - After conversion, a warning was noted; system detected that the CERS 500 suppressor had been paired with the AG22 and AS22 columns
 - Operator was able to remove the cation suppressor and replace it with the anion suppressor before beginning work



Case Study: Jorge – Operates with Strict Preventative Maintenance Guidelines

- Jorge is a lab manager who manages up to four IC systems for monitoring steam water
 - He also manages up to 3 analytical chemists who operate the IC systems
- The IC systems are used primarily for anion analysis in steam water
- They cannot afford to be down for more than a few hours as it will affect downstream processes
- The regulated method requires that consumables are tracked and replaced on a regular basis as part of routine preventative maintenance
- Jorge is very sensitive to maximizing instrument up-time, but is also concerned about the budget for replacement of consumables



Operating with Strict Preventative Maintenance Guidelines – Challenges

- Log-book entries may not be effective
 - Responsible chemist may not reliably update log-book
- Consumables may be switched out before needed
 - During an audit any discrepancies require an immediate change of consumable if it is not possible to determine the consumable's exact age
- Risk is that perfectly good consumables are often being replaced unnecessarily
- Jorge is concerned that he is paying for consumables without being certain that the current replacement schedule is effective



Operates with Strict Preventative Maintenance Guidelines – Solution



- The lab just received a new Dionex Integrion RFIC HPIC system
- With the new consumables tracking feature, Jorge has switched to a digital log-book
 - The Dionex Integrion IC system automatically logs the installation date on each consumable without human intervention
- Logging mistakes are now eliminated and preventative maintenance can be scheduled more accurately



Case Study: Danielle – Minimizing Downtime

- Danielle runs a food testing laboratory that performs QC analysis on fruit juices and carbonated beverages
- Danielle runs one IC system for organic acid analysis
- The lab is very focused on maintaining a low cost of operation per sample while meeting its management's requirements for fast and accurate results
- Danielle is very sensitive to maximizing instrument up-time, but is also concerned about her budget for replacement of consumables



Minimizing Downtime – Challenges

- Consumables may be replaced unnecessarily only to find that it was a different consumable responsible for the issue
- On multiple occasions, Danielle suspected that the background conductivity was abnormally high
 - Unfortunately, she couldn't remember how long the system had been like that
 - She also couldn't remember what the conductivity was the last time the suppressor or CR-ATC was replaced
 - Sometimes when she notices this, the suppressor or CR-ATC fails soon afterward



Minimizing Downtime – A Solution

- The lab just received a new Dionex Integrion RFIC HPIC system to replace their aging IC
- 12 months after setting up the instrument, Danielle suspected that the background conductivity was abnormally high
 - Danielle queried the background trend data and determined that the background was indeed elevated
 - She determined that the suppressor voltage spiked at the same time that the background changed
 - The CR-TC current was unchanged
- Danielle ordered a replacement suppressor
 - The replacement suppressor restored the system to full working condition within just a few hours, thus reducing down-time



Consumables Tracking – Benefits

- A level of intelligence built into the system
- Logging and monitoring of individual consumables over time
- Trend plotting of key parameters over the lifetime of a consumable
- Monitoring of the age and usage of individual consumables



Remote Services Monitoring

- Unity Labs' Remote Services Monitoring is an add-on feature for Dionex Integrion IC systems
 - Enhances the Consumables Device Monitoring feature
- Tracks and stores instrument diagnostic data on a cloud server
 - Includes Consumables Device Monitoring parameters
 - Provides a mechanism for field service engineers to monitor key consumable parameters
- Parameters are monitored on an hourly basis
 - Data trends include a finer level of detail, making diagnostics more accurate



Remote Services Monitoring

- Provides data to your field service engineer
- Provides an alert service if any instrument parameter exceeds a predefined limit
- Remote Services and Monitoring is a value added feature and is included during the instrument warranty period (3 years)
 - The service is also available with Unity Lab Services “Critical Support” and “Essential Support” Service Plans for coverage beyond the warranty period



Remote Services Monitoring – Benefits

- Provides an alert service if any instrument parameter exceeds a predefined limit
- Generates trend plots of instrument parameters
- Allows a qualified Unity Lab Services technician to view instrument status



Complete Your Analytical Workflow

- High Purity Reagents and Chemicals

- NaOH, low in carbonate 50% solution
- High purity reagents
- High purity water

- Laboratory Equipment

- Gloves and safety glasses
- Syringe filters
- Vials



Please join me in the **Dionex Integrion IC** system section of our booth and I will address any further comments and questions.

