



ThermoFisher
S C I E N T I F I C

iQuan: Best Practices For Peptide Quantitation On a Triple Quadrupole Mass Spectrometer

The world leader in serving science

Application outline for peptide quantitation

- **(1) Background and Workflows**
- **(2) Peptide selection and standards**
- **(3) Collision energy (CE) optimization**
- **(4) Liquid chromatography (LC) gradient optimization**
- **(5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection**
- **(6) Experimental set up for peptide quantitation – an example workflow**
- **(7) Results**
- **(8) Data processing with Thermo Scientific™ TraceFinder 4.1™**

Application outline for peptide quantitation

- **(1) Background and Workflows**
- (2) Peptide selection and standards
- (3) Collision energy (CE) optimization
- (4) Liquid chromatography (LC) gradient optimization
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- (6) Experimental set up for peptide quantitation – an example workflow
- (7) Results
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™


Protein Quantitation with Confidence

- **Goal:** Addressing protein and peptide quantitation with the same ease as that of traditional small molecule quantitation
- **Challenges:**
 - Transitioning from traditional workflows for small molecules
 - Developing robust, sensitive, reproducible methods regardless of expertise
 - Addressing sample preparation complexity
 - Reducing cost/sample
- **Confidence from Sample Preparation, LC-MS, to data analysis and reporting**



Sample Prep with SMART Digest Kit: Why?

- Ensuring reproducibility is a critical challenge
 - Existing in-solution trypsin digestion of proteins is
 - Time consuming
 - Requires multiple steps (protein assay, denaturation, reduction and alkylation)
 - Labor resulting in increased chances of errors, and lack of reproducibility
- SMART Digest Kit
 - Highly reproducible
 - Quick and easy to use
 - Detergent free
 - Less prone to chemically-induced PTMs
 - Autolysis-free
 - Amenable to automation




<https://www.thermofisher.com/order/catalog/product/60109-101>

Four basic components of Peptide Quantitation


Challenge:

- Quantifying with confidence multiple targets
- Complex matrices
- Limited sample volume



How Do We Get Started?

Optimization Process



For Targets Without Empirical MS (Discovery data)...

Do the math!

- For the Top-Class QqQ available today, 10's-100's amol can be detected

Questions

- How much total sample do you need to load to get this range (or higher) of your target peptides?
- What do you need to do to enrich your sample to get ~100 amol on column?

Sample Enrichment

1 ug neat plasma (no enrichment)

Plasma Depletion (protein level)

Orthogonal Chromatography (peptide level)

Antibody Enrichment (peptide level)

Equivalent Sample Volume

~ 15 nL

~220 nL per inj.

~500 - 1500 nL per inj.

~100-1000 μ L per inj.

Application outline for peptide quantitation

- (1) Background and Workflows
- **(2) Peptide selection and standards**
- (3) Collision energy (CE) optimization
- (4) Liquid chromatography (LC) gradient optimization
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- (6) Experimental set up for peptide quantitation – an example workflow
- (7) Results
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™

How Many Peptide Candidates per Protein is Ideal?

- Considerations:
 - Can you get multiple peptides from each protein?
 - Do you need to look at site-specific PTMs?
 - What do you do if you can only detect 1-2 peptides per protein?



nature
biotechnology

Perspective

Nature Biotechnology **24**, 971 - 983 (2006)
Published online: 9 August 2006 | doi:10.1038/nbt1235

Protein biomarker discovery and validation:
the long and uncertain path to clinical utility

Nader Rifai¹, Michael A Gillette² & Steven A Carr²

Better biomarkers are urgently needed to improve diagnosis, guide molecularly targeted therapy and monitor activity and therapeutic response across a wide spectrum of disease. Proteomics methods based on mass spectrometry hold special promise for the discovery of novel biomarkers that might form the foundation for new clinical blood tests, but to date their contribution to the diagnostic armamentarium has been disappointing. This is due in part to the lack of a coherent pipeline connecting marker discovery with well-established methods for validation. Advances in methods and technology now enable construction of a comprehensive biomarker pipeline from six essential process components: candidate discovery, qualification, verification, research assay optimization, biomarker validation and commercialization. Better understanding of the overall process of biomarker discovery and validation and of the challenges and strategies inherent in each phase should improve experimental study design, in turn increasing the efficiency of biomarker development and facilitating the delivery and deployment of novel clinical tests.



HHS Public Access

Author manuscript

Nat Methods. Author manuscript; available in PMC 2014 August 01.

Published in final edited form as:

Nat Methods. 2014 February ; 11(2): 149–155. doi:10.1038/nmeth.2763.

Demonstrating the feasibility of large-scale development of standardized assays to quantify human proteins

Jacob J. Kennedy¹, Susan E. Abbatiello², Kyunggon Kim³, Ping Yan¹, Jeffrey R. Whiteaker¹, Chenwei Lin¹, Jun Seok Kim⁴, Yuzheng Zhang¹, Xianlong Wang¹, Richard G. Ivey¹, Lei Zhao¹, Hophil Min³, Youngju Lee⁴, Myeong-Hee Yu⁴, Eun Gyeong Yang⁴, Cheolju Lee⁴, Pei Wang¹, Henry Rodriguez⁵, Youngsoo Kim³, Steven A. Carr², and Amanda G. Paulovich¹

¹Fred Hutchinson Cancer Research Center, 1100 Fairview Ave. N., Seattle, WA 98109

²Broad Institute of MIT and Harvard, 7 Cambridge Center, Cambridge, MA 02142



³Department of Biomedical Engineering, Seoul National University College of Medicine, 28 Yongon-Dong, Seoul 110-799 Republic of Korea

⁴Center for Theragnosis, Korea Institute of Science and Technology, Seoul 136-791, Republic of Korea

⁵Office of Cancer Clinical Proteomics Research, National Cancer Institute, National Institutes of Health, Bethesda, MD 20892


Peptide Selection Criteria

Differentially Expressed Proteins from DDA



External Standards

Sample



SMART Digest Kit




LC/MS

Peak
Integration



Target Protein
Concentration



Standard



SMART Digest Kit




LC/MS


Peak
Integration



Standard 1







Stable Isotope Labeled (SIL) Internal Standards



<https://www.thermofisher.com/br/en/home/life-science/protein-biology/peptides-proteins/custom-peptide-synthesis-services/peptides-targeted-quantitation.html>

Advantages of SIL Standards

- Helps with ID due to same RT and fragment ratio pattern
- Serves as Internal Standard for Sample Prep
 - Each SIL spiked in at same concentration, every sample should have similar SIL Peak Area
 - Using SIL's for each peptide can improve precision (reproducibility) of quantitation



Questions



Application outline for peptide quantitation


- (1) Background and Workflows
- (2) Peptide selection and standards
- **(3) Collision energy (CE) optimization**
- (4) Liquid chromatography (LC) gradient optimization
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- (6) Experimental set up for peptide quantitation – an example workflow
- (7) Results
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™

Two methods for CE optimization

- Large number of targets
 - “Automated”
 - Specific softwares (Skyline)
 - LC multiple targets/run
- Few targets
 - Manual Optimization
 - Direct infusion
 - One target at a time

Collision Energy (CE) Optimization

- On a triple quadrupole MS, each product ion is made separately by CID, and transferred to Q3 and the detector




- CE can be optimized for each Q3 product ion individually
- Multiple product ions should be monitored in complex matrices
 - Confirms peptide ID and provides a means to determine if there are interferences

<http://www.nature.com/nmeth/journal/v9/n6/full/nmeth.2015.html>

CE optimization for large number of precursors

- Efficiently performed utilizing the LC system and software such as Skyline
- Start with >5 fmol on-column for nanoflow (above LOQ)
- Select fragment ions to target
- If >10 peptides require CE optimization, first find RTs for the LC method, then schedule CE optimization runs
- Most vendors provide a starting “CE regression equation” for their platform



https://skyline.gs.washington.edu/labkey/wiki/home/software/Skyline/page.view?name=tutorial_optimize_ce

Start with Recommended CE Regression Equations

- TSQ Quantiva™ has the following recommended CE regression equations for peptides at CID gas pressure of 1.5 mTorr.
- 2+ precursors: $CE = m/z * 0.0339 + 2.3597$
- 3+ precursors: $CE = m/z * 0.0295 + 1.5123$
- CE regression equations can be adjusted based on your peptide set.

Edit Collision Energy Equation

Name:

OK Cancel

Regression parameters:

	Charge	Slope	Intercept
▶	2	0.0339	2.3597
	3	0.0295	1.5123
*			

Optimization: _____

Step size: Step count:

Use Results Show Graph...

LC-SRM-MS Methods Can Be Generated and Imported

- Skyline generates the CE step gradient for each transition
- Transition lists can be imported into the TSQ Method Editor
- Automatic Method Exports from Skyline are available in TSQ v 2.1 software
- Product ion m/z is stepped in 0.01 m/z to associate steps in CE
- Large numbers of peptides (>10) may require scheduled methods for CE optimization

Method Editor

Method Timeline

Method Duration (min): 34


Experiment 1

SRM Table						
Compound	Start Time (min)	End Time (min)	Polarity	Precursor (m/z)	Product (m/z)	Collision Energy (V)
1 GGPFSDSYR(+2).-4	12.48	15.48	Positive	493.2	212.063	11.1
2 GGPFSDSYR(+2).-3	12.48	15.48	Positive	493.2	212.073	13.1
3 GGPFSDSYR(+2).-2	12.48	15.48	Positive	493.2	212.083	15.1
4 GGPFSDSYR(+2).-1	12.48	15.48	Positive	493.2	212.093	17.1
5 GGPFSDSYR(+2)	12.48	15.48	Positive	493.2	212.103	19.1
6 GGPFSDSYR(+2).1	12.48	15.48	Positive	493.2	212.113	21.1
7 GGPFSDSYR(+2).2	12.48	15.48	Positive	493.2	212.123	23.1
8 GGPFSDSYR(+2).3	12.48	15.48	Positive	493.2	212.133	25.1
9 GGPFSDSYR(+2).4	12.48	15.48	Positive	493.2	212.143	27.1
10 NLSVEDAAR(+2).-4	10.33	13.33	Positive	487.7	317.153	10.9
11 NLSVEDAAR(+2).-3	10.33	13.33	Positive	487.7	317.163	12.9
12 NLSVEDAAR(+2).-2	10.33	13.33	Positive	487.7	317.173	14.9
13 NLSVEDAAR(+2).-1	10.33	13.33	Positive	487.7	317.183	16.9
14 NLSVEDAAR(+2)	10.33	13.33	Positive	487.7	317.193	18.9
15 NLSVEDAAR(+2).1	10.33	13.33	Positive	487.7	317.203	20.9
16 NLSVEDAAR(+2).2	10.33	13.33	Positive	487.7	317.213	22.9
17 NLSVEDAAR(+2).3	10.33	13.33	Positive	487.7	317.223	24.9
18 NLSVEDAAR(+2).4	10.33	13.33	Positive	487.7	317.233	26.9
19 VSFELFADK (heavy)(+2).-4	19.34	22.34	Positive	532.2	341.151	12.3
20 VSFELFADK (heavy)(+2).-3	19.34	22.34	Positive	532.2	341.161	14.3
21 VSFELFADK (heavy)(+2).-2	19.34	22.34	Positive	532.2	341.171	16.3
22 VSFELFADK (heavy)(+2).-1	19.34	22.34	Positive	532.2	341.181	18.3
23 VSFELFADK (heavy)(+2)	19.34	22.34	Positive	532.2	341.191	20.3
24 VSFELFADK (heavy)(+2).1	19.34	22.34	Positive	532.2	341.201	22.3
25 VSFELFADK (heavy)(+2).2	19.34	22.34	Positive	532.2	341.211	24.3
26 VSFELFADK (heavy)(+2).3	19.34	22.34	Positive	532.2		


CE Step Gradient for On-Line Optimization

- Import data for acquired CE optimization runs back into Skyline
- Choose for “precursor” or “product” optimization (optimal CEs may be different)
- Save as a CE Optimization Library, and the settings will be exported for all methods with these optimized peptides.

Precursor Level Optimization



Product Level Optimization



Optimal CE for Each Fragment

Automated CE optimization complete using LC

The screenshot displays the Thermo Xcalibur Instrument Setup Method Editor interface. The main window is titled "Method Editor" and shows the "Scan Parameters" tab. The "Method Timeline" at the top indicates a duration of 40 minutes with a scan labeled "SRM" starting at 6.7 minutes and ending at 40 minutes. The "SRM Table" is the central focus, listing 19 scan entries with columns for Compound, Retention Time (min), RT Window (min), Polarity, Precursor (m/z), Product (m/z), and Collision Energy. A green box highlights the "Collision Energy" column, showing values ranging from 12.9 to 26.9. To the right, the "SRM Properties" panel is visible, containing various parameters such as Chromatographic Peak Width (sec), Use Cycle Time, Cycle Time (sec), Use Calibrated RF Lens, Q1 Resolution (FWHM), Q3 Resolution (FWHM), CID Gas (mTorr), Source Fragmentation (V), Use Chromatographic Filter, Use Retention Time Reference, Display Retention Time, Use Quan Ion, and Show Visualization. A "Copy Experiment Time" button is located at the bottom of this panel. The interface also includes a sidebar with scan type buttons (SRM, Full Scan Q1, Full Scan Q3, Product Ion Scan, Precursor Ion Scan, Neutral Loss Scan, SIM Q1, SIM Q3, QED, CO) and a status bar at the bottom left indicating "Ready".

	Compound	Retention Time (min)	RT Window (min)	Polarity	Precursor (m/z)	Product (m/z)	Collision Energy
1	TAENFR(+2)	9.1	2	Positive	369.182	322.187	13
2	TAENFR(+2)	9.1	2	Positive	369.182	436.23	18.9
3	TAENFR(+2)	9.1	2	Positive	369.182	565.273	12.9
4	TAENFR(+2)	9.1	2	Positive	369.182	636.31	14.9
5	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	332.196	13
6	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	446.239	18.9
7	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	575.281	12.9
8	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	646.318	14.9
9	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	466.277	26.9
10	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	563.33	22.9
11	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	660.383	16.9
12	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	731.42	18.9
13	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	802.457	18.9
14	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	476.286	26.9
15	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	573.338	22.9
16	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	670.391	16.9
17	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	741.428	18.9
18	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	812.465	18.9
19	EVPTI SPOCK(+3)	11.7	2	Positive	492.227	493.23	13.4

Questions



Manual Tune Page Optimization

Select precursor m/z

Quad resolution

Optimize RF lens
(best ion transmission
for each precursor)

Choose CID gas

Pressure

Collision energy range

CE step number

Click optimize to start

TSQ Quantiva Tune Application 2.1.1514

Thermo TSQ QUANTIVA SCIENTIFIC

ION SOURCE DEFINE SCAN CALIBRATION

Scan Optimization

Sample Injection Mode Syringe

Mass Input Option Formula (add adduct to formula) m/z

m/z value 493.768

Compound Name S5AAPPPPR(heavy)

Q1 Resolution (FWHM) 0.7

Charge State 2

Source Fragmentation 0

Precursor - Optimize RF Lens

Adjust Precursor Mass

Product

CID Gas (mTorr) 1.5

Product Input Option Unknown Product Ions Known Product Ions

Collision Energy Start (V) 10

Collision Energy End (V) 35

Use Collision Energy Step

Collision Energy Steps 10

Adjust Product Mass

Product Mass

	Product Mass
1	476.285
2	573.338
3	670.391
4	741.428
5	812.465

Optimize

- Infusion-based approach
- Samples must be more concentrated and have some organic solvent for better ESI performance
- MS ramps voltages
- XICs are plotted to determine optimal settings
- PDF report is generated
- Optimized settings can be pasted into Method Editor

Optimization Without Known Fragment Ions

Select Unknown Product Ions

Top N
(how many fragment ions to optimize)

Low mass exclusion
(doesn't select ions below this m/z)

Minor mass losses can also be excluded.

Click on Edit Mass List, select masses to exclude.

The screenshot displays the 'Product Input Options' section of a software interface. It includes radio buttons for 'Unknown Product Ions' (selected) and 'Known Product Ions'. Below are input fields for 'Collision Energy Start (V)' (5), 'Collision Energy End (V)' (55), 'Collision Energy Steps' (5), and checkboxes for 'Use Collision Energy Step', 'Adjust Product Mass', and 'Exclude Loss Masses'. The 'top N' field is set to 10, and 'Low Mass Exclusion' is set to 400. An 'Edit Mass List' button is visible next to the 'Exclude Loss Masses' checkbox. A status bar at the bottom indicates 'Optimization in progress...' and 'Optimization completed successfully.'. A 'Setup Mass Loss List' dialog box is open, showing a list of 'Excludes' with checkboxes: 'Water (18)', 'Ammonia (17)', and 'Methyl (15)' are checked, while two 'User Define' entries are unchecked. The dialog also has 'OK' and 'Cancel' buttons.

Example Optimization Report



Compound Optimization Report




Date & Time: Monday, April 11, 2016 11:45 AM
 Instrument Model: TSQ Quantiva
 Instrument Serial: TQH-Q1-0188
 Software Version: 2.0.1292.15

Source Parameters		Compound Optimization Input	
Source Type	H-ESI	Source Fragmentation (V)	0
Spray Voltage (V)	3500	Step Collision Energy Value	5
Sheath Gas (Arb)	10	Collision Energy Start (V)	5
Aux Gas (Arb)	0	Collision Energy End (V)	55
Sweep Gas (Arb)	0	Charge State	2
Ion Transfer Tube Temp (°C)	350	Step Collision Energy	Yes
Vaporizer Temp (°C)	25	Unknown/known Products	Known
		Exclude Loss Masses	None
		Compound Name	compound
		Product	Yes
		CID Gas (mTorr)	1.5
		Adjust Precursor Mass	Yes
		m/z Value	801.411
		Q1 Resolution	0.7
		Adjust Product Mass	Yes
		Number of Products	5
		Optimize RF Lens	Yes

Compound Optimization Results						
Compound Name	Precursor m/z	Product m/z	Collision Energy (V)	RF Lens (V)	Intensity	Source Fragmentation (V)
compound	801.533	772.315	28.202	117.506	5182.305	0
compound	801.533	928.45	28.051	117.506	2476.773	0
compound	801.533	1015.491	29.466	117.506	1815.119	0
compound	801.533	1072.531	29.77	117.506	7327.675	0
compound	801.533	1171.53	27.747	117.506	2064.839	0



Breakdown Curve for 801.533 to 928.45 transition



Serial Number: TQH-Q1-0188

Signature:
 Monday, April 11, 2016 11:45 AM
 1 of 5

Serial Number: TQH-Q1-0188


Signature:
 Monday, April 11, 2016 11:45 AM
 2 of 5

Application outline for peptide quantitation


- (1) Background and Workflows
- (2) Peptide selection and standards
- (3) Collision energy (CE) optimization
- (4) **Liquid chromatography (LC) gradient optimization**
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- (6) Experimental set up for peptide quantitation – an example workflow
- (7) Results
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™

Ion ratio confirmation in matrix background


Neat Standard




PRTC-60mingrad #18130-18865 RT: 23.46-24.40 AV: 49 NL: 2.01E5
F: + c NSI SRM ms2 422.432 [559.309-559.311, 674.329-674.331, 731.350-731.352]



Standard in Matrix




1fmole-PRTC2 #3262-3507 RT: 24.07-24.47 AV: 61 NL: 7.21E3
F: + c NSI SRM ms2 422.432 [559.309-559.311, 674.329-674.331, 731.350-731.352]




If CE optimization done neat, check signal in the matrix

- Ensure none of the target peaks have interference from the matrix background (peak splitting, shoulders).
- Modify the gradient to help shift interferences

Short Gradient




Long Gradient



Optimization of LC gradient

Goal: Maximize the time between targets and minimize the effect of the background matrix, while eluting in the shortest time possible.



Gradient Length	Chromatographic Resolution	Sensitivity	Reproducibility	Sample Throughput
Short	↓	↓	↓	↑
Long	↑	↑	↑	↓

Application outline for peptide quantitation

- (1) Background and Workflows
- (2) Peptide selection and standards
- (3) Collision energy (CE) optimization
- (4) Liquid chromatography (LC) gradient optimization
- **(5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection**
- (6) Experimental set up for peptide quantitation – an example workflow
- (7) Results
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™

Method Editor for Thermo Scientific™ TSQ Quantiva™

Method Editor for Thermo Scientific™ TSQ Quantiva™

Method Editor | Global Parameters | Scan Parameters | Summary

Method Timeline

Method Duration (min): 40

Method Timeline: 6.7, 13.3, 20, 26.7, 33.3, 40

Experiment 1

Scans: SRM, Full Scan Q1, Full Scan Q3, Product Ion Scan, Precursor Ion Scan, Neutral Loss Scan, SIM Q1, SIM Q3, QED, CD

SRM Table							
	Compound	Retention Time (min)	RT Window (min)	Polarity	Precursor (m/z)	Product (m/z)	Collision Energy
1	TAENFR(+2)	9.1	2	Positive	369.182	322.187	13
2	TAENFR(+2)	9.1	2	Positive	369.182	436.23	18.9
3	TAENFR(+2)	9.1	2	Positive	369.182	565.273	12.9
4	TAENFR(+2)	9.1	2	Positive	369.182	636.31	14.9
5	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	332.196	13
6	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	446.239	18.9
7	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	575.281	12.9
8	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	646.318	14.9
9	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	466.277	26.9
10	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	563.33	22.9
11	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	660.383	16.9
12	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	731.42	18.9
13	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	802.457	18.9
14	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	476.286	26.9
15	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	573.338	22.9
16	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	670.391	16.9
17	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	741.428	18.9
18	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	812.465	18.9
19	EXPT SCOOK(+3)	11.7	2	Positive	402.227	403.23	13.4

SRM Properties

- Chromatographic Peak Width (sec): 10
- Use Cycle Time:
- Cycle Time (sec): 1
- Use Calibrated RF Lens:
- Q1 Resolution (FWHM): 0.7
- Q3 Resolution (FWHM): 0.7
- CID Gas (mTorr): 1.5
- Source Fragmentation (V): 0
- Use Chromatographic Filter:
- Use Retention Time Reference:
- Display Retention Time:
- Use Quan Ion:
- Show Visualization:

Copy Experiment Time

Ready

Method Editor for Thermo Scientific™ TSQ Quantiva™

Method Editor | Global Parameters | Scan Parameters | Summary

Method Timeline

Method Duration: 1 | 6.7 | 13.3 | 20 | 26.7 | 33.3 | 40 | New

Compound	Retention Time (min)	RT Window (min)	Polarity	Precursor (m/z)	Product (m/z)	Collision Energy	
1	TAENFR(+2)	9.1	2	Positive	369.182	322.187	13
2	TAENFR(+2)	9.1	2	Positive	369.182	436.23	18.9
3	TAENFR(+2)	9.1	2	Positive	369.182	565.273	12.9
4	TAENFR(+2)	9.1	2	Positive	369.182	636.31	14.9
5	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	332.196	13
6	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	446.239	18.9
7	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	575.281	12.9
8	TAENFR (heavy)(+2)	9.1	2	Positive	374.186	646.318	14.9
9	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	466.277	26.9
10	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	563.33	22.9
11	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	660.383	16.9
12	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	731.42	18.9
13	SSAAPPPPPR(+2)	11.23	2	Positive	488.764	802.457	18.9
14	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	476.286	26.9
15	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	573.338	22.9
16	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	670.391	16.9
17	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	741.428	18.9
18	SSAAPPPPPR (heavy)(+2)	11.23	2	Positive	493.768	812.465	18.9
19	EXPT SCOOK(+3)	11.7	2	Positive	402.227	403.23	13.4

Experiment 1

SRM Table

SRM Properties

- Chromatographic Peak Width (sec): 10
- Use Cycle Time:
- Cycle Time (sec): 1
- Use Calibrated RF Lens:
- Q1 Resolution (FWHM): 0.7
- Q3 Resolution (FWHM): 0.7
- CID Gas (mTorr): 1.5
- Source Fragmentation (V): 0
- Use Chromatographic Filter:
- Use Retention Time Reference:
- Display Retention Time:
- Use Quan Ion:
- Show Visualization:

Copy Experiment Time

Ready

Parameter Selection for Targeted Peptide Quantitation

Set chromatographic peak width

Use cycle time for best sampling frequency

Set cycle time to get 10-15 points across the peak

If RF Lens was not optimized for each target, use calibrated value

Q1 and Q3 Resolution can be individually set for targets in the SRM table, or globally here

Choose the CID gas pressure you optimized with

Chromatographic Filter provides "on the fly" signal averaging

The screenshot shows the 'SRM Properties' configuration window. The parameters and their values are as follows:


Parameter	Value
Chromatographic Peak Width (sec)	10
Use Cycle Time	<input checked="" type="checkbox"/>
Cycle Time (sec)	1
Use Calibrated RF Lens	<input checked="" type="checkbox"/>
Q1 Resolution (FWHM)	0.7
Q3 Resolution (FWHM)	0.7
CID Gas (mTorr)	1.5
Source Fragmentation (V)	0
Use Chromatographic Filter	<input checked="" type="checkbox"/>
Use Retention Time Reference	<input type="checkbox"/>
Display Retention Time	<input checked="" type="checkbox"/>
Use Quan Ion	<input type="checkbox"/>
Show Visualization	<input type="checkbox"/>

Additional features highlighted in green boxes include 'Use Retention Time Reference', 'Use Quan Ion', and 'Show Visualization'. A 'Copy Experiment Time' button is located at the bottom right of the window.

**Features
new to
v2.1
software**

Cycle time determination


Use the peak widths from narrowest target peaks (at lower concentrations).



$$\frac{\text{peak width}}{\text{cycle time}} = \text{Scans across the peak}$$

$$9.6 \text{ seconds} / 0.64 \text{ seconds} = 15 \text{ scans}$$


Quadrupole Unit Resolution (0.7 FWHM)





Note: 0.7 FWHM is equivalent to 1.0 mass width at base of peak (m/z scale)

The Effect of Q1 Resolution on Intensity

- PRTC peptide mix spiked into HeLa and analyzed with the same gradient and method, only Q1 resolution was changed




Isobaric Discrimination – Effect of Increasing Resolution



SRM vs HSRM

- 200 amol on column GISNEGQNASIK[HeavyK]
- While the signal drops with Q1 setting of 0.2, the background and chemical noise are drastically reduced



SRM


Q1 = 0.7

H-SRM

Q1 = 0.2

Parameters to consider- Dwell Time

- Increasing Dwell Time improves S/N by increasing ion statistics
- Reproducibility of signal is also improved
- Dwell Times of 10 msec or greater provide high quality data for low abundant analytes
- TSQ Quantiva™ can detect signal with dwell times as low as 1 msec
- RT scheduling will increase dwell times vs unscheduled SRM methods



Dwell time can be variable using cycle time

$$\frac{\text{cycle time}}{\# \text{ product ions}} = \text{dwell time per fragment}$$

$$0.5/10 = 0.05 \text{ seconds}$$

$$0.5/50 = 0.01 \text{ seconds}$$

- Minimum dwell time should be 10 milliseconds for low abundance targets.
- Shorter dwell times have poor ion statistics and can make quantitation more difficult.
- We want to have optimal chromatographic resolution of targets and minimal RT overlap.

Determining Cycle Time when “setting” Fixed Dwell

Method: 1 ms Fixed Dwell

Transition List: 10 Peptide, 8 transitions each (80 Transitions total)

```
PRTC_10Peptides_3Transitions_lmsFixedDwell_lfmol_060414_1#48820 RT: 15.52
Total Ion Current:      10479.21
Scan Low Mass:         628.29
Scan High Mass:        1171.56
Scan Start Time (min): 15.52
Scan Number:           48820
Base Peak Intensity:   3288.78
Base Peak Mass:        1072.49
Scan Mode:             + c NSI SRM ms2 801.411 [628.291-628.293, 685.312-685.314, 772.344-772.346, 871.413-871.415, 928.434-928.436, 1015.466-1015.468, 1072.488-1072.490, 1171.556-1171.558]
```

Scan Header

```
TSQ Quantiva Data:
=====
Elapsed Scan Time
Average Scan by I
Micro Scan Count:
TSQ Quantiva Data:
=====
Elapsed Scan Time (sec):    0.018
Average Scan by Inst:      No
Micro Scan Count:          1
```

- Elapsed scan time of 18 ms for 8 transitions
- $18 \text{ ms} / 8 \text{ transitions} = 2.25 \text{ msec}$ dwell time per transition
- With a setting of 1 ms dwell, an interscan delay of $\sim 1.25 \text{ ms}$ is calculated
- Elapsed scan time of $18 \text{ ms} \times 10 \text{ peptides} = \sim 180 \text{ ms}$ Cycle Time
- *Note: The interscan delay will be dynamic and specific to the mass jump*

Application outline for peptide quantitation

- (1) Background and Workflows
- (2) Peptide selection and standards
- (3) Collision energy (CE) optimization
- (4) Liquid chromatography (LC) gradient optimization
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- **(6) Experimental set up for peptide quantitation – an example workflow**
- (7) Results
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™

Experimental Design Example for Peptide Quantitation

Goal: To quantify >360 peptides in a single QqQ method


- **Sample Details**

- HeLa Lysate Digest (Pierce, 88328)
 - 0.5 ug/uL
- Retention Time Calibration Peptides (Pierce, 88320)
 - Response curve from 25 amol-100 fmol/uL
 - Light versions spiked at a fixed 10 fmol/uL
- 6 x 5 LC-MS/MS Peptide Reference Mix (Promega, V7495)
 - Isoforms for each peptide ranged: 20 amol – 200 fmol/uL
- MS Qual/Quant QC Mix (Sigma, MSQC1-1VL)
 - 14 light and heavy peptides at various L:H ratios
 - Peptide concentrations ranged from 160 amol/uL – 300 fmol/uL


- **Discovery Experiment to Select Target Peptides**
 - Top 20 data dependent acquisition on a QE HF
 - 30K Resolution
 - Max fill time: 20 msec
 - AGC target: $1e^6$

Peptide Selection and Spectral Library Generation

Top 20 dd-MS2 Screen QE HF



Database Search



Peptide Selection

- Select peptides from search engine results
- Import list into Skyline Software


- Generate Spectral Library
- Select and refine peptide list
- Export targeted list of transitions



- Balance RT windows and number of transitions
- Adjust MS parameters for best cycle time
- Optimize CE (optional)




tSRM Method TSQ Quantiva



LC Conditions

- Trap-Elute configuration (Easy nanoLC 1000)
 - Trap column: Acclaim PepMap 100, 100 μm x 2cm, 5 μm beads
 - Analytical column: EasySpray PepMap, 75 μm x 25 cm, 2 μm beads, 50°C.
- Gradient
 - A = 2% ACN/0.1% formic;
 - B = 90% ACN/0.1% formic



- Total run time = 60 minutes

Pre-column equilibration

Volume: μl

Flow: $\mu\text{l} / \text{min}$

Max. pressure: Bar

Analytical column equilibration

Volume: μl

Flow: $\mu\text{l} / \text{min}$

Max. pressure: Bar

Sample pickup

Volume: μl (Max. is "loop size - 2 μl ")

Flow: $\mu\text{l} / \text{min}$

Sample loading

Volume: μl


Flow: $\mu\text{l} / \text{min}$

Max. pressure: Bar

Targeted SRM for 361 peptides (1842 transitions)

- Cycle time = 1 sec
- CID gas pressure = 1.5 mTorr
- Q1/Q3 Resolution: 0.7/0.7
- RT windows: 2 minutes
- Chrom Filter = 8
- Collision Energy = optimized by transition
- Dwell time range = 3.5 – 125 msec
- Total method time = 1 hour

RT: 8.95 - 36.67



NL:
7.79E7
TIC MS
071616_Qu
antiva_100
_fmoI_034


Questions




Application outline for peptide quantitation

- (1) Background and Workflows
- (2) Peptide selection and standards
- (3) Collision energy (CE) optimization
- (4) Liquid chromatography (LC) gradient optimization
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- (6) Experimental set up for peptide quantitation – an example workflow
- **(7) Results**
- (8) Data processing with Thermo Scientific™ TraceFinder 4.1™


Reproducibility – Peak Areas



Reproducibility – Peak Areas




Reproducibility – Retention Times





Reproducibility – Retention Times

- Eight XICs overlaid, show very good RT reproducibility
- SIL peptide is plotted
- Time range for shown data is 48 hours
- RT windows of 2 minutes were used, but could reduce to 1.5 or 1 min.




TSQ Quantiva™ : Consistent Sampling Frequency

PRTC peptide ELGQSGVDTYLQTK




Most congested part of chromatogram

Linear Response While Monitoring 1842 Transitions



Linear Response While Monitoring 1842 Transitions



Application outline for peptide quantitation

- (1) Background and Workflows
- (2) Peptide selection and standards
- (3) Collision energy (CE) optimization
- (4) Liquid chromatography (LC) gradient optimization
- (5) Thermo Scientific™ TSQ Quantiva™ method editor and parameter selection
- (6) Experimental set up for peptide quantitation – an example workflow
- (7) Results
- (8) **Data processing with Thermo Scientific™ TraceFinder 4.1™**

Summary

- Remember to follow key points when establishing your assay.
 - Careful selection of peptide targets, give yourself multiple options.
 - Select m/z fragments above precursor and optimize CE's for each.
 - Optimize chromatography, ensure no interference with fragment ions.
 - Utilize SIL internal standards if possible.
 - Establish standards for endogenous protein concentration.
 - Run "Scheduled" with RT windows when possible to maximize dwell times.
 - Ensure you have enough scans across the peak.
 - Utilize ion ratio confirmation in your processing method.
 - Quantitate utilizing fragment with best signal/noise.

Acknowledgements

Susan Abbatiello
Kristi D. Akervik
Brad Groppe
Nick Molinaro
Kent Seeley
Tara Schroeder
Alan Atkins
Katie Southwick
Detlef Schumann

Parts – Reversed Phase Columns

Part Number	Description	Column Dimensions (ID x L)	Packing Phase	Bead Diameter (µm)
ES800	EASY-Spray	75 µm x 15 cm	PepMap C18	3
ES801	EASY-Spray	50 µm x 15 cm	PepMap C18	2
ES802	EASY-Spray	75 µm x 25 cm	PepMap C18	2
ES803	EASY-Spray	75 µm x 50 cm	PepMap C18	2

- http://www.proxeon.com/productrange/easy_spray/intro/index.html

Parts – Reversed Phase Columns

Part Number	Description	Analytical Column Dimensions (ID x L)	Packing Phase	Bead Diameter (μm); Pore Size (\AA)
164738	Acclaim	75 μm x 15 cm	Acclaim PepMap 100 C18	3; 100
164739	Acclaim	75 μm x 50 cm	Acclaim PepMap 100 C18	3; 100
164940	Acclaim	75 μm x 15 cm	Acclaim PepMap 100 C18	2; 100
164941	Acclaim	75 μm x 25 cm	Acclaim PepMap 100 C18	2; 100
164942	Acclaim	75 μm x 50 cm	Acclaim PepMap 100 C18	2; 100
164568	Acclaim	75 μm x 15 cm	Acclaim PepMap 100 C18	3; 100
164569	Acclaim	75 μm x 25 cm	Acclaim PepMap 100 C18	3; 100

- <https://www.thermofisher.com/order/catalog/product/160321>

Parts – Reversed Phase Trap Columns

Part Number	Description	Trap Column Dimensions (ID x L)	Packing Phase	Bead Diameter (μm)
164535	Acclaim	75 μm x 2 cm	PepMap 100 C18	3
164564-CMD also known as 164564	Acclaim	100 μm x 2 cm	PepMap 100 C18	5

- <https://www.thermofisher.com/order/catalog/product/160321>

Parts – Metal Needle Kits and Ion Source Accessories

Part Number	Description
00950-00954	32-Gauge metal needle kit for high flow LC flow rates between 5 μ L/min to 400 μ L/min
97144-20040	34-Gauge metal needle kit for low flow LC flow rates between 500nL/min to 10 μ L/min
00106-10511	Tubing, 0.075mm ID x 0.193mm OD fused-silica capillary for low flow applications instead of metal needle
ES232	EASY-Spray Emitter Positioning Tool
ES235	EASY-Spray Emitter Wash Cap
106868-0064	Cleaving Stone, 1" x 1" x 1/32"