



ThermoFisher
S C I E N T I F I C

Conquer the Challenges of Small Molecule Analysis


Dr. Richard Jack

The world leader in serving science

Background


What are PFAS?

- PFASs are Per- and PolyFluorinated Alky Substances. Exclusively anthropogenic.
- Structures contain a hydrophobic perfluoroalkyl backbone and a hydrophilic end group
- Include a diverse range of compounds with a variety of chain lengths and end groups




Perfluorooctanoic acid

- PFOA
- Teflon®



8:2 Fluorotelomer sulfonate

- 8:2 FTS



Industrial Uses

- PFAS are used in a variety of applications because of their chemical and physical properties. These include:
 - Industrial polymers (Teflon® - PFOA)
 - Stain repellants (Scotch Guard® - PFOS)
 - Aqueous film forming foams (AFFF) – fire fighting applications

Sources

- Can be found anywhere at differing (generally lower) concentrations,
- Areas of elevated concentration and concern are:
 - Airports
 - Run-off from incidents of fire
 - Landfill leachate
 - WWTP effluent

PFAA Drinking - , Wastewater and Soil Clean up > Non-targeted and targeted



- **FEDERAL LEVEL ACTIONS**

- Water and Soil Method Validations
- EPA Off. Of Ground Water and Drinking Water
- EPA Office of Water
- Approx. 600 military bases in US have PFC contamination

- **STATE LEVEL ACTIONS**

- PFOS and PFOA in WW regulated in 7 States
- **California** - “Expressed support for including the **broader panel** of perfluoroalkyl and polyfluoroalkyl substances (PFASs).”



- **North Carolina** – New Legislation specifically focused on PFC/PFAA monitoring




SECTION 7.(d) The sum of eight million dollars (\$8,000,000) ...
(ii) ...of essential scientific instruments,
(iii) ...sample collection and analysis, training

In May 2017, Administrator Scott Pruitt established a task force to restore EPA's Superfund program to its rightful place at the center of the Agency's core mission to protect health and the environment.



Dr. P. Lee Ferguson
Assoc. Prof. of Civil and Environmental Engineering
Duke University



EPA is Proposal to Regulate PFAS in DW



U.S. ENVIRONMENTAL PROTECTION AGENCY
NEWS RELEASE
WWW.EPA.GOV/NEWSROOM

- Announced in Feb. 2019
- First ever comprehensive nationwide Action Plan to help states address concerns.
- Proposing a Federal Maximum Contaminant Level
- Begin the process to **propose** a regulation
- ***This doesn't mean PFAS are regulated!***

EPA's PFAS Action Plan: A Summary of Key Actions




EPA's PFAS Action Plan outlines concrete steps the agency is taking to address PFAS and to protect public health.

EPA's Per- and Polyfluoroalkyl Substances (PFAS) Action Plan:

- Demonstrates the agency's critical national leadership by providing both short-term solutions and long-term strategies to address this important issue.
- Provides a multi-media, multi-program, national research and risk communication plan to address this emerging environmental challenge.
- Responds to the extensive public input the agency has received over the past year during the PFAS National Leadership Summit, multiple community engagements, and through the public docket.

EPA is taking a proactive, cross-agency approach to addressing PFAS. The key actions EPA is taking to help provide the necessary tools to assist states, tribes, and communities in addressing PFAS are summarized below.



- Drinking Water
- Clean up
- Monitoring
- Research – *Water and Toxicology*
- Enforcement

Profiling

e.g.
500+
cpds

What's in my sample?
What knowledgebase do
I leverage to find out?

Screening

e.g.
100+
cpds

Are these compounds
really in my sample and
how much?


Targeted Quantitation

e.g.
20+
cpds

Am I absolutely
confident I should focus
on these compounds in
my sample?

Thermo Scientific™ Small Molecule Analysis Solutions

Profiling/Screening



Thermo Scientific™ Orbitrap™ Hybrid Mass Spectrometers

Targeted Quantitation



Thermo Scientific™ TSQ Triple Quadrupole Mass Spectrometers


Data Analysis



Thermo Scientific™ Compound Discoverer™ with mzCloud™ mass spectral library

mzCloud is a trademark of HighChem LLC, Slovakia

Orbitrap Mass Analyzer Technology for Screening of Unknowns




- High and ultra-high resolution enables discrimination between ions of interest and interfering ions in the very low and low m/z range
- Mass accuracy – the superior resolution enables accurate mass assignments with sub 1-ppm mass accuracy to eliminate false positives
- Retrospective analysis – enables investigation of new analytes in the same samples because it collects data on all analytes in the sample!

Learn More: www.thermofisher.com/orbitrap

Enhanced Resolution Using Orbitrap Technology

Resolution

$$R = \frac{m}{\Delta m}$$



- Quadrupole MS $R = \frac{m}{\Delta m} = \frac{500}{0.6} = 833$

- Orbitrap MS $R = \frac{m}{\Delta m} = \frac{500}{0.0035} = 140000$



Mass Accuracy

$$\Delta m / z = \frac{m_{meas} - m_{true}}{m_{true}} \cdot 10^6$$

- Quadrupole MS $\Delta m / z = \frac{5001 - 5000}{500} \cdot 10^6 = 200ppm$

- Orbitrap MS TOF MS $\Delta m / z = \frac{50010314 - 50010214}{50010314} \cdot 10^6 = 2ppm$

Resolution: 10k, 30k, 50k, 100k



Compounding Insights to Match Analysis Requirements

Compound Discoverer 3.0 Software for Small Molecule Unknown Identification

Efficiently extract high-confidence insights from information-rich small molecule HRAM data


Serves as a hub to seamlessly connect users to the tools they need to analyze productively and confidently



mzCloud Mass Spectral Library

Rank search more effectively with industry leading online spectra fragmentation library

mzCloud is a trademark of HighChem LLC, Slovakia




mzLogic Algorithm


From 1000's of candidates and hours of work to fast automated logical analysis

Learn More: www.thermofisher.com/compounddiscoverer


Identifying Unknowns – Spectral Libraries and Compound Databases




← Online spectral library
> 2.8 million spectra




← Chemical structure database with >480
data sources, 59 million structures



← Local spectral libraries



← Local databases



HRAM MS/MS and MSⁿ

HCD and CID fragmentation


Multiple Energy Levels

100% Professionally Curated


>8000 Compounds

>1M Fragment Structures


Ideally



Spectral Library Search



But what if ...



Spectral Library Search

NO reference data found in spectral library

Database Search
(HMDB, KEGG, BioCyc, PubChem)

441 hits
(based on formula)

MS conditions:

Spray voltage 3800V, Capillary temperature 295°C, sheath gas 32 au, Aux gas 7 au, S-lens RF level 55

Full MS/data dependent MS²

System templates

- Full MS - SIM
- AIF
- Full MS / AIF
- Full MS / dd-MS² (TopN)**
- Targeted-SIM
- Targeted-MS²
- Targeted-SIM / dd-MS²
- Full MS / AIF / NL dd-MS²

- full MS:resolution **140,000 FWHM**
scan range **100-1500 m/z**
- ddMS²:resolution **35,000 FWHM**
Inclusion list **on**
isolation **1 m/z**
exclusion **5 s**
NCE **10,40,200**

List of 146 PFAS target analytes


| | Mass (m/z) | Formula [M] | Formula type | Species | CS [z] | Polarity |
|----|------------|-------------|------------------|---------|--------|----------|
| 1 | 176.97731 | C3HF3O3 | Chemical formula | M-H | 1 | Negative |
| 2 | 112.9807 | C3HF3O2 | Chemical formula | M-H | 1 | Negative |
| 3 | 148.9826 | CF3SO3H | Chemical formula | M-H | 1 | Negative |
| 4 | 162.9823 | C3F5SO3H | Chemical formula | M-H | 1 | Negative |
| 5 | 176.9773 | C3HF3O3 | Chemical formula | M-H | 1 | Negative |
| 6 | 197.9804 | C3HF5SO2NH | Chemical formula | M-H | 1 | Negative |
| 7 | 198.9804 | C2HF5O3 | Chemical formula | M-H | 1 | Negative |
| 8 | 212.9792 | C4HF7O2 | Chemical formula | M-H | 1 | Negative |
| 9 | 228.97411 | C4HF7O3 | Chemical formula | M-H | 1 | Negative |
| 10 | 241.9816 | C4H9F7N3O3H | Chemical formula | M-H | 1 | Negative |
| 11 | 242.00723 | C6H5F7O3 | Chemical formula | | 1 | Negative |
| 12 | 244.9803 | C4HF7O4 | Chemical formula | M-H | 1 | Negative |
| 13 | 247.9816 | C3HF7NSO2H | Chemical formula | M-H | 1 | Negative |
| 14 | 247.9823 | C3HF7SO2NH | Chemical formula | M-H | 1 | Negative |
| 15 | 248.9461 | C3HF7SO3 | Chemical formula | M-H | 1 | Negative |
| 16 | 252.97601 | C3HF9O2 | Chemical formula | M-H | 1 | Negative |
| 17 | 264.0002 | C6H9F9O | Chemical formula | | 1 | Negative |
| 18 | 276.9411 | C4F7SO4H | Chemical formula | M-H | 1 | Negative |
| 19 | 278.97050 | C3HF9O3 | Chemical formula | M-H | 1 | Negative |
| 20 | 282.9481 | C4HF9O2 | Chemical formula | M-H | 1 | Negative |
| 21 | 291.9884 | C5H9F7NSO3H | Chemical formula | M-H | 1 | Negative |
| 22 | 297.9584 | C4HF11NSO2H | Chemical formula | M-H | 1 | Negative |
| 23 | 298.94298 | C4F9O3 | Chemical formula | | 1 | Negative |
| 24 | 310.9807 | C3HF9O5 | Chemical formula | M-H | 1 | Negative |
| 25 | 310.98078 | C3HF9O5 | Chemical formula | M-H | 1 | Negative |
| 26 | 312.97281 | C6H7ClO2 | Chemical formula | M-H | 1 | Negative |
| 27 | 325.96647 | C8H4F8O3 | Chemical formula | M-H | 1 | Negative |




146





Full Scan MS Data: Screening



Full Scan MS/MS Data:




GenX Gets Reported As HFPO Anion.



Finding Unknown PFAA's > Conc. than EPA Targeted List

**Continued: Showing
>1E+04<2.6E+06 area Count.**




Eaton Analytical

44 Total suspects/hits from a single sample (NC).

23

Targeted Quantitation Pain Points For Every Analytical Laboratory



NEW Thermo Scientific TSQ Triple Quadrupole Mass Spectrometers

Confident quantitation for any compound, any matrix, any user

Thermo Scientific™
TSQ Fortis™ Triple
Quadrupole MS



New level of confidence

Thermo Scientific™
TSQ Quantis™ Triple
Quadrupole MS



Everyday excellence

Thermo Scientific™
TSQ Altis™ Triple
Quadrupole MS




For what's now and what's next

Learn More: www.thermofisher.com/confidentquan

Sample Preparation

□ EPA Draft SW-846 Method 8327



* All standards were obtained from Wellington Laboratories

**Acrodisc GxF/0.2 µm GHP membrane syringe-driven filters were washed twice with LC-MS grade methanol (2x 10 mL) and acetonitrile (2x 10 mL)


*** Silanized-amber glass autosampler vials sealed with polypropylene caps were free of contaminants and interferences



Thermo Scientific™ Vanquish™ Flex Binary UHPLC System fitted with PFC-free kit and interfaced with the TSO Altis Triple Quadrupole Mass Spectrometer

Thermo Scientific TSQ Altis MS – Calibration Curves: Range 5 – 200 ppt

| Compounds | | LV 1 – 5 ppt % Deviation | LV2 - 10 ppt % Deviation | LV3 – 20 ppt % Deviation | RT (min) | R ² |
|-----------|----------|-----------------------------|-----------------------------|-----------------------------|----------|----------------|
| N-1 | PFTreA | 15% | 7% | -2% | 14.83 | 0.9906 |
| N-2 | PFTriA | 2% | 2% | 4% | 14.63 | 0.9966 |
| N-3 | PFDoA | -8% | 4% | 7% | 14.3 | 0.9989 |
| N-4 | PFUnA | -3% | 4% | -1% | 13.72 | 0.9996 |
| N-5 | PFDA | 12% | 1% | -5% | 13.03 | 0.9979 |
| N-6 | PFNA | 3% | 6% | -2% | 12.21 | 0.9983 |
| N-7 | PFOA | 13% | -1% | -2% | 11.22 | 0.9972 |
| N-8 | PFHpA | 1% | 6% | -2% | 9.91 | 0.9991 |
| N-9 | PFHxA | 8% | 0% | -1% | 7.94 | 0.9983 |
| N-10 | PFPeA | -16% | 4% | 14% | 4.98 | 0.9984 |
| N-11 | PFBA | 0% | -1% | 3% | 2.68 | 0.9993 |
| N-12 | PFDS | 4% | -14% | 11% | 13.7 | 0.9939 |
| N-13 | PFNS | 0% | -6% | 15% | 13.04 | 0.9976 |
| N-14 | PFOS | -7% | -4% | 9% | 12.24 | 0.9981 |
| N-15 | PFHpS | 16% | -7% | -9% | 11.3 | 0.9979 |
| N-16 | PFHxS | 13% | -5% | 0% | 10.11 | 0.9985 |
| N-17 | PFPeS | 5% | 4% | -3% | 8.42 | 0.9991 |
| N-18 | PFBS | 2% | 1% | -2% | 5.73 | 0.9995 |
| N-19 | PFOSA | 10% | 5% | -4% | 13.66 | 0.9931 |
| N-20 | FtS 8:2 | 6% | -3% | -1% | 13 | 0.9997 |
| N-21 | FtS 6:2 | 7% | 3% | -7% | 11.12 | 0.9977 |
| N-22 | FtS 4:2 | 23% | -10% | -9% | 7.66 | 0.9976 |
| N-23 | NEtFOSAA | 4% | -13% | 9% | 14.04 | 0.9985 |
| N-24 | NMeFOSAA | -10% | 8% | 6% | 13.64 | 0.9993 |



Overlay of all PFC compounds analyzed in this method

Calibration curves

- ✓ Linearity over the range 5 – 200 ppt
(2-fold dilution not taken into consideration)
- ✓ R² > 0.99 for all compounds
- ✓ Deviation < 20%


60 EPA water samples

| Ground Water | Reagent Water | Surface Water | Waste Water |
|--|--|--|--|
| <ul style="list-style-type: none">• 5x 5 mL native analytes• 5x 5 mL low spike• 5x 5 mL high spike | <ul style="list-style-type: none">• 5x 5 mL native analytes• 5x 5 mL low spike• 5x 5 mL high spike | <ul style="list-style-type: none">• 5x 5 mL native analytes• 5x 5 mL low spike• 5x 5 mL high spike | <ul style="list-style-type: none">• 5x 5 mL native analytes• 5x 5 mL low spike• 5x 5 mL high spike |

- ✓ Field samples divided into 3 batches of 20 samples
- ✓ 2 Method blanks, 2 LLOQs levels (10 and 20 ppt) and 2 Lab controls (LCS) were prepared for each batch

Perfluorotri- and Tetradecanoic Acids (PFTriDA and PFTeDA): Solubility Issues

Variable recoveries obtained due to their low solubility in water
(highest spike concentration was diluted in 50:50 methanol: water)



ICAL – LV6 (80 ppt)

BATCH 1 – LCS (80 ppt)

BATCH 2 – LCS (80 ppt)

BATCH 3 – LCS (80 ppt)



EPA Field Samples Results

**RSD < 20% for majority of the compounds among different water matrices.
All of the Reagent Water RSD <20%**

| N = 5 | GROUND WATER | | | | REAGENT WATER | | | | SURFACE WATER | | | | WASTE WATER | | | |
|--------------|-----------------|-------|------------------|-------|-----------------|-------|------------------|-------|-----------------|-------|------------------|-------|-----------------|-------|------------------|-------|
| | Low level spike | | High level spike | | Low level spike | | High level spike | | Low level spike | | High level spike | | Low level spike | | High level spike | |
| | Average (ppt) | RSD % | Average (ppt) | RSD % | Average (ppt) | RSD % | Average (ppt) | RSD % | Average (ppt) | RSD % | Average (ppt) | RSD % | Average (ppt) | RSD % | Average (ppt) | RSD % |
| N1_PFTeDA | 26.13 | 12% | 69.5 | 12% | 24.43 | 10% | 71.84 | 6% | 26.78 | 18% | 82.63 | 13% | 23.92 | 11% | 77.91 | 4% |
| N2-PFTrDA | 22.81 | 11% | 65.76 | 12% | 22.29 | 9% | 74.45 | 6% | 23.8 | 17% | 76.87 | 15% | 22.24 | 8% | 77.47 | 7% |
| N3-PFDoA | 21.22 | 8% | 64.08 | 12% | 20.01 | 9% | 73.25 | 9% | 22.17 | 19% | 71.55 | 14% | 20.54 | 13% | 75.8 | 9% |
| N4_PFUdA | 22.05 | 11% | 65.41 | 7% | 21.23 | 12% | 72.08 | 9% | 23.19 | 21% | 69.81 | 11% | 21.74 | 12% | 74.95 | 9% |
| N5_PFDA | 22.98 | 9% | 64.83 | 8% | 21.63 | 11% | 72.83 | 8% | 23.55 | 17% | 69.98 | 11% | 23.16 | 3% | 76.89 | 7% |
| N6-PFNA | 22.29 | 8% | 66.83 | 8% | 21.21 | 9% | 73.4 | 7% | 23.15 | 16% | 70.65 | 10% | 22.7 | 9% | 75.64 | 7% |
| N7-PFOA | 22.89 | 12% | 65.88 | 10% | 21.26 | 9% | 71.72 | 9% | 24.15 | 15% | 69.87 | 9% | 27.08 | 8% | 79.79 | 7% |
| N8-PFHpA | 23.34 | 11% | 65.96 | 11% | 21.89 | 10% | 72.23 | 8% | 25.51 | 18% | 69.48 | 11% | 26.35 | 12% | 77.48 | 9% |
| N9-PFHxA | 22.33 | 12% | 64.9 | 10% | 21.05 | 8% | 71.47 | 8% | 23.36 | 16% | 69.43 | 11% | 40.48 | 7% | 93.5 | 8% |
| N10_PFPeA | 27.41 | 13% | 63.44 | 11% | 23.42 | 12% | 71.5 | 10% | 32.24 | 23% | 67.87 | 9% | 38.32 | 3% | 84.01 | 10% |
| N11_PFBA | 18.99 | 24% | 64.11 | 4% | 21.04 | 10% | 70.09 | 6% | 19.59 | 16% | 63.82 | 6% | 14.01 | 45% | 66.61 | 6% |
| N12-PFDS | 20.79 | 14% | 58.96 | 20% | 18.2 | 14% | 60.87 | 12% | 22.02 | 24% | 56.67 | 13% | 19.43 | 24% | 69.12 | 14% |
| N13-PFNS | 20.32 | 27% | 56.17 | 16% | 19.67 | 8% | 59.42 | 10% | 20.08 | 17% | 57.1 | 12% | 18.19 | 20% | 63.55 | 9% |
| N14-PFOS | 22.57 | 20% | 60.07 | 8% | 21.26 | 16% | 63.04 | 9% | 24.71 | 14% | 62.52 | 10% | 24.29 | 10% | 75.43 | 13% |
| N15-PFHpS | 20.81 | 12% | 61.61 | 9% | 19.53 | 11% | 67.95 | 10% | 21.7 | 17% | 63.27 | 10% | 19.43 | 15% | 71.37 | 10% |
| N16-PFHxS | 19.69 | 15% | 59.56 | 10% | 18.38 | 10% | 63.56 | 9% | 20.89 | 14% | 61.66 | 7% | 19.94 | 11% | 67.58 | 9% |
| N17-PFPeS | 20.32 | 11% | 62.4 | 9% | 19.35 | 6% | 65.83 | 9% | 21.52 | 18% | 63.48 | 8% | 19.87 | 10% | 68.3 | 9% |
| N18-PFBS | 22.55 | 13% | 61.07 | 10% | 19.98 | 11% | 63.39 | 10% | 21.34 | 26% | 61.46 | 9% | 24.88 | 15% | 52.03 | 36% |
| N19-PFOSA | 20.51 | 15% | 60.19 | 12% | 18 | 10% | 60.89 | 12% | 22.6 | 17% | 63.44 | 10% | 19.66 | 10% | 67.62 | 9% |
| N20_FtS8_2 | 22.18 | 7% | 63.6 | 6% | 21.55 | 10% | 69.86 | 8% | 20.12 | 18% | 66.8 | 8% | 21.51 | 8% | 71.42 | 7% |
| N21-FtS6_2 | 20.64 | 13% | 64.29 | 11% | 21.46 | 17% | 69.2 | 7% | 29.54 | 37% | 80.12 | 13% | 33.15 | 25% | 82.55 | 11% |
| N22-FtS4_2 | 21.11 | 6% | 62.04 | 5% | 20.4 | 11% | 69.78 | 6% | 21.02 | 15% | 65.49 | 11% | 19.75 | 7% | 69.17 | 8% |
| N23_NeFOSAA | 23.66 | 12% | 63.56 | 11% | 22.96 | 18% | 70.55 | 7% | 23.72 | 25% | 70.56 | 14% | 22.16 | 16% | 76.7 | 8% |
| N24_NMeFOSAA | 23.88 | 9% | 65.42 | 9% | 21.8 | 8% | 73.34 | 7% | 23.97 | 17% | 71.28 | 8% | 22.9 | 10% | 73.38 | 7% |

- ✓ **TSQ Altis Triple Quadrupole MS** is able to quantitate at 5 ppt the list of PFCs listed in EPA draft SW-846 method 8327 using direct injection (2-fold dilution not taken into consideration)
- ✓ **TSQ Altis Triple Quadrupole MS** is able to quantitate certain PFCs 5 times lower than the LLOQ reported by EPA
- ✓ **Thermo Scientific™ Accucore™ RP-MS LC Columns** provide similar performance as the Waters Acquity CSH Phenyl-hexyl described in EPA draft SW-846 method 8327
- ✓ Retention time stability was very good inter-batch
- ✓ Large injection volumes overload the analytical columns and chromatographic peak fronting is observed. Reduced injection volumes maintain sensitivity and improve assay robustness.
- ✓ Longest chain perfluorocarboxylic acids showed high variability within batches mainly due their low solubility in water (higher spike concentration)
- ✓ EPA field samples showed RSDs below 20% for most of the compounds among different water matrices

Thank You for Listening




Questions

For more information visit:
<http://thermofisher.com/quantitation>

Thank You

Please return our survey to receive a drink ticket for our daily networking event where you can continue discussions with our experts!



Visit Booth #2632

Join us Tuesday through
Thursday from 3:30 – 5:00 p.m.
to collaborate with our technical experts

Redeem this ticket
for a complimentary beverage!

ThermoFisher
SCIENTIFIC