


Thermo Scientific Orbitrap Exploris 4.3 Instrument Control Software (OES 4.3 ICSW) –

Overview

February 2024


 The world leader in serving science



Thermo Scientific Orbitrap Exploris MS Portfolio – one ICSW



Thermo Scientific™ Orbitrap Exploris™ MX
Mass Detector



Thermo Scientific™ Orbitrap Exploris™ 120
Mass Spectrometer



Thermo Scientific™ Orbitrap Exploris™ 240
Mass Spectrometer



Thermo Scientific™ Orbitrap Exploris™ 480
Mass Spectrometer



Thermo Scientific™ Orbitrap Exploris™ GC
Mass Spectrometer




For Research Use Only. Not for use in diagnostic procedures.



Thermo Scientific™ Orbitrap Exploris™ GC 240
Mass Spectrometer


Comparison of Orbitrap Exploris Portfolio



Orbitrap Exploris MX



Orbitrap Exploris 120



Orbitrap Exploris 240



Orbitrap Exploris 480

	Orbitrap Exploris MX	Orbitrap Exploris 120	Orbitrap Exploris 240	Orbitrap Exploris 480
Max Resolution (FWHM) @ m/z 200	180,000	120,000	240,000	480,000
Mass range	40 – 3,000 (8,000 *)	40 – 3,000	40 – 6,000 (8,000 *)	
Precursor ion selection	n/a	≤ 2,500		
Sensitivity		S/N 250 @ 200 fg reserpine (tSIM)		S/N 150 @ 50 fg reserpine (tSIM)
MSMS scan rate (Hz)	22 Hz (Full Scan)	22 Hz		40 Hz
Mass accuracy - external	< 3 ppm RMS drift over 24 hours			
Mass accuracy w/ EASY-IC - internal	< 1 ppm over 5 days			
Spectral multiplexing	n/a	20		
Polarity switching * : one cycle equals (pos./switch/neg./switch)	60 k Full Scan* < 700 ms (equals > 1.4 Hz)	60 k Full Scan* < 700 ms (equals > 1.4 Hz) 60 k tSIM Scan* < 600 ms (equals > 1.6 Hz)		
Calibration	One-click calibration with FlexMix and dedicated calibration probe - with harmonization and improved user experience across all TNG platforms (TSQs, Hybrids, Tribbrids)			
One-Point Mass Calibration	One-Point (Self) Mass Calibration achieves < 3 ppm RMS drift over at least 4 weeks			
Scan modes	<ul style="list-style-type: none"> Full Scan 	<ul style="list-style-type: none"> Full Scan ddMS2 (Top1-4) tSIM (targeted mass list) ddMS2 (Top1-4) Full Scan ddMS2 (targeted list) (Top1-4) 	<ul style="list-style-type: none"> Full Scan ddMS2 (topN) Full Scan ddSIM tSIM (targeted mass list) ddMS2 Full Scan ddMS2 (targeted mass list) With options for <ul style="list-style-type: none"> 'Number of Scans' (= TopN) 'Cycle Time' 	<ul style="list-style-type: none"> Full Scan ddMS2 (topN) Full Scan ddSIM tSIM (targeted mass list) ddMS2 Full Scan ddMS2 (targeted mass list) With options for <ul style="list-style-type: none"> 'Number of Scans' (= TopN) 'Cycle Time' 'Scans per Outcome' (branching)
Full MS AIF t-SIM DIA MS2				
combinable within in one single experiment, such as:	In addition, multiple experiments can be created combining various Full Scan experiments	In addition, up to 5 experiments can be created combining the above listed scan types	In addition, multiple experiments can be created combining the above listed scan types	In addition, multiple experiments can be created combining the above listed scan types
Advanced acquisitions	APD	AcquireX (chargeable option)	AcquireX, APD, AcquireX AB TMT @ 45k resolution setting	<ul style="list-style-type: none"> 16 msec transient (7,500 min resolution) System Templates supporting BoxCar and SureQuant approaches TurboTMT with TMT reagents up to 18-plex

Release Notes

Thermo Scientific Orbitrap Exploris Series 4.3 Instrument Control Software Release Notes							
Version	Version No.	Orbitrap Exploris 480	Orbitrap Exploris 240	Orbitrap Exploris 120	Orbitrap Exploris GC	Orbitrap Exploris GC 240	Orbitrap Exploris MX
1.0	1.0.77.7	✓	—	—	—	—	—
1.1	1.1.117.22	✓	—	—	—	—	—
1.1 SP1	1.1.117.26	✓	—	—	—	—	—
2.0	2.0.182.18	✓	✓	✓	—	—	—
2.0 SP1	2.0.182.25	✓	✓	✓	—	—	—
2.0 SP2	2.0.182.35	✓	✓	✓	—	—	—
3.0	3.0.261.13	✓	✓	✓	✓	✓	—
3.1	3.1.279.9	✓	✓	✓	✓	✓	—
4.0	4.0.309.27	✓	✓	✓	✓	✓	✓
4.0 SP1	4.0.309.28	✓	✓	✓	✓	✓	✓
4.1	4.1.335.19	✓	✓	✓	✓	✓	✓
4.2	4.2.362.16	✓	✓	✓	✓	✓	✓
4.2 SP1	4.2.362.21	✓	✓	✓	✓	✓	✓
4.2 SP2	4.2.362.26	✓	✓	✓	✓	✓	✓
4.2 SP3	4.2.362.36	✓	✓	✓	✓	✓	✓
4.2 SP4	4.2.362.42	✓	✓	✓	✓	✓	✓
4.3	4.3.458.15	✓	✓	✓	✓	✓	✓



Source: Release Notes for OES 4.3 ICSW

Thermo Scientific Orbitrap Exploris Series 4.3 Instrument Control Software Release Notes

Installation Notes

Supported Target Systems

Thermo Scientific Orbitrap Exploris 120 mass spectrometer
Thermo Scientific Orbitrap Exploris 240 mass spectrometer
Thermo Scientific Orbitrap Exploris 480 mass spectrometer
Thermo Scientific Orbitrap Exploris MX mass spectrometer
Thermo Scientific Orbitrap Exploris GC mass spectrometer
Thermo Scientific Orbitrap Exploris GC 240 mass spectrometer

System Requirements

The minimum hardware and software configurations required for the Orbitrap Exploris Series 4.3 Instrument Control Software operation are as follows:

System	Requirements
PC	3.0 GHz Quad Core Intel™ Processor 32 GB RAM 512 GB SSD Hard Drive Display Monitor Resolution of 1920 × 1080 Two Network Interface Cards (NIC), 1000 MBit/s
Software	Microsoft™ Windows™ 10 Enterprise 2016 LTSC, 2019 LTSC or 2021 LTSC Thermo Scientific Xcalibur 4.7

Tip The Orbitrap Exploris Series 4.3 Instrument Control Software was only tested within the delivered composition.

Note:
Xcalibur 4.7 software applies
Foundation 3.1 SP9.

Source: Release Notes for OES 4.3 ICSW

List of New Features and Improvements in OES 4.3 ICSW

New Features

General

- Operating Manual, Pre-Installation Requirements Guide, and Software manuals are updated
- Additional resolution settings are also accessible for OE 120 and OE MX (11.25k, 22.5k, and 90k) in Tune and Method Editor

Method Editor

- Display of absolute AGC target
- MaxIT dynamic for FSddMS2 experiments
- Lock Mass injection per experiment

OE GC:

- ME: User Defined LM correction - RunStart mode
- ME: Production Ions Scan - Product Ion Window Modes
- Tune: improved Manual Calibration tables separated by polarity
- Tune: Improved Leak Check

Improvements

General

- Scheduled one-point mass calibration (via Tune Preferences)
 - Addition of a Calendar view

Tune

- System calibration
- Tune Preferences: Hotlink function is now limited to the Ion Source settings

Method Editor

- EASY-IC “On/Off” available for all experiments in the timeline
- EASY-IC “On/Off” available for Lock Mass Correction Mode: Timed (for all experiments in the timeline)
- Method Execution: the instrument now preserves the Ions Source settings and polarity after the sample acquisition has finished.

Diagnostics

- All MALDI related functionalities are now applicable to OT Exploris MX as well

Display of the Absolute AGC Value

AGC Target	Custom	★
Normalized AGC Target (%)	100	★
Absolute AGC Value	1.000e6	★

- If the user chooses AGC Target „Custom“, the absolute AGC Value will be displayed underneath.
- The absolute AGC target value is a conversion of the AGC percentage target set in the method into absolute values. This read-only value will be automatically updated according to the scan type and the target percentage.

Isolation width for DDA/tMS2 for OE240 is extended

Targeted MS² Scan Properties [Show Favorites](#)

tMS²

Multiplex Ions ★

Isolation Window (m/z) ★

Offset ★

Energy Type ★

Specify the Isolation Window for the precursor ions
Range: 0.4-2000
[Learn more...](#)

- In tMS2 and FSddMS2 experiments the isolation width has been extended to 2000 Da

Additional Resolution Settings

➔ now also for OE120 and MX



Added Resolution to Provide More Flexibility in Short LC Gradient Analysis

Tune

thermo scientific
Orbitrap Exploris 480

ION SOURCE DEFINE SCAN CALIBRATION

Scan Type: Full Scan

Orbitrap Resolution: 15000

Scan Range (m/z): 7500

RF Lens (%): 11250

AGC Target: 22500

Maximum Injection Time: 30000

Time (ms): 45000

Microscans: 90000

Source Fragmentation: 120000

Use EASY-IC™: 180000

Method Editor

Full Scan Properties Show All

Orbitrap Resolution: 120000

RF Lens (%): 7500

Polarity: 11250

15000

22500

30000

45000

60000

90000

120000

180000

240000

480000

OE120

Full Scan Properties Show All

Orbitrap Resolution: 120000

Scan Range (m/z): 15000

RF Lens (%): 22500

30000

Polarity: 45000

60000

90000

120000

OE MX

Full Scan Properties Show All

Orbitrap Resolution: 120000

Scan Range (m/z): 15000

RF Lens (%): 22500

30000

Polarity: 45000

60000

90000

120000

180000

- Reminder
 - General concept → Only when the full calibration is passed calibration values are stored
 - If a procedure does not pass, repeating the „system calibration“ will resume from a checkpoint before the not passed procedure
 - Update of the date for the next recommended System calibration when ALL procedures have passed
- **NEW**: In negative mode the calibration values get saved even if the ICS HV calibration has not passed. The recommended date will not be updated until also the ICS HV calibration has passed.

„Use EASY-IC“ On/Off

Improvement: for **EACH** experiment

Targeted MS ² Scan Properties		Show All
Isolation Window (m/z)	2	
Collision Energy Type	Normalized	▼
HCD Collision Energy (%)	30	
Orbitrap Resolution	15000	▼
Scan Range Mode	Auto	▼
RF Lens (%)	70	
Polarity	Positive	▼
Use EASY-IC™	On	▼

- If On is selected, it provides an internal reference mass that is used for mass correction during a run (Lock Mass Correction: EASY-IC™, Mode ≠ RunStart)
- Now available for each experiment in the timeline – here displayed for a tMS2 experiment
- This provides the possibility in mixed experiments to perform Scan-to-Scan EASY-IC only in the FullScan experiment and not in an e.g. additional tMS2 experiment
- If no lock mass is found in one scan, the system will apply the last successful locking information to this scan. Time duration of last locking and lock mass correction are provided in the scan header of the individual scan.

„Use EASY-IC“ On/Off - Timed


Improvement: also for Timed EASY-IC “On/Off“ can be chosen

Targeted MS ² Scan Properties		Show All
Isolation Window (m/z)	2	
Collision Energy Type	Normalized	▼
HCD Collision Energy (%)	30	
Orbitrap Resolution	15000	▼
Scan Range Mode	Auto	▼
RF Lens (%)	70	
Polarity	Positive	▼
Use EASY-IC™	On	▼


- If “On” is selected, it provides an internal reference mass that is used for mass correction during the defined time window
- If no lock mass is found in one scan, the system will apply the last successful locking information to this scan. Time duration of last locking and lock mass correction are provided in the scan header of the individual scan.

Lock Mass Injection parameter moved to scan properties

≤ OES 4.2



OES 4.3



- Lock Mass Injection function available for every experiment in the timeline
- Whereas a FullScan might not need the custom lock mass to be injected (no time penalty) the lock mass could be injected for tMS2 experiments.
- Please Note: EASY-IC™ injects the lock mass to every chosen experiment (see previous slide)

- Additionally we renamed „Internal Mass Calibration“ to „Lock Mass Correction“ as the old naming was confusing for many customers.

MaxIT Mode „Dynamic“ for FSddMS2*

No additional properties are exposed when selecting dynamic Max IT but method execution is affected

Number of Scans

- Max IT for every scan in the experiment is scaled by $(1+TopN)/(1+SumCand)$
- TopN = user-defined max. number of dd scans, SumCand = sum of actually available candidates from master scan,
- Increment of 1 to account for the master scan itself
- Scaling factor must be >1 to take effect.
- **Example:**
with TopN=100 but only 50 available candidates in a cycle, the expected max. IT values that result from parallel acquisition of scans would be scaled by a factor of ~ 2
- \rightarrow for a Resolution of 15k (ddMS2) the maxIT will be extended from 32 ms to $\sim 2 \times 32 = \sim 64$ ms (minus overheads)

Cycle time

- MaxIT for every scan is determined as $\max(\text{“Expected_max_IT”}, \text{Cycle Time} / \text{SumCand})$
- With “Expected_max_IT” = maxIT value pertaining to parallel acquisition.
- **Example:**
Cycle Time: 0.6 sec; 15k Res – 32 ms transient length
 $\rightarrow \sim 18$ MS2 spectra (theoretical) with a maxIT of 32 ms
If only 9 candidates available \rightarrow maxIT extended by a factor of ~ 2

*previously only available for targeted experiments.

One-Point Self Calibration

Improved Capabilities – Calendar View

The image shows a software interface for configuring self-calibration. On the left is a dialog box titled "Mass Self-Calibration Options". It contains the following elements:

- A header bar with the title "Mass Self-Calibration Options".
- Text: "System Self-Calibration is enabled. To abort or change options, please uncheck the system self calibration check box below".
- A checked checkbox: "Run One-Point Mass Self-Calibration".
- A section titled "Schedule Self-Calibration" containing:
 - "Cadence" with radio buttons for "Daily" and "Weekly" (selected).
 - "Time" with a dropdown menu currently showing "12:00 AM".
 - "Day" with checkboxes for "Monday", "Tuesday", "Wednesday", "Thursday", "Friday", "Saturday", and "Sunday".

On the right is a vertical time selection dropdown menu showing a list of times from 12:00 AM to 6:00 AM in 30-minute increments. The list is scrollable, with 12:00 AM at the top and 6:00 AM at the bottom.

- The One-point Self-calibration can be set up in a more detailed way
- E.g. it can be set to be performed every Monday, Wednesday and Friday.
- Also, there are more options for the timing.


Updates to Manuals for Orbitrap Exploris Series

Pre-Installation Requirements Guide

Operator Manuals

Model specific **Software Manuals** and **online help**

 The world leader in serving science




Pre-Installation Requirements Guide And Operating Manual


Pre-Installation Requirements Guide



Operating Manual



Software Manuals for these Orbitrap Exploris models



Thermo Scientific™
Orbitrap Exploris™ MX Mass Detector




Thermo Scientific™ Orbitrap Exploris™ 120
Mass Spectrometer



Thermo Scientific™ Orbitrap Exploris™ 240
Mass Spectrometer



Thermo Scientific™ Orbitrap Exploris™ 480
Mass Spectrometer





Updated Software Manuals and updated online help are part of the delivered ISO-Image and installed upon the installation of OES 4.3 ICSW



Display Online Help – of Tune – and – of Method Editor - via Fct F1 Key

Display online help – of Tune - via Fct F1 key

or via the gear wheel symbol in Tune – located to the right (top) corner in Tune



Display online help – of Method Editor - via Fct F1 key






- Working under Chromeleon CDS * software

* CDS = Chromatography Data System

 The world leader in serving science



OES 4.3 ICSW and Chromeleon CDS Software

LC-MS data acquisition under Chromeleon

- **OES 4.3 driver** validated for use with Chromeleon CDS 7.2.10 MUh software and Chromeleon CDS 7.3.2 MUb
- Defect fixes:
The previous issue that the OT Exploris model difference in Chromeleon is not recognized is fixed
- The Chromeleon Driver Compatibility matrix is updated when new combinations of software versions are tested. For more information and to view the compatibility matrix, sign on to <https://support.thermoinformatics.com/downloads/default.aspx>, and then select **Chromeleon > Chromeleon > Related Drivers > Driver Compatibility Matrix**.