

**Trace analysis of  
polychlorinated dibenzo-p-  
dioxins/furans (PCDD/F) in  
food and feed using GC-MS/MS**

**Dr. Nicholas A. Warner**


Product Specialist for GC and GC/MS

# Dioxins

## Organochlorines


Synthesis byproduct

Incineration byproduct




2,3,7,8-TCDD

(most toxic known chemical)



2,3,7,8-TCDF

Bio-accumulative



Toxic



## POPs regulated by Stockholm Convention

- |   |
|---|
| Aldrin  |
| Chlordane                                       |
| DDT   |
| Dieldrin  |
| Endrin  |
| Heptachlor                                      |
| Hexachlorobenzene                               |
| Mirex   |
| Toxaphane                                       |
| Polychlorinated biphenyls (PCB)                 |
| <b>Polychlorinated dibenzo-p-dioxins (PCDD)</b> |
| <b>Polychlorinated dibenzofurans (PCDF)</b>     |

The Dirty Dozen



# Adoption of European safety guidelines


**2001** - EU Scientific Committee on Food - risk assessment of PCDD/Fs and dioxin-like PCBs in food

- Weekly intake of 14 pg WHO TEQ/kg body weight limit established
- Maximum and action limits for PCDD/F within food and feed set highlighting the need for constant monitoring



## Globalization of food market

- No one jurisdiction is risk free from exposure
- Monitoring of food and feed exports and imports need to avoid future accidental exposure



# Instrumentation requirements

Meeting the needs of compliance

- High resolution magnetic sector mass spectrometry (HRMS)
  - Gold standard in Dioxin analysis (Worldwide compliance)
- **2014** – EU allows GC-MS/MS for dioxin analysis in food and feed



Thermo Scientific™  
TSQ™ 9610 GC-MS/MS



Thermo Scientific™  
DFS™ GC-HRMS



**Project Goal**  
Compare performance between platforms in meeting Dioxin analysis criteria in food and feed

**Internal evaluation**  
➤ Thermo Fisher Scientific

**External evaluation**  
➤ EURL (Freiburg, Germany)  
➤ Wageningen Food Safety Research Institute (Netherlands)

# Injection and separation methodology

Thermo Scientific iConnect PTV injector



Thermo Scientific TG-Dioxin column (60 m)



Injection		Oven temperature program	
Injector	iConnect PTV	Initial temperature (°C)	140
Liner	PTV liner with concentric baffle (P/N 453T2845-UI)	Hold time (min)	2
Injection mode	Large volume	Rate 1 (°C·min <sup>-1</sup> )	25
Injection volume (µL)	5	Temperature 1 (°C)	250
Split flow (mL·min <sup>-1</sup> )	100	Rate 2 (°C·min <sup>-1</sup> )	2.5
Initial injector temperature (°C)	75	Temperature 2 (°C)	260
Injection time (min)	0.4	Hold time (min)	5
Transfer rate (°C·s <sup>-1</sup> )	2.5	Rate 3 (°C·min <sup>-1</sup> )	2.5
Final temperature (°C)	300	Temperature 3 (°C)	285
Cleaning phase rate (°C·s <sup>-1</sup> )	14.5	Rate 4 (°C·min <sup>-1</sup> )	10
Cleaning phase temperature (°C)	330	Temperature 4 (°C)	320
Cleaning phase hold time (min) / flow rate (mL·min <sup>-1</sup> )	5 / 200	Hold time (min)	15
Carrier gas flow rate (mL·min <sup>-1</sup> )	1.2		



## Trace analysis of polychlorinated dibenzo-p-dioxins/dibenzofurans using GC-MS/MS in accordance with EU Regulations 2017/644 and 2017/771 for food and feed

### Authors

Nicholas Warner<sup>1</sup>, Adam Ladak<sup>2</sup>, Daniel Kutscher<sup>1</sup>, Alexander Schächtele<sup>3</sup>, and Kerstin Krätschmer<sup>4</sup>

<sup>1</sup>Thermo Fisher Scientific, Germany

<sup>2</sup>Thermo Fisher Scientific, UK

<sup>3</sup>European Union Reference Laboratory (EURL) for Halogenated POPs in Feed and Food, Freiburg, Germany

<sup>4</sup>Wageningen Food Safety Research, Part of Wageningen University & Research, Wageningen, The Netherlands

### Goal

To demonstrate the suitability of the Thermo Scientific™ TSQ™ 9610 triple quadrupole GC-MS/MS system with Thermo Scientific™ Chromelion™ Chromatography Data System (CDS) software for the routine and regulatory compliance testing of polychlorinated dibenzo-p-dioxins and polychlorinated dibenzofurans (PCDD/Fs) in food and feed samples in accordance with Commission Regulations (EU) 2017/644 and 2017/771.

### Introduction


Polychlorinated dibenzo-p-dioxins/dibenzofurans (PCDD/F), or dioxins as they are commonly referred to, are classified as persistent organic pollutants (POPs) under the Stockholm Convention. Due to their chemical stability and high fat solubility, these chemicals can accumulate within food chains and pose exposure/health risks to humans through consumption of food items (i.e., dairy, meat, and fish). Regulations are in place to monitor food and feed for the presence of dioxins to protect the population.<sup>1</sup> Current maximum levels allowable for PCDD/Fs in food and feedstuffs are at the pg·g<sup>-1</sup> concentration range due to the toxicity risk they pose.<sup>2,3</sup> However, the European Commission has announced plans to lower the current maximum levels in certain feedstuffs in 2024 based on updated risk assessment.<sup>4</sup> As of 2014, a change of EU regulations permitted the use of gas chromatography-triple quadrupole mass

### Keywords

GC-MS/MS, time-resolved SRM, polychlorinated dibenzo-p-dioxins, polychlorinated dibenzofurans, sensitivity, Regulation (EU) 2017/644, Regulation (EU) 2017/771, food, feed

# Differences in mass selectivity


## DFS



- Double focusing mass analyzer
  - Magnetic and Electrostatic focusing
  - High mass resolution data obtained
    - $R = 10,000$  (10% valley definition)
    - $R = 20,000$  FWHM
  - Robust performance from source design

vs.

## GC-MS/MS



- Mass separation in alternating electric field
  - AC/DC applied across quadrupoles
  - Unit (low) mass resolution technique
- MS/MS selectivity
  - Monitor compound selective fragmentation
  - Fast scanning capabilities
    - -> 800 simultaneous transitions possible



# LOQ determination

## Calibration based approach for GC-MS/MS analysis<sup>1</sup>


- Defined as the lowest level of analyte detected and meets performance criteria:

- Within the specified retention time window (for all monitored ions).
- ***Ion ratio intensities***  $\leq 15\%$  of theoretical (or calculated from standards) values
- Deviation from the average ***relative response factor***  $\leq 30\%$  for all calibration points
- ***The LOQ calculated by taking the lowest point of the calibration curve and correcting for the final sample volume, sample intake weight, and associated internal standard recovery.***


LOQ determined experimentally

Standard at LOQ

0.020pg



16pg



$$\text{Sample LOQ (pg/g)} = \sum_{n=PCDD/F}^{17} \text{Min Conc}_n(\text{pg}/\mu\text{L}) * \left( \frac{\text{Sample volume}(\mu\text{L})}{\text{Sample weight(g)} * \text{Recovery } l(\%)} \right)$$

1. Wenzl, T., Haedrich, J., Schaechtele, A., Robouch, P., Stroka, J., Guidance Document on the Estimation of LOD and LOQ for Measurements in the Field of Contaminants in Feed and Food; EUR 28099, Publications Office of the European Union, Luxembourg, 2016, ISBN 978-92-79-61768-3; doi:10.2787/8931u

# Thermo Scientific™ Chromeleon Chromatography Data System (CDS) 7.3.2

## Dioxin eWorkflow

The screenshot displays the Chromeleon Console #2 interface. On the left, a file browser shows the 'Dioxin analyzer initial setup' folder. The 'eWorkflows' panel is open, showing a list of workflows under 'ChromeleonLocal', with 'Dioxin Analyzer - v1.3' selected. A red arrow points from the file browser to the eWorkflows panel. Below the eWorkflows panel, a table lists various workflow items with columns for Name, Date Modified, and Type.

Name	Date Modified	Type
eur 28099 en_lod loq guidance document.pdf	3/28/2019 3:50:19 PM +01:00	Attachment
PCB review	5/29/2019 12:54:44 PM -05:00	View Settings
PCBs	3/29/2019 2:32:45 PM +00:00	Report Template
PCBs	3/29/2019 2:50:50 PM +00:00	Processing Method
PCBs - PTV(bk) - AI-AS1610 - 1µL	5/31/2024 4:17:24 PM +02:00	Instrument Method
PCBs - PTV(bk) - TriplusRSHSmart - 1µL	5/31/2024 4:20:13 PM +02:00	Instrument Method
PCBs - PTV(fr) - AI-AS1610 - 1µL	5/31/2024 4:17:40 PM +02:00	Instrument Method
PCBs - PTV(fr) - TriplusRSHSmart - 1µL	5/31/2024 4:20:01 PM +02:00	Instrument Method
PCBs - SSL(bk) - AI-AS1610 - 1µL	5/31/2024 4:17:47 PM +02:00	Instrument Method
PCBs - SSL(bk) - TriplusRSHSmart - 1µL	6/3/2024 9:54:57 AM +02:00	Instrument Method
PCBs - SSL(fr) - AI-AS1610 - 1µL	5/31/2024 4:19:53 PM +02:00	Instrument Method
PCBs - SSL(fr) - TriplusRSHSmart - 1µL	6/3/2024 9:55:06 AM +02:00	Instrument Method
PCDDF	4/25/2019 11:45:45 PM -05:00	Report Template
PCDDF review	5/29/2019 12:51:55 PM -05:00	View Settings
PCDDFs	3/29/2019 3:07:36 PM +00:00	Processing Method
PCDDFs - PTV(bk) - AI-AS1610 - 1µL	5/31/2024 4:18:10 PM +02:00	Instrument Method
PCDDFs - PTV(bk) - TriplusRSHSmart - 1µL	5/31/2024 4:19:21 PM +02:00	Instrument Method
PCDDFs - PTV(fr) - AI-AS1610 - 1µL	5/31/2024 4:18:30 PM +02:00	Instrument Method
PCDDFs - PTV(fr) - TriplusRSHSmart - 1µL	5/31/2024 4:19:15 PM +02:00	Instrument Method
PCDDFs - SSL(bk) - AI-AS1610 - 1µL	5/31/2024 4:18:37 PM +02:00	Instrument Method
PCDDFs - SSL(bk) - TriplusRSHSmart - 1µL	5/31/2024 4:19:08 PM +02:00	Instrument Method
PCDDFs - SSL(fr) - AI-AS1610 - 1µL	5/31/2024 4:18:47 PM +02:00	Instrument Method
PCDDFs - SSL(fr) - TriplusRSH - 1µL	4/25/2019 7:22:55 PM -05:00	Instrument Method
TSQ ion ratio intensity calculator PCB.xlsx	3/4/2019 10:16:46 AM +01:00	Attachment
TSQ ion ratio intensity calculator PCDD.xlsx	3/5/2019 11:44:09 AM +01:00	Attachment
TSQ ion ratio intensity calculator PCDF.xlsx	3/1/2019 12:17:38 PM +01:00	Attachment
WHO-TEQ_Upperbound calculator - standards.xlsx	3/28/2019 3:42:37 PM +01:00	Attachment


Automated full method setup for PCBs and PCDD/Fs including:

- Instrument methods for TriPlus RSH Smart and AI-AS 1610 autosamplers
- Isotopic dilution quantification methods
- Reporting layouts/templates
- Supporting compliance documentation




# Meeting performance criteria

Analysis repeatability at trace levels






**Assessment of the dioxin sensitivity using TCDD step standard ranging from 2 to 100 fg·µL<sup>-1</sup> with a 5 µL injection and reproducibility over 10 injections.**

Peak number	Compound	Concentration (fg On-Column)	Repeatability (%RSD) n = 10
1	1368-TCDD	10	10.3
2	1379-TCDD	20	7.3
3	1378-TCDD	50	8.7
4	1478-TCDD	125	6.6
5	1234-TCDD	250	6.0
6	2378-TCDD	500	6.6





# Meeting performance criteria

## Ion ratio performance at trace levels



# Meeting performance criteria

## Response factor variance and accuracy




	A	B	H	J
1	Peak Name	Ret. Time		
2		min		
3	First Injection	First Injection	EURL external check standard	EURL analysis LOQ/4 -4
5	2378-TCDF	19.292	Pass	Pass
6	2378-TCDD	19.852	Pass	Pass
7	12378-PeCDF	23.319	Pass	Pass
8	23478-PeCDF	24.704	Pass	Pass




# MS/MS selectivity in complex matrices

Minimal matrix impact at sub –pg/g concentrations





TCDF in laying hen egg  
using GC – MS/MS




TCDF in broiler fat  
using GC – MS/MS




TCDF in mixed animal fat  
using GC – MS/MS




13C-2378-TCDF IS




13C-2378-TCDF IS




13C-2378-TCDF IS



# EURL sample analysis (Freiburg, Germany)



# Wageningen Food Safety Research (Netherlands) analysis





# Conclusions

- Femtogram level sensitivity achievable by TSQ 9610 GC-MS/MS in food and feed samples
- Check standard performance within ion ratio thresholds for EU regulatory compliance at femtogram levels
- Good data agreement between TSQ 9610 and DFS through internal and external evaluation highlights accuracy and precision performance of the TSQ 9610
- Easy implementation with the Chromeleon 7.3.2 CDS Dioxin analyzer eWorkflow



# Thank you

Any Questions

