

Hydraulic Fracturing Wastewaters

Monika Verma, Carl Fisher, Richard Jack
Thermo Fisher Scientific, Sunnyvale, CA, USA

Executive Summary

Wastewater produced from hydraulic fracturing (fracking) contains salts (anions and cations), metals, and radioisotopes that have been mobilized from the hydrocarbon-rich shale layer. When high-pressure is used to crack the shale, wastewater flows to the surface where it is subject to different fates. Fracking wastewater can be disposed of into injection wells if no further use is required. If the wastewater will be used in subsequent fracking events, the water is first desalinated prior to mixing with additives. Lastly, wastewater can be sent to publicly owned treatment works (POTWs) to be treated prior to surface water discharge. Due to health and environmental safety concerns, anions, cations and organic acids found in fracking wastewaters that will not be directly disposed of should be monitored. Reagent-Free™ Ion Chromatography (RFIC™) systems can be used to determine the concentrations of anions, cations, and organic acids that are commonly found in fracking waters.

Keywords

Fracking, flowback water, produced water, Reagent-Free Ion Chromatography, publicly owned treatment works (POTW)




Figure 1. Fracking overview.

Introduction

Hydraulic fracturing (or fracking) extracts hard-to-reach natural gas and oil from deep underground deposits by drilling a well into bedrock (primarily shale) and then injecting fluid under high-pressure causing cracks to form. These cracks release trapped hydrocarbons that can be captured. Fracking fluid contains approximately 90% water which can be sourced from ground water, surface water, or treated wastewater.¹ The remainder of the fluid includes sand and chemical additives such as friction reducers, anti-bacterial agents, and corrosion inhibitors.² Sand is used to prop open cracks, facilitating the flow of hydrocarbons.

Upon relieving the injection pressure, the fracking fluid's direction of flow reverses which returns fracking wastewater to the surface. The percent of fracking wastewater that returns can range from 10–70%, depending on the nature of rock layers present.

Fracking Wastewater

There are two types of fracking wastewater: flowback and produced. Flowback water refers to the fluid that returns to the surface upon completion of a hydraulic fracturing event prior to production. Produced water refers to the fluid that returns to the surface once the well starts producing natural gas and oil. Both types of water can contain liquid that was resident in the shale layer (formation brine); however flowback water is primarily composed of fracking fluid and its associated chemicals.

Figure 2 depicts a typical hydraulic fracturing workflow. Fracking wastewater is either disposed of by pumping into injection wells, partially treated by desalination, or discharged to surface water post-treatment. The above disposal and transfer processes do run the risk of unintentional spills during transport which can impact water sources in the immediate vicinity.




Figure 2. Hydraulic fracturing workflow.

Reuse of Fracking Wastewater

Due to the scarcity and cost of water, fracking wastewaters are reused for additional fracking events in the same proximity. Prior to reuse, these waters are desalinated to remove salts, which can inhibit the effectiveness of additives and are known to cause scaling. For example, cations such as calcium, barium, and strontium may cause scale formation in water pumps, pipes, etc., reducing extraction efficiency. To counter this, additional anti-scaling agents would need to be added to the fracking fluid prior to reuse.³ Occasionally organic acids such as formate and acetate are added to control pH. These organic acids are sources of carbon for bacterial growth. Bacterial growth in hydraulic fracturing wastewaters can result in the production of hydrogen sulfide which is very toxic and may cause corrosion. Biocides are added to prevent bacterial growth.⁴

Discharge of Fracking Wastewater into Surface Water

Another option is the treatment of fracking wastewaters onsite or at POTWs prior to discharge. If not properly treated, these wastewaters will have a negative impact on downstream drinking water sources. For example, significant concentrations of chloride and bromide can be found in hydraulic fracturing wastewaters. If ozonation is used during the drinking water purification process, chloride and bromide react with ozone to produce chlorite and bromate. The National Primary Drinking Water Regulations regulate both bromate and chlorite, as the former is carcinogenic and the latter can lead to nervous system damage.⁵ Sulfate, which is regulated under the National Secondary Drinking Water Regulations, can make drinking water unpalatable and have a laxative effect.⁶

Ion Analysis of Fracking Wastewater

The anions, cations, and organic acids mentioned above can be easily determined in hydraulic fracturing wastewater using Thermo Scientific™ Dionex™ RFIC systems which electrolytically create the precise amount of eluent required for ion chromatography (IC) applications. This capability eliminates the need to handle the acids and bases required for the preparation of IC eluents, allowing chromatographers to run a full range of gradient and isocratic separations more effectively than with hand-made eluents. These systems allow for a simpler and more reliable way to deliver superior results while simultaneously saving time and labor. Furthermore, because the instrument pump seals and pistons only come into contact with deionized water, overall pump maintenance is significantly reduced.⁷ Figure 3 a, b shows examples of cation, anion, and organic acid determinations in flowback water using RFIC systems.




Figure 3a. Determination of cations in fracking flowback water using capillary IC.⁸




Figure 3b. Determination of inorganic anions and organic acids in fracking flowback.

References

1. Study of the Potential Impacts of Hydraulic Fracturing on Drinking Water Resources PROGRESS REPORT; EPA/601/R-12/011; US Environmental Protection Agency Office of Research and Development: Washington, DC, December 2012;
2. FracFocus Chemical Disclosure Registry, What Chemicals are Used Page. <http://fracfocus.org/chemical-use/what-chemicals-are-used> (accessed May 6, 2014).
3. Vidic, R.D.; Brantley, S.L.; Vandenbossche, J.M.; Yoxtheimer, D.; Abad, J.D. Impact of Shale Gas Development on Regional Water Quality. *Science*. 2013, 340, 1235009.
4. Thermo Fisher Scientific. Whitepaper 70693: The Importance of Anion and Organic Acid Determinations in Fracking Wastewater by Ion Chromatography. Sunnyvale, CA 2013. [Online] <http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/CMD%20Documents/WP-70693-Anion-Organic-Acid-Determinations-Fracking-Wastewater-WP70693-E.pdf> (accessed May 6, 2014)
5. U.S. Environmental Protection Agency. Basic Information about Disinfection Byproducts in Drinking Water: Total Trihalomethanes, Haloacetic Acids, Bromate, and Chlorite. <http://water.epa.gov/drink/contaminants/basicinformation/disinfectionbyproducts.cfm#> (accessed May 6, 2014).
6. U.S. Environmental Protection Agency. Secondary Drinking Water Regulations: Guidance for Nuisance Chemicals. <http://water.epa.gov/drink/contaminants/secondarystandards.cfm> (accessed May 6, 2014)
7. Thermo Fisher Scientific. Thermo Scientific Dionex Reagent-Free Ion Chromatography (RFIC) System Capabilities Brochure. Sunnyvale, CA 2013. [Online] <http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/CMD%20Documents/29368-BR-70667-Reagent-Free-Ion-Chromatography-RFIC-System-Capabilities.pdf> (accessed May 6, 2014).
8. Thermo Fisher Scientific. Thermo Scientific Application Note 1094, Determination of Cations in Hydraulic Fracturing Flowback Water from the Marcellus Shale, AN71085_E, Sunnyvale, CA, 2014.

www.thermoscientific.com/chromatography

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific, Sunnyvale, CA
USA is ISO 9001:2008 Certified.

Australia +61 3 9757 4300	Denmark +45 70 23 62 60	Japan +81 6 6885 1213	Singapore +65 6289 1190
Austria +43 810 282 206	Finland +358 9 3291 0200	Korea +82 2 3420 8600	Sweden +46 8 556 468 00
Belgium +32 53 73 42 41	France +33 1 60 92 48 00	Latin America +1 561 688 8700	Switzerland +41 61 716 77 00
Brazil +55 11 3731 5140	Germany +49 6103 408 1014	Netherlands +31 76 579 55 55	Taiwan +886 2 8751 6655
Canada +1 800 530 8447	India +91 22 6742 9494	New Zealand +64 9 980 6700	UK/Ireland +44 1442 233555
China 800 810 5118 (free call domestic) 400 650 5118	Italy +39 02 950 591	Norway +46 8 556 468 00	USA +1 800 532 4752

WP71116-EN 0714S

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand