

New Mass Spectrometry Tools to Transform Metabolomics and Lipidomics

July.3.2013

Ken Miller Vice President of Marketing, Life Sciences Mass Spectrometry

The world leader in serving science

Omics & the Central Dogma of Biology



Addressing the Metabolomics Challenges



Developing SOPs for major sample types



Innovating LC/MS to improve separation & detection





Developing streamlined, workflow-oriented software

Building and connecting to new and existing database



Orbitrap Delivers the Highest Quality Data

Better data, Better results!

- Highest Resolution
- Best Mass Accuracy
- Widest Dynamic Range
- Excellent Sensitivity





LCMS Metabolomics Offerings

Quan / Qual analysis of endogenous molecules of low molecular mass (<1200-1500Da)



Thermo Fisher

Orbitrap Fusion Tribrid Mass Spectrometer



Unmatched Analytical Performance

Revolutionary performance Unprecedented usability Exceptional versatility



Orbitrap Fusion Tribrid Mass Spectrometer





Flexible, Intuitive, Easy to Learn Software

thod Timeline					
tethod Duration (mins) 140 # 233	46.7 70 93.3	116.7 140	- Q + Nex · · · Delete		
Experiment 1	Fime Range 0 - 140 mins				
Save as Template			ddMSn Properties		
ystem Templates	MS		MSn Level	2	
Metabolomics •	Monoisotopic Precursor		Isolation Mode	Quadrupole	
	Charge State		lsol. Window (m/z)	2	
Froteomics	Charge State		Use Isolation m/z Offset		
nall Molecules	Dynamic Exclusion		Activation Type	ETD	
	-0		Use calibrated charge dependent ETD parameters		
istom Templates			Detector Type	lon Trap	
ty Experiments	Decisions		Ion Trap Scan Rate	Rapid	
Others 3			Scan Range Mode	Auto: m/z Normal	
			First Mass (m/z)	120	
	ddMS2 OT HCD		AGC Target	1e4	
	Product Ion Trigger		Minimum AGC Target	[1e4	
			Maximum Injection Time (ms)	35	
			Microscans	1	
	ddMS2 IT ETD		Data Type	Profile	
			Polarity	Positive	
				- parti	

The Thermo Scientific[™] Orbitrap Fusion[™] Method Editor allows for drag and drop method building.

Dynamic Scan Management for maximum data quality



Orbitrap Fusion Tribrid MS



Scan rate FTMS ²	15 Hz
Scan rate ITMS ²	20 Hz
Max resolution	450K at m/z 200
Quad isolation	down to 0.4 amu
lon trap isolation	down to 0.2 amu
Mass Accuracy	1 ppm
Dissociation	CID, HCD, ETD
MSn	Up to MS^{10} in ion trap or Orbitrap
Analyzers	Q, FTMS, ITMS
Detectors	Ion Trap, Orbitrap
Compact	1186 x 674 x 650 mm (w, d, h)

Unmatched Analytical Performance

450,000 Resolution to remove spectral interferences

CID/HCD/ETD detected by the **Ion Trap or Orbitrap** at any level of MSⁿ for maximum experimental flexibility

Powered by the new **Dynamic Scan Management** architecture that lets every user easily get the maximum information from every sample

Mass Accuracy Matters!

Accurate mass improves the ability to confidently identify unknown compounds ppm error = (mass error / MW) x 10⁶

EXAMPLE:

• Reserpine $(C_{33}H_{40}N_2O_9)$ has a protonated ion at 609.28066

Mass Error (ppm)	# of Possible Chemical Formulae (using only C, H, N & O)
165 (0.1 daltons)	209
10	13
5	7
3	4
2	2
1	2
0.5 (0.0003 daltons)	1



Orbitrap Fusion Tribrid Mass Spectrometer



Challenge

Generate absolutely most accurate data to facilitate compound identification.



Orbitrap Fusion MS

EASY IC (internal calibration): new ETD source allows use of internal calibration for sub-ppm mass accuracies. On/off toggle in SW. No ion suppression by IC.



Reagent/Internal Calibrant Source





Data Type	Profile	•
Polarity	Positive	•
Use internal calibration		
Source Fragmentation		















Easy-IC Providing Sub-ppm Mass Accuracy in MS and MS/MS

Orbitrap Fusion MS operating @120K for FS, 15k for HCD MS2, with UHPLC. Sample, Urine. Analyte: L-Tryptophan, C₁₁H₁₂N₂O₂ [M+H]⁺=205.09715



MS

HCD MS/MS

ASMS 2013 poster **MP045** by Junhua Wang *et al.*



Orbitrap Fusion Tribrid Mass Spectrometer



Challenge

Separation of distinct isobars in complex mixtures



Orbitrap Fusion MS

Orbitrap Fusion can achieve up to 450,000 resolution.



See the Fine Isotope Pattern with Ultra-high Resolution

L-Methionine $C_5H_{11}NO_2S$



ASMS 2013 poster MP045 by Junhua Wang et al.



Resolving the Most Complex Sample

Nano-infusion of NIST SRM1950 plasma extract, Negative mode



Thermo Fishel S C I E N T I F I C

Orbitrap Fusion Tribrid Mass Spectrometer



Challenge

Achieve maximum coverage of the metabolome/lipidome.

Many molecules require special fragmentation conditions to reveal their structural information.



Orbitrap Fusion MS

Instrument architecture allows for extremely fast MS/MS and full flexibility of MSⁿ experiments for structural analysis including novel workflows such as HCD³ and MS² HCD followed by MS³ CID.



Method 1: The "Super" Q Exactive MS Mode



MS OT, R=60K

Full MS MS² Top1 MS² Top2 MS² Top3 MS² Top4 MS² Top5





MS² Top1 MS² Top2 MS² Top3 MS² Top4 MS² Top5



Method 3: MSⁿ Enables *de novo* Compound ID



Sheldon, Mistrik, Croley, JAm Soc Mass Spectrom, 2009, 20: 370-376



Method 3: HR/AM Ion Trees for Unknown Identification



De novo structural elucidation: capture substructure with MSⁿ ion tree



HCD MSⁿ is Unique to Orbitrap Fusion MS





Structural Homology Identified by MSⁿ Ion Tree Search





The Industry's Leading Portfolio of MS Solutions





SCIENTIFIC

Transform Your Science