



**ThermoFisher**  
S C I E N T I F I C

# Creation, Implementation and Transferability of a Glucose Unit (GU) Library for HILIC-FLD Analysis of APTS-Labelled Released N-Glycans


Phil Widdowson, Ph.D., European Pharma & Biopharma Applications Development Scientist  
Thermo Fisher Scientific, Runcorn, UK

The world leader in serving science

- *N*-linked glycosylation
  - Importance of glycans in biotherapeutics
  - Regulations
- Analysis of released *N*-glycans
  - Applied Biosystems™ GlycanAssure™ HyPerformance APTS kit
- Glucose Unit (GU) libraries
- Construction of APTS GU library
- Transferability of APTS GU library
  - eWorkflow

# N-linked Glycans

- ~50% of biopharmaceuticals are glycosylated
  - N-glycan attached to Asn residues
    - Consensus sequence N-X-S/T
- Glycosylation can have an impact on
  - Immunogenicity
  - Biological activity, drug clearance and half-life
  - Stability
- Therefore is a Critical Quality Attribute (CQA)
  - ICH Q6B
- Microheterogeneity affected by manufacturing process
  - Both a manufacturing and analytical challenge




Reviewed in:

Liu, *J Pharm Sci.* 2015; 104: 1866-1884

Abes & Teillaud, *Pharmaceuticals* 2010; 3:146-57


ADCC: Antibody-Dependent Cell-mediated Cytotoxicity

CDC: Complement-Dependent Cytotoxicity




**ICH (Q6B)** recommends 6 test approaches for characterization and confirmation of biological products:

- Amino acid sequence
- Amino acid composition
- Terminal amino acid sequence
- Peptide map
- Sulfhydryl group(s) and disulfide bridges
- **Carbohydrate structure**
  - *“For glycoproteins, the carbohydrate content and structure (neutral sugars, amino sugars, and sialic acids) is determined.”*




**Early development (High throughput clone screening) → QC (CQA monitoring)**

# Analysis of N-linked Glycans




# GlycanAssure HyPerformance Kit – Denaturation and Deglycosylation




- Denatures and solubilizes protein
- Enables rapid and complete deglycosylation

# GlycanAssure HyPerformance Kit – APTS Labelling




# GlycanAssure HyPerformance Kit – Sample Clean-up







# Overview of Glycan Labelling Using GlycanAssure



- External retention calibration approach
- Uses labelled dextran (poly-glucose) ladder standard
  - Generates retention time calibration curve
  - Converts chromatographic time to GU value
- Critical for glycan library search robustness
  - Minimizes chromatographic variation associated with HILIC
    - System-to-system
    - Day-to-day
    - Site-to-site
- Requires method standardization



**APTS-labelled dextran ladder (GU 3-12)**



**Retention time/GU calibration curve**


# Importance of GU Calibration

- Human IgG – different
  - Analysts
  - Mobile phase
  - Sample preparation
  - Instruments
  - Columns...
- Retention times differ by ~1 min. across the run
- Dextran ladder calibration applied to correct this difference

For G0F peak (n=12)

RSD of RT = 3.5%


RSD of GU = **0.9%**



Human IgG analysed across different sites

# Method Parameters


Mobile phase A	100 mM ammonium formate, pH 4.4
Mobile phase B	Acetonitrile
Column temperature	50°C
Injection volume	2 µL
Excitation wavelength	455 nm
Emission wavelength	500 nm
Sensitivity	7
Lamp mode	High power
Data collection rate	10 Hz



# Overview of Library Construction


- Biopharmaceutical relevant samples
  - NIST mAb
  - RNase B
  - Human IgG
- Individually labelled, *N*-glycan standards
  - Supplement the sample glycan pools
  - Aid and confirm glycan identifications
- Two analysts; Two different sites; Two column lots
  - Different UHPLC instruments, mobile phase preparations, sample preparations

# Typical Sample Chromatograms Used in Library Construction



NISTmAb, RNase B and hulG representative samples

n=9 overlay of hulG performed in the same analysis




# Glycans Incorporated into APTS GU Library

- APTS GU library contains 29 glycans
  - Biopharmaceutical relevant samples
  - Individually purified standards
- Average standard deviation (StDEV) for all glycans = 0.06
  - Robust and reproducible
  - Important for instrument-to-instrument and site-to-site transferability
  - Confident glycan assignments

Glycan	Average GU	StDEV
G0N	3.768	0.037
G0	3.897	0.057
G0N'	3.930	0.041
G0FN	4.030	0.031
G0FN'	4.171	0.042
G0F	4.193	0.060
G0FB	4.655	0.071
G1	4.777	0.065
G1'	4.891	0.050
G1Fa	4.945	0.061
G1Fb	5.173	0.059
G1FB	5.339	0.074
M5	5.510	0.040
G1FB'	5.589	0.051
G2	5.641	0.041
G2F	5.902	0.060
G2FB	6.157	0.074
M6	6.298	0.042
G1S1	6.470	0.066
G1FS1	6.729	0.081
M7	7.134	0.023
A1	7.181	0.078
A1F	7.447	0.082
A1FB	7.671	0.093
M8	8.039	0.046
M9	8.723	0.051
A2	8.755	0.098
A2F	8.975	0.101
A2FB	9.104	0.110

# GU Annotated Chromatograms

Biotherapeutic samples annotated with library GU values





# Preparation of eWorkflow

The screenshot displays the eWorkflow Editor software interface, which is used for configuring and managing laboratory workflows. The interface is divided into several panels:

- Chromleon Console:** Located at the top left, it shows a navigation menu with options like 'Folder', 'Sequence', 'Instrument', 'Processing', 'Report Template', 'Injection', 'Spectral Library', 'View Settings', 'eWorkflow', and 'Electronic Signature'.
- eWorkflow Editor - new\_eWorkflow\*:** This is the main workspace, which is split into multiple tabs: 'eWorkflow General', 'Sequence General', and 'Sequence Layout'.
  - General Tab:** Contains settings for the workflow type (currently 'None'), tags, and associated items (including 'Vanquash\_DAD').
  - Sequence General Tab:** Shows default settings for report templates, view settings, and channels. It also includes options for auto-reporting, notification settings, and electronic signatures.
  - Sequence Layout Tab:** Displays a table of sequence items, including headers, sample blocks, brackets, and footers.
- Sequence Preview:** A detailed view of the sequence items, showing columns for Chromatogram, Name, Type, Level, Spike Group, Position, Volume, Instrument Method, Processing Method, Status, and Inject Time. The preview shows 12 items, including blanks, samples, and calibration standards.

- APTS GU library has been created
  - GlycanAssure HyPerformance Kit
- 29 biopharmaceutically relevant glycans incorporated
- Use of GU library increases robustness and confidence in glycan assignments
  - Without the requirement for mass spectrometry
- Library in the process of internal/external validation
- Creation of an eWorkflow
  - Seamless transfer of methods between laboratories/sites