


IonPac® Cryptand C1 Concentrator Column



The IonPac Cryptand C1 Concentrator Column is designed primarily for the analysis of trace perchlorate in drinking water. A method using the Cryptand C1 Concentrator Column can quantify 140 ng/L (ppt) of perchlorate in a background of total dissolved solids (chloride, carbonate, and sulfate) at combined concentrations as high as 3000 mg/L. The Cryptand C1 Concentrator Column is the specified concentrator column for sample preconcentration in U.S. EPA Method 314.1.

The Cryptand C1 is an adjustable-capacity concentrator column containing a macroporous, 17.5-µm resin that has been grafted with the macrocyclic 2,2,2 cryptand compound. This column has approximately 30 µeq/col of cryptand capacity available for use. The functional capacity depends on the eluent concentration and type of cation bound within the

cage of the cryptand molecule. The unique ability to adjust the capacity of the column from a fixed amount to nearly zero capacity can be achieved by changing either the eluent concentration or type of cation. This unique feature makes the Cryptand C1 Concentrator Column a powerful tool for the determination of trace perchlorate in drinking water and high-ionic-strength water.

Now sold under the
Thermo Scientific brand

Thermo
SCIENTIFIC



Figure 4 illustrates the system flow path for the determination of trace perchlorate according to U.S. EPA Method 314.1. The sample is first loaded onto the Cryptand C1 Concentrator Column manually or with an autosampler. Most environmental water samples contain dissolved solids in the form of chloride, carbonate, and sulfate at concentrations from 10 mg/L (ppm) up to as high as 1000 mg/L (ppm) each. The majority of the dissolved solids are rinsed off the Cryptand C1 Concentrator Column with a rinse step of 10 mM sodium hydroxide. Because perchlorate is highly retained relative to the dissolved solids, most if not all the perchlorate remains on the Cryptand C1 Column. The Cryptand C1 Concentrator Column is then switched in-line with the guard/analytical column set and the eluent concentration is stepped down to 0.5 mM sodium hydroxide. At 0.5 mM sodium hydroxide, the Cryptand C1 Concentrator Column has nearly zero capacity, and the perchlorate will elute off and onto the head of the guard/analytical column set, along with any residual dissolved solids.

The analytical column set (AS16 or AS20) has a fixed capacity, and the perchlorate concentrates at the head of the guard/analytical column set and is not eluted due to the low hydroxide eluent concentration (0.5 mM). A step change to 65 mM sodium hydroxide is used to elute the perchlorate in ~20 min. The perchlorate peak is well resolved from the residual dissolved solids. A final step change to 100 mM sodium hydroxide is used to elute off any residual ions that might be present in certain water samples. This step also ensures that the Cryptand C1 Concentrator Column is converted completely to the sodium form prior to the next analysis. The Cryptand C1




Figure 3. Determination of trace perchlorate in drinking water using the AS20 column following concentration on a Cryptand C1 Concentrator Column.




Figure 4. System flow path for trace perchlorate in drinking water.

Concentrator Column is then switched out of line from the guard/analytical column set and is ready for loading a new sample.

Low- to sub-µg/L (ppb) levels of perchlorate can easily be quantified using the AS16 and AS20 columns and a 2-mL sample injection, as shown in Figures 2 and 3.

SPECIFICATIONS

Dimensions: IonPac Cryptand C1

Concentrator Column:

4 × 35 mm

Maximum Operating Pressure:

4000 psi

Mobile Phase Compatibility:

pH 0–14; 0–100% HPLC solvents

Substrate Characteristics:

Bead Diameter (µm): 17.5

Pore Size: 100 Å

Cross-Linking (%DVB): 55%

Ion-Exchange Group:

Grafted Cryptand (with bound metal cation)

Functional Group Characteristics:

Variable hydrophobicity (eluent dependent)

Capacity:

4 × 35 mm concentrator column:
30 µeq (sodium form)

Column Construction:

Polyetheretherketone (PEEK) with 10-32 threaded ferrule-style end fittings. All components are nonmetallic.

ORDERING INFORMATION

To order in the U.S., call 1-800-346-6390, or contact the Dionex Regional Office nearest you. Outside the U.S., order through your local Dionex office or distributor. Refer to the following part number.

IonPac Cryptand C1 Concentrator Column

4 × 35 mm P/N 062893

System Requirements for EGC-NaOH Eluent Generation

For NaOH eluents generated using the EGC-NaOH cartridge, use with ICS-2000, ICS-2500, or ICS-3000 RFIC systems equipped with eluent generation. The eluent generator is used to automatically produce sodium hydroxide gradients from deionized water. The EGC-NaOH cartridge can only be controlled using Chromeleon 6.7 and subsequent Chromeleon releases.

ASRS, AutoSuppression, and IonPac are registered trademarks of Dionex Corporation.

Passion. Power. Productivity.



Dionex Corporation

1228 Titan Way
P.O. Box 3603
Sunnyvale, CA
94088-3603
(408) 737-0700

North America

U.S. (847) 295-7500
Canada (905) 844-9650

South America

Brazil (55) 11 3731 5140

Europe

Austria (43) 1 616 51 25 Belgium (32) 3 353 4294 Denmark (45) 36 36 90 90
France (33) 1 39 30 01 10 Germany (49) 6126 991 0 Italy (39) 02 51 62 1267
The Netherlands (31) 20 683 9768 Switzerland (41) 62 205 9966
U.K. and Ireland (44) 1276 691722

Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 28475235
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190

www.dionex.com



LPN 1678-01 PDF 10/06
©2006 Dionex Corporation