

# Thermo Scientific TriPlus RSH Autosampler

The Thermo Scientific™ TriPlus™ RSH Autosampler offers exceptional precision, flexibility, and productivity in robotic, sample-handling solutions. Compatible with Thermo Scientific GC and GC/MS systems, the autosampler sets new standards in automation and provides advanced liquid-handling cycles that enable automated functionality beyond traditional liquid, headspace, and solid-phase microextraction (SPME) injections.

## Integrated sampling system

The TriPlus RSH Autosampler expands automated capabilities beyond sample injections to advanced sample-handling cycles. Automating the whole workflow, from sample preparation to injection, increases precision and reproducibility and enables unprecedented flexibility and productivity in GC and GC-MS sampling.

## Ultimate productivity

Designed for expanded productivity, this integrated sampling system accommodates large sample capacities supporting full, unattended 24/7 operation. As an example, a maximum of 972, 2 mL vials combined with multiple 100 mL wash/waste bottles enable weekend long unattended operations – a goal not attainable with most other sampling systems currently on the market.



## Unmatched flexibility

The TriPlus RSH robotic sample-handling system offers a multitude of configurations including, liquid, headspace and solid-phase microextraction. In addition to other capabilities, this sampling system delivers the precision you demand for achieving high quality results.

Scalable capabilities enable expanded GC and GC-MS application ranges and the best matching of techniques to sample types.

## Seamless operation

Accurate automation of a multitude of traditionally manual tasks is enabled by the smart Automatic Tool Change capability (ATC). The ATC feature enables the user to set up a sequence using up to six different syringes, automatically loaded by the autosampler to accurately perform dilutions, standard additions, calibrations set-up, and sample injections. The ability to exchange syringes for different tasks enables accurate and highly precise sample-handling steps in a single, unattended sequence prior to automated sample injection.

## TriPlus RSH configurations

The TriPlus RSH Autosampler is available in several base configurations for:

- Automated liquid sampling and injection
- Static headspace sampling and injection
- SPME Arrow sampling and injection
- Sample preparation routine workflows

All the configurations exist with standard or extended X-arm length. Additional upgrade kits, tools and accessories are available to extend the sample vial capacity or to add other sampling capabilities and transform any configuration into multi-technique platforms.

## Features and technical specifications

<b>Description</b>	<ul style="list-style-type: none"><li>• XYZ robotic sample-handling apparatus with ATC capability, allowing different injection techniques executed within the same sample sequence in a fully unattended way</li><li>• Up to six different tools can be managed at the same time, using two ATC stations</li><li>• Using optional accessories and dedicated programmed workflows, the system is also capable of automating the most common sample preparation steps and deliver the highest level of sample-handling flexibility</li></ul>
<b>Communication</b>	<ul style="list-style-type: none"><li>• Two independent LAN ports</li></ul>
<b>Local user interface</b>	<ul style="list-style-type: none"><li>• LED status indicators</li><li>• Optional control panel with 4 keys, round knob and graphical LCD display</li></ul>
<b>Instrument control</b>	<ul style="list-style-type: none"><li>• Local controller for direct access to instrument configuration and movements (optional)</li><li>• Thermo Scientific chromatography data systems integrated with Virtual Terminal software to completely mimic the local controller</li></ul>
<b>Teaching functions</b>	<ul style="list-style-type: none"><li>• Manual without using tools or external devices</li></ul>
<b>Injector compatibility</b>	<ul style="list-style-type: none"><li>• Compatible with on-column (COC), programmable temperature vaporizing (PTV), packed (PKD), purged packed (PPKD), split-splitless (SSL) injectors</li></ul>
<b>High throughput configuration</b>	<ul style="list-style-type: none"><li>• Dual GC set-up with Double Pro and Confirmation modes: single TriPlus RSH serving two independent GC or GC/MS systems, for liquid, HS or SPME sample injection or for sample preparation, using the same or two different software systems. Two different sampling methods can be used (for liquid injection only)</li><li>• Rapid Mode starts the syringe washing cycle during the current GC cooling phase</li></ul>
<b>Barcode reader</b>	<ul style="list-style-type: none"><li>• Two active laser scanners for all standard vials using 1-dimension barcodes in a horizontal orientation</li></ul>
<b>Vortexer</b>	<ul style="list-style-type: none"><li>• Intensive sample mixing with an agitation speed up to 2000 rpm. Compatible with 0.5, 0.7, 2, 5, 10, or 20 mL vials</li></ul>
<b>Incubator/agitator</b>	<ul style="list-style-type: none"><li>• Capacity 6 × 20 mL vials (compatible with 2 and 10 mL vials with adapters)</li><li>• 30–200 °C temperature range</li><li>• 250–750 rpm agitation speed</li></ul>
<b>Temperature-controlled vial trays</b>	<ul style="list-style-type: none"><li>• Heated and cooled trays expand the range of applications from sample injection to sample/standard preparation</li></ul>
<b>Dilutor</b>	<ul style="list-style-type: none"><li>• Single- or multi-solvent dispensing tool for single- or multi-solvent dilution</li><li>• External solvent reservoirs for increased solvent capacity and longer unattended operations</li></ul>
<b>Mounting kits</b>	<ul style="list-style-type: none"><li>• Thermo Scientific™ TRACE™ 1300 Series GC, TRACE GC Ultra and FOCUS GC</li><li>• Mounting kits for major GCs on the market also available</li></ul>
<b>Sampling techniques</b>	<ul style="list-style-type: none"><li>• Liquid, Static Headspace, SPME, SPME Arrow, ITEX-DHS</li></ul>
<b>Advanced sample prep</b>	<ul style="list-style-type: none"><li>• Dedicated PrepCycles are available to perform routine sample handling workflows</li><li>• Tools such as Vortexer, Incubator/Agitator or Dilutor as well as multiple large volume syringes, large solvent station, <math>\mu</math>SPE option and the ATC can be used to automate routine sample preparation procedures, such as standard dilution, standard addition, sequential dilution, derivatization and sample clean-up</li></ul>
<b>Sampling Workflow Editor software</b>	<ul style="list-style-type: none"><li>• This software tool allows the user to easily program custom sample preparation workflows through an intuitive drag-and-drop visual programming interface</li></ul>

## Liquid sampling

<b>Vial volumes</b>	<ul style="list-style-type: none"><li>• 300 µL fixed insert vials, 0.5, 0.7, 2, 2.5, 10, and 20 mL vials. 96/384 Microtiter or Deep Well plates with Automatic Foil Cutter to pierce alumina or plastic foils prior the needle penetration</li></ul>
<b>Bottom sensing for vials</b>	<ul style="list-style-type: none"><li>• Capable of liquid injection starting from small-volume samples. Capability to inject from samples as low as 5 µL into a vial. Possibility of performing up to three 1 µL injections from a 5 µL sample, depending on vial type</li></ul>
<b>Height from vial bottom</b>	<ul style="list-style-type: none"><li>• User selectable between 0.1 and 32 mm in 0.1 mm steps</li></ul>
<b>Injection speed for liquid samples</b>	<ul style="list-style-type: none"><li>• Selectable from 0.1 µL/sec up to 2000 µL/sec and fully programmable</li></ul>
<b>Sample capacity</b>	<ul style="list-style-type: none"><li>• Depending on autosampler configuration<ul style="list-style-type: none"><li>– Up to 4608 well plates or 6912 well plates with extended X-arm</li><li>– Up to 840 × 0.5/0.7 mL vials or 1260 × 0.5/0.7 mL vials with the extended X-arm</li><li>– Up to 648 × 2 mL sample vials or 972 × 2 mL vials with the extended X-arm</li><li>– Up to 240 × 10 mL or 20 mL vials or 360 × 10 mL or 20 mL vials with the extended X-arm</li></ul></li></ul>
<b>Syringes</b>	<ul style="list-style-type: none"><li>• Capable of handling liquid volumes in the range 0.1 µL – 10 mL</li><li>• Capable of using 0.5 µL, 1.0 µL, 5 µL, 10 µL (standard), 25 µL, 50 µL, 100, µL, 250 µL, 500 µL, 1000 µL, 10000 µL syringes for sample injection and/or volume transfer</li><li>• Needle lengths: 57 mm or 85 mm</li></ul>
<b>Syringe cleaning</b>	<ul style="list-style-type: none"><li>• Wash stations for up to 4 different solvents for a total of 40 mL (standard)</li><li>• Optional large washing stations for up to 3 × 100 mL solvent bottles</li><li>• Waste: 1 × 10 mL or drain to external waste bottle</li><li>• Possibility to install multiple solvent stations to expand solvent and waste volumes</li></ul>
<b>Injection volume</b>	<ul style="list-style-type: none"><li>• Range from 0.1 to 10,000 µL in 0.1 µL steps up to 100 µL, and 1 µL steps between 100 µL and 10 mL</li></ul>
<b>Liquid injection modes</b>	<ul style="list-style-type: none"><li>• 8 fully customizable method-specific preset menus available:<ul style="list-style-type: none"><li>– Basic enrichment</li><li>– Enrichment needle solvent wash</li><li>– Internal standard double</li><li>– Internal standard post</li><li>– Needle solvent wash</li><li>– Solvent flush double</li><li>– Solvent flush post</li></ul></li></ul>
<b>Optional tools</b>	<ul style="list-style-type: none"><li>• Peltier-controlled drawer for well plates, 300 µL fixed insert vials, 2 and 10 mL vials. Temperature selectable between 4 and 40 °C</li><li>• Cooled tray holders for well plates, 300 µL fixed insert vials, 2, 10 and 20 mL vials. Temperature selectable between 4 and 70 °C. Requires external circulator bath</li><li>• Large Solvent Station – 3 × 100 mL</li><li>• Large Volume Wash Station – 2 × 100 mL and drain</li><li>• Fast Syringe Washing module with two solvents – 2 × 1000 mL</li></ul>
<b>Typical liquid injection repeatability</b>	<ul style="list-style-type: none"><li>• &lt;0.3 RSD % obtained under standard Thermo Scientific instrument conditions</li></ul>

## Headspace sampling

<b>Vial volumes</b>	<ul style="list-style-type: none"><li>• Compatible with 2, 10 and 20 mL vials</li></ul>
<b>Syringe sizes</b>	<ul style="list-style-type: none"><li>• Gastight 1, 2, 5, and 5 mL</li></ul>
<b>Needle lengths</b>	<ul style="list-style-type: none"><li>• 65 mm, compatible with every injector port</li></ul>
<b>Sample capacity</b>	<ul style="list-style-type: none"><li>• Depending on autosampler configuration<ul style="list-style-type: none"><li>– Up to 180 × 10 or 20 mL vials or 300 × 10 or 20 mL vials with the extended X-arm</li></ul></li></ul>
<b>Injection volume range</b>	<ul style="list-style-type: none"><li>• 0.1 to 5 mL in 0.1 mL steps, depending on syringe</li></ul>
<b>Injection speed</b>	<ul style="list-style-type: none"><li>• 1 to 100 mL/min, in 1 mL/min increments</li></ul>
<b>Syringe temperature</b>	<ul style="list-style-type: none"><li>• OFF or 40 °C to 150 °C in 1 °C steps</li></ul>
<b>Incubation oven</b>	<ul style="list-style-type: none"><li>• Capacity 6 × 20 mL vials (compatible with 2 and 10 mL vials with adapters)</li><li>• 30–200 °C temperature range, in 1 °C steps</li><li>• 250–750 rpm agitation speed</li></ul>
<b>Incubation time</b>	<ul style="list-style-type: none"><li>• 0.1 to 600.0 min in 0.1 min increments</li></ul>
<b>Syringe flush capability</b>	<ul style="list-style-type: none"><li>• With inert gas</li></ul>
<b>Solvent syringe cleaning</b>	<ul style="list-style-type: none"><li>• Optional washing stations (4 × 10 mL or 2 × 100 mL vials)</li></ul>
<b>Multiple Headspace Extraction (MHE)</b>	<ul style="list-style-type: none"><li>• Yes (optional accessory)</li></ul>
<b>Enrichment sampling</b>	<ul style="list-style-type: none"><li>• Yes with optional kit for cold trap</li></ul>
<b>Optional tools</b>	<ul style="list-style-type: none"><li>• Peltier-cooled tray holder for 300 µL fixed insert vials, 2 and 10 mL vials; temperature selectable between 4 °C and 40 °C</li><li>• Cooled tray holders for 300 µL fixed insert vials, 2, 10 and 20 mL vials; requires external circulator bath; temperature selectable between 4 °C and 70 °C</li></ul>
<b>Typical headspace injection repeatability</b>	<ul style="list-style-type: none"><li>• &lt;0.7 RSD % under Thermo Scientific standard conditions</li></ul>

## Solid-phase microextraction (SPME)

<b>Tool</b>	<ul style="list-style-type: none"><li>• Includes fiber holder for standard SPME fibers, 23-gauge</li><li>• Compatible with SSL and PTV injectors</li></ul>
<b>Vial volumes</b>	<ul style="list-style-type: none"><li>• Compatible with 2, 10 and 20 mL vials</li></ul>
<b>Sample capacity</b>	<ul style="list-style-type: none"><li>• Depending on autosampler configuration<ul style="list-style-type: none"><li>– Up to 648 × 2 mL sample vials or 972 × 2 mL vials with the extended X-arm</li><li>– Up to 180 × 10 or 20 mL vials or 300 × 10 or 20 mL vials with the extended X-arm</li></ul></li></ul>
<b>Incubation oven</b>	<ul style="list-style-type: none"><li>• Capacity 6 × 20 mL vials (compatible with 2 and 10 mL vials with adapters)</li><li>• 30–200 °C temperature range, in 1 °C steps</li><li>• 250–750 rpm agitation speed</li></ul>
<b>Incubation time</b>	<ul style="list-style-type: none"><li>• 0.1 to 600.0 min in 0.1 min increments</li></ul>
<b>Vial penetration depth</b>	<ul style="list-style-type: none"><li>• Standard or custom between 20 mm and 70 mm, suitable for headspace or direct immersion (DI) extraction</li></ul>
<b>Fiber conditioning station</b>	<ul style="list-style-type: none"><li>• Optional, 2-ports, 40–350 °C, inert gas purged</li><li>• Suitable for both SPME and SPME Arrow fibers</li></ul>
<b>Fiber types</b>	<ul style="list-style-type: none"><li>• 10 mm fiber length</li><li>• PDMS (7, 30, 100 µm), Polyacrylate (85 µm), Carbon WR/PDMS (95 µm), DVB/PDMS (65 µm), DVB/Carbon WR/PDMS (50-30 µm)</li></ul>

## SPME Arrow

<b>Tool</b>	<ul style="list-style-type: none"><li>• Includes fiber holder for SPME Arrow fibers, 1.1 and 1.5 mm o.d.</li><li>• Compatible with SSL injector</li></ul>
<b>Vial volumes</b>	<ul style="list-style-type: none"><li>• Compatible with 10 and 20 mL vials</li></ul>
<b>Sample capacity</b>	<ul style="list-style-type: none"><li>• Depending on autosampler configuration<ul style="list-style-type: none"><li>– Up to 180 × 10 or 20 mL vials or 300 × 10 or 20 mL vials with the extended X-arm</li></ul></li></ul>
<b>Incubation oven</b>	<ul style="list-style-type: none"><li>• Capacity 6 × 20 mL vials (compatible with 2 and 10 mL vials with adapters)</li><li>• 30–200 °C temperature range, in 1 °C steps</li><li>• 250–750 rpm agitation speed</li></ul>
<b>Incubation time</b>	<ul style="list-style-type: none"><li>• 0.1 to 600.0 min in 0.1 min increments</li></ul>
<b>Vial penetration depth</b>	<ul style="list-style-type: none"><li>• Standard or custom between 20 mm and 70 mm, suitable for headspace or direct immersion (DI) extraction</li></ul>
<b>Heatex-Stirrer</b>	<ul style="list-style-type: none"><li>• For intensive heating and stirring during the extraction step, 40–200 °C, 0–1600 rpm</li></ul>
<b>Fiber conditioning station</b>	<ul style="list-style-type: none"><li>• Optional, 2-ports, 40–350 °C, inert gas purged. Suitable for both SPME and SPME Arrow fibers</li></ul>
<b>Fiber types</b>	<ul style="list-style-type: none"><li>• 20 mm fiber length</li><li>• PDMS (100 µm/1.1 mm o.d., 250 µm/1.5 mm o.d.), Polyacrylate (100 µm/1.1 mm o.d.)</li><li>• Carbon WR/PDMS (120 µm/1.1 mm o.d.), DVB/PDMS (120 µm/1.1 mm o.d.), DVB/Carbon WR/PDMS (120 µm/1.1 mm o.d.)</li></ul>

## In-Tube Extraction Dynamic Headspace (ITEX-DHS) option

<b>Tool</b>	<ul style="list-style-type: none"><li>• Includes sampling gas-tight syringe, focusing trap, built-in trap heating and cooling fan and trap cleaning capability</li></ul>
<b>Temperatures</b>	<ul style="list-style-type: none"><li>• Trap 30–350 °C</li><li>• Syringe 40–150 °C</li></ul>
<b>Extraction parameters</b>	<ul style="list-style-type: none"><li>• Flow rate 10–1000 µL/s, stroke cycles 0-1000, volume 0-1300 µL, incubation time up to 600 min, water removal step</li></ul>
<b>Vial volumes</b>	<ul style="list-style-type: none"><li>• Compatible with 20 mL vials</li></ul>
<b>Vial penetration depth</b>	<ul style="list-style-type: none"><li>• Standard or custom between 10 and 35 mm</li></ul>
<b>Sample capacity</b>	<ul style="list-style-type: none"><li>• Depending on autosampler configuration<ul style="list-style-type: none"><li>– Up to 180 × 10 or 20 mL vials or 300 × 10 or 20 mL vials with the extended X-arm</li></ul></li></ul>
<b>Incubation oven</b>	<ul style="list-style-type: none"><li>• Capacity 6 × 20 mL vials (compatible with 2 and 10 mL vials with adapters)</li><li>• 30–200 °C temperature range, in 1 °C steps</li><li>• 250–750 rpm agitation speed</li></ul>
<b>Incubation time</b>	<ul style="list-style-type: none"><li>• 0.1 to 600.0 min in 0.1 min increments</li></ul>
<b>Traps</b>	<ul style="list-style-type: none"><li>• Tenax TA 80/100 mesh as standard, other single- or multi-layer microtraps available for volatile and semi-volatile compounds enrichment</li></ul>

## Dilutor tool

<b>Tool</b>	<ul style="list-style-type: none"><li>• Available for single- or multi-solvent dispensing (up to four), compatible with standard X-arm length only</li></ul>
<b>Dispensing syringe size</b>	<ul style="list-style-type: none"><li>• 100 µL, 1 mL, 5 mL, 10 mL</li></ul>
<b>Dispensed volume</b>	<ul style="list-style-type: none"><li>• From 10 µL to 1000 µL, depending on the syringe volume</li></ul>
<b>Dilution cycle</b>	<ul style="list-style-type: none"><li>• Dedicated PreCycles for in-batch single- or multi-solvent addition, optional pre- and post-washing step, optional mixing step (requires Vortexer module)</li></ul>

**Micro Solid-Phase Extraction ( $\mu$ SPE)**

<b><math>\mu</math>SPE tool kit</b>	<ul style="list-style-type: none"> <li>Includes hardware for <math>\mu</math>SPE handling, script, standard operation procedure and quick installation guide</li> </ul>
<b>Sample capacity</b>	<ul style="list-style-type: none"> <li>54 (standard)–108 <math>\times</math> 2-mL sample vials, elution vials and <math>\mu</math>SPE cartridges</li> </ul>
<b><math>\mu</math>SPE syringe volume</b>	<ul style="list-style-type: none"> <li>1000 <math>\mu</math>L for conditioning/elution solvent and raw sample</li> </ul>
<b>Liquid syringe volume</b>	<ul style="list-style-type: none"> <li>10 and 25 <math>\mu</math>L for ISTD/protectant addition and clean sample injection</li> </ul>
<b>Elution speed</b>	<ul style="list-style-type: none"> <li>Optimized at 2 <math>\mu</math>L/s</li> </ul>
<b>Elution solvent capacity</b>	<ul style="list-style-type: none"> <li>3 <math>\times</math> 100 mL</li> </ul>
<b>Washing solvent capacity</b>	<ul style="list-style-type: none"> <li>2 <math>\times</math> 1000 mL</li> </ul>
<b><math>\mu</math>SPE workflow</b>	<ul style="list-style-type: none"> <li>Optimized QuEChERS extracts clean-up workflow with internal standard and analyte protectant addition (optional), <math>\mu</math>SPE cartridges conditioning (optional) and on-line GC injection</li> </ul>

Find out more at [thermofisher.com/triplusrsh](https://thermofisher.com/triplusrsh)