


# Thermo Scientific Dionex IonPac AS22 and AS22-Fast Anion-Exchange Columns

The Thermo Scientific™ Dionex™ IonPac™ AS22 column is a carbonate based anion-exchange column designed for the determination of inorganic anions and low-molecular-weight organic acids including fluoride, acetate, chloride, nitrite, bromide, nitrate, phosphate, and sulfate. The Dionex IonPac AS22 column is used with isocratic carbonate/bicarbonate eluents and suppressed conductivity detection. The Dionex IonPac AS22 column can be used in combination with Dionex Eluent Generators and the Electrolytic pH Modifier (EPM), which automatically produce potassium carbonate/bicarbonate eluents from water. The Dionex IonPac AS22-Fast column is ideal for fast IC using eluent generation or eluent regeneration.

Analysis of a Municipal Drinking Water Sample Using the  
Dionex IonPac AS22-Fast 4 × 150 mm Column



## Column Performance

The Dionex IonPac AS22 column meets the performance requirements specified in U.S. Environmental Protection Agency (EPA) Method 300.0 (A). The common inorganic anions are easily separated in a variety of sample matrices, including drinking water, wastewater, process streams, and scrubber solutions. The Dionex IonPac AS22 column selectivity provides excellent retention of fluoride from the water dip and resolution of fluoride, acetate, and formate. Solvent compatibility permits easy column cleanup after the analysis of samples with hydrophobic components. The Dionex IonPac AS22 column is available in 4 mm, 2 mm and 0.4 mm formats.

## Fast IC

The Dionex IonPac AS22-Fast column is ideal for Fast IC as it is designed to have sufficient capacity to maintain resolution even in a short column format. Fast separations are achieved on any Thermo Scientific Dionex system at higher flow rates. In a short column format, backpressures produced at higher flow rates are reduced while allowing overall shorter run times. This allows for the determination of anions with high resolution even in drinking, surface, groundwater, and wastewater matrices in under 5 min. Laboratories can achieve higher productivity and increased throughput. The Dionex IonPac AS22-Fast column is available in 4 mm and 2 mm formats.

## Analyze Inorganic Anions in Diverse Sample Matrices

- Source water and drinking water
- Municipal and industrial wastewater
- Industrial cooling water
- Hazardous waste extracts and dump site leachates
- Acid rain
- Foods and beverages
- Anionic counterions in pharmaceutical preparations and synthetic peptides

- Polymers such as polyols and polysulfonates
- Scrubber solutions

## Superior Chromatographic Performance

- Fast isocratic separation of the common inorganic anions in simple sample matrices in 8 min.
- Isocratic separation of common inorganic anions in complex sample matrices in 12 min.
- High speed separation of the common inorganic anions in less than 5 min using the Dionex IonPac AS22-Fast Column.
- Carbonate peak well resolved from common inorganic anions.
- Meets performance requirements specified in EPA Method 300.0 (A).
- Ideal alternative for Dionex IonPac AS4A-SC, AS12A, AS14, and AS14A inorganic anion applications.
- Simplified Reagent-Free™ ion chromatography (RFIC™) System operation provided by Dionex Eluent Generators and EPM, which require only a deionized water source to produce potassium carbonate/bicarbonate eluent.
- Ideal column for use with Eluent Regeneration, enabling simple, non-stop operation for up to 28 days.

- Simple, accurate eluent preparation with the Dionex IonPac AS22 Eluent Concentrate—just dilute in deionized water and start operation.
- Eluent suppression using the Thermo Scientific Dionex AERS 500 Anion Electrolytically Regenerated Suppressor or Thermo Scientific™ Dionex™ ACES™ Anion Capillary Electrolytic Suppressor technology provides RFIC operation with low backgrounds and enhanced analyte sensitivity.
- High capacity: 210 µeq/col. (4 × 250 mm column).
- Operate at ambient or elevated temperatures. Column selectivity is optimized for a 30 °C operating temperature to ensure reproducible retention times.
- Compatible with organic solvents to enhance analyte solubility, modify column selectivity, or for effective column cleanup.

## High Efficiency Particle Structure

The Dionex IonPac AS22 column was developed using a unique polymer bonding technology. The stationary phase consists of a novel hyper-branched anion-exchange condensation polymer, electrostatically attached to the surface of a wide-pore polymeric substrate. The substrate is surface-sulfonated in exactly the same manner as is common in Dionex latex coated anion-exchange materials. However, in this anion-exchange resin, alternating treatments of epoxy monomer and amines produce a coating which is grown directly off the substrate surface as illustrated in Figure 1. Resin capacity is controlled through the number of alternating coating cycles. The Dionex IonPac AS22 column uses a high-capacity resin (210 µeq/4 mm column) with optimized selectivity for the common inorganic anions in a variety of sample matrices.




Figure 1. Structure of a Dionex IonPac AS22 column particle.

## Economical Capillary Operation

The Dionex IonPac AS22 column is available in 0.4 mm i.d. for capillary operation offering the advantage of reduced operating costs.

- Ideal for limited sample volumes due to higher mass sensitivity.
- One hundred fold reduction in eluent consumption and waste disposal.
- 4 mm applications can be directly transferred to the 0.4 mm format by reducing flow rate by one hundred fold.

Figure 2 is an example of the separation of inorganic anions and acetate using the Dionex IonPac AS22 capillary column.

Column:	Dionex IonPac AG22/AS22 (0.4 × 250 mm)	Peaks:	1. Fluoride	1.25 mg/L (ppm)
Eluent:	4.5 mM sodium carbonate/ 1.4 mM sodium bicarbonate		2. Acetate	5.0
Temperature:	30 °C		3. Chloride	2.5
Flow Rate:	12 µL/min		4. Nitrite	3.75
Inj. Volume:	0.4 µL		5. Bromide	6.25
Detection:	Suppressed conductivity, Dionex ACES 300 suppressor, 4 mm, AutoSuppression recycle mode		6. Nitrate	6.25
			7. Phosphate	10.0
			8. Sulfate	7.5




Figure 2. Separation of inorganic anions and acetate using the Dionex IonPac AG22/AS22 capillary column.

## Isocratic Separation of Common Inorganic Anions

The Dionex IonPac AS22 column provides excellent separation of fluoride, chloride, nitrite, bromide, nitrate, phosphate and sulfate using an isocratic carbonate/bicarbonate eluent and suppressed conductivity detection. Using a 4.5 mM carbonate/1.4 mM bicarbonate eluent, the common inorganic anions plus acetate can be resolved in approximately 12 min as shown in Figure 3.

Column:	Dionex IonPac AG22, AS22, 4 mm	Peaks:	1. Fluoride	5 mg/L (ppm)
Eluent:	4.5 mM sodium carbonate/ 1.4 mM sodium bicarbonate		2. Acetate	20
Temperature:	30 °C		3. Chloride	10
Flow Rate:	1.2 mL/min		4. Nitrite	15
Inj. Volume:	10 µL		5. Bromide	25
Detection:	Suppressed conductivity, Dionex ASRS 300 suppressor, 4 mm, AutoSuppression recycle mode		6. Nitrate	25
			7. Phosphate	40
			8. Sulfate	30




Figure 3. Separation of common inorganic anions plus acetate on a Dionex IonPac AS22 column using a 4.5 mM carbonate/1.4 mM bicarbonate eluent.

## Inorganic Anions in Drinking Water

The Dionex IonPac AS22 column is ideal for compliance monitoring of drinking water and wastewater. The column meets or exceeds the performance requirements of U.S EPA Method 300.0 (A). Common inorganic anions are separated in approximately 12 min in drinking water as shown in Figure 4. The analysis of a drinking water sample using the Dionex IonPac AS22 capillary column is shown in Figure 5.

Column:	Dionex IonPac AG22, AS22, 4 mm	Peaks:	1. Fluoride	0.84 mg/L (ppm)	7. Bromide	0.02
Eluent:	4.5 mM sodium carbonate/ 1.4 mM sodium bicarbonate		2. Formate	0.03	8. Nitrate	0.89
Temperature:	30 °C		3. Chloride	15.59	9. Carbonate	NQ
Flow Rate:	1.2 mL/min		4. Nitrite	0.01	10. Phosphate	0.22
Inj. Volume:	25 µL		5. Unknown	NQ	11. Sulfate	20.45
Detection:	Suppressed conductivity, Dionex ASRS 300 suppressor, 4 mm, AutoSuppression recycle mode		6. Chlorate	0.18		
					NQ: Not Quantified	




Figure 4. Determination of inorganic anions in a municipal drinking water sample using the Dionex IonPac AS22 column.

Column:	Dionex IonPac AG22/AS22 (0.4 × 250 mm)	Peaks:	1. Fluoride	2.06 mg/L (ppm)	7. Nitrate	0.86
Eluent:	4.5 mM sodium carbonate/ 1.4 mM sodium bicarbonate		2. Formate	NQ	8. Phosphate	0.10
Temperature:	30 °C		3. Chloride	16.41	9. Sulfate	21.67
Flow Rate:	0.012 mL/min		4. Nitrite	0.04	10. Unknown	NQ
Inj. Volume:	0.4 µL		5. Unknown	NQ		
Detection:	Suppressed conductivity, Dionex ACES 300 suppressor, 4 mm, AutoSuppression recycle mode				NQ: Not Quantified	




Figure 5. Analysis of municipal drinking water using the Dionex IonPac AS22 capillary column.

The Dionex IonPac AS22-Fast column is suitable for separations requiring higher flow rates for the fast analysis of inorganic anions in simple sample matrices as shown in Figure 6. The Dionex IonPac AS22-Fast column (4 × 150 mm) was operated with the same eluent at 1.2 and 2.0 mL/min. Using a higher flow rate in a shorter column format, the overall run time was reduced to 4.5 min with optimal resolution of the common inorganic anions.

The use of a Dionex IonPac AS22-Fast 2 × 150 mm column format allows greater linear velocities and reduced eluent consumption compared to a 4 mm format. Figure 7 shows the excellent resolution for the common inorganic anions in less than 4 min.

Column:	Dionex IonPac AS22-Fast (4 × 150 mm)	Peaks:	1. Fluoride	5.0 mg/L (ppm)
Eluent:	4.5 mM sodium carbonate/ 1.4 mM sodium bicarbonate		2. Chloride	10.0
Flow Rate:	A. 1.2 mL/min B. 2.0 mL/min		3. Nitrite	15.0
Inj. Volume:	10 µL		4. Bromide	25.0
Temperature:	30 °C		5. Nitrate	25.0
Detection:	Suppressed conductivity, Dionex ASRS 300 suppressor, 4 mm, AutoSuppression recycle mode		6. Phosphate	40.0
			7. Sulfate	30.0




Figure 6. Separation of common inorganic anions using the Dionex IonPac AS22-Fast 4 × 150 mm column.

Column:	Dionex IonPac AS22-Fast (2 × 150 mm)	Peaks:	1. Fluoride	5.0 mg/L (ppm)
Eluent:	4.5 mM sodium carbonate/ 1.4 mM Sodium bicarbonate		2. Chloride	10.0
Flow Rate:	A. 0.3 mL/min B. 0.5 mL/min C. 0.63 mL/min		3. Nitrite	15.0
Inj. Volume:	2.5 µL		4. Bromide	25.0
Temperature:	30 °C		5. Nitrate	25.0
Detection:	Suppressed conductivity, Dionex ASRS 300 suppressor, 2 mm, AutoSuppression recycle mode		6. Phosphate	40.0
			7. Sulfate	30.0




Figure 7. Separation of common inorganic anions using the Dionex IonPac AS22-Fast 2 × 150 mm column.

Figure 8 illustrates the analysis of a municipal drinking water sample using the Dionex IonPac AS22-Fast ( $4 \times 150$  mm) column. At 2.0 mL/min, the anions can be separated in less than 5 min. The bottom trace shows an enlarged image of the separation. The anions are well resolved even in the presence of disparate anion concentrations including fluoride/formate and chloride/nitrite.

### Extended Application Capabilities

The unique selectivity and high capacity of the Dionex IonPac AS22 column makes it ideal for methods development of specialized anion applications. Figure 9 shows the separation of a variety of environmental anions including inorganic anions, oxyhalides, oxyanions, and organic acids using an isocratic carbonate/bicarbonate eluent. These 18 analytes are easily separated in less than 30 min.




Figure 8. Analysis of a municipal drinking water sample using the Dionex IonPac AS22-Fast  $4 \times 150$  mm column.

Column:	Dionex IonPac AG22, AS22, 4 mm	Peaks:	1. Fluoride	3 mg/L (ppm)	10. Phosphate	40
Eluent:	4.8 mM sodium carbonate/ 1.0 mM sodium bicarbonate		2. Acetate	20	11. Sulfate	30
Flow Rate:	1.5 mL/min		3. Formate	10	12. Selenate	30
Temperature:	30 $^{\circ}$ C		4. Bromate	20	13. Iodide	40
Inj. Volume:	10 $\mu$ L		5. Chloride	5	14. Arsenate	30
Detection:	Suppressed conductivity, Dionex ASRS 300 suppressor, 4 mm, AutoSuppression recycle mode				15. Thiocyanate	40
			6. Nitrite	15	16. Perchlorate	40
			7. Chlorate	25	17. Thiosulfate	40
			8. Bromide	25	18. Chromate	40
			9. Nitrate	25		




Figure 9. Anion separation including inorganic anions, organic acids, oxyanions, and oxyhalides using the Dionex IonPac AS22 column.

The Dionex IonPac AS22 column is also ideal for the determination of anions in food and beverage samples containing high concentrations of carbonate. Carbonated bottled water contains up to 500 mg/L (ppm) of carbonate. The unique selectivity of the Dionex IonPac AS22 column positions carbonate well away from the common inorganic anions and does not interfere with quantification of the analytes of interest. Figure 10 shows the determination of inorganic anions in carbonated water using the Dionex IonPac AS22 column. Figure 11 shows the analysis of carbonated water using the Dionex IonPac AS22 capillary column. The sample was sonicated for 15 min to remove some of the carbonate.

The Dionex IonPac AS22 column can be used to evaluate the mass balance of drugs and synthetic peptide preparations. Figure 12 illustrates the use of the Dionex IonPac AS22 column to determine the anionic counterion amount and type. The common inorganic anions plus trifluoroacetate can be determined in less than 14 min.




Figure 10. Determination of inorganic anions in carbonated water using the Dionex IonPac AS22 column.





Figure 11. Analysis of the common inorganic anions in carbonated water using the Dionex IonPac AS22 capillary column.





Figure 12. Analysis of the common inorganic anions and TFA using the Dionex IonPac AS22 column.

## System Requirements

The Dionex IonPac AS22 Capillary Column is recommended for use with the Thermo Scientific™ Dionex™ ICS-5000+ or ICS-4000 Capillary Reagent-Free™ HPIC™ system. The Dionex IonPac AS22 and AS22-Fast Analytical Columns are recommended for use with Dionex ICS-2100 or ICS-5000+ RFIC system equipped with an eluent generator and EPM. The Dionex IonPac AS22 Analytical Column can also be used with older Dionex IC systems equipped with an eluent generator of a Thermo Scientific Dionex RFC-30 Reagent-Free Controller. The eluent generator is used to automatically produce carbonate and bicarbonate eluents from deionized water.

### Suppressor Recommendation

For optimum ease-of-use and economy, the Dionex IonPac AS22 and AS22-Fast columns should be used with the Dionex AERS 500 Anion Electrolytically Regenerated Suppressor or the Dionex ACES-300 Anion Capillary Electrolytic Suppressor.

### Concentrator Columns

For concentrator work with a 2 mm or 4 mm Dionex IonPac AS22 column, use the: Dionex IonPac AG22 guard column; Ultratrace Anion Concentrator Columns (Dionex IonPac UTAC-ULP1, UTAC-XLP1, UTAC-ULP2, or UTAC-XLP2) or Trace Anion Concentrator Column (Dionex IonPac TAC-ULP1) when a single piston pump such as the Thermo Scientific Dionex AXP Auxiliary Pump (pulse damper required) is used for sample delivery. In addition to the concentrator columns listed above, use the Dionex IonPac UTAC-LP1, UTAC-LP2 or TAC-LP1 when the sample is delivered using a syringe or a low-pressure AutoSampler, e.g., the Thermo Scientific Dionex AS-DV AutoSampler.

For concentrator work with a 0.4 mm capillary column, use the Dionex IonPac AG22 Capillary Guard Column or the Thermo Scientific™ Dionex™ IonSwift™ MAC-100 concentrator column.

SPECIFICATIONS	
<b>Dimension</b>	Dionex IonPac AS22 Analytical Column (2 × 250 mm), (4 × 250 mm) Dionex IonPac AS22 Capillary Column (0.4 × 250 mm) Dionex IonPac AS22-Fast Analytical Column (2 × 150 mm), (4 × 250 mm) Dionex IonPac AG22 Guard Column (2 × 50 mm), (4 × 50 mm) Dionex IonPac Capillary Guard Column (0.4 × 50 mm) Dionex IonPac AG22-Fast Guard Column (2 × 30 mm), (4 × 30 mm)
<b>Maximum Operating Pressure</b>	3000 psi
<b>Mobile Phase Compatibility</b>	pH 0–14; 0–100% HPLC solvents
<b>Substrate Characteristics</b>	
<b>Analytical Column</b>	Supermacroporous resin
Bead Diameter (µm)	6.5
Pore Size Å	2000
Cross-Linking (%DVB)	55
<b>Guard Column</b>	Microporous resin
Bead Diameter (µm)	11
Pore Size Å	<1
Cross-Linking (%DVB)	55%
<b>Ion-Exchange Functional Group</b>	
Alkanol quaternary ammonium ion	
<b>Functional Group Characteristics</b>	
Hydrophobicity Ultralow	
<b>Capacity</b>	
2.1 ueq	(0.4 x 250 mm capillary column)
0.06 ueq	(90.4 x 50 mm capillary guard column)
52.5 µeq	(2 × 250 mm analytical column)
1.5 µeq	(2 × 50 mm guard column)
210 µeq	(4 × 250 mm analytical column)
6 µeq	(4 × 50 mm guard column)
126 µeq	(4 × 150 mm analytical column)
4.0 µeq	(4 × 30 mm guard column)
31.5 µeq	(2 × 150 mm analytical column)
1.0 µeq	(2 × 30 mm guard column)
<b>Column Construction</b>	PEEK™ with 10–32 threaded ferrule style end fittings. All components are nonmetallic.



## Ordering Information

For more information or to place an order, contact the Thermo Scientific Dionex Products office nearest you or your local distributor. Phone numbers and addresses for worldwide subsidiaries can be found in the About Us section of [www.thermoscientific.com/dionex](http://www.thermoscientific.com/dionex).

Product Description	Part Number
Dionex IonPac AS22 Capillary Column (0.4 × 250 mm)	079057
Dionex IonPac AG22 Capillary Guard Column (0.4 × 50 mm)	079058
Dionex IonPac AS22 Analytical Column (4 × 250 mm)	064141
Dionex IonPac AG22 Guard Column (4 × 50 mm)	064139
Dionex IonPac AS22 Analytical Column (2 × 250 mm)	064137
Dionex IonPac AG22 Guard Column (2 × 50 mm)	064135
Dionex IonPac AS22-Fast Analytical Column (4 × 150 mm)	079936
Dionex IonPac AG22-Fast Guard Column (4 × 30 mm)	072784
Dionex IonPac AS22-Fast Analytical Column (2 × 150 mm)	079937
Dionex IonPac AG22-Fast Guard Column (2 × 30 mm)	072785
Dionex IonPac AMC-1 Anion MicroConcentrator (2 × 15 mm)	051760
Dionex IonPac TAC-2 Trace Anion Concentrator (3 × 35 mm)	043101
Dionex IonPac TAC-LP1 Low Pressure Trace Anion Concentrator (4 × 35 mm)	046026
Dionex IonPac TAC-ULP1 Ultra Low Pressure Trace Anion Concentrator (5 × 23 mm)	061400
Dionex IonPac UTAC-LP1 Ultra Trace Anion Concentrator–Low Pressure (4 × 35 mm)	063079
Dionex IonPac UTAC-ULP1 Ultra Trace Anion Concentrator–Ultra Low Pressure (5 × 23 mm)	063475
Dionex IonPac UTAC-XLP1 Ultra Trace Anion Concentrator–Extremely Low Pressure (6 × 16 mm)	063459
Dionex IonPac UTAC-LP2 Ultra Trace Anion Concentrator- Low Pressure (4 × 35 mm)	079917
Dionex IonPac UTAC-ULP2 Ultra Trace Anion Concentrator- Ultra Low Pressure (5 × 23 mm)	079918
Dionex IonPac UTAC-XLP2 Ultra Trace Anion Concentrator- Extremely Low Pressure (6 × 16 mm)	072781
Dionex IonPac AS22 Sodium Carbonate/Bicarbonate Eluent Concentrate (250 mL of 100X concentrate)	063965

## [www.thermoscientific.com/dionex](http://www.thermoscientific.com/dionex)

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. PEEK is a trademark of Victrex PLC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific, Sunnyvale, CA  
USA is ISO 9001:2008 Certified.

<b>Africa</b> +43 1 333 50 34 0	<b>Denmark</b> +45 70 23 62 60	<b>Japan</b> +81 6 6885 1213	<b>Russia/CIS</b> +43 1 333 50 34 0
<b>Australia</b> +61 3 9757 4300	<b>Europe-Other</b> +43 1 333 50 34 0	<b>Korea</b> +82 2 3420 8600	<b>Singapore</b> +65 6289 1190
<b>Austria</b> +43 810 282 206	<b>Finland</b> +358 9 3291 0200	<b>Latin America</b> +1 561 688 8700	<b>Sweden</b> +46 8 556 468 00
<b>Belgium</b> +32 53 73 42 41	<b>France</b> +33 1 60 92 48 00	<b>Middle East</b> +43 1 333 50 34 0	<b>Switzerland</b> +41 61 716 77 00
<b>Brazil</b> +55 11 3731 5140	<b>Germany</b> +49 6103 408 1014	<b>Netherlands</b> +31 76 579 55 55	<b>Taiwan</b> +886 2 8751 6655
<b>Canada</b> +1 800 530 8447	<b>India</b> +91 22 6742 9494	<b>New Zealand</b> +64 9 980 6700	<b>UK/Ireland</b> +44 1442 233555
<b>China</b> 800 810 5118 (free call domestic) 400 650 5118	<b>Italy</b> +39 02 950 591	<b>Norway</b> +46 8 556 468 00	<b>USA</b> +1 800 532 4752

# Thermo

SCIENTIFIC

A Thermo Fisher Scientific Brand