

Thermo Scientific Delta Ray

Isotope Ratio Infrared Spectrometer

Thermo Scientific™ Delta Ray™ is a high performance, mid-infrared laser based, isotope ratio infrared spectrometer. It offers simultaneous determination of $\delta^{13}\text{C}$ and $\delta^{18}\text{O}$ in CO_2 at ambient concentrations with a precision as low as 0.05‰.

- $\delta^{13}\text{C}$ and $\delta^{18}\text{O}$ in CO_2
- Precision as low as 0.05‰
- Mid-infrared laser based isotope ratio infrared spectrometer
- Portable, field deployable
- Universal Reference Interface – URI
- Up to 1 Hz data
- 200 ppm to 100% CO_2 with dilution
- Global support



Precise, Verifiable Isotope Ratio Measurements

The Delta Ray analyzer comes with a universal reference interface (URI) that provides fully automated referencing and calibration for verifiable measurements and long term confidence. Smart referencing adjusts the reference gas concentration to the sample to achieve superior performance. With the optional dilution box, samples with CO_2 concentration from 200 ppm to 100% can be analyzed.

Isotope Ratios To-Go

The modular design, low weight and built-in referencing make the Delta Ray the analyzer of choice for demanding field applications such as greenhouse gas monitoring, ecology, plant science, carbon sequestration and storage research, or volcanic monitoring. No need to fill vials and wait for the analysis back in the laboratory. Instead feature rich data is acquired autonomously 24 hours a day, 7 days a week, with a time resolution as low as 1 second directly where the sample originates. A quickstart guide gets you going in no time.

Real-Time Results

Powered by the Thermo Scientific Qtegra™ Intelligent Scientific Data Solution™ (ISDS), the data acquisition is driven by powerful workflows. Select sample ports, drive external valves, or synchronize your measurements with external triggers all with a few clicks. Qtegra software then acquires the data into a Qtegra LabBook without user interaction. The fully referenced and calibrated results can be inspected immediately and are easily exported to a spreadsheet or your favorite data analysis package. Qtegra software runs on the on-board Windows 7 computer, which in connection with an internet connection allows you to control your Delta Ray analyzer from anywhere in the world.

Thermo
SCIENTIFIC

Simple, Robust, Powerful Mid-Infrared Technology

At the core of the analyzer is a difference frequency generation (DFG) laser that operates in the mid-infrared, where the absorption of CO₂ and its isotopologues are so strong, that 5 m path length are sufficient to achieve a precision as low as 0.05‰. It is simple, robust and powerful.

Sample gas is dried in the analyzer to prevent any interaction of water and CO₂ (isotope exchange) with the added benefit of providing dry mole fraction concentration data. In the measurement cell, where the laser analyzes the gas, pressure and temperature are precisely controlled.

Analytical Performance

| Performance Specifications, in Air | δ ¹³ C | δ ¹⁸ O | CO ₂ Concentration |
|--|--|-------------------|-------------------------------|
| Precision, 60 sec averaging, 30 replicates | 0.15‰ 1 SD | 0.2‰ 1 SD | 70 ppb |
| External reproducibility*, 300 sec averaging, 10 replicates | 0.07‰ 1 SD | 0.1‰ 1 SD | |
| Reproducibility, 1 h resolution, 24 h | 0.05‰ 1 SD | 0.05‰ 1 SD | |
| Operating range | 200 - 3,500 ppm | | |
| Guaranteed specifications range without dilution (100% CO ₂ with optional dilution) | 380 - 1,500 ppm | | |
| Measurement interval | 1, 10, 60 sec | | |
| Cell exchange time (response time) | 35 sec @ 80 sccm sample flow | | |
| Measurement technique | Tunable laser direct absorption in the mid infrared (4.3 μm) TDLAS | | |

*demonstrated with Zero enrichment template: sample and reference are the same gas, 5 min sample average, 5 min reference average, repeat at least 10 times.

Operating Conditions

| | |
|----------------------|---|
| Temperature range | 10 - 35°C |
| Temperature gradient | 0.2°C/min |
| Humidity | 10 - 80% R.H. up to 31°C 10 - 50% R.H. at 35°C non condensing |
| Size | 19" rack mount chassis Analyzer: 588 x 424 x 219 mm (D x W x H) URI: 588 x 424 x 131 mm (D x W x H) |
| Weight | Analyzer: 25 kg URI: 12 kg |
| Power consumption | 100 - 240 V, 50 - 60 Hz Analyzer: 150 W typical / 500 W startup at 25°C URI: < 160 W |
| Sample temperature | -10 to 45°C |
| Sample pressure | Sensor input 700 - 1,200 mbar |
| Sample flow | Sensor input 80 sccm @ 1,000 mbar abs. |
| Pump | Internal, oil free |
| Qtegra ISDS licenses | Each Delta Ray system comes with a Qtegra license for installation on one (1) instrument and up to three (3) additional desktop installations for data review and manipulation. Additional desktop licenses can be purchased as an option |

Connections

| | |
|------------|--|
| Electrical | Video: DVI port (incl. VGA) 4 USB ports 1 ethernet 1 Gbit 4 digital outs for 12 V or 24 V valves max 3.3 W each 2 digital in for trigger (TTL) 3 analog outs 0-5 V (0.5 mV resolution) 1 analog in 0-5 V (1 mV resolution) |
| Gas | Sample A (1 - 12 bar relative) Sample B (700 - 1,200 mbar (abs.)), passing through dryer. Carrier out: for external dilution User provided calibration gas: 2 x pure CO ₂ with known δ ¹³ C and δ ¹⁸ O, synthetic air (CO ₂ < 0.2 ppm) |
| Fittings | 1/16" Swagelok, 10 - 32 coned port |

Transport

| | |
|-------------|----------------------------|
| Vibration | MIL-STD 810 G |
| Temperature | -33 to 63°C non condensing |

www.thermoscientific.com/DeltaRay

© 2013 Thermo Fisher Scientific Inc. All rights reserved. Microsoft and Windows are trademarks of Microsoft Corp. Swagelok is a trademark of Swagelok Company. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Thermo Fisher Scientific (Bremen) GmbH
Management System Registered to
ISO 9001:2008

| | | | |
|-------------------------------------|---|--------------------------------------|-------------------------------------|
| Africa-Other +27 11 570 1840 | Europe-Other +43 1 333 50 34 0 | Japan +81 45 453 9100 | Spain +34 914 845 965 |
| Australia +61 3 9757 4300 | Finland/Norway/Sweden +46 8 556 468 00 | Latin America +1 561 688 8700 | Switzerland +41 61 716 77 00 |
| Austria +43 1 333 50 34 0 | France +33 1 60 92 48 00 | Middle East +43 1 333 50 34 0 | UK +44 1442 233555 |
| Belgium +32 53 73 42 41 | Germany +49 6103 408 1014 | Netherlands +31 76 579 55 55 | USA +1 800 532 4752 |
| Canada +1 800 530 8447 | India +91 22 6742 9434 | New Zealand +64 9 980 6700 | |
| China +86 10 8419 3588 | Italy +39 02 950 591 | Russia/CIS +43 1 333 50 34 0 | |
| Denmark +45 70 23 62 60 | | South Africa +27 11 570 1840 | |

PS30277_E 11/13G

Thermo

SCIENTIFIC

Part of Thermo Fisher Scientific