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References

Dionex OnGuard II A

1. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Determination of Sialic Acids Using UHPLC with Fluorescence Detection. Application Notes 266; LPN 2662. Sunnyvale, CA, 2011.
2. Thermo Fisher Scientific. Determination of Trace Sodium in Cranberry Powder. Application Note 1003; LPN 3080. Sunnyvale, CA, 2012.
3. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Application Note 87, Determination of Sugar Alcohols in Confections and Fruit Juices by HPAEC with Pulsed Amperometric Detection.
4. van den Hoop, Marc A. G. T.; Cleven, Rob F. M. J.; van Staden, Johannes J.; Neele, Jos. "Analysis of Fluoride in Rain Water—Comparison of Capillary Electrophoresis with Ion Chromatography and Ion-Selective Electrode Potentiometry". *J. Chromatogr., A*. **1996**, 739, 241–248.
5. Rohrer, Jeffrey; Thayer, Jim; Avdalovic, Nebojsa; Weitzhandler, Michael. "HPAEC-PAD Analysis of Monoclonal Antibody Glycosylation". *Techniques in Protein Chemistry VI*, Academic Press, **1995**, 65–73.

Dionex OnGuard II Ag

1. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Determination of Perchlorate in Drinking Water by Ion Chromatography . Application Update 145; LPN 1540. Sunnyvale, CA, 2004.
2. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Determination of Perchlorate in Vegetation Samples Using Accelerated Solvent Extraction and Ion Chromatography. Application Note 356; LPN 1830. Sunnyvale, CA, 2006.
3. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Application Note 101, Trace Level Determination of Bromate in Ozonated Drinking Water Using Ion Chromatography.
4. Slingsby, R. W.; Pohl, C. A. "Approaches to Sample Preparation for IC-Sulfate Precipitation on Barium-Form Ion Exchangers". *J. Chromatogr., A*. **1996**, 739, 49–55.
5. Charles, L.; Pepin, D. "Analysis of Oxyhalides in Water by Ion Chromatography—Ionspray Mass Spectrometry". *J. Chromatogr., A*. **1998**, 804, 105–111.
6. Sarzanini, Corrado; Bruzzoniti, Maria Concetta; Mentasti, Edoardo. "Preconcentration and Separation of Halocetic Acids by Ion Chromatography". *J. Chromatogr., A*. **1999**, 850, 197–211.

7. Wicks, Richard J.; Moran, Mary Ann; Pittman, Laura J.; Hodson, Robert E. "Carbohydrate Signatures of Aquatic Macrophytes and Their Dissolved Degradation Products as Determined by a Sensitive High-Performance Ion Chromatography Method". *Applied and Environmental Microbiology*. Nov. 1991, 3135–3143.
8. Shotyk, William; Immenhauser-Potthast, Ina; Vogel, Hubert A. "Determination of Nitrate, Phosphate and Organically Bound Phosphorus in Coral Skeletons by Ion Chromatography". *J. Chromatogr., A*. **1995**, 706, 209–213.
9. Adam, Luke C.; Gordon, Gilbert. "Direct and Sequential Potentiometric Determination of Hypochlorite, Chlorite and Chlorate Ions when Hypochlorite Ion is Present in Large Excess". *Anal. Chem.* **1995**, 67, 535–540.
10. Joyce, Robert J.; Dhillon, Harpreet S. "Trace Level Determination of Bromate in Ozonated Drinking Water Using Ion Chromatography". *J. Chromatogr., A*. **1994**, 671, 165–171.
11. Stratford, Michael R. L.; Dennis, Madeleine; Cochrane, Ralph; Parkins, Charles S.; Everett, Steven A. "The Role of Nitric Oxide in Cancer, Improved Methods for Measurement of Nitrite and Nitrate by High-Performance Ion Chromatography". *J. Chromatogr., A*. **1997**, 770, 151–155.
12. Weinberg, Howard. "Pre-Concentration Techniques for Bromate Analysis in Ozonated Waters". *J. Chromatogr., A*. **1994**, 671, 141–149.

Dionex OnGuard II Ba

1. Dionex Corporation (Now Part of Thermo Fisher Scientific). Determination of Perchlorate in Drinking Water by Ion Chromatography. Application Update 145; LPN 1540. Sunnyvale, CA, 2004.
2. Dionex Corporation (Now Part of Thermo Fisher Scientific). Determination of Perchlorate in Vegetation Samples Using Accelerated Solvent Extraction. Application Note 356; LPN 1830. Sunnyvale, CA, 2006.
3. Accelerated Solvent Extraction and Ion Chromatography. Application Note 356; LPN 1830. Sunnyvale, CA, 2006.
4. Dionex Corporation (Now Part of Thermo Fisher Scientific). Application Note 101, Trace Level Determination of Bromate in Ozonated Drinking Water Using Ion Chromatography.
5. Slingsby, R. W.; Pohl, C. A. "Approaches to Sample Preparation for IC-Sulfate Precipitation on Barium-Form Ion Exchangers". *J. Chromatogr., A*. **1996**, 739, 49–55.
6. Charles, L.; Pepin, D. "Analysis of Oxyhalides in Water by Ion Chromatography-Ionspray Mass Spectrometry". *J. Chromatogr., A*. **1998**, 804, 105–111.

Dionex OnGuard II H

1. Dionex Corporation (Now Part of Thermo Fisher Scientific). Determination of Chelating Agents in Drinking Water and Wastewater Samples. Application Note 268; LPN 2664. Sunnyvale, CA, 2011
2. Dionex Corporation (Now Part of Thermo Fisher Scientific). Application Note 101, Trace Level Determination of Bromate in Ozonated Drinking Water Using Ion Chromatography.
3. Slingsby, R. W.; Pohl, C. A. "Approaches to Sample Preparation for IC-Sulfate Precipitation on Barium-Form Ion Exchangers". *J. Chromatogr., A*. **1996**, 739, 49–55.
4. Charles, L.; Pepin, D. "Analysis of Oxyhalides in Water by Ion Chromatography-Ionspray Mass Spectrometry". *J. Chromatogr., A*. **1998**, 804, 105–111.
5. Kuo, Cing-Yuan. "Improved Application of Ion Chromatographic Determination of Carboxylic Acid in Ozonated Drinking Water". *J. Chromatogr., A*. **1998**, 804, 265–272.
6. Chen, Sau Soon; Spiro, Michael. "Rose-hip Tea: Equilibrium and Kinetic Study of Mineral Ion Extraction". *Food Chemistry* **1993**, 48, 47–50.
7. Grinnell, Brian W.; Hermann, Robert B.; Yan, S. Betty. "Human Protein C Inhibits Selectin-Mediated Cell Adhesion: Role of Unique Fucosylated Oligosaccharide". *Glycobiology* **1994**, 4 (2), 221–115

Dionex OnGuard II M

1. Slingsby, R. W. "Sample Preparation for Ion Chromatography—New Developments, Considerations and Options". International Ion Chromatography Symposium 1999, Paper #33.
2. Slingsby, R. W. "When It's Not Just "Dilute and Shoot": Developments in Sample Preparation for Ion Chromatography". Pittsburgh Conference 2000, Paper #145.

Dionex OnGuard II P

1. Dionex Corporation (Now Part of Thermo Fisher Scientific). Determination of Hexavalent Chromium in Dyes. Application Brief 107; LPN 2323. Sunnyvale, CA, 2009.
2. Shotyk, William. "Ion Chromatography of Organic-Rich Natural Waters from Peat Lands, I". *J. Chromatogr.* **1993**, 640, 309–322.
3. Van Nifterik, L.; Xu, J.; Laurent, J. L.; Mathieu, J. "Analysis of Cellulose and Kraft Pulp Ozonolysis Products by Anion-Exchange Chromatography with Pulsed Amperometric Detection". *J. Chromatogr.* **1993**, 640, 335–343.

4. Steinmann, P.; Shoty, William. "Ion Chromatography of Organic-Rich Natural Waters from Peat Lands, III. Improvements for Measuring Anions and Cations". *J. Chromatogr., A*. **1995**, 706, 281–286.
5. Steinmann, P.; Shoty, William. "Ion Chromatography of Organic-Rich Natural Waters from Peat Lands, IV. Dissolved Free Sulfide and Acid-Volatile Sulfur". *J. Chromatogr., A*. **1995**, 706, 287–292.
6. Steinmann, P.; Shoty, William. "Ion Chromatography of Organic-Rich Natural Waters from Peat Lands, V. Fe²⁺ and Fe³⁺". *J. Chromatogr., A*. **1995**, 706, 293–299.
7. Wicks, Richard J.; Moran, Mary Ann; Pittman, Laura J.; Hodson, Robert E. "Carbohydrate Signatures of Aquatic Macrophytes and Their Dissolved Degradation Products as Determined by a Sensitive High-Performance Ion Chromatography Method". *Appl. and Environ. Microbiol.* Nov. 1991, 3135–3143.
8. Blank, Robert R.; Allen, Fay; Young, James A. "Extractable Anions in Soils Following Wildfire in a Sagebrush-Grass Community". *Soil Sci. Soc. Am. J.* **1994**, 58, 564–570.
7. Chen, Jian. "Determination of Organic Acids in Industrial Streams by Ion Chromatography After Solid-Phase Extraction". *J. Chromatogr., A*. **1996**, 739, 273–280.
8. Kadnar, Rainer. "Determination of Alkali and Alkaline Earth Metals in Oilfield Waters by Ion Chromatography". *J. Chromatogr., A*. **1998**, 804, 217–221.
9. Wicks, Richard J.; Moran, Mary Ann; Pittman, Laura J.; Hodson, Robert E. "Carbohydrate Signatures of Aquatic Macrophytes and Their Dissolved Degradation Products as Determined by a Sensitive High-Performance Ion Chromatography Method". *Appl. and Environ. Microbiol.* Nov. 1991, 3135–3143.

Dionex OnGuard II RP

1. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Determination of Perchlorate in Vegetation Samples Using Accelerated Solvent Extraction and Ion Chromatography. Application Note 356; LPN 1830. Sunnyvale, CA, 2006.
2. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Higher Resolution Separation of Organic Acids and Common Inorganic Anions in Wine. Application Notes 273; LPN 2727. Sunnyvale, CA, 2011.
3. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Determination of Inorganic Counterions in Pharmaceutical Drugs Using Capillary IC. Application Brief 136; LPN 2917. Sunnyvale, CA, 2011.
4. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Application Note 140, Determination of Iodide in Urine.
5. Dionex Corporation (Now Part of Thermo Fisher Scientific) . Application Note 37, Determination of Iodide in Milk Products.
6. Kadnar, P.; Rieder, J. "Determination of Anions in Oilfield Waters by Ion Chromatography", *J. Chromatogr., A*. **1995**, 706, 301–305.

Dionex OnGuard II Na

Dionex Corporation (Now Part of Thermo Fisher Scientific) . HPAE-PAD Determination of Infant Formula Sialic Acids. Application Notes 253; LPN 2561. Sunnyvale, CA, 2011.

Dionex OnGuard II Ag/H

Dionex Corporation (Now Part of Thermo Fisher Scientific) . Determination of Phytic Acid in Soybeans and Black Sesame Seeds. Application Note 295; LPN 3028. Sunnyvale, CA, 2012.

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© 2015-2017, 2020-2022 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

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