

Upgrade Your Maps


Fast, confident and more reliable peptide mapping.



Smarter Navigation...

Peptide mapping is a core analytic in biotherapeutic development. Yet it is slow, involves painstaking sample preparation, and is prone to errors.

Let us guide you through the rough terrain of peptide mapping and upgrade you with a solution that's better, wherever on the map you are.



We're ready to upgrade your peptide mapping capabilities. Thermo Scientific[™] SMART Digest[™] kits give you more confident and reproducible digestion. Thermo Scientific[™] Vanguish[™] Flex UHPLC speeds up your mapping, whilst giving you better retention time stability. Our range of sub-2µm and solid core particle LC columns including the Thermo Scientific[™] Acclaim[™] and Accucore[™] ranges ensure fast, reliable separations. Thermo Scientific[™] BioPharma Finder[™] software is easier to use, and ensures you get the greater coverage you demand. Powerful, innovative new solutions for simply better peptide mapping.

...then let us be your guide to simply better mapping.


- Do you need greater reproducibility?
 - Do you want faster peptide navigation?
 - How about greater sequence coverage?
 - A simple, easy to use solution?
 - More **confidence** in your peptide assignation?



bigestion

Sample preparation for peptide mapping, such as by tryptic digestion, is laborious, difficult, and prone to error. Not only that, it's not always consistent - different

analysts can produce different results. You need confidence in your digestion, that's where SMART Digest Kits come in.




Easier to Use


SMART Digest Kits are an innovative way of rapidly and reproducibly digesting proteins into their constituent peptides, ready for characterization. The process is simplicity itself: Firstly add the protein and a buffer to the SMART Digest tube. Then incubate and shake, typically for about an hour. You can then optionally centrifuge or filter, and you're done!

Proteins are heat-denatured for digestion, so denaturing agents and reduction/alkylation steps are not needed. It's a high throughput technique that can be readily automated. The kit comes as 12 x 8 stripe digest tubes in PCR format or in bulk format for use with 96 well plates. Optionally a Thermo Scientific[™] SOLAµ[™] Solid Phase Extraction (SPE) Plate can be included for use with more complex plasma-based sample sets. Reproducibility is typically significantly better than that of traditional tryptic protein digestion and often with increased sensitivity.

SMART Digest Kits







" Amazingly stable David Goodsell, Protein Data Bank

Better Reproducibility

Overlay of 13 consecutive chromatographic runs of a peptide sample separated on an analytical Thermo Scientific[™] Acclaim[™] 120 C18 column (RP18, 2.2 µm, 2.1 x 250 mm) and prepared from a mAb, digested with the SMART Digest kit.

Better Sensitivity

Measurement of serum Thyroglobulin after tryptic digestion of serum samples.

- SMART Digest: 25% plasma, 3.5 h digestion
- In-solution digest: 20% plasma, R/A, 4 + 16 h digestion

For in-solution protocol, see for example: Clarke et al. (2012), J. Investigative Medicine, 60(8)

Separation

The characterization of biotherapeutics places rigorous demands on liquid chromatography systems. Biocompatibility, reproducibility and high performance are necessary. However, it's also important for an LC system to offer flexibility and ease of use. Vanquish Flex UHPLC has all of this and more.

What about the Choice of Column?


For perfect peptide mapping, the LC column needs to produce excellent resolution and reliable reproducibility. The Acclaim RSLC 120 C18 column is ideal for the job, and is our recommended LC column of choice for peptide mapping applications.

Faster Mapping

Total ion chromatograms obtained from peptide mapping experiments of rituximab applying gradient lengths from 30 to 5 minutes.

More Confidence

Faster Mapping



Vanquish Flex UHPLC

Easy Transfer from LC-MS to LC-UV


Peptide mapping is typically done using mass spectrometric detection. However, for certain batch testing and QA/QC applications, analysts are increasingly turning to LC-UV, a simpler and faster technique. This relies on LC performance that is high enough to provide data of sufficient quality.

In the data shown right, mass spectrometric data is compared to LC-UV data. The Vanquish Flex has a highly stable flow delivery as well as sample pre-compression, and this results in excellent retention time precision and increased confidence in the peak assignment. The data comparison clearly shows good confirmation between the LC-MS data (upper trace) and the LC-UV data (lower trace). Vanquish Flex is therefore extremely well suited to being used in a routine LC-UV peptide mapping in a quality control environment.



Overlaid chromatograms of the total ion current (TIC, upper trace), measured using LC-MS and the UV trace at 214 nm (measured using Vanquish Flex UHPLC, lower trace) of a SMART Digest Kit digested rituximab sample.







- Extremely high pressure flow path to ensure high peak capacity and peak area precision
- Easy integration with MS
- Multiple detector options


^d Mass Informatics

The collation, interpretation and deconvolution of peptide mapping mass spectrometry data is as important as the analysis. You need a software package that allows the identification of all peptides, gives confirmation of the amino acid sequence, identifies the site and type of any modifications and maps the disulfide bonds. Welcome to BioPharma Finder.

What about Mass Spectrometry?

Thermo Scientific[™] Q Exactive[™] Hybrid Quadrupole-Orbitrap[™] mass spectrometers are ideal for peptide mapping applications where maximum characterization is required and full sequence information needed. It is extremely accurate, typically identifying component compounds to within plus minus 3ppm. Furthermore, it is linear across four and a half orders of magnitude, proving ideal for complex and varied peptide mixtures.





MS/MS spectra of the glyco-peptide aa 290-302 (TKPREEQFN*STFR, *=AG20), as measured using Vanquish Flex UHPLC - Q Exactive HF MS and detected with BioPharma Finder.

BioPharma Finder is our new

integrated software for both peptide mapping and intact analysis. It's easy to use, with an interactive wizard for editing process methods and allowing for easy visualization of results. It also has an automated protein coverage map with chromatographic peak shading.


The user can readily produce a mapping modification summary. This clearly describes any detected modifications and gives details about type, peptide, sequence, confidence value and percent abundance. This information is ideal for the analyst when determining any modifications to the peptide sequence.

BioPharma Finder Software

100% coverage peptide mapping and intact protein characterization in one package

The ability to perform Sequence Variant Searching, disulfide link analysis, de novo sequencing with quantitation of post-translational modifications assures you wont miss anything. Multiple interactive plot display options mean you can see what you want, when you need it.






BioPharma Finder is also an excellent tool for batchto-batch comparability, essential in routine peptide mapping applications. Interactive visualization tools such as chromatographic shading afford increased efficiency in data mining and interpretation.



Minism MS Signd = 100000 Das File = Speed_20150731_Cal20P3SPE_NISTspaked_HCD_Top20_125_4al_03 raw

Data File = Spaceu_aver Protesse = Trypsin	Number of MS Peaks MS Peak Area		Semence Coverage	Abundance (mol
	Number of MS Peaks	MS Peak Area	100.0%	45.53%
Proteins	464	26.7%	100.000	54,47%
1: Light Chain	1149	57.8%	100.0%	J4.4774
2: Heavy Chain Heidenfiled	4122	15.5%		

MS Signal Threshold = 6.3e+001





Integrated workflows for intact protein

- to peptide mapping analysis
- Interactive visualization tools
- Superior deconvolution and peptide
- identification algorithms

Smarter Navigation

Greater

Coverage

It's time for an upgrade. The peptide mapping workflow is ready for new technologies, new solutions, and better answers. Let us guide you to a world of better, more confident, and more reproducible peptide mapping.

> This workflow has made my mapping much easier, and frees my lab staff for

more valuable tasks





Ease of Use SMART Digest Kits

Reproducibility Acclaim RSLC 120 C18



Bf

Speed Vanquish Flex UHPLC

Confidence Q Exactive Plus Mass Spectrometer

Coverage **BioPharma Finder Software**





variants gives me the confidence "



Find out about technologies that can upgrade your maps **www.thermofisher.com/peptidemaps**

f facebook.com/ChromatographySolutions У twitter.com/ChromSolutions ♠ analyteguru.com **Stay Ahead!** Scan the QR code or follow the URL to watch the latest webinars on Peptide Mapping.



thermofisher.com

©2016 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

