

Thermo Scientific TC/EA



*High Temperature Conversion  
Elemental Analyzer*

## Thermo Scientific TC/EA High Temperature Conversion Elemental Analyzer

For many years the on-line isotope ratio analysis of organic bulk samples for  $^{13}\text{C}/^{12}\text{C}$ ,  $^{15}\text{N}/^{14}\text{N}$  and  $^{34}\text{S}/^{32}\text{S}$  has been performed rapidly, easily and precisely using a conventional Elemental Analyzer – continuous flow isotope ratio mass spectrometers (IRMS). Measurement of  $^{18}\text{O}/^{16}\text{O}$  and D/H ratios in organic and inorganic matter had been restricted to off-line sample preparation in the past. The Thermo Scientific high temperature conversion elemental analyzer (TC/EA) is the technological breakthrough offering the benefits of continuous flow IRMS to on-line oxygen and hydrogen isotope ratio analysis of solid and liquid bulk samples.

### Direct Analysis of $^{18}\text{O}/^{16}\text{O}$ and D/H in Organic and Inorganic Samples

- Automated analysis with high sample throughput
- Only small sample amounts required
- Low blanks and negligible memory
- Oxygen isotope ratios from organic compounds, water and selected inorganic compounds
- Hydrogen isotope ratios from organic compounds and water
- Nitrogen and oxygen isotope ratios from nitrates



# Principle of Operation

Quantitative high temperature conversion, also referred to as pyrolysis, is a new technique in which oxygen present in a compound is converted to CO, and hydrogen contained in a compound is converted to H<sub>2</sub>. The process is rapid and quantitative in a reducing environment at high temperatures, typically exceeding 1400 °C.


The Thermo Scientific TC/EA preparation system and Thermo Scientific ConFlo IV universal interface are outlined in the schematic diagram. This combination facilitates the conversion reaction, separation of reaction gases, transfer into the IRMS and referencing against standard gases.

The reactor consists of a glassy carbon tube with glassy carbon filling, ensuring that neither sample nor reaction gases can get into contact with oxygen containing surfaces (e.g. Al<sub>2</sub>O<sub>3</sub>) while at high temperatures. Only this technology enables memory-free conversion reactions with no restrictions on compound type.


The reaction gases are separated in an isothermal gas chromatograph, which is also part of the TC/EA. The gases are admitted to the IRMS via the ConFlo IV universal interface. The ConFlo IV allows automatic sample gas dilution and generation of reference gas pulses, enabling individual referencing of each sample gas peak. The reference gas pulses can be automatically adjusted in signal intensity which also allows fully automated H<sup>3+</sup> factor determination at any time within a sample sequence.

- No fractionation
- Lowest memory
- No restrictions on organic samples
- Usable for selected inorganic samples
- Referencing of each peak using reference gas

## Schematic of the TC/EA, ConFlo IV and IRMS



## Principle of High Temperature Conversion



# Specifications and Installation Requirements

The Thermo Scientific TC/EA can be connected to any current Thermo Scientific Isotope Ratio MS equipped for continuous flow applications. If hydrogen isotope ratios are to be analyzed, the IRMS must be equipped with an energy filter to suppress  $^4\text{He}^+$  ions on the DH collector (Thermo Scientific MAT 253, DELTA V Plus, and DELTA V Advantage IRMS).

## Gases

High purity helium (99.999% or better).

Reference gases ( $\text{CO}$  and  $\text{H}_2$ ) with pressure regulators.

To use the TC/EA with  $\text{CO}$  and  $\text{H}_2$  reference gas, the laboratory must be equipped with  $\text{CO}$  and  $\text{H}_2$  detectors.

## Power

230 V, single phase, 8 A.

## Dimensions and Weight

45 x 70 x 50 cm (W x D x H), 59 kg.

### External Precision for Isotope Ratios, H, O, (n = 10), $\delta$ -Isotope

Benzoic acid at natural abundance for H and O

	$^{18}\text{O}/^{16}\text{O}$	D/H
50 $\mu\text{g}$ O, 25 $\mu\text{g}$ H	0.4 ‰	3 ‰

Water for H and O\*

	$^{18}\text{O}/^{16}\text{O}$	D/H
0.5 $\mu\text{L}$	0.2 ‰	2 ‰

\* requires liquid injector



# Applications

The Thermo Scientific TC/EA high temperature elemental analyzer can be used for simultaneous hydrogen and oxygen isotope ratio determination of all organic compounds. Due to its very high maximum operation temperature restrictions to individual sample classes do not exist for the analysis of hydrogen isotope ratios the IRMS must be capable of  $^4\text{He}^+$  suppression at the DH collector (Thermo Scientific MAT 253, DELTA V Plus, DELTA V Advantage). Selected inorganic compounds can also be analyzed, such as nitrates (N and O), phosphates (O) and sulfates (O).

One important application is the isotope ratio analysis of hydrogen and oxygen in water samples. For this application an autosampler for liquid samples can be added to the system, replacing the standard autosampler for solid samples. Typical sample amounts for water are  $< 0.5 \mu\text{L}$ . The analysis of smaller water samples depends on the quality of syringe transfer and injection. Both isotope ratios can be determined within one run, making this setup perfectly suitable for e.g. high throughput doubly-labeled water analysis or screening of water resources.

Typical measurement time is about 5 minutes per sample for oxygen and less than 3 minutes for hydrogen isotope ratio determination. A determination of both isotope ratios in one run takes about 6 minutes.

## Options


### Liquid Injection Kit

Isotope analysis of water and liquids by high temperature conversion requires rapid transfer of the injected liquid into the reaction zone. The  $\text{H}_2\text{O}$  kit includes a special insert for the reactor, a septum-equipped injection head, a microliter syringe, screw top vials and spare parts. A special autosampler for liquids, e.g. Thermo Scientific AS3000, is recommended achieving high throughput and repeatedly reproducible injection for highest precision.


### Bottom Feed Connector

The bottom feed connector (BFC) is the most recent improvement for high temperature conversion elemental analysis. By feeding the carrier gas flow from the bottom instead of the top, as in the conventional design, maintenance intervals become significantly less. Channeling of the complete flow through the glassy carbon tube reduces the effects of ash deposition which makes cleaning procedures occasional (less frequent). The BFC increases the analytical capacity of the high temperature conversion by at least a factor of five for solid samples and improves peak shapes. The design of the Thermo Scientific BFC allows the use of the original reactor material. No change of dimensions of the tubes is required.

### Flow Path with Conventional and Bottom Feed Connector



## One Injection - Two Isotope Ratios



### Isotope Ratios of Single Elements

O	→	CO	→	$\delta^{18}\text{O}$ , wt% O
H	→	$\text{H}_2$	→	$\delta\text{D}$ , wt% H

### Isotope Ratios from Two Elements

H + O	→	$\text{H}_2 + \text{CO}_2$	→	$\delta\text{D} + \delta^{18}\text{O}$ , wt% H + wt% O
-------	---	----------------------------	---	--

Solids: organic and inorganic samples

Water\*:  $\text{H}_2\text{O} + \text{C} \rightarrow \text{CO} + \text{H}_2$

\* requires liquid injector

## Laboratory Solutions Backed by Worldwide Service and Support

Tap our expertise throughout the life of your instrument. Thermo Scientific Services extends its support throughout our worldwide network of highly trained and certified engineers who are experts in laboratory technologies and applications. Put our team of experts to work for you in a range of disciplines – from system installation, training and technical support, to complete asset management and regulatory compliance consulting. Improve your productivity and lower the cost of instrument ownership through our product support services. Maximize uptime while eliminating the uncontrollable cost of unplanned maintenance and repairs. When it's time to enhance your system, we also offer certified parts and a range of accessories and consumables suited to your application.

To learn more about our products and comprehensive service offerings, visit us at [www.thermoscientific.com](http://www.thermoscientific.com).



In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

### **Africa-Other**

+27 11 570 1840 • [analyze.sa@thermo.com](mailto:analyze.sa@thermo.com)

### **Australia**

+61 3 9757 4300 • [analyze.au@thermo.com](mailto:analyze.au@thermo.com)

### **Austria**

+43 1 333 50 34 0 • [analyze.at@thermo.com](mailto:analyze.at@thermo.com)

### **Belgium**

+32 53 73 42 41 • [analyze.be@thermo.com](mailto:analyze.be@thermo.com)

### **Canada**

+1 800 530 8447 • [analyze.ca@thermo.com](mailto:analyze.ca@thermo.com)

### **China**

+86 10 8419 3588 • [analyze.cn@thermo.com](mailto:analyze.cn@thermo.com)

### **Denmark**

+45 70 23 62 60 • [analyze.dk@thermo.com](mailto:analyze.dk@thermo.com)

### **Europe - Other**

+43 1 333 50 34 0 • [analyze.emea@thermo.com](mailto:analyze.emea@thermo.com)

### **Finland / Norway / Sweden**

+46 8 556 468 00 • [analyze.se@thermo.com](mailto:analyze.se@thermo.com)

### **France**

+33 1 60 92 48 00 • [analyze.fr@thermo.com](mailto:analyze.fr@thermo.com)

### **Germany**

+49 6103 408 1014 • [analyze.de@thermo.com](mailto:analyze.de@thermo.com)

### **India**

+91 22 6742 9434 • [analyze.in@thermo.com](mailto:analyze.in@thermo.com)

### **Italy**

+39 02 950 591 • [analyze.it@thermo.com](mailto:analyze.it@thermo.com)

### **Japan**

+81 45 453 9100 • [analyze.jp@thermo.com](mailto:analyze.jp@thermo.com)

### **Latin America**

+1 561 688 8700 • [analyze.la@thermo.com](mailto:analyze.la@thermo.com)

### **Middle East**

+43 1 333 50 34 0 • [analyze.emea@thermo.com](mailto:analyze.emea@thermo.com)

### **Netherlands**

+43 1 333 50 34 0 • [analyze.nl@thermo.com](mailto:analyze.nl@thermo.com)

### **New Zealand**

+64 9 980 6700 • [analyze.au@thermo.com](mailto:analyze.au@thermo.com)

### **Russia/CIS**

+27 11 570 1840 • [analyze.emea@thermo.com](mailto:analyze.emea@thermo.com)

### **South Africa**

+27 11 570 1840 • [analyze.sa@thermo.com](mailto:analyze.sa@thermo.com)

### **Spain**

+34 914 845 965 • [analyze.es@thermo.com](mailto:analyze.es@thermo.com)

### **Switzerland**

+41 61 716 77 00 • [analyze.ch@thermo.com](mailto:analyze.ch@thermo.com)

### **UK**

+44 1442 233555 • [analyze.uk@thermo.com](mailto:analyze.uk@thermo.com)

### **USA**

+1 800 532 4752 • [analyze.us@thermo.com](mailto:analyze.us@thermo.com)

[www.thermoscientific.com](http://www.thermoscientific.com)

Thermo Fisher Scientific (Bremen) GmbH  
Management System Registered to ISO 9001:2008

©2004, 2008, 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

BR30034\_E 09/10C

**Thermo**  
SCIENTIFIC