Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





The European Water Framework Directive

Thermo Scientific Environmental Solutions Reference Guide

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Introduction to EU WFD 2013/39/EU

Introduction

- The EU Water Framework Directive (WFD) was introduced in August 2013, amending the EU directives 2000/60/EC, 2008/105/EC, 2000/60/EC, and 2008/105/EC.
- The European Water Framework Directive is a directive that commits all member states to actively control for a long list of environmental contaminants in all water bodes of the various member states.
- It lays down a strategy against water pollution, to be applied in all European Union member states.
- It involves identifying priority substances and monitoring different classes of contaminants; and it includes the first watch list (used for future prioritization exercises).
- Member States have the flexibility to apply an EQS (Environmental Quality Standard) for an alternative matrix or, where relevant, an alternative biota taxon, for example sub-phylum Crustacea, paraphylum "fish", class Cephalopoda or class Bivalvia (mussels and clams).
- Compound levels are expressed as Environmental Quality Standard or (EQS) and the annual average (AA) and the Maximum Allowable Concentration (MAC) is given.
- The EU WFD encourages the development of novel monitoring methods such as passive sampling and other tools



Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis
- Reference List





UK Specific Regulations: Introduction to CIP 2

- The UK Chemical Investigations Programme (CIP 2) looks into the occurance, sources, and removal of trace substances in waste water treatment facility effluent.
- This regulation helps to establish priorities for remediative action to ensure surface waters meet new Environmental Quality Standards (EQS).
- The CIP 1 program was managed by UK Water Industry Research (UKWIR) and implemented from 2010-2013.
- The CIP 2 program is a follow-up program of sampling and analysis to be implemented between 2014 and 2020.
- The primary objective of CIP 2 is to identify and characterize sites where EQS levels are breached.
- In the program, 70 priority substances were determined from 162 sewage treatment works (STW) effluents.
- Additionally, 11 pharmaceuticals were identified as priority monitoring candidates. It is important to note that EQS is defined for only three pharmaceutical compounds: Diclofenac, E2, and EE2. All substances selected for monitoring analysis were detected previously in waste water effluent samples.
- The determined environmental concentrations of many priority substances in effluent exceeded EQS.

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Chapter 1: Index of Compounds and Techniques

Priority hazardous substances	GC-MS	LC-MS	ICP-MS	Priority hazardous substances	GC-MS	LC-MS	ICP-MS	Priority hazardous substances	GC-MS	LC-MS	ICP-MS
Anthracene	\checkmark			Perfluorooctane sulfonic acid and its derivatives (PFOS)				Naphthalene			
Brominated diphenylethers	\checkmark			Quinoxyfen	\checkmark			Nickel and its compounds			
Cadmium and its compounds				Dioxins and Dioxin-like compounds	\checkmark			Octylphenols	\checkmark		
C ₁₀₋₁₃ Chloroalkanes	\checkmark			Hexabromocyclododecane (HBCDD)				Pentachlorophenol	\checkmark		
Di-(2-ethylhexyl)phthalate (DEHP)	\checkmark			Heptachlor and heptachlor epoxide	\checkmark			Simazine	\checkmark	\checkmark	
Endosulfan	\checkmark			Alachlor	\checkmark			Trichlorobenzenes	\checkmark		
Hexachlorobenzene (HCB)	\checkmark			Atrazine	\checkmark			Trichloromethane (chloroform)			
Hexachlorobutadiene (HCBD)	\checkmark			Benzene	\checkmark			Aclonifen			
Hexachlorocyclohexane	\checkmark			Chlorfenvinphos	\checkmark			Bifenox			
Mercury and its compounds				Chlorpyrifos (ethyl)	\checkmark			Cybutryne	\checkmark		
Nonylphenols	\checkmark			1,2-dichloroethane	\checkmark			Cypermethrin	\checkmark		
Pentachlorobenzene	\checkmark			Dichloromethane	\checkmark			Dichlorvos	\checkmark		
Polyaromatic hydrocarbons (PAHs)	\checkmark			Diuron				Terbutryn	\checkmark	\checkmark	
Tributyltin compounds	\checkmark			Fluoranthene				Diclofenac (watchlist)		\checkmark	
Trifluralin				Isoproturon				17-beta-estradiol (watchlist)			
Dicofol	\checkmark			Lead and its compounds				17-alpha-ethinylestradiol (watch list)	\checkmark	\checkmark	

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids

ICP-MS compounds


- Highlights & configuration
- Metal analysis


Reference List



Chapter 2: Sample Preparation Techniques

General sample preparation workflow





Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Automated Sample Preparation Solutions

The Thermo Scientific[™] Dionex[™] AutoTrace[™] 280 Solid-Phase Extraction instrument saves time, solvent, and labor ensuring high reproducibility and productivity for analytical laboratories. The instrument can process up to six samples in 2–3 hours with only 15 minutes of operator involvement. The Dionex AutoTrace 280 instrument uses powerful pumps (no check valves) and proven constant-flow technology to efficiently process even the most difficult samples.

Save time, save solvent

With Dionex AutoTrace and Thermo Scientific[™] Dionex[™] ASE[™] Accelerated Solvent Extraction systems, laboratories can effectively automate the solvent extraction process for liquid and solid matrices.



Learn more about the accelarared solvent extraction systems

Download the accelerated solvent extraction environmental summary application notebook

Introduction to WFD

- Introduction to UK CIP 2 List of compounds and analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Solid Phase Extraction (SPE)

LC-MS amenable compounds

Loading	 Material: Thermo Scientific[™] HyperSep[™] Retain PEP cartridge pH 2
Wash	• 10% MeOH/H ₂ O
Elution	 ACN (0.1%FA) – ACN – ACN (0.1% NH₄OH) EtOAc (PFCs, HBCDD, Phenols, Triclosan, Triazoles)
Evaporation	• Final Volume (0.5 mL)
Reconstitution	Addition of 4.5 mL of 0.1% FA
Analysis	 10 μL direct injection 1-5 mL Injection Thermo Scientific[™] EQuan Max[™] LC-MS system

Introduction to WFD

Introduction to UK CIP 2 List of compounds and

analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids

ICP-MS compounds

- Highlights & configuration
- Metal analysis

Reference List



Liquid Liquid Extraction (LLE)

GC-MS amenable compounds

Classical off-line sample preparation: Liquid-liquid extraction or LLE

To 100 mL of surface water, 10 g of NaCl was added. The extraction solvent was 10 mL dichloromethane with 1 % iso-octane. The sample was shaken vigorously for 20 minutes. The organic phase was removed and dried with $NaSO_4$.

The extract was evaporated to 1 mL under a low flow of nitrogen and transferred to a 2 mL (12x32 mm) GC vial for injection.

No extra clean up was performed.

Iso-octane acted as a keeper to retain the most volatile compounds, such as trichlorobenzene and naphthalene. The starting temperature of the GC was increased to 65 °C, saving cooling time and decreasing the run-to-run time.



Introduction to WFD

Introduction to UK CIP 2 List of compounds and

- analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS


GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Solid Phase Extraction (SPE) Procedure

GC-MS amenable pesticides

EPA method 508 – Analysis of chlorinated pesticides, herbicides and organohalides

Using 10 g 75 mL Thermo Scientific™ HyperSep™ C18 SPE cartridge (P/N 60108-703)

Sample preparation

1 L water collected Add MgCl₂ (final conc. 10 mL/L)

Condition HyperSep C18 SPE cartridge

1 x 1 mL 1:1 EtAC/CH $_3$ Cl $_2$ 1 x 10 mL CH $_3$ OH 1 x 10 mL H $_2$ O

Apply sample

Add 1 x 5 mL CH_3OH to sample Mix Take 50 µL of sample and mix

Load sample at 1 to 2 mL/minute

Elute

Insert fresh collection tubes into manifold 1 x 10 mL EtAc 1 x 10 mL CH₃Cl₂ 1 x 3 mL EtAc/CH₃Cl₂

Evaporation

Evaporate cluates to 0.8 mL under a gentle stream of nitrogen in a heated water bath 40 °C Add internal standard Adjust volume to 1 mL

Analysis

Add 1 to 2 µL onto GC

Recommended GC column	Part number
Thermo Scientific™ TraceGOLD™ TG-OCP I 30 m x 0.25 m x 0.25 µm	26078-1420



Introduction to WFD

- Introduction to UK CIP 2 List of compounds and analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Sample preparation

1 x 3 mL CH₂Cl

1 x 3 mL CH OH

Apply sample

1 x 3 mL 0.05N HCL

Collect 1 L of H_oO

Adjust pH to 2 with 6N HCI

Condition HyperSep Retain PEP SPE cartridge

NOTE: *Do not allow the cartridge to dry out*

Load 1 L of water sample at 20 mL/minute

Dry column for 10 to 15 minutes

SPE Procedure: Phenols by GC-MS

GC-MS determination of phenols in drinking water

Using 500 mg 6 mL HyperSep Retain PEP SPE cartridge (P/N 60107-206)

Wash column

1 x 10 mL H₂0 Dry column (1 minute at > 10 °Hg)

Elute

1 x 10 mL CH CI x 3 mL CH CI CH_{3} CI Concentrate the extract to 0.9 mL in water bath (40 °C) under a gentle stream of nitrogen

Analysis

Adjust final volume to 1.0 mL with CH_Cl Analyze the extract with using GC-MS $^{\rm S}$

Recommended GC column	Part number
TraceGOLD TG 5MS 30 m x 0.25 m x 0.25 µm	26098-1420



Introduction to WFD

- Introduction to UK CIP 2 List of compounds and analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis
- Reference List



Automated Solvent Extraction (ASE) Procedure

Dioxins in biota by GC-MS

Environmental contaminants in fish and egg samples

Equipment

GC-FCD

HPI C

Dionex ASE 200 accelerated solvent extractor* with ASE solvent controller Choose either 11 mL stainless steel extraction cells (P/N 049560) or 22 mL stainless steel extraction cells (P/N 049561) or 33 mL stainless steel extraction cells (P/N 049562) Cellulose filters (P/N 049458) Collection vials, 40 mL (P/N 048783) or 60 mL (P/N 048784) Dionex SE 500 solvent evaporation system (P/N 063221) Analytical balance (to read to nearest 0.0001 g or better) Tissue homogenizer (Buchi B-400 or equivalent) Freeze drier (for PCB extraction) Centrifuge (for organotin extraction) Mechanical shaker (for organotin extraction) Dionex ASE 150 and 350 systems can be used for equivalent results Analysis GC GC-MS

Pressure:	1500 psi
Temperature:	175 °C
Solvent:	100% Toluene
Static time:	10 min
Static cycles:	2
Flush volume:	60%
Purge time:	60 sec
Static time:	1 or 2 min**
Flush volume:	60%
Purge time:	60 sec
Cycles:	2
Total time:	12 min
Total solvent:	20 mL

Extraction conditions

* Petroleum ether and hexane were found to be equivalent as extraction solvents for fat in meat.

**When extracting more than 1 g of a high-fat sample, a 2 min static time may be beneficial.

	ľ	1	E			n	n		D	
S	С	I	Е	Ν	Т	I	F	I	С	

Download the accelerated solvent extraction environmental summary application notebook

Results

PCDDs/PCDFs in fish tissue samples (ng/kg or ppt)using ASE.

Compound	Soxhlet	ASE	Certified
2,3,7,8-TCDD	7.6	7.6	6.6
1,2,3,4,8-PCDD	4.3	4.3	4.4
1,2,3,4,7,8-HCDD	1.4	1.4	1.9
2,3,4,7,8-TCDF	13.4	12.6	11.9
1,2,3,7,8-PCDF	5.4	5.1	5.0
1,2,3,4,7,8-HCDF	12.5	12.2	12.2
OCDD	12.4	6.4	6.3
Total TEQ	21.4	21.1	21.0

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Sample Preparation

ICP-MS amenable components

Aqueous samples

Aqueous environmental samples that contain less than 0.2% (m/v) dissolved solids at the point at which the sample enters the plasma can be measured directly with ICP-MS, following filtration of the samples if they contain suspended particles that could block the nebulizer of the instrument. For samples containing more than 0.2% (m/v) dissolved solids at the point of entry to the plasma, dilution is required either of the sample itself, by on-line or off-line liquid dilution, or by dilution of the nebulized sample aerosol using an additional flow of argon gas. Liquid dilution of samples is most commonly achieved using diluted nitric acid (typically 1 to 2% (v/v)), with addition of Au (at 0.2 to 1 ppm) to stabilize Hg to minimize memory effects with this element.

Solids and slurry samples

These samples generally require digestion with concentrated nitric acid followed by dilution with water and then filtration to remove undigested material. Digestion can be achieved using hot plate or microwave-based methods.

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids

ICP-MS compounds

- Highlights & configuration
- Metal analysis

Reference List





Chapter 3: GC-MS Compounds

GC-MS amenable compounds

The section below describes the GC amenable contaminants for the surface water directive, containing the following compound classes:

- Volatile organic compounds analyzed with Headspace GC-MS in SIM mode
- Polybrominated diphenyl ethers using GC-MS/MS
- Organotin with derivatization with GC-MS/MS
- Chlorinated alkanes using GC-MS/MS
- Dioxins and PCB using GC-MS/MS
- The remaining pesticides, PAHs, phthalate, and phenol compounds with GC-MS/MS

The groups of compounds above are divided according to analytical workflows and required detection limits. Measurement in SIM mode requires monitoring three ions for each compound for which the ion ratio deviation is monitored. Measurements in MS/MS mode require two transitions per compound for which the ion ratio will be monitored according to the European Directive 2002/657/EC.

GC-MS Highlights

Table of Contents

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



The Thermo Scientific ${}^{\scriptscriptstyle\rm TM}$ TriPlus ${}^{\scriptscriptstyle\rm TM}$ RSH autosampler offers:

- Automated sample preparation capabilities
- Unattended switching from headspace to liquid injection to SPME
- Support for dual GC injection configuration

The Thermo Scientific ${}^{\scriptscriptstyle\rm TM}$ TRACE ${}^{\scriptscriptstyle\rm TM}$ 1300 series GC offers:

- Instant connect modularity for easy maintenance, removal, and module replacement
- Fast and reliable heating and cooling of the oven for rapid analysis

The Thermo Scientific ${}^{\scriptscriptstyle \rm T}$ TSQ ${}^{\scriptscriptstyle \rm T}$ series triple quadrupole and ISQ ${}^{\scriptscriptstyle \rm T}$ single quad series offers:

- The Thermo Scientific[™] Extractabrite[™] removable ion source, for removal without venting the mass spec
- Enhanced Velocity Optics (EVO), allowing the optimal number of transition scans without comprimising data quality for large sample compound lists.
- User-friendly software, providing tools for automated method development and integration of compound names, retention times, and transitions between the instrument and processing methods







Introduction to WFD

- Introduction to UK CIP 2
- List of compounds and
- analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Compound group	Autosampler	TRACE GC 1300/1310	MS analyzer (acquisition mode)
VOC	Headspace	SSL or PTV	ISQ or TSQ Series (EI-Timed SIM)
BDE	LV option	SSL or PTV	TSQ Series (EI-SRM or NCI-SIM)
Organotin compounds	Standard	SSL or PTV	TSQ Series (EI-SRM)
Polychlorinated alkanes	Standard	SSL or PTV	TSQ Series (EI-SRM)
Dioxins/Dioxin-like compounds	Standard	SSL or PTV	TSQ Series (EI-SRM) Thermal Scientific™ DFS™ High Resolution GC/MS
Pesticides	Standard	SSL or PTV	TSQ Series (EI-SRM)
PAHs	Standard	SSL or PTV	TSQ Series (EI-SRM)
Phenols	Standard	SSL or PTV	TSQ Series (EI-SRM)
Phthalates	Standard	SSL or PTV	TSQ Series (EI-SRM)

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- **ICP-MS** compounds
- Highlights & configuration
- Metal analysis
- Reference List



Volatile Organic Compounds (HS) GC-MS

Sample preparation

- 10 mL of surface water
- 2 g NaCl

Headspace conditions

60 °C, 20 min., 1.5 mL injection; 1/40 split

Column


TraceGOLD TG-VMS 20 m x 0.18 mm ID x 1.00 µm (P/N 26080-4950)


Instrument methodology


GC method: 30 °C (4 min); 18 °C/min; 100 °C (0 min); 40 °C/min;

230 °C (3 min). Injector: 200 °C ISQ LT MS: Timed SIM mode

Source temperature: 250 °C







Chromatogram of volatile compounds monitored in Timed SIM mode with monochloro ethene eluting first.



Learn more about analyzing Volatile Organic Compounds with GC-MS

Download the poster presentation on GC-MS amenable EU WFD compounds

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- **ICP-MS** compounds
- Highlights & configuration
- Metal analysis
- Reference List





Polybrominated Diphenylethers by GC-MS/MS

Methodology

GC method

120 °C (1 min); 20 °C/min; 320 °C (5 min)

PTV method


80 °C (0.1 min); 10 °C/sec; 300 °C (6 min); 14.5 °C/sec; 340 °C (20 min); 2 min splitless

TSQ Quantum Ultra system Timed SRM mode and/or Cineg with NH₄


Source temperature 260 °C Column

TraceGOLD TG-5HT column (5% phenyl film) 15 m x 0.25 mm x 0.10 µm (P/N 26095-0350)

Liner Siltek Baffled (P/N 453T21210)



Compound	IDL LOQ in pg in El SRM	IDL LOQ in pg in NCI SIM
BDE28	<0.2	<0.1
BDE47	<0.2	<0.1
BDE99	<0.5	<0.1
BDE100	<0.5	<0.1
BDE154	<0.5	<0.1
BDF153	< 0.5	<0.1



Learn more about PBDE analysis with GC-MS/MS

Download the poster presentation on GC-MS amenable EU WFD compounds

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS


GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Transfer line 300 °C

PTV method

Column

Temperature program: 45 °C (Hold 2 min)

Ramp 55 °C/min – 175 °C – Ramp 35 °C – 300 °C (Hold 2 min);

Injector Temperature: 50 °C – Spitless injection 0.1 min

PTV Clean: 350 °C - 11 min - Clean flow 50 mL/min

TSQ Quantum system in El SRM Source Temp. 250 °C

PTV Transfer: 8 °C/sec - 280 °C (Hold 1 min)

TraceGOLD TG-5HT column (5% phenyl film) of

30 m x 0.25 m x 0.25 µm (P/N 26095-1420)

Liner Siltek Baffled (P/N 453T21210)


Organotins by GC-MS/MS

Sample preparation

400 mL water sample Adjust pH to 5 Ethylation by adding a 2% w/v sodium tetraethyl borate; solution in 0.1M NaOH; Extraction with pentane; Evaporate to 400 μL; 3 μL injection volume



Compound	IDL in µg/I in El SRM
Tributyl tin	0.00007





Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration

ENTIFIC

- Metal analysis

Reference List



Polychlorinated Alkanes by GC-MS/MS

Methodology

GC method: 100 °C (1 min); 40 °C/min; 320 °C (3 min). Injector PTV: 60 °C (0.1 min); 14.5 °C/sec; 280 °C (1 min); 2 μ l, 1 min splitless.


TSQ Series Timed SRM mode


Source temperature 280 °C

Column


TraceGOLD TG-5SilMS 20 m x 0.18 mm x.18µm (P/N 26096-5780)

Liner Siltek Baffled (P/N 453T21210)





SCOP AT 5 19 1 SCOP 16 5 1



Download the poster presentation on GC-MS amenable EU WFD compounds

Introduction to WFD

Introduction to UK CIP 2 List of compounds and

- analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Dioxins in Biota by GC-MS/MS

niection No

WHO-PCDO/F-TEQ ub

GC method

Initial 100 °C, hold 2 min, ramp 25.0 °C /min - 250 °C, Ramp 2.5 °C/ min, 285 °C , ramp 10 °C/ min to 330 °C Hold 5.0 min. transfer line: 280 °C

SSL method

Injector temperature 260 °C, splitless injection 2 min

TSQ 8000 Evo System in El SRM mode


Source temp: 300 °C

Ionization: EI, 40 eV


Column

TraceGOLD TG-5SiIMS 60 m \times 0.25 mm l.D. \times 0.25 μm (P/N 26096-1540)

Liner SSL single taper (P/N 453A2342)



2378TCDD and 13C2378 TCDD from left to right: at 95 fg in fishmeal; 98 fg in animal fat and 28 fg in egg fat sample.



Results in fish meal compared with GC-HRMS

PCDD/Fs limit of quantitation						
Compound	pg/µL					
2378-TCDF	0.01					
2378-TCDD	0.01					
12378-PeCDF	0.02					
23478-PeCDF	0.02					
12378-PeCDD	0.02					
123478-HxCDF	0.04					
123678-HxCDF	0.04					
234678-HxCDF	0.03					
123478-HxCDD	0.06					
123678-HxCDD	0.03					
123789-HxCDD	0.04					
123789-HxCDF	0.04					
1234678-HpCDF	0.02					
1234678-HpCDD	0.05					
1234789-HpCDF	0.02					
OCDD	0.05					
OCDF	0.03					

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS


GC-MS compounds


- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- **ICP-MS** compounds
- Highlights & configuration
- Metal analysis

Reference List





Pesticides, PAH, and Phenols by GC-MS/MS

Methodology

GC method: 100 °C (1 min); 40 °C/min; 320 °C (3 min). Injector PTV: 60 °C (0.1 min); 14.5 °C/sec; 280 °C (1 min); 2 µl, 1 min splitless.


TSQ series system Timed SRM mode

Source temperature 280 °C


Column


TraceGOLD TG-5SiIMS 20 mm x 0.18 mm x 0.18 µm (P/N 26096-5780)

Liner Siltek Baffled (P/N 453T21210)



Octylphenol at 0.4 pg in river water, cal curve, r²= 0.9923





Benzo(a)pyrene at 0.4 pg in river water; cal curve, r²=0.9991



Introduction to WFD

Introduction to UK CIP 2 List of compounds and

analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Combined Pesticides, PAH, and Phenols: Results

	2013/39/EU		IDL shashes		20	IDL absolute	
Compounds	AA-EQS inland surface water required limit		IDL absolute amount injected in pg	Compounds	AA-EQS inland surface water	EQS/3 required limit	amount injected in pg
	(µg/L)	(µg/L)			(µg/L)	(µg/L)	injootou in pg
Alachlor	0.3	0.1	0.3	1,2,3-TCB			0.05
Atrazine	0.6	0.2	1.5	1,2,4-TCB			0.1
Chlorfenvinphos	0.1	0.034	4	1,3,5-TCB			0.1
Chlorpyrifos (ethyl)	0.03	0.01	0.3	Trifluralin	0.03	0.01	0.1
Cyclodiene	∑=0.01	0.0034		Dicofol	0.0013	0.0004	0.5
Aldrin	0.0025	0.0008	0.1	Quinoxyfen	0.15	0.05	0.2
Dieldrin	0.0025	0.0008	0.1	Aclonifen	0.12	0.04	2
Endrin	0.0025	0.0008	0.2	Bifenox	0.012	0.004	2
Isodrin	0.0025	0.0008	1.5	Cybutryne	0.0025	0.00083	0.5
DDT total	0.025	0.0083		Cypermethrin	0.00008	0.00003	2
pp DDT	0.01	0.0033	0.1	Dichlorvos	0.0006	0.0002	0.3
op DDT	0.025	0.0083	0.1	Heptachlor	0.000002	0.00000067	0.5
pp DDE	0.025	0.0083	0.15	Heptachlorepoxide (Cis and Trans)	0.000002	0.00000067	2
pp DDD	0.025	0.0083	0.5	Terbutryn	0.065	0.022	0.3
op DDE	0.025	0.0083	0.1	Anthracene	0.1	0.034	0.15
op DDD	0.025	0.0083	0.1	Benzo(a)pyrene	0.00017	0.000057	0.25
Endosulfan a	0.005	0.0017	1.5	Benzo(b)fluoranthene	0.00017	0.000057	0.25
Endosulfan b	0.005	0.0017	2	Benzo(k)fluoranthene	0.00017	0.000057	0.25
Hexachlorobenzene	(MAC) 0.05	0.017	0.08	Benzo(ghi)perylene	0.00017	0.000057	0.5
Hexachlorobutadiene	(MAC) 0.6	0.2	0.02	Indeno(123cd)pyrene	0.00017	0.000057	0.5
Hexachlorocyclohexane	0.02	0.007		Fluoranthene	0.0063	0.0021	0.25
a-HCH			0.2	Naphtalene	2	0.67	0.25
b-HCH			0.2	Pentachlorobenzene	0.007	0.0023	0.15
g-HCH			0.2	Nonylphenol (4-nonylphenol)	0.3	0.1	0.2
d-HCH			0.5	Octylphenol	0.1 0.034		0.00
Simazine	1	0.33	1.2	(4-1.1'. 3.3'-tetramethylbutyl)-phenol			0.08
Trichlorobenzenes	0.4	0.133		Pentachlorophenol	0.4	0.14	2

Download the poster presentation on GC-MS amenable EU WFD compounds

Learn more about analyzing alkylphenoles with GC-MS/MS

Learn more about analyzing pesticides with GC-MS/MS

Learn more about analyzing pyrethroids in water with GC-MS

 (Δ)

Introduction to WFD

- Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





LC-MS System Highlights

The TSQ Quantiva triple quadrupole mass spectrometer exceeds the most stringent analytical requirements for quantitative performance with attogram-level sensitivity, unprecedented usability and exceptional robustness.

- Unprecedented Quantitative Performance Thermo Scientific[™] active ion management (AIM[™]) technology — electrodynamic ion funnel, ion beam guide with neutral blocker, Thermo Scientific[™] HyperQuad[™] quadrupole mass filter, and active collision cell — enables attogram-level sensitivity.
- Ultrafast selected-reaction monitoring (SRM) of 500 SRM/s, with up to 30,000 definable SRMs, enables quantification of more compounds in less time.
- Easy, robust, and reliable intuitive drag-and-drop method editor software with application templates simplifies method development and operation.
- Thermo Scientific[™] Ion Max NG[™] ion source makes all gases and voltages automatically on installation for ease of use, while allowing flexible spray position for ultimate performance.



Introduction to WFD

- Introduction to UK CIP 2 List of compounds and analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Recommended LC-MS Configuration

Compound group	LC system	Injection volume (µL)	MS analyzer (Acquisition mode)
Herbicides	Thermo Scientific™ Dionex™UltiMate™ 3000 EquanMAX Plus	10 (offline) 100 (Large Vol Inj) 1000 (online)	TSQ Endura (SRM)
Insecticides	Insecticides UltiMate 3000 EquanMAX Plus		TSQ Endura (SRM)
Fungicides	UltiMate 3000 EquanMAX Plus	10 (offline) 100 (Large Vol Inj) 1000 (online)	TSQ Endura (SRM)
Pharmaceuticals	UltiMate 3000 EquanMAX Plus	10 (offline) 2000 (online)	TSQ Quantiva (SRM)
Endocrine disruptors hormones	UltiMate 3000 EquanMAX Plus	50 (offline) 2000 (online)	TSQ Quantiva (SRM)
Perfluorooctane sulfonic acid (PFOS)	UltiMate 3000 EquanMAX Plus	10 (offline) 2000 (online)	TSQ Quantiva (SRM)

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Pesticides by LC-MS/MS

Confirming Ions

PTT IL 40 | WE can taget study

.

LC Method

UltiMate 3000 HPLC system


Column Thermo Scientific[™] Accucore[™] C18 100 x 2.1 mm, 2.6 µm (P/N 17126-102130) at 30 °C Mobile phase: (A) 2mM ammonium acetate (B) MeOH Gradient: Start at 10% B (Hold 1 min), ramp to 90% B in 5min. Ramp to 100% B in 1 min (Hold 1 min). Ramp to 10% B in 0.10 min (Hold 4 min)


Flow rate 0.4 mL/min

Injection Volume 10 µL

MS Quantiva ESI (SRM)


Capillary temperature (°C): 270; vaporizer temperature (°C): 400 sheath gas pressure (Arb): 50; aux gas pressure (Arb): 20 sweep gas pressure (Arb): 1





Isoproturon at LOD; showing two transitions

Quinoxyfen at LOD; showing two transitions



₹ ×

4 1/1 1

6.40 (WE real level Broad

RT 6.40 AA 104953.03 AH 10051.10 SN 39.19

Quan Peak

Learn how this customer benefits from the Thermo Scientific LC-MS/MS systems

Diuron at LOD; showing two transitions

Introduction to WFD

- Introduction to UK CIP 2 List of compounds and analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Pesticides by LC-MS/MS

Results

	2013/39	/E11	Required LOD	River	Effluent	Influent
	2013/39	/EU	(pg on column)	LOD	LOD	LOD
Compounds	AA-EQS inland surface water (µg/L)	EQS/3		(pg on column)	(pg on column)	(pg on column)
		(µg/L)				
Atrazine	0.6	0.2	2	0.01	0.01	0.01
Chlorfenvinphos	0.1	0.034	0.34	0.01	0.01	0.01
Diuron	0.2	0.067	0.67	0.1	0.1	0.1
Isoproturon	0.3	0.1	1	0.01	0.01	0.01
Simazine	1	0.34	3.4	0.03	0.03	0.01
Quinoxyfen	0.15	0.05	0.5	0.01	0.01	0.01
Terbutryn	0.065	0.022	0.22	0.05	0.01	0.1

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds


- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Pharmaceuticals by LC-MS/MS

LC Method

UltiMate 3000 HPLC system

Column Accucore C18 100 x 2.1 mm, 2.6 μm (P/N: 17126-102130) at 30 °C


Mobile phase (A) 2mM ammonium acetate (B) MeOH Gradient: Start at 10% B (Hold 1 min), ramp to 90% B in 5 min. Ramp to 100% B in 1 min (Hold 1 min). Ramp to 10% B in 0.10 min (Hold 4 min)


Flow rate 0.4 mL/min


Injection volume 10 µL

MS Quantiva ESI (SRM)

Capillary temperature (°C): 270; vaporizer temperature (°C): 400 Sheath gas pressure (Arb): 50; aux gas pressure (Arb): 20 Sweep gas pressure (Arb): 1







Atorvastatin

Erthromycin

Download application note on analysis of basic and acidic pharmaceutical compounds

Introduction to WFD

- Introduction to UK CIP 2 List of compounds and
- analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Pharmaceuticals by LC-MS/MS

	CIP	2	River LOD	Effluent	Influent
Compounds	Required LOD (µg/L)	LOD In pg	(pg on column)	(pg on column)	(pg on column
Diclofenac (WFD)	0.01	0.1	0.03	0.1	0.03
Ibuprofen	0.01	0.1	1	0.8	0.8
Atorvastatin	0.01	0.1	0.05	0.05	0.05
ortho-hydroxyatorvastatin	0.01	0.1	0.05	0.1	0.1
para-hydroxyatorvastatin	0.01	0.1	0.05	0.1	0.1
Propanolol	0.01	0.1	0.4	0.2	0.2
Atenolol	0.01	0.1	0.2	0.15	0.1
Amoxicillin	Not on draft - I	March 2014	1	1	1
Erythromycin	0.1	1	0.6	0.15	0.15
Norerythromycin	0.1	1	Standard not available		
Azithromycin	0.005	0.05	2	tbd	0.4
Clarithromycin	0.01	0.1	0.1	0.1	0.03
Ciprofloxacin	0.01	0.1	1	tbd	tbd
Ranitidine	0.1	1	0.3	0.1	0.1
Carbamazepine	0.1	1	0.01	0.1	0.1
10,11-epoxycarbamazepine	0.1	1	0.01	0.1	0.04
Sertraline	0.01	0.1	0.5	1	0.1
Norsertraline	0.01	0.1	0.5	1	0.1
Fluoxetine	0.01	0.1	0.15	0.1	0.1

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Phenols and Triazoles with LC-MS/MS

LC Method

UltiMate 3000 HPLC system

Column Accucore C18 100 x 2.1 mm, 2.6 μm (P/N 17126-102130) at 30 °C

Mobile phase (A) 2mM ammonium acetate (B) MeOH


Gradient Start at 10% B (Hold 1 min), ramp to 90% B in 5min. Ramp to 100% B in 1 min (Hold 1 min). Ramp to 10% B in 0.10 min (Hold 4 min);


Flow rate 0.4 mL/min

Injection volume 10 μL

MS Quantiva ESI (SRM)

Capillary temperature (°C): 270; vaporizer temperature (°C): 400 Sheath gas pressure (Arb): 50; aux gas pressure (Arb): 20 Sweep gas pressure (Arb): 1





Tris(isopropylphenyl)phosphate at LOD; showing two transitions

		CIP 2					
Compounds	Required LOD (µg/L)	Required LOD pg on column	River LOD pg on column	Effluent pg on column	Influent pg on column		
Nonylphenol/4-nonylphenol	0.1	0.4	GC-MS	GC-MS	GC-MS		
Octylphenols	0.1	0.14	GC-MS	GC-MS	GC-MS		
4-nonylphenol triethoxylate	0.1	0.4	0.2	0.1	0.1		
Tris(isopropylphenyl)	0.01	-	0.1	0.1	0.1		
Phosphate	0.01		0.1	0.1	0.1		
Benzotriazole	0.002	-	0.1	-	0.1		

Benzo triazole at LOD; showing two transitions

Introduction to WFD

Introduction to UK CIP 2 List of compounds and

analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids

ICP-MS compounds

- Highlights & configuration
- Metal analysis

Reference List



Perfluorooctanesulfonic Acid (PFOS) by LC-MS/MS

LC Method

UltiMate 3000 HPLC system

Column Accucore C18 100 x 2.1 mm, 2.6 μm (P/N 17126-102130) at 30 °C

Mobile phase (A) 2 mM ammonium acetate (B) MeOH Gradient: Start at 10% B (Hold 1 min), ramp to 90% B in 5min. Ramp to 100% B in 1 min (Hold 1 min). Ramp to 10% B in 0.10 min (Hold 4 min);

Flow rate 0.4 mL/min


Injection volume 10 μL


MS Quantiva ESI (SRM)

Capillary temperature (°C): 270; vaporizer temperature (°C): 400 Sheath gas pressure (Arb): 50; aux gas pressure (Arb): 20 Sweep gas pressure (Arb): 1

To reduce the contamination coming from the system we recommend to replace parts of the UHPLC system and introduce a trap column before the injector: Thermo ScientificTM Hypersil GOLDTM column $50 \times 2.1 \text{ mm} \times 3 \text{ µm} (P/N 25003-052130)$

		2013/39/EU		LODs
Compounds	AA-EQS inland surface water (µg/L)	AA-EQS Other surface water (µg/L)	Biota µg/kg	(µg/L)
PFOS	0.00065	0.00013	9.1	0.1
PFOA	-	-	-	0.1





PFOS linear and branched at 5 μ g/l; showing two transitions

Learn how this customer benefits from the Thermo Scientific LC-MS/MS

Introduction to WFD

Introduction to UK CIP 2 List of compounds and

analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Hexabromocyclododecane (HBCDD) by LC-MS/MS

Results

LC Method

UltiMate 3000 HPLC system

Column Accucore C18 100 x 2.1 mm, 2.6 μm (P/N 1726-102130) at 35 °C

Mobile phase (A) 2mM ammonium acetate B) MeOH Gradient: Start at 0% B, ramp to 100% B in 10 min (Hold 2 min). Ramp to 0% B in 0.10 min (Hold 3.3 min)


Flow rate 0.3 mL/min

Injection volume 10 μL

MS Quantiva ESI (SRM)


Capillary temperature (°C): 280; vaporizer temperature (°C): 120 Sheath gas pressure (Arb): 50; aux gas pressure (Arb): 15 sweep gas pressure (Arb): 1

	2013/39/EU			LODs	
Compounds	AA-EQS inland surface water (µg/L)	AA-EQS other surface water (µg/L)	Biota µg/kg	(required}	achieved
HBCDD	0.0016	0.0008	167	2000	2000



Gradient elution using 2mM ammonium acetate and MeOH Accucore C18 - 100 x 2.1mm, 2.6 μm

HBCDD: α , β , γ isomer separation



Thermo scientific

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Steroid Analysis by LC-MS/MS

LC Method

UltiMate 3000 HPLC system

Column Hypersil GOLD 50 x 2.1 mm, 1.9 μm (P/N 25002-052130) at 30°C

Mobile phase (A) 0.02 % ammonium hydroxyde in H_20 (B) 0.02 % ammonium hydroxyde in MeOH

Flow rate 0.4 mL/min

Injection volume 50 μL

MS Quantiva ESI (SRM)

Capillary temperature (°C): 320; vaporizer temperature (°C): 350 Sheath gas pressure (Arb): 50; aux gas pressure (Arb): 40 sweep gas pressure (Arb): 1

40 10 10 2.96 2.96 5.76 6.38 6.85 7.69 8.82 9.77 100 2.99 100 1.91 3.64 4.23 4.90 5.76 6.38 6.85 7.69 8.82 9.77 100 2.99 1.91 3.63 4.20 5.05 5.63 6.22 6.76 7.71 8.58 9.43 100 0.74 1.91 2.63 3.63 4.30 5.05 5.63 6.22 6.76 7.71 8.58 9.43 100 0.93 1.85 2.32 3.68 4.53 5.35 6.14 6.78 7.44 8.05 8.55 8.98 9.66 100 0.93 1.85 2.32 3.68 4.53 5.35 6.14 6.78 7.44 8.05 8.55 8.98 9.66 100 0.23 .041 3.64 4.46 5.10 5.84 6.45 7.12 7.90 8.77 9.66	RT: 0.00 - 10.00	↓	Blank River	Water					
80 River Water Sample # 1 spiked at 0.35 ng/L 90 0.74 1.91 90 3.63 4.30 90 3.63 4.30 90 3.01 90 8 91 3.63 92 0.93 93 1.85 93 1.85 93 3.68 94 4.53 94 5.95 94 6.14 95 6.14 96 8.55 98 9.66 90 2.99 90 1.07 91 1.55 92 0.155 94 5.10 95 6.14 96 8.72 96	20-1 1 WWWWWWWW	2.96	4.23 4.9	0 5.76	6.00 6.0	- 7	69	0.02	9.77
80 River Water Sample # 1 spiked at 0.35 ng/L 90 0.74 1.91 90 3.63 4.30 90 3.63 4.30 90 3.01 90 8 91 3.63 92 0.93 93 1.85 93 1.85 93 3.68 94 4.53 94 5.95 94 6.14 95 6.14 96 8.55 98 9.66 90 2.99 90 1.07 91 1.55 92 0.155 94 5.10 95 6.14 96 8.72 96	0-h	2.99	Mr. M.M.	M-M	0.30 0.0		<u>^</u>	0.02	
60 River Water Sample # 1 spiked at 0.35 ng/L 0.74 1.91 0 3.63 0 3.63 0 3.63 0 3.63 0 3.63 0 3.63 0 3.63 0 3.61 0 3.63 0 3.63 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.01 0 3.03 0 3.03 0 3.03 0 3.04 0 3.05 0 3.04 0 3.05 0 3.04 0 3.04 0 5.04									
River Water Sample # 1 spiked at 0.35 ng/L 20 0.74 1.91 3.63 4.30 5.05 5.63 6.22 6.76 7.71 8.58 9.43 00 3.01 3.01 3.01 3.01 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00 8.00									
20 0.74 1.91 2.63 3.63 4.30 5.05 5.63 6.22 6.76 7.71 8.58 9.43 3.01 3.01 3.01 8.59 9.43 0.00			River Water	Sample # 1 spil	ked at 0.35 ng/L				
00 3.01 80 60 40 0.93 0.93 1.85 0.23 0.93 1.07 2.99 80 60 40 0.101 1.07 200 1.07 201 1.07 202 1.07 203 1.07 204 1.07 205 1.07 204 1.07 205 204 1.07 205 204 1.07 205 204 1.07 205 206 1.07 207 208 804 209 804 807 807 807 807 807 807 808 808 809 807 808 808 809 809 800 801 802 803									
100 3.01 80 60 40 0.93 0.93 1.85 0.23 1.85 0.23 1.85 2.99 80 80 60 4.0 0.107 2.99 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 80 81 81 82 82 82 82 82 82 82 82 82 8341	20 Martin Martin 20	my hanna	4.30 5	5.63	6.22 6.76	7	.71	8.58	9.43
60- 40- 20- 0.33 1.85 2.32 0- 0.23 M ^M M M M M M M M M M M M M M M M M M		3.01							
River Water Sample # 2 spiked at 0.35 ng/L 20 0.93 1.85 2.32 3.68 4.53 5.35 6.14 6.78 7.44 8.05 8.55 8.98 9.66 100 2.99 2.99 2.99 80 60 60 40 2.99 8.05 8.55 8.98 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66 9.66	80 -								
40 - 0.93 0.23 MM M M M M M M M M M M M M M M M M M	60 -								
200 0.23 MM/mmmmmmmmmmmmmmmmmmmmmmmmmmmmmmmmmmm	40 -		River Water	Sample # 2 spil	ked at 0.35 ng/L				
00 2.99 80 60 40 1.07 20 0.21 0.07 3.41 3.41 2.04 4.6 5.10 5.64 6.45 7.12 7.00 8.77 9.65	20 0.93 1.85 2.32	2.60							
100 80 60 40 20 02 1.07 20 02 1.55 2.10 341 204 465 510 584 645 712 700 877 966	0=	day have a	4.53	5.35	6.14 6.78	7.44	8.05	8.55 8.98	9.66
60- 40- 20- 20- 20- 20- 20- 20- 20- 20- 20- 2		2.55							
40 - 1.07 20 - 1.07 20 - 0.21 M M M M M M M M M M M M M M M M M M M	80 -								
20 - 1.07 20 - 0.21 $M - 1.55$ 2.10 3.41 204 446 510 584 845 712 700 877 966	60 -								
20 - 0.21 Martin 155 2.10 341 204 446 5.10 5.04 6.45 7.12 7.00 8.77 9.66	40 - 1.07	1	River Water	Sample # 3 spik	ed at 0.35 ng/L				
$0 - 0.31$ 3° 3.94 4.40 5.10 5.84 6.45 7.12 7.90 8.77 9.66	20	3.41 0		F 40 - F 0 4	0.45	7.40	=	0.77	0.00
		ma land	<u> </u>	1	1	1		· · · ·	
0 1 2 3 4 5 6 7 8 9 1 Minutes	0 1 2	3			6	7	8	9	10

	CIP 2	COM(2011)876	2013/39/EU	LOD influent/ effluent
Compounds	Required LODs			(ng/L)
	Influent/effluent	AA-EQS inland surface water (ng/L)	AA-EQS inland surface water	
	(ng/L)			
17-alpha-ethinylestradiol (EE2)	0.03	0.35	Watch list	100
17-beta-estradiol (E2)	0.3	4	Watch list	10
Estrone	1	_	-	15

Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique Sample preparation

- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



The Thermo Scientific[™] iCAP[™] Q ICP-MS system has been developed with groundbreaking technology to enable advanced high-performance analysis combined with total reliability and ultra-flexibility.

- Ergonomically designed quadrupole ICP-MS system with smallest footprint
- Innovative RAPID (Right Angle Positive Ion Deflection) lens technology (90° ion optics) for separation of ions and neutrals
- Innovative interface with skimmer cone featuring unique, user-replaceable inserts mounted behind the cone tip to minimize memory effects
- The only ICP-MS system to include proprietary QCell technology combining proven He KED interference reduction
- High-performance quadrupole analyzer pumped by a novel split flow turbo pump backed by a single rotary pump
- New, simultaneous analog/PC detector with real time multi-channel analyzer electronics provides >9 orders of dynamic range suitable for both steady state and transient signal analysis
- The combination of an advanced inlet system, the RAPID lens and QCell provide enhanced performance for outstanding signal to background
- New, exceptionally robust RF generator with fast dynamic frequency impedance matching for high plasma stability without need for a shield





Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List





Peristaltic pump speed Pump


	tubing for both
Carrier	(sample) and internal
	standard
Nebulizer	PFA-ST
Interface	cones Nickel
RF Power	1550 W
Cool gas flow	14 L/min
Auxiliary gas flow	0.8 L/min
Nebulizer gas flow	0.97 L/min
Number of sweeps	20
Dwell Times	0.001 - 0.02 ms
Data points per isotope	1
Replicates per analysis	3
Signal stabilization time	15 s
FAST sample loop	0.5 mL
FAST uptake time	3 s
FAST rinse time (dual stations)	1+3 s
Total analysis time (sample to sample)	56 s

40 rpm

tubing orange/green

Metal Analysis by ICP-MS: Results

Applications



Fast analysis for all compounds is achieved

Element	Lowest required AA- EQS (in µg/L)	Method detection limit shown in AN43127 (in μg/L)
Cd	≤ 0.08	0.017
Pb	1.2	0.007
Hg	0.07	0.011
Ni	4	0.014

Introduction to WFD

Introduction to UK CIP 2

- List of compounds and
- analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis

Reference List



Heinz Singer et al. Determination of biocides and pesticides by on-line solid phase extraction coupled with mass spectrometry and their behavior in wastewater and surface water published at Environmental Pollution – 2010 (158) 3054-3064

References

Liza Viglino et al. On-line solid phase extraction and liquid chromatography/tandem mass spectrometry to quantify pharmaceuticals. pesticides and some metabolites in wastewaters. drinking and surface waters published at Journal of Environmental Monitoring – 2008 (10) 482-489

Laszlo Hollosi et al. Analytical sensitivity of direct injection LC-MS methods for the compounds of EU Water Frame Work directive presented at SETAC Glasgow 2013

Marina Gorga et al. Multi-residue analytical method for the determination of endocrine disruptors and related compounds in river and waste water using dual column liquid chromatography switching system coupled to mass spectrometry published at Journal of Chromatography A – 2013 (1295) 57-66

Rebeca López-Serna et al. Direct Analysis of pharmaceuticals, their metabolites and transformation products in environmental waters using on-line TurboFlow chromatography-liquid chromatography-tandem mass spectrometry published at Journal of Chromatography A – 2012 (1252) 115-129

Marta Llorca et al. Analysis of perfluoroalkyl substances in waters from Germany and Spain published at Science of the Total Environment – 2012 (431) 139-150

Salima Idder et al. Quantifitative on-line pre-concentration-liquid chromatography coupled with tandem mass spectrometry method for the determination of pharmaceutical compounds in water published at Analytica Chimica Acta-2013 (805) 107-115

Introduction to WFD

- Introduction to UK CIP 2
- List of compounds and
- analytical technique
- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS
- GC-MS compounds
- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols
- LC-MS compounds
- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- ICP-MS compounds
- Highlights & configuration
- Metal analysis
- Reference List



Sara Valsecchi et al. Optimization of an on-line SPE-UHPLC/MS/MS method for the analysis of perfluoroalkylated acid in Italian drinking and surface waters presented at SETAC Glasgow 2013

Technical note 148 : Automated extraction and determination of human hormones in drinking water using solid-phase extraction and HPLC with UV detection

Application note 52389 - Consolidated GC-MS/MS analysis of OCPs, PAHs, and PCBs in environmental samples

Application note 51899 – Latest advances in the analysis of volatile organic compounds by single quadrupole GC-MS

Technical note 10319 – Simplifying complex multi-residue pesticide methodology in GC-MS/MS

Technical note 52099 – The determination of organotins in water using triple quadrupole GC-MS/MS

Application note 30098 – DFS-analysis of brominated flame retardants with high resolution GC/MS

Application note 43098 – Speciation analysis of Cr (III) and Cr (IV) in drinking waters using anion exchange chromatography coupled to the Thermo Scientific iCAP Q ICP-MS

Application note 40849 - Lead in natural waters by graphite furnace atomic absorption using EPA Method 200.9

Application note 40851 – Arsenic in natural waters by graphite furnace atomic absorption using EPA Method 200.9

Application note 40853 - Environmental series - US EPA Method 200.7 using the iCAP 6500 Duo ICP



Introduction to WFD

Introduction to UK CIP 2 List of compounds and analytical technique

- Sample preparation
- General sample preparation
- Automated solutions
- LC-MS compounds
- GC-MS compounds
- ICP-MS

GC-MS compounds

- Highlights & configuration
- Volatile Organics
- PBDE
- Organotins
- Chlorinated alkanes
- Dioxins
- Pesticides, PAH, phenols

LC-MS compounds

- Highlights & configuration
- Pesticides
- Pharmaceuticals
- Phenoles and triazoles
- PFOS
- HBCDD
- Steroids
- **ICP-MS** compounds
- Highlights & configuration
- Metal analysis

Reference List



www.thermoscientific.com/water-analysis

©2015 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0 Australia +61 3 9757 4300 Austria +43 810 282 206 Belgium +32 53 73 42 41 Brazil +55 11 3731 5140 Canada +1 800 530 8447 China 800 810 5118 (free call domestic) Italy +39 02 950 591 400 650 5118

Denmark +45 70 23 62 60 Europe-Other +43 1 333 50 34 0 Korea +82 2 3420 8600 Finland +358 10 3292 200 France +33 1 60 92 48 00 Germany +49 6103 408 1014 India +91 22 6742 9494

Japan +81 6 6885 1213 Latin America +1 561 688 8700 Middle East +43 1 333 50 34 0 Netherlands +31 76 579 55 55 **New Zealand** +64 9 980 6700 Norway +46 8 556 468 00



Russia/CIS +43 1 333 50 34 0 **Singapore** +65 6289 1190 Sweden +46 8 556 468 00 Switzerland +41 61 716 77 00 **Taiwan** +886 2 8751 6655 UK/Ireland +44 1442 233555 **USA** +1 800 532 4752



A Thermo Fisher Scientific Brand