


# Scale up

your structural  
analysis



Orbitrap Ascend  
Structural Biology Tribrid  
mass spectrometer



# Scale up insight from complex native structures

## Unlock the intricacies of complex molecular structures

Built to meet the demands of structural biology, the Thermo Scientific™ Orbitrap™ Ascend Structural Biology Tribid™ mass spectrometer adds innovations to scale up characterization of complex structures, including intact proteoforms, post-translational modifications (PTMs) and complex native structures, as well as peptide-level analyses such as hydrogen exchange MS and cross-linked peptides identification. Analyze challenging species up to  $m/z$  16,000 with the Native MS option. With these capabilities and more, access unprecedented experimental throughput, versatility and usability to meet tomorrow's challenges while achieving more insight today.



### Native proteomics

Analyze complex mixtures, cover higher mass range and perform orthogonal fragmentation techniques



### Quantitation

Achieve accurate, proteome-wide quantitation for **structural characterization**



### Gas-phase fractionation

Completely characterize protein complexes and proteoforms




### Peptide-level structural analysis

Choose from the **widest selection** of methods for **identification and quantification of cross-linked peptides**, HDX-DIA



### Ligand binding analysis


Characterize complex samples with wide dynamic range and analyze labile compounds



Orbitrap Ascend Structural Biology Tribid mass spectrometer


# Innovative technology addresses structural analysis needs

The Orbitrap Ascend Structural Biology Tribrid mass spectrometer features innovations to scale up your native structural characterization with throughput, versatility and ease-of-use.



**OPTIONS** IC | ETD | PTCR | Native MS\* | UVPD | FAIMS Pro Duo interface

\*New on this platform



# Native complex

# Identify more species across a broad $m/z$ range in complex mixtures with DIA-PTCR

## Scale-up complex sample characterization

The Native MS option extends the quadrupole isolation range to  $m/z$  2,000–8,000 from  $m/z$  50–2,000, and isolation width to  $m/z$  5–3,000 from  $m/z$  0.4–1,200. The proton transfer charge reduction (PTCR) ion source option generates perfluoroperhydrophenanthrene (PFPP) ions for subsequent gas-phase, ion-ion reactions. When narrower data independent acquisition (DIA) windows are used for PTCR analysis of complex native mixtures, the PTCR spectra generated are simplified and sensitivity is increased. In addition, extending the mass range detected enables elucidation of hidden peaks in the charge state envelope. For the analysis of membrane proteins in nanodisc-based experiments, researchers can use this capability to enhance structural insight, such as the **stoichiometry of different lipid classes** and bound proteins.




**“I am very impressed by the high-mass quadrupole isolation, which enables exciting new DIA-PTCR experiments to characterize complex samples.”**


Michael T. Marty, PhD  
Associate Professor of Chemistry & Biochemistry  
University of Arizona

## DIA-PTCR of membrane proteins in nanodiscs


### (A) Full scan




### (B) PTCR on $m/z$ 7000





### (C) PTCR on $m/z$ 7300



### (D) PTCR on $m/z$ 7600



### (E) Sliding window deconvolution of DIA-PTCR spectrum



Orbitrap Ascend Structural Biology Tribid mass spectrometer analysis of nanodiscs using extended quadrupole isolation and DIA-PTCR. Deconvolution of DIA-PTCR spectra produces a molecular weight profile of approximately 140–175 kDa. The spacing between adjacent deconvoluted peaks at 700–750 Da fall into the molecular weight range associated with nanodiscs lipids.

Data courtesy of Associate Professor Michael T. Marty, University of Arizona.

# Identify native membrane proteins using PTCR and top-down techniques

# Native top-down

## Scale up to realize the potential of native proteomics


With the Orbitrap Ascend Structural Biology Tribrid mass spectrometer, it's now possible to realize the potential of native proteomics, particularly for difficult-to-analyze membrane proteins. Using native protein sample preparation techniques, complexes can be directly injected into the mass spectrometer and the proteins elucidated using their precursor and fragmentation patterns. At the intact protein level, PTCR simplifies the spectra produced from complex samples. For top-down analysis, difficult-to-break proteins can be fragmented using multiple fragmentation options, generating enough sequence coverage for high-confidence identification of proteins within complex samples.





**“I’m delighted that we now have an Orbitrap Ascend Structural Biology Tribrid mass spectrometer because it has transformed our ability to study complicated systems. It gives us the additional flexibility to break the intact molecule into its components and really identify what it is, and also what’s there. Often, in these very large assemblies, there’s a small molecule hidden. That’s what we want to find. Is there a metabolite? Is there a drug that’s changing the properties of this protein? That’s very important for me.”**

Professor Dame Carol Robinson DBE FRS FMedSci FRSC  
Director of the Kavli Institute for Nanoscience Discovery  
University of Oxford

## Fragments mapped on to structure of $\beta$ 1AR



5 m/z width quadrupole isolation of precursor m/z 3152.92 and HCD CE95



Data courtesy of Corinne Lutomski, Jack Bennett and Tarick El-Baba, Professor Dame Carol Robinson's lab, University of Oxford, and Ildir Liko, OMass Therapeutics.



## Scale up cross-link identification

The presence of mono-links complicates the process of cross-linked peptide identification by MS even after the enrichment. By analyzing intensity ratios of specific peaks with enrich-able cross-linkers, partial differentiation between the cross-links and mono-links in samples can be achieved. Real-Time Library Search triggering of high-resolution MS<sup>2</sup> scans based on the intensity ratios significantly enhances cross-linked precursors identification, predominantly for unenriched samples and short chromatography gradients. The approach particularly improves productivity for higher throughput and resource-constrained biological projects.



**“The unique feature of Real-Time Library Search on the Orbitrap Ascend Structural Biology Tribrid mass spectrometer enables on-the-fly selection of cross-linked precursors for MS<sup>2</sup> sequencing. This is really exciting and has always been what we hoped to achieve over the last years—selective targeting and sequencing of cross-linked peptides.”**

Fan Liu, PhD  
 Professor of Biochemistry  
 Leibniz Research Institute for Molecular Pharmacology



Sixty-four human proteins were cross-linked with PhoX (p/n A52286). Cross-linked peptides were enriched and spiked into HeLa digest at a 1:20 ratio and analyzed on the Orbitrap Ascend Structural Biology Tribrid mass spectrometer with (+) and without (-) Real-Time Library Search. Cross-links were identified using Thermo Scientific™ Proteome Discoverer™ software revision 3.1 with XlinkX node.

*Data courtesy of Max Ruwolt and Fan Liu, Leibniz-Research Institute for Molecular Pharmacology.*

# Perform disulfide mapping and quantitation with EThcD and the XlinkX node

# Disulfide mapping

## Scale up disulfide bridge assessment certainty


Quality control assessment of therapeutic recombinant proteins requires accurate measurement of critical attributes such as sequence fidelity, proper folding and PTMs. Errors can lead to diminished bioactivity and, in the context of therapeutic proteins, an elevated risk for immunogenicity. Though analytically challenging, assessment of disulfide bridges—the covalent bonds linking two Cysteine residues in a protein—is essential because they determine the correct folding and stability of proteins and thus have a major influence on their efficacy.

When equipped with the Thermo Scientific™ FAIMS Pro Duo interface and electron transfer dissociation (ETD) with supplemental activation (EThcD), the Orbitrap Ascend Structural Biology Tribrid mass spectrometer offers an efficient workflow for disulfide mapping. Microwave Assisted Acid Hydrolysis (MAAH) sample preparation effectively digests the protein of interest, generating massive peptide redundancy. Because MAAH is prone to background noise, the FAIMS Pro Duo interface uses differential ion mobility to remove interferences prior to MS analysis, bringing the desired signals into focus. EThcD generates highly informative spectra with ions diagnostic of disulfide-bridged peptides. XlinkX, a node in Proteome Discoverer software, enables rapid, high confidence identification of disulfide bridges at very low false discovery rates and quantifies the occupancy correctly.



“By combining novel sample preparation and advances in data analysis with the amazing FAIMS and EThcD capabilities of the Orbitrap Ascend Structural Biology Tribrid mass spectrometer, we are obtaining information on disulfide bridges at unprecedented detail in less than one hour. This transforms high performance QC of recombinant protein products.”

Richard Scheltema, PhD  
Professor of Structural Proteomics  
University of Liverpool



Trastuzumab is digested using MAAH. The subsequent background noise produced by MAAH is removed prior to MS analysis using the orthogonal selectivity provided by the FAIMS Pro Duo interface, enhancing signal-to-noise ratios for the ions of interest, including ions diagnostic of disulfide bridges.


*Data courtesy of Professor Richard Scheltema, University of Liverpool, and Søren Heissel, The Rockefeller University.*

# Proteoform detection

## Access higher mass proteoforms with PTCR fractionation of the MS<sup>1</sup>

### Scale up to measure larger proteoforms by top-down proteomics

Top-down proteomics enables scientists to identify and structurally characterize PTMs and protein isoforms—known collectively as proteoforms—as they occur in the cell. By combining PTCR and the Native MS option, which extends instrument mass range to  $m/z$  16,000, proteoform characterization was expanded to a new level, exploring larger proteins above the 30 kDa range by simplifying complex MS<sup>1</sup> spectra using targeted isolation windows. The same  $m/z$  range is again isolated and fragmented for proteoform identification and primary structure characterization using Thermo Scientific™ ProSightPD™ software.



HEK293 cells were grown to 80–100% confluency before harvesting by trypsinization. Proteins were fractionated by molecular weight using 10% T gel. Fractions 2–14 were analyzed using tPRCT methodology.

*Sample courtesy of Assistant Professor Luca Fornelli, University of Oklahoma.*



# Experience more high-quality results with less hassle using automated, remote and schedulable system checks and calibrations


## Scale up convenience and ease-of-use

The Auto-Ready ion source is a fully integrated, standard, easy-to-use feature that increases laboratory productivity with automated, remote and schedulable system checks and internal calibrations. Because there is no need to remove the source (HESI, nESI or high-field asymmetric waveform ion mobility spectrometry [FAIMS]), there are no experimental setup interruptions required to perform internal calibrations. The user can automate the calibration to start at a scheduled time—for example, every week—when there are no experiments planned to run on the instrument. The calibration can run completely remotely, regardless of the nature of the last experiment. Because the calibration can be scheduled to occur regularly and automatically without interrupting vital work, users can expect to maintain mass spectrometer performance, improve data consistency and achieve more accurate and precise quantitation.




## Auto-Ready ion source

Separate ion transfer tube



Dedicated emitter



Robust delivery system



## Automated weekly calibrations

Status

Self-check is scheduled to run every Wednesday at 12:00 AM in Check, Calibrate if Required mode.

Polarity (+)    Polarity (+/-)

Orbitrap Mass

- Recommended Calibration: 5/2/2021

System

- Recommended Calibration: 5/9/2021
- Optional Calibrations

FlexMix Volume    Full (>= 70%)

Calibration

Mode: Check, Calibrate if required

Polarity: Positive and Negative

Type: Orbitrap Mass & System

Optional Calibrations

Easy-IC:

System self-check completed successfully at 01:59 PM on Feb 19

Passed

## Resources and support



### Services Central—All your service information at your fingertips

Spend less time searching for support and more time focusing on your important work. This online platform has what you need to easily manage your instruments and equipment.

Learn more at [thermofisher.com/servicescentral](https://thermofisher.com/servicescentral)



### Technical and online support

Helping you keep your instruments running at peak performance is our goal. Whether you're looking for an instrument manual or spare parts, want to submit a repair request or check on the status of your warranty or service contract, we have every support option you're looking for.

Learn more at [thermofisher.com/technicalresources](https://thermofisher.com/technicalresources)



### Protect your investments with expert lab services

Unity™ Lab Services provides a single source for integrated lab service, support and supply management. Our customized service offerings and world-class service experts have the flexibility and experience to address your laboratory's needs. We provide a complete portfolio of services and support solutions designed to help you improve productivity, reduce total cost of ownership and help ensure performance throughout your laboratory.

Learn more at [unitylabservices.com](https://unitylabservices.com)

## More resources

Thermo Scientific™ Ardia™ Platform  
[thermofisher.com/ArdiaProteomics](https://thermofisher.com/ArdiaProteomics)

Thermo Scientific™ BioPharma Finder™ software  
[thermofisher.com/biopharmafinder](https://thermofisher.com/biopharmafinder)

Thermo Scientific OptiMSe workflow  
[thermofisher.com/optimse](https://thermofisher.com/optimse)

Thermo Scientific Proteome Discoverer software  
[thermofisher.com/ProteomeDiscoverer](https://thermofisher.com/ProteomeDiscoverer)

 Learn more at [thermofisher.com/OrbitrapAscendStructuralBiology](https://thermofisher.com/OrbitrapAscendStructuralBiology)

**General Laboratory Equipment – Not For Diagnostic Procedures.** © 2024 Thermo Fisher Scientific Inc. All rights reserved. OMASS is a trademark of OMass Therapeutics Limited. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. **24-085-1056 | BR002862-EN 0524C**