


thermo scientific



## Celebrate simplicity

Thermo Scientific Dionex Easion  
Ion Chromatography System

**ThermoFisher**  
SCIENTIFIC



## Worth celebrating: A small, simple IC system for everyday ion analysis

The new Thermo Scientific™ Dionex™ Easion™ Ion Chromatography (IC) system is anything but “routine.” It’s a perfect fit for any lab doing everyday water analysis or other typical IC applications. Small in size but big in value, the system’s advanced technology, innovative design and intuitive Thermo Scientific™ Chromeleon™ Chromatography Data System (CDS) make it easy even for new operators to get first-rate results for either anions or cations.

# Simply better value




Dionex Easion IC system

- **Everything you need to get up and running in about 90 minutes**—including columns, suppressor, eluent concentrate, and regenerant concentrate—no additional pumps or equipment required
- **Membrane-based chemical suppression** ensures maximum uptime and reduces baseline drift to minimize sample-to-sample variability
- **Space-saving footprint** fits anywhere there's a bit of bench space
- **Simple operation and maintenance**—no extra items to change or maintain
- **Preconfigured analysis kits** take the guesswork out of system configuration for anions or cations
- **Industry-leading quality and reliability** give you confidence in your purchase decision—and assure your satisfaction for years to come



**“Simply Intelligent” Chromeleon CDS** speeds up daily operations and helps you achieve consistent, high quality analyses. An intuitive, friendly interface shortens and smooths the path to results and makes it easy to train new users.

# Simple, trouble-free design



## Here's how easy it is to get high quality, reproducible results



### Step 1:

#### Prepare your eluent


Add 10 mL concentrate to the eluent bottle, fill to 1 L with high purity water, and mix



### Step 2:

#### Add your regenerant


Add concentrate to the regenerant bottle, dilute with high purity water, and mix



### Step 3:

#### Connect your regenerant

Attach the tubing from the regenerant bottle to the "Regen In" on the suppressor, then attach the tubing from the "Regen Out" to waste




System equilibrated. Ready for regeneration.

### Step 4:

#### Equilibrate your system


Start the pump and flow eluent through system until baseline shows low noise and drift



### Step 5:

#### Plug in your sample information

Enter the sample sequence into Chromeleon CDS



### Step 6:

#### Inject your standards and samples

You can use either the manual syringe port or an autosampler



### Step 7:

#### Process your data


Generate first-class results using Chromeleon CDS

# Simply superior resolution and reproducibility

If your lab has a small number of samples for routine water analysis—or any other standard IC application, for that matter—you'll find the Dionex Easion IC system is a great fit. Industry-leading chemical suppression technology allows quick, hassle-free online regeneration, and ensures optimum performance for both anions and cations.


## Inorganic anions

The Dionex Easion IC system shows excellent resolution for common anions found in drinking water.




## Ten replicate injections of a seven-ion standard

Despite the system's compact size and operational simplicity, there's no compromise in the quality and consistency of results.




## Inorganic anions in a drinking water sample

The technology of the Dionex Easion IC system reduces baseline noise, giving you the sensitivity needed to identify and quantify all ions of interest.




The versatile Dionex Easion IC system delivers superior results for both cations and anions.



### Separation of six common cations in less than 15 min

The Thermo Scientific™ Dionex™ IonPac™ CS12A analytical column eases analysis and increases time efficiency.



### Common inorganic cations and ammonium

Are you especially interested in analyzing ammonium in your samples? The Thermo Scientific™ Dionex™ IonPac™ CS16 analytical column can be used on the Dionex Easion IC system for increasing resolution between sodium and ammonium; which is very helpful when sodium is present in high amounts.



# The simpler path to IC success

Here's one more reason to celebrate:  
Right out of the box, you're ready to run!

Processing a lot of samples? Add a Thermo Scientific™ Dionex™ AS-DV Autosampler for walk-away convenience, and you'll also enjoy the benefits of built-in sample filtration.



## Ordering information

Description	Part number
Dionex Easion IC system	22138-60101
Thermo Scientific™ Dionex™ Easion™ Anion Regenerant and Eluent Concentrates	057555 & 063965
Thermo Scientific™ Dionex™ Easion™ Cation Regenerant and Eluent Concentrates	057556 & 080388
Chromeleon CDS, version 7.2	7100.0108-ICSP
Dionex AS-DV Autosampler*	068907
Thermo Scientific™ Dionex™ Easion™ Vials with Filter Cap*	038141
Microsoft® Windows® 10 PC*	7200.0850-ICSP

\*Optional items

Find out more at [thermofisher.com/easion](https://thermofisher.com/easion)

**For Research Use Only. Not for use in diagnostic procedures.** © 2020-2024 Thermo Fisher Scientific Inc. All rights reserved. Microsoft and Windows are registered trademarks of Microsoft Corporation. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all locations. Please consult your local sales representative for details. **BR73688-EN 0124M**

**ThermoFisher**  
SCIENTIFIC