


Food safety

Fully automated method implementation for compliant analysis of polychlorinated dibenzodioxins/furans by GC-MS/MS

Polychlorinated dibenzo-*p*-dioxins/furans (PCDDs/PCDFs), commonly known as dioxins, are a group of persistent organic pollutants (POPs) that have been the focus of environmental scrutiny and regulation for decades. Due to their extreme toxicity, the European Union (EU) has set strict thresholds on allowable limits of dioxins within food and feed ($\mu\text{g}\cdot\text{g}^{-1}$ concentration) to protect public health. This requires constant monitoring of food and feed items imported into the EU from jurisdictions with different threshold limits or less stringent regulatory requirements.^{1,2}



Performance-based methodology for dioxin analysis is applied within the EU to ensure accurate reporting of dioxin levels within food and feed. Results are evaluated according to strict quality control and reporting criteria. For analytical testing laboratories, this often means a significant time investment for adoption and implementation of the methodology for dioxin analysis, as well as a requirement for highly skilled operators. Even for laboratories experienced with GC-MS-based test procedures, the addition of dioxin testing to the overall portfolio is demanding. To simplify this process, a dedicated eWorkflow™ template is available within the Thermo Scientific™ Chromeleon™ Chromatography Data System. This eWorkflow template is available for the Thermo Scientific™ TSQ™ 9610 triple quadrupole GC-MS/MS system and is free to download from the Thermo Scientific™ AppsLab™ Library of Analytical Applications ([Dioxin eWorkflow](#)).

Figure 1. Chemical structure of polychlorinated dibenzodioxins/furans (PCDD/F)




Figure 2. Overview of Dioxin eWorkflow template and its contents in Chromeleon CDS

Processing methods within the Dioxin eWorkflow template use isotopic dilution for quantitation, allowing users to perform data acquisition and real-time quantitation within the same software. Spiking levels of both syringe and internal standard can be defined per injection, permitting samples spiked with different levels of internal or syringe standard concentration to be analyzed in a single sequence. This can be seen in Figure 3 where internal standard concentration level within the sequence is defined at different levels for calibration standards (i.e., PCDD/F CS1), check standards (LOQ, LOQ/2, LOQ/4), and blank/sample injections (Sample).

A summary of the sample results is presented in an interactive table (Figure 4). The information provided by default includes the following parameters:

- Internal standard recovery – recovery assessment from analyte extraction/cleanup process.
- Targeted and observed ion ratio confirmation – validates that the detected compound meets the theoretical ratio between quantification and qualifier ions.

- Limit of quantification (LOQ) – choice of lowest calibration standard or user defined limit.
- Lower, middle, and upper bound concentration – concentrations reported when all results below LOQ are replaced with a zero value (lower bound), half the value of the LOQ (middle bound), or the value of the LOQ (upper bound).
- Deviation upper/lower bound – difference between upper and lower bound limits.

A dedicated sample LOQ is calculated for each component for every sample injection, accounting for the analyte extraction recovery (or internal standard recovery), injection volume, and sample weight. In addition, lower, middle, and upper bound values of the World Health Organization PCDD/F Toxic Equivalence (WHO-PCDD/F-TEQ)¹ are automatically calculated together with percent deviation between upper and lower bound limits, which is required for reporting purposes.

| Name | Type | *CM7:IntStd_Level |
|--|----------------------|-------------------|
| Dioxin analysis CS1 | Calibration Standard | PCDD/F CS1 |
| Dioxin analysis CS2 | Calibration Standard | PCDD/F CS1 |
| Dioxin analysis CS3 | Calibration Standard | PCDD/F CS1 |
| Dioxin analysis CS4 | Calibration Standard | PCDD/F CS1 |
| Dioxin analysis CS5 | Calibration Standard | PCDD/F CS1 |
| Dioxin analysis LOQ/4 | Check Standard | PCDD/F LOQ/4 |
| Dioxin analysis LOQ/2 | Check Standard | PCDD/F LOQ/2 |
| Dioxin analysis LOQ | Check Standard | PCDD/F LOQ |
| Dioxin analysis RBW - Reagent blank | Blank | Sample |
| Dioxin analysis QK1 - Mixed animal fat | Unknown | Sample |
| Dioxin analysis QK3 - Egg fat | Unknown | Sample |
| Dioxin analysis QK6 - Fish oil | Unknown | Sample |

Internal standard concentration is defined per injection, allowing analysis of samples/standards containing varying levels of internal standard

Figure 3. Interactive injection sequence table

| Peakname | Ret.Time | ISTD recovery | Target IR | Obs IR | IR deviation | TEF | Sample LOQ | <LOQ? | Lowerbound | Middlebound | Upperbound | SUM Upperbound | Deviation Upper/Lowerbound |
|---------------|----------|---------------|-----------|--------|--------------|--------|------------|-------|-----------------------|-----------------------|-----------------------|-----------------------|----------------------------|
| | min | % | % | % | % | | pg/g | | WHO-PCDD/F-TEQ (pg/g) | WHO-PCDD/F-TEQ (pg/g) | WHO-PCDD/F-TEQ (pg/g) | WHO-PCDD/F-TEQ (pg/g) | |
| 2378-TCDF | 19.32 | 85.26 | 93.1 | 94.7 | 1.7 | 0.1 | 0.044 | | 0.046 | 0.046 | 0.046 | | |
| 2378-TCDD | 19.90 | 74.29 | 96.0 | 90.4 | -5.9 | 1 | 0.050 | | 0.070 | 0.070 | 0.070 | | |
| 12378-PeCDF | 23.36 | 80.30 | 79.7 | 88.8 | 11.3 | 0.03 | 0.232 | <LOQ | 0.000 | 0.003 | 0.007 | | |
| 23478-PeCDF | 24.75 | 91.77 | 77.5 | 79.2 | 2.2 | 0.3 | 0.203 | | 0.072 | 0.072 | 0.072 | | |
| 12378-PeCDD | 25.00 | 90.59 | 78.5 | 76.6 | -2.4 | 1 | 0.206 | | 0.347 | 0.347 | 0.347 | | |
| 123478-HxCDF | 28.14 | 88.71 | 64.2 | 61.6 | -4.0 | 0.1 | 0.210 | <LOQ | 0.000 | 0.011 | 0.021 | | |
| 123678-HxCDF | 28.25 | 91.33 | 62.8 | 61.5 | -2.0 | 0.1 | 0.204 | <LOQ | 0.000 | 0.010 | 0.020 | | |
| 234678-HxCDF | 28.93 | 87.78 | 64.5 | 52.5 | -18.7 | 0.1 | 0.213 | <LOQ | 0.000 | 0.011 | 0.021 | | |
| 123478-HxCDD | 29.04 | 80.16 | 63.9 | 68.3 | 6.8 | 0.1 | 0.233 | | 0.059 | 0.059 | 0.059 | | |
| 123678-HxCDD | 29.12 | 86.13 | 64.9 | 61.9 | -4.6 | 0.1 | 0.217 | | 0.123 | 0.123 | 0.123 | | |
| 123789-HxCDD | 29.43 | 86.13 | 64.2 | 66.0 | 2.8 | 0.1 | 0.217 | | 0.066 | 0.066 | 0.066 | | |
| 123789-HxCDF | 29.79 | 99.20 | 63.5 | 70.4 | 10.8 | 0.1 | 0.188 | <LOQ | 0.000 | 0.009 | 0.019 | | |
| 1234678-HpCDF | 31.43 | 90.29 | 79.8 | 84.6 | 6.1 | 0.01 | 0.207 | | 0.003 | 0.003 | 0.003 | | |
| 1234678-HpCDD | 32.86 | 90.88 | 79.6 | 79.2 | -0.6 | 0.01 | 0.205 | | 0.091 | 0.091 | 0.091 | | |
| 1234789-HpCDF | 33.59 | 111.19 | 79.6 | 92.5 | 16.3 | 0.01 | 0.168 | <LOQ | 0.000 | 0.001 | 0.002 | | |
| OCDD | 37.45 | 106.88 | 96.0 | 93.9 | -2.2 | 0.0003 | 0.349 | | 0.002 | 0.002 | 0.002 | | |
| OCDF | 37.68 | 106.88 | 93.3 | 86.9 | -6.9 | 0.0003 | 0.349 | <LOQ | 0.000 | 0.000 | 0.000 | 0.968 | 9.3% |

Upper/lower bound deviation should not exceed 20%

Full overview of Internal standard recovery and ion ratio performance

Sample LOQ automatically corrected for injection volume, sample weight, and internal standard recovery

Figure 4. Sample summary table provided with the Dioxin eWorkflow template within Chromeleon CDS

The Dioxin eWorkflow template allows analytical testing laboratories to easily implement compliant analytical methodology for dioxin analysis. It delivers validated data acquisition/processing methods and reporting templates for dioxin analysis, helping remove time-consuming steps around method development and data processing.

Included within the Dioxin eWorkflow template:

- Injection sequence templates with pre-optimized analytical parameters for method injection and data acquisition.
- Incorporation of isotopic dilution quantitation, removing the need for additional software or script writing to perform data calculations.
- Interactive tables and charts to facilitate quick data evaluation for compliance purposes.
- Fully customizable data evaluation and report layout templates based on the analyst's needs.

References

1. Commission Regulation (EU) 2017/644 (2017). Laying down methods of sampling and analysis for the control of levels of dioxins, dioxin-like PCBs and non-dioxin-like PCBs in certain foodstuffs and repealing Regulation (EU) No 589/2014.
2. Commission Regulation (EU) 2017/771 (2017). Amending Regulation (EC) No 152/2009 as regards the methods for the determination of the levels of dioxins and polychlorinated biphenyls.

 Learn more at thermofisher.com/triple-quadrupole-gc-ms