Product Spotlight



Focus on your samples and instruments

Introduction

When discussing sample analysis, it is normally expected that this also includes the injection of blanks, calibration standards, check standards and so on, all combined in a series of injections, a sequence.

However, this is not always the case. For example, in the chemical or petrochemical industries, the analysis is typically driven by the need to get a fast result for an individual sample without additional injections. The main chromatography instrument used is a gas chromatograph (GC) and it is typical for one operator to use several of these instruments with each GC dedicated to a specific analysis. A calibration curve for a specific method is only required periodically, for example once every six months so running a daily blank and a check standard is sufficient to confirm the validity of the calibration.

All of this drives a more sample and instrument-centric approach to the analysis, rather than a sequence-centric view. To support these specific workflows required by (petro)chemical laboratories, Thermo Scientific[™] Chromeleon[™] Chromatography Data System (CDS) provides a comprehensive set of tools that streamline and simplify operation.



SCIENTIFIC

Instruments Overview

As mentioned before, one analyst in the laboratory can be responsible for multiple instruments. It is key that they can instantly see the status of all these instruments at a glance. Not all instruments have fully automated injection systems and may need manual interaction, for example loading a gas cylinder which was used to take the sample. Monitoring multiple instruments for their status and readiness to perform this manual action should be possible in one view.

With the Chromeleon Instruments Overview (Figure 1) the status of all the analyst's favorite instruments are visible at once, clearly indicating which instruments are running (including run time and the operator), available for a new sample analysis or ready for a new injection. Analysts can also control the instrument queue directly from this overview.

With additional privileges, system administrators can select a range of visible instruments or have an overview of all instruments in the network (Figure 2), this can quickly identify issues, for example when an Instrument Controller is offline.

C Back C Creat	te File Edit Vie							
Instruments	**	Launch eWorkflow	• 🗐 View Instrument 🧠	View Queue 🔳 Sto	p Queue 👂 Start Queue 👻 🍸 Filtering 🚪	Grouping 🗳 Detach View		
🔕 🖌	Filter V	Instrument Name	Instrument Status	Controller	Sequence	Injection	Queued By	Queue Status
		07_GC	Run: 0.17 of 0.24 [min]	NLBRE-6R2RV02	Aromatics%202021-09-20%2015-03	15 of 31 (Sample 10)	Charlie	Pending sequences: 1
Overview		🗐 08_GC	ldle	NLBRE-6R2RV02				
NLBRE-6R2RV02		🚍 09_GC	ldle	NLBRE-6R2RV02	DHA%202021-09-20%2016-40		Charlie	Cancelled sequences: 1
607_GC	*	■ 10_GC	ldle	NLBRE-6R2RV02				
08_GC	****	🗐 11_GC	ldle	NLBRE-6R2RV02				
09_GC	×	12_GC	Idle	NLBRE-6R2RV02				
		9. 						
Data								

Figure 1. Favorite Instruments Overview

nstruments	« DLaunch eWorkflow	• 🗐 View Instrument 🧠	View Queue 🔳 Sto	op Queue 🗼 Start Queue 🕞 🍞 Filtering 🚍 Groupin	ng 🖳 Detach View		
Fiter S	/ Instrument Name	Instrument Status	Controller	Sequence	Injection	Queued By	Queue Status
	07_GC	Run: 0.17 of 0.24 [min]	NLBRE-6R2RV02	Aromatics %202021-09-20%2015-03	15 of 31 (Sample 10)	Charlie	Pending sequences: 1
Overview	🗐 08_GC	Idle	NLBRE-6R2RV02				
NLBRE-6R2RV02	09_GC	Idle	NLBRE-6R2RV02	DHA%202021-09-20%2016-40		Charlie	Cancelled sequences:
07 GC	🗐 10_GC	Idle	NLBRE-6R2RV02				
08_GC	🗐 11_GC	Idle	NLBRE-6R2RV02				
09_GC	12_GC	Idle	NLBRE-6R2RV02				
10_GC	💓 GC-1	Idle	NLBRE-44DT0N2				
12 GC	🗐 GC-2	Idle	NLBRE-44DT0N2				
NLBRE-44DTON2	ICS-6000-LT	Idle	NLBRE-44DT0N2				
GC-1	Trace-GC-LT	Idle	NLBRE-44DT0N2				
GC-2	UltiMate3000_RS-LT	ldle	NLBRE-44DT0N2	2021-09-20%2014-53-09%20Alkylphenone%20Linearity			
ICS-6000-LT	Vanquish-LT	Idle	NLBRE-44DTON2				
Trace-GC-LT UttiMate3000_RS-LT	02_ICS-3000_Dual_1	Remote Connecting	NLBRE-BC093G2				
Vanguish-LT	03_ICS-3000_Dual_2	Remote Connecting	NLBRE-BC093G2				
NLBRE-BC093G2	04_ICS-2100	Remote Connecting	NLBRE-BC093G2				
02_ICS-3000_Dual_1	Vanquish_Flex	Remote Connecting	NLBRE-BC093G2				
03_ICS-3000_Dual_2 04_ICS-2100 Vanquish_Rex							
Data y eWorkflows							

Figure 2. Full Instruments Overview

Instant sample run

From the Instrument Overview samples can be started instantly. Once a suitable and available instrument is selected for the required analysis, the run can be started in only two clicks.

- Select the analysis method to run (Figure 3)
- 2. Click Finish to start the analysis (Figure 4)

Chromeleon software's built-in eWorkflow[™] procedures provide a unique framework to support this ultra-streamlined operation. They ensure that whichever instrument is selected, the correct methods and files are used for the analysis every time. This eliminates errors, simplifies operation and minimizes training requirements.

Back Creat	File Edit	view Tool	s Help						
nstruments		OLaun	ch eWorkflow •	🖓 🗐 View Instrument 🧠 View	Queue 🔳 Sto	op Queue 🗼 Start Queue 🕣 🍸 Filtering 📒	Grouping 🖳 Detach View		
9 🙀	Filter S		Aromatics on	NLBRE-CM01		Sequence	Injection	Queued By	Queue Status
		(OHA on NLBR	E-CM01		Aromatics %202021-09-20%2015-03	15 of 31 (Sample 10)	Charlie	Pending sequences:
Overview				Gasoline on NLBRE-CM0					
NLBRE-6R2RV02			RGA on NLBRE-CM01			DHA%202021-09-20%2016-40		Charlie	Cancelled sequences
	2		SimDist on NI		/02				
09_GC	2	12.0			BRE-6R2RV02				
10_GC	÷	120	<i></i>	NUI NLI					
11_GC	****	· · · · · ·			Unit-Unit 1992			1	
I 12.60	*]	

Figure 3. Select Method



Figure 4. Start Run

Instrument control

The online controls for a specific instrument can be directly accessed from the Overview pane (Figure 5A). With a single click, an ePanel opens (Figure 5B) displaying the real time instrument status along with controls for direct operation.

Direct data access

In order to quickly and easily find data that was just generated, results can be directly accessed without leaving the Overview pane. One click opens the queue tab of the instrument ePanel (Figure 6A) where the recent sequence can be found. The sequence is displayed, including MiniPlots[™] of the actual chromatograms (Figure 6B) and even custom columns, which can include results, for instant information. A double click on the sequence will open the data in the Chromatography Studio for detailed viewing or further processing if required.

Conclusion

By providing intelligent tools to enable operators to work with multiple instruments, start analyses quickly and review results immediately, Chromeleon CDS supports the workflows in (petro)chemical laboratories, while delivering ease of use and minimizing training requirements.



Figure 5. View Instrument



Figure 6. View Queue for Direct Data Access

thermo scientific



Join the Chromeleon CDS community to future proof your investment: **thermofisher.com/chromeleonsupport**

Like Charlie Chromeleon on Facebook to follow his travels and get important updates on chromatography software! facebook.com/CharlieLovesChromatography

Visit AppsLab Library for online access to applications for GC, IC, LC, MS and more. **thermofisher.com/AppsLab**

Find out more at thermofisher.com/chromeleon

© 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. **SP73778-EN 0122M**

