Thermo Fisher

Next-generation sequencing

Complete next-generation sequencing solution for pharmacogenomics research Ion AmpliSeq Pharmacogenomics Research Panel

The Ion AmpliSeq[™] Pharmacogenomics (PGx) Research Panel is ideally suited to clinical research studies investigating how genomic variations in drug metabolizing enzymes (DMEs) impact drug efficacy and safety. The panel is part of a complete end-to-end next-generation sequencing (NGS) research solution for genotyping single-nucleotide polymorphisms (SNPs), as well as insertion/deletion (indel) and copy number variation (CNV) analysis in 40 known DMEs (Figure 1). The panel focuses on 136 well-documented SNPs and indel variants, and captures *CYP2D6* copy number variation at both the gene and exon 9 levels, enabling the screening of broad selection of haplotypes including *CYP2D6*36* alleles.

The Ion AmpliSeq Pharmacogenomics Research Panel utilizes Ion AmpliSeq[™] chemistry to offer:

- Customizable content, allowing addition or removal of gene targets using Ion AmpliSeq Designer
- Outstanding performance on buccal swab samples from only 10 ng of input
- High coverage and unbiased amplification, resulting in robust performance and minimizing sample reruns
- Genotyping and copy number analysis in a single workflow
- Sample multiplexing flexibility enabling processing of 8–96 samples in a single batch
- Integrated downstream analysis solution

Ion AmpliSeq Pharmacogenomics Research Panel			
Total amplicons	119		
Genes covered	40		
SNV/indels detected	136		
CYP2D6 CNV events*	2		

* CNV detection at gene level and exon 9(*36) level.

ABCB1	CYP2C8	F2	MTHFR	VKORC1
ABCG2	CYP2C9	F5	OPRD1	HLA-A*3101
ADRA2A	CYP2D6	GABRA6	OPRK1	HLA-B*5701
ANKK1	CYP3A4	GABRP	OPRM1	HLA-B*1502
APOE	CYP3A5	GRIK4	SLCO1B1	
COMT	DBH	HTR2A	TPMT	
CYP1A2	DPYD	HTR2C	UGT1A1	
CYP2B6	DRD1	ITGB3	UGT2B15	
CYP2C19	DRD4	KIF6	UGT2B7	

Figure 1. The Ion AmpliSeq Pharmacogenomics Research Panel analyzes genetic variants including SNPs, indels, and CNVs associated with 40 DMEs. Panel design is customizable, enabling flexible editing to an existing gene target list.

ion torrent

Complete sample-to-results workflow solution

The Ion AmpliSeq Pharmacogenomics Research Panel is part of a complete NGS workflow solution, enabling rapid and flexible analysis of buccal swab samples in 2 days, from start to finish (Figure 2). Each step in the workflow can be automated with minimal hands-on time. A flexible system together with a broad spectrum of throughput needs enables parallel processing of 48–384 samples in a single workflow—reducing data turnaround time significantly when project size is variable. An integrated data analysis solution completes the process, providing report-ready haplotype answers.

Our support team is ready to help new and existing customers set up and maintain all components of the system, enabling you to fully focus on your projects.



Exceptional sample success

The low per-sample input requirement of only 10 ng, robust compatibility of Ion AmpliSeq chemistry with buccal swab– derived genomic DNA, and broad gene coverage enable lower assay failure rate and sample reruns (Figure 3). This can mean the difference between missing deadlines and delivering research results on time. There are 9 additional panels included in the panel that can generate unique sample-specific genetic barcodes useful for sample identification purposes. Table 1 shows the recommended numbers of samples for Ion 510, 520, and 530 Chips to achieve >500x average coverage per amplicon using estimated throughputs of 3 million, 6 million, and 20 million reads, respectively.

Table 1. Sample multiplexing recommendations for each type of lon chip.

Chip	Multiplexed samples
lon 510	48 samples
lon 520	96 samples
lon 530	384 samples

Α

Genotyping across 136 genetic variants				
Cell line samples	400			
Buccal swab samples	432			
Genotype accuracy	>99.9%			
Genotype call rate	>99.9%			

В

CYP2D6 gene and exon 9 CNVs	
Buccal swab samples	432
CYP2D6 gene-level calls	429
CYP2D6 gene-level no calls	3
CYP2D6 gene-level false calls	0
CYP2D6 gene-level call rate	99.3%
CYP2D6 exon 9-level calls	415
CYP2D6 exon 9-level no calls	17
CYP2D6 exon 9-level false calls	0
CYP2D6 exon 9-level call rate	96.1%

Figure 3. Genotyping and CNV data. (A) Genotyping results from over 800 cell line and buccal swab samples demonstrate >99.9% accuracy with known genotypes. **(B)** Focusing on the 432 buccal swab samples, *CYP2D6* CNVs were called in 99.3% and 96.1% of samples at the gene and exon 9 levels, respectively, with no false positives.

Primer pools for the Ion AmpliSeq Pharmacogenomics Research Panel can be ordered through the Ion AmpliSeq Designer at <u>ampliseq.com</u>.

Ordering information

Product	Quantity	Cat. No.		
Sample DNA extraction using the KingFisher Apex system				
MagMAX DNA Multi-Sample Ultra Kit	500 preps	A25597		
	2,500 preps	A25598		
Sequencing library construction				
Ion AmpliSeq Library Kit 2.0	8 reactions	4475345		
	96 reactions	4480441		
	384 reactions	4480442		
Dual Barcode Kit 1-96	1 kit	A39360		
Ion Library Equalizer Kit	96 reactions	4482298		
Templating and sequencing using the Ion Chef and Ion GeneStudio S5 systems				
lon 510 & lon 520 & lon 530 Kit – Chef	8 reactions	A34461		
Ion 510 Chip Kit	8 chips	A34292		
Ion 520 Chip Kit	8 chips	A27762		
Ion 530 Chip Kit	8 chips	A27764		

Learn more at thermofisher.com/pharmacogenomics

ion torrent