

Ask your pathologist about the Oncomine Dx Target Test

A paradigm change in testing for targeted therapies in NSCLC, cholangiocarcinoma, and thyroid cancer

The Ion Torrent™ Oncomine™ Dx Target Test is the first targeted next-generation sequencing (NGS) *in vitro* diagnostic test for non-small cell lung cancer (NSCLC), cholangiocarcinoma (CC), and thyroid cancer (TC), simultaneously delivering multiple biomarker results for multiple targeted therapies from one sample within 4 days.

Did you know:


- Many biopsy samples are so small that they cannot be analyzed by some NGS tests, especially panels containing hundreds of genes, leading to tissue exhaustion
- It can take several weeks to get results with alternative NGS tests, potentially delaying a treatment decision

Choosing the right NGS test can make a difference for your patient

The Oncomine Dx Target Test is an FDA-approved, NGS companion diagnostic (CDx) test that can:

- Identify patients who are candidates for multiple therapies— one test indicated as a CDx device to aid in selecting NSCLC, CC, TC, and medullary thyroid cancer (MTC) patients for treatment with targeted therapies
- Accept small samples (10 ng DNA and 10 ng RNA) so that more patients can potentially access targeted therapies
- Generate results in a laboratory within 4 days, enabling faster treatment decisions

This test is reimbursed by Medicare and over 40 commercial payers, covering more than 200 million US enrollees.



Cancer type	Gene	Targeted therapies
NSCLC	<i>BRAF</i>	TAFINLAR® (dabrafenib) in combination with MEKINIST® (trametinib)
	<i>EGFR</i> L858R and exon 19 deletions	IRESSA® (gefitinib)
	<i>EGFR</i> exon 20 insertions	EXKIVITY™ (mobocertinib) RYBREVANT™ (amivantamab-vmjw)
	<i>ERBB2/HER2</i> activating mutations (SNVs and exon 20 insertions)	ENHERTU® (fam-trastuzumab deruxtecan-nxki)
	<i>RET</i>	GAVRETO™ (pralsetinib) RETEVMO® (selpercatinib)
	<i>ROS1</i>	XALKORI® (crizotinib)
CC	<i>IDH1</i>	TIBSOVO® (ivosidenib)
MTC	<i>RET</i> mutations (SNVs, MNVs, and deletions)	RETEVMO® (selpercatinib)
TC	<i>RET</i> fusions	RETEVMO® (selpercatinib)

Figure 1. List of genes for therapeutic use.

Oncomine Dx Target Test—performance

Concordance with FDA-approved or validated reference methods based on FISH, PCR, Sanger sequencing, or NGS was established for all CDx biomarkers included in the test:

- 100% overall percent agreement (OPA), positive percent agreement (PPA) and negative percent agreement (NPA) for *BRAF* mutation, *EGFR* exon 20 insertions, and *ROS1* fusions
- 99% OPA, PPA, and NPA for *EGFR* exon 19 deletions and L858R
- 99% OPA, 100% PPA, and 99% NPA for *ERBB2/HER2* activating mutations (SNVs and exon 20 insertions)
- 92% OPA, 91% PPA, and 92% NPA in the first study and 95% OPA, 92% PPA, and 97% NPA in the second study for *RET* fusions in NSCLC
- 98% OPA, 99% PPA, and 97% NPA for *IDH1* mutations
- 99% OPA, 100% PPA, and 98% NPA for *RET* mutations in MTC and 100% OPA, PPA, and NPA for *RET* fusions in TC

Oncomine Dx Target Test—report

Sequence variations to indicate therapeutic use for NSCLC (for illustrative purposes only; *EGFR*, *BRAF*, *ERBB2/HER2*, *ROS1*, and *RET* are mutually exclusive)

DNA sequence variants						
Gene	Display name	Amino acid change	Nucleotide change	Test result	Hotspot ID	Associated therapy
<i>EGFR</i>	<i>EGFR</i> L858R	p.Leu858Arg	c.2573T>G	Positive	COSM6224	IRESSA® (gefitinib)
<i>EGFR</i>	<i>EGFR</i> exon 20 insertions	p.Ala767_Ser768 insSerValAsp	c.2311_2312ins GCGTGGACA	Positive	COSM13428	EXKIVITY™ (mobocertinib) RYBREVANT™ (amivantamab-vmjw)
<i>BRAF</i>	<i>BRAF</i> V600E	p.Val600Glu	c.1799T>A	Positive	COSM476	TAFINLAR® + MEKINIST® (dabrafenib in combination with trametinib)
<i>ERBB2</i>	<i>ERBB2</i> exon 20 insertions	p.Gly776delinsLeuCys	c.2326_2326delGinsCTTT	Positive	COSM12554	ENHERTU® (fam-trastuzumab deruxtecan-nxki)
Gene fusions (RNA)						
Gene	Display name			Test result	Associated therapy	
<i>ROS1</i>	<i>ROS1</i> fusions			Positive	XALKORI® (crizotinib)	
<i>RET</i>	<i>RET</i> fusions			Positive	GAVRETO™ (pralsetinib) RETEVMO® (selpercatinib)	

Sequence variations to indicate therapeutic use for CC (for illustrative purposes only)

DNA sequence variants						
Gene	Display name	Amino acid change	Nucleotide change	Test result	Hotspot ID	Associated therapy
<i>IDH1</i>	<i>IDH1</i> R132G	p.Arg132Gly	c.394C>G	Positive	COSM28749	TIBSOVO® (ivosidenib)
<i>RET</i>	<i>RET</i> A883T	p.Ala883Thr	c.2647G>A	Positive	COSM100081	RETEVMO® (selpercatinib)

Figure 3. Example of the Oncomine Dx Target Test report format. The report includes a section with results of the validated biomarkers and information about relevant treatment indication, as well as a separate section with the other biomarkers not validated for treatment selection (not shown).

If your pathology laboratory does not perform the Oncomine Dx Target Test, you can send samples to one of these reference laboratories.

Reference lab	Telephone number	Website
NeoGenomics Laboratories, Inc.	866-776-5907	neogenomics.com
OncoCyte	615-639-0710	oncocyte.com

Find out more at thermofisher.com/oncomine-dxtarget

For In Vitro Diagnostic Use. © 2017–2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. XALKORI is a registered trademark of Pfizer Inc. TAFINLAR and MEKINIST are registered trademarks of Novartis. IRESSA is a registered trademark of AstraZeneca. GAVRETO is a trademark of Blueprint Medicines Corporation. TIBSOVO is a registered trademark of Agios Pharmaceuticals, Inc. EXKIVITY is a trademark of Takeda Pharmaceuticals. RYBREVANT is a trademark of Johnson & Johnson. ENHERTU is a registered trademark of Daiichi Sankyo Company, Limited. RETEVMO is trademark of Eli Lilly and Company. **COL119232 0123**

Oncomine Dx Target Test—content

Gene targets for NSCLC				
Gene targets for therapeutic use				
<i>BRAF</i> : V600E	<i>EGFR</i> : L858R, exon 19 deletions, and exon 20 insertions	<i>ERBB2/HER2</i> activating mutations (SNVs and exon 20 insertions)	<i>RET</i> : fusions	<i>ROS1</i> : fusions
Analytically validated gene targets				
<i>KRAS</i>	<i>MET</i> *	<i>PIK3CA</i>		
Additional gene targets**				
<i>AKT1</i> <i>ALK</i> *	<i>ERBB2</i> <i>ERBB3</i>	<i>HRAS</i> <i>KIT</i>	<i>MTOR</i> <i>NRAS</i>	<i>RET</i>
<i>CDK4</i> <i>DDR2</i>	<i>FGFR2</i> <i>FGFR3</i>	<i>MAP2K1</i> <i>MAP2K2</i>	<i>PDGFRA</i> <i>RAF1</i>	<i>ROS1</i>
Gene targets for CC		Gene targets for TC	Gene targets for MTC	
Gene targets for therapeutic use		<i>IDH1</i> : mutations	<i>RET</i> : fusions	
			<i>RET</i> : fusions	

Figure 2. Oncomine Dx Target Test complete gene list.

* The test reports fusion/translocation variants for *ROS1* and *RET* only. This test only reports mutations for *ALK* and *MET*.

** Performance for the additional gene target variants has been validated based on a representative method.