



Comprehensive liquid biopsy analysis with only one tube of blood

Relevant insights from DNA, RNA, and all classes of variants

The Ion Torrent™ Oncomine™ Pan-Cancer Cell-Free Assay is a targeted next-generation sequencing (NGS) assay that delivers current, relevant insights for oncogenes associated with more than 18 different cancer types. Through simultaneous multibiomarker analysis of DNA and RNA from just one tube of blood, results are now only three days away.

- **Multiple biomarkers from one blood sample**—this complementary, noninvasive option to traditional tissue biopsies only requires one tube of blood for detection of variants across 52 oncogenes that are associated with 18 different cancer types
- **One assay, two programs**—two distinct programs are available for initial evaluation of the Oncomine Pan-Cancer Cell-Free Assay, or for potential adoption into your menu; these programs are available through Life Technologies Clinical Services Lab (lifelabdx.com or toll-free 1.888.734.8588)
- **One streamlined workflow**—go from blood sample to variant report in three days
- **Established performance**—tissue orthogonal concordance was conducted on plasma samples with corresponding solid tumor molecular characterization; reported concordance for SNVs and CNVs at >99%; key confirmed alterations observed in clinical samples include *EGFR* p.E746_A750del, *EGFR* p.T790M, *EGFR* p.G719S, *KRAS* G12/13, *NRAS* codon 61, *ERBB2* amp, *EGFR* amp




Table 1. Gene targets included in Oncomine Pan-Cancer Cell-Free Assay.

Hotspot genes			Tumor suppressor genes	Copy number genes	Gene fusions
AKT1 ALK AR ARAF BRAF CHEK2 CTNNB1 DDR2 EGFR ERBB2 ERBB3 ESR1 FGFR1	FGFR2 FGFR3 FGFR4 FLT3 GNA11 GNAQ GNAS HRAS IDH1 IDH2 KIT KRAS MAP2K1 MAP2K2	MET MTOR NRAS NTRK1 NTRK3 PDGFRA PIK3CA RAF1 RET ROS1 SF3B1 SMAD4 SMO	APC FBXW7 PTEN TP53	CCND1 CCND2 CCND3 CDK4 CDK6 EGFR ERBB2 FGFR1 FGFR2 FGFR3 MET MYC	ALK BRAF ERG ETV1 FGFR1 FGFR2 FGFR3 MET NTRK1 NTRK3 RET ROS1

Table 2. 18 cancer types included in Oncomine Pan-Cancer Cell-Free Assay.

Cancer types	Top 10 genes ranked by mutation frequency
Bladder	TP53, PIK3CA, FGFR3, HRAS, ERBB2, KRAS, CTNNB1, BRAF, NRAS, FBXW7
Brain and CNS	TP53, IDH1, PIK3CA, EGFR, CHEK2, ALK, CTNNB1, BRAF, PTEN, PDGFRA
Breast	PIK3CA, TP53, ERBB2, PTEN, SF3B1, AKT1, ERBB3, ESR1, KRAS, FGFR2
Cervical	PIK3CA, FBXW7, TP53, KRAS, ERBB2, PTEN, ERBB3, MTOR, CTNNB1, SMAD4
Colorectal	PIK3CA, FBXW7, TP53, KRAS, ERBB2, PTEN, ERBB3, MTOR, CTNNB1, SMAD4
Endometrial	PIK3CA, TP53, CTNNB1, PTEN, KRAS, FGFR2, FBXW7, MTOR, NRAS, ERBB2
Esophageal	TP53, PIK3CA, SMAD4, FBXW7, KRAS, ERBB2, APC, CTNNB1, PTEN, SMO
Gastric	TP53, PIK3CA, KRAS, FBXW7, ERBB3, ERBB2, SMAD4, CTNNB1, APC, MAP2K1
Head and neck	TP53, PIK3CA, HRAS, PTEN, FBXW7, RET, KRAS, FGFR3, BRAF, ERBB2
Kidney	TP53, MTOR, CHEK2, PIK3CA, PTEN, MET, FGFR3, EGFR, KRAS, SF3B1
Liver	CTNNB1, TP53, PIK3CA, KRAS, PTEN, KIT, IDH1, GNAS, APC, FGFR2
Lung	TP53, KRAS, EGFR, PIK3CA, BRAF, NRAS, PTEN, FBXW7, APC, CTNNB1
Melanoma	BRAF, NRAS, TP53, MAP2K1, CTNNB1, GNA11, PTEN, IDH1, KIT, GNAQ
Ovarian	TP53, PIK3CA, KRAS, EGFR, CTNNB1, CHEK2, ERBB2, MET, FBXW7, NRAS
Pancreatic	KRAS, TP53, SMAD4, APC, GNAS, CTNNB1, SF3B1, PIK3CA, BRAF, FGFR1
Prostate	TP53, PIK3CA, CTNNB1, AR, CHEK2, APC, PTEN, IDH1, AKT1, BRAF
Sarcoma	TP53, IDH1, NRAS, PIK3CA, KRAS, FGFR4, ERBB2, IDH2, HRAS, CTNNB1
Thyroid	BRAF, NRAS, HRAS, RET, TP53, KRAS, PIK3CA, AKT1, GNAS, CCND1

Next-generation sequencing (NGS) workflow—from sample to report



A comprehensive liquid biopsy NGS workflow for streamlined detection and analysis of variants from 52 oncogenes (Table 1) that are associated with 18 different cancer types (Table 2). Receive a report in three days with relevant insights that may provide guidance regarding future treatment decisions.

(NGS) workflow—optimized for operational efficiency from sample to report, the NGS workflow for the

OncoPrint Pan-Cancer Cell Free Assay consists of three key steps, enabling you to go from blood sample to report in three days. During sample preparation, cell-free nucleic acids are extracted, enriched, and amplified. These amplicon-based libraries are then assembled overnight before targeted resequencing. Our integrated informatics solution then takes you from variant caller to a finished report that provides contextual insights about sample-specific variants and their use with respect to labels, guidelines, and current global clinical trials.

Strong concordance with tissue samples enables confidence in your liquid biopsy results

Tissue orthogonal concordance was evaluated on plasma samples with corresponding solid tumor molecular characterization. Reported concordance for SNVs and CNVs were shown to be >99%. Key confirmed alterations observed in samples include *EGFR* p.E746_A750del, *EGFR* p.T790M, *EGFR* p.G719S, *KRAS* G12/13, *NRAS* codon 61, *ERBB2* amp, and *EGFR* amp.



**EXAMPLE
HEALTH
SYSTEM**



1234567890

Example Health System
123 Street
City, ST 12345 USA
Tel (123) 123-1234
email@examplehealth.com
www.examplehealth.com

Tracking Number: 00-123456789
Case Number: 00-123456789
Date: 08 Oct 2018
1 of 81

Sample Information

Year of Birth: 1966
Gender: Female
Smoking Status: Never Smoker
Case Number: 00-123456789

Primary Tumor Site:
Sample Type:
Sample ID:
Sample Collected:

Colon
FFPE
00-123456789
02/01/2018

Sample Cancer Type: Colorectal Cancer

Clinically Significant Biomarkers

■ Indicated
■ Contraindicated

Genomic Alteration	Relevant Therapies (On this cancer type)	Relevant Therapies (On other cancer type)	Clinical Trials
BRAF V600E BRAF proto-oncogene, serine/threonine kinase Test ID: Assay Frequency: 27.96%	cetuximab + vemurafenib + ■ chemotherapy panitumumab + vemurafenib + chemotherapy	binimetinib + encorafenib ¹ cobimetinib + vemurafenib ² dabrafenib + trametinib ³ trametinib ⁴ dabrafenib ⁵ vemurafenib ⁶ nivolumab BRAF inhibitor + MEK inhibitor	27
NF1 R243Q* neurofibromin 1 Test ID: Assay Frequency: 38.54%	None	■ erastinib	6
PKC3CA R89Q phosphatidylinositol-3-OH kinase class III isoform p115alpha subunit alpha Test ID: Assay Frequency: 22.15%	None	None	10
PTEN R173C phosphatidylinositol-3-OH kinase, class II Test ID: Assay Frequency: 30.33%	None	None	9
NOTCH1 P2415del notch 1 Test ID: Assay Frequency: 0.94%	None	None	2
Tumor Mutational Burden All 99 nucleotide sequences	None	None	2

Sources included in relevant therapies: PD4, NOTCH, BRAF, PTEN

© 2018 Thermo Fisher Scientific Inc. All rights reserved.

Document: The data presented here is a result of the analysis of published drug sources, but may not be exhaustive. The data version is 2018-04-05

Sample report from OncoPrint Reporter links variant data to current labels, guidelines, and clinical trials. Information is provided in an easy-to-consume format.

Performance data of the Oncomine Pan-Cancer Cell-Free Assay

Validation item	SNV/indel	CNV	Fusion
Analytical sensitivity	<ul style="list-style-type: none">>99.9% at 0.5% allele fraction80% at 0.1% allele fraction	>99% at >1.35 fold amplification	>99% at 0.4% fusion fraction
Analytical specificity	>99%	>99%	>99%
Analytical accuracy	>99%	>99%	>99%
Precision within run	98%	>99%	>99%
Precision across runs	99%	>99%	>99%
Tissue concordance	>99% for informative variants	>99% for informative variants	N/A

Two programs to get you started

If your laboratory is interested in a trial evaluation of your samples using the Oncomine Pan-Cancer Cell-Free Assay, sequencing and analysis is available through the Proof-of-Principle program.

If your laboratory is interested in validating and incorporating the Oncomine Pan-Cancer Cell-Free Assay into your menu, analytical validation support is also available.

Both programs are supported and executed by our CLIA-certified, CAP-accredited clinical services lab facility, Life Technologies Clinical Services Lab (LTCSL). LTCSL can help pharma and other customers find comprehensive, accurate answers. They specialize in developing oncology molecular assays for pharmaceutical clients with small, challenging samples.

Life Technologies Clinical Services Lab offers a comprehensive end-to-end oncology workflow, including histology, pathology, and NGS. They also provide custom assay development and analytical validation services to pharma customers using multiple technology platforms (qPCR, Sanger sequencing, and NGS), custom design, and existing analytically validated panels. For more information or to get started with their programs, please visit lifelabdx.com or call toll free 1.888.734.8588.

Proof-of-Principle (POP) program

The POP program is designed for customers interested in evaluating the Oncomine Pan-Cancer Cell-Free Assay with their own research samples. A small batch will be run on our NGS workflow, and a report with optional raw data will be provided back to you.

Sample requirements for the POP program

- The specimen must be NSCLC
- Plasma obtained from two 10 mL K₂EDTA tubes (lavender top) at a minimum
- Whole blood should be processed within 6 hours post-collection
- Whole blood samples should be spun for 10 minutes at 1,000–2,000 x g in a refrigerated centrifuge

Plasma should be separated from whole blood after centrifugation is completed, taking care not to disturb the buffy coat.

Analytical validation support

This program is designed for customers interested in validating and incorporating the Oncomine Pan-Cancer Cell-Free Assay into their menu. During the validation period, services will be provided to the customer to help minimize or potentially eliminate downtime as a result of the analytical validation efforts.

For more information, visit thermofisher.com/cfna-assays

ThermoFisher
SCIENTIFIC